

# Package ‘AlpsNMR’

March 29, 2021

**Type** Package

**Title** Automated spectral Processing System for NMR

**Version** 3.0.5

**Encoding** UTF-8

**Description** Reads Bruker NMR data directories both zipped and unzipped.

It provides automated and efficient signal processing for untargeted NMR metabolomics.

It is able to interpolate the samples, detect outliers, exclude regions, normalize, detect peaks, align the spectra, integrate peaks, manage metadata and visualize the spectra.

After spectra processing, it can apply multivariate analysis on extracted data.

Efficient plotting with 1-D data is also available.

Basic reading of 1D ACD/Labs exported JDX samples is also available.

**License** file LICENSE

**LazyData** FALSE

**Depends** R (>= 4.0), dplyr (>= 0.7.5), future (>= 1.10.0), magrittr (>= 1.5)

**Imports** utils, graphics, stats, grDevices, signal (>= 0.7-6), assertthat (>= 0.2.0), rlang (>= 0.3.0.1), stringr (>= 1.3.1), tibble (>= 1.3.4), tidyr (>= 1.0.0), readxl (>= 1.1.0), plyr (>= 1.8.4), purrr (>= 0.2.5), glue (>= 1.2.0), reshape2 (>= 1.4.3), GGally (>= 1.4.0), mixOmics (>= 6.3.2), matrixStats (>= 0.54.0), writexl (>= 1.0), fs (>= 1.2.6), rmarkdown (>= 1.10), speaq (>= 2.4.0), htmltools (>= 0.3.6), ggrepel (>= 0.8.0), pcaPP (>= 1.9-73), furrr (>= 0.1.0), ggplot2 (>= 3.1.0), baseline (>= 1.2-1), zip (>= 2.0.4), tidyselect (>= 0.2.5), BiocParallel, SummarizedExperiment, S4Vectors

**Suggests** DT (>= 0.5), testthat (>= 2.0.0), plotly (>= 4.7.1), ChemoSpec, knitr

**biocViews** Software, Preprocessing, Visualization, Classification, Cheminformatics, Metabolomics, DataImport

**RoxygenNote** 7.1.1

**Roxygen** list(markdown = TRUE)

**VignetteBuilder** knitr

**git\_url** <https://git.bioconductor.org/packages/AlpsNMR>

**git\_branch** RELEASE\_3\_12

**git\_last\_commit** c83bb59

**git\_last\_commit\_date** 2021-01-20

**Date/Publication** 2021-03-29

**Author** Ivan Montoliu Roura [aut],

Sergio Oller Moreno [aut] (<<https://orcid.org/0000-0002-8994-1549>>),

Francisco Madrid Gambin [aut] (<<https://orcid.org/0000-0001-9333-0014>>),

Luis Fernandez [aut, cre] (<<https://orcid.org/0000-0001-9790-6287>>),

Héctor Gracia Cabrera [aut],

Nestlé Institute of Health Sciences [cph],

Institute for Bioengineering of Catalonia [cph]

**Maintainer** Luis Fernandez <[lfernandez@ibecbarcelona.eu](mailto:lfernandez@ibecbarcelona.eu)>

## R topics documented:

|   |    |
|---|----|
| AlpsNMR-package . . . . .                 | 4  |
| AUC_model . . . . .                       | 5  |
| bp_kfold_VIP_analysis . . . . .           | 5  |
| bp_VIP_analysis . . . . .                 | 7  |
| computes_peak_width_ppm . . . . .         | 9  |
| confusion_matrix . . . . .                | 10 |
| files_to_rDolphin . . . . .               | 10 |
| file_lister . . . . .                     | 12 |
| filter.nmr_dataset_family . . . . .       | 13 |
| format.nmr_dataset . . . . .              | 13 |
| format.nmr_dataset_1D . . . . .           | 14 |
| format.nmr_dataset_peak_table . . . . .   | 15 |
| hmdb . . . . .                            | 16 |
| HMDB_blood . . . . .                      | 16 |
| HMDB_cell . . . . .                       | 16 |
| HMDB_urine . . . . .                      | 17 |
| is.nmr_dataset . . . . .                  | 17 |
| is.nmr_dataset_1D . . . . .               | 18 |
| is.nmr_dataset_peak_table . . . . .       | 18 |
| load_and_save_functions . . . . .         | 19 |
| models_stability_plot_bootstrap . . . . . | 20 |
| models_stability_plot_plsda . . . . .     | 22 |
| model_VIP . . . . .                       | 23 |
| MUVR_model_plot . . . . .                 | 23 |
| new_nmr_dataset . . . . .                 | 24 |
| new_nmr_dataset_1D . . . . .              | 25 |
| new_nmr_dataset_peak_table . . . . .      | 26 |
| nmr_align_find_ref . . . . .              | 27 |
| nmr_baseline_removal . . . . .            | 28 |
| nmr_baseline_threshold . . . . .          | 29 |
| nmr_batman . . . . .                      | 30 |
| nmr_batman_options . . . . .              | 31 |
| nmr_data . . . . .                        | 33 |
| nmr_dataset . . . . .                     | 34 |
| nmr_dataset_1D . . . . .                  | 34 |
| nmr_dataset_family . . . . .              | 35 |

|  |    |
|--|----|
| nmr_dataset_peak_table . . . . .                         | 35 |
| nmr_dataset_peak_table_to_SummarizedExperiment . . . . . | 35 |
| nmr_data_1r_to_SummarizedExperiment . . . . .            | 36 |
| nmr_data_analysis . . . . .                              | 37 |
| nmr_data_analysis_method . . . . .                       | 39 |
| nmr_exclude_region . . . . .                             | 40 |
| nmr_export_data_1r . . . . .                             | 41 |
| nmr_identify_regions_blood . . . . .                     | 41 |
| nmr_identify_regions_cell . . . . .                      | 42 |
| nmr_identify_regions_urine . . . . .                     | 43 |
| nmr_integrate_regions . . . . .                          | 44 |
| nmr_interpolate_1D . . . . .                             | 45 |
| nmr_meta_add . . . . .                                   | 46 |
| nmr_meta_export . . . . .                                | 48 |
| nmr_meta_get . . . . .                                   | 49 |
| nmr_meta_get_column . . . . .                            | 50 |
| nmr_normalize . . . . .                                  | 51 |
| nmr_pca_build_model . . . . .                            | 52 |
| nmr_pca_outliers . . . . .                               | 54 |
| nmr_pca_outliers_filter . . . . .                        | 55 |
| nmr_pca_outliers_plot . . . . .                          | 56 |
| nmr_pca_outliers_robust . . . . .                        | 57 |
| nmr_pca_plots . . . . .                                  | 58 |
| nmr_ppm_resolution . . . . .                             | 59 |
| nmr_read_bruker_fid . . . . .                            | 60 |
| nmr_read_samples . . . . .                               | 61 |
| nmr_zip_bruker_samples . . . . .                         | 62 |
| Parameters_blood . . . . .                               | 63 |
| Parameters_cell . . . . .                                | 63 |
| Parameters_urine . . . . .                               | 63 |
| permutation_test_model . . . . .                         | 64 |
| permutation_test_plot . . . . .                          | 65 |
| Pipelines . . . . .                                      | 67 |
| plot.nmr_dataset_1D . . . . .                            | 71 |
| plot_bootstrap_multimodel . . . . .                      | 72 |
| plot_interactive . . . . .                               | 73 |
| plot_plsda_multimodel . . . . .                          | 74 |
| plot_plsda_samples . . . . .                             | 75 |
| plot_vip_scores . . . . .                                | 76 |
| plot_webgl . . . . .                                     | 78 |
| plsda_auroc_vip_compare . . . . .                        | 79 |
| plsda_auroc_vip_method . . . . .                         | 80 |
| ppm_resolution . . . . .                                 | 81 |
| ppm_VIP_vector . . . . .                                 | 81 |
| print.nmr_dataset . . . . .                              | 82 |
| print.nmr_dataset_1D . . . . .                           | 83 |
| print.nmr_dataset_peak_table . . . . .                   | 84 |
| p_value_perm . . . . .                                   | 85 |
| random_subsampling . . . . .                             | 86 |
| rdCV_PLS_RF . . . . .                                    | 87 |
| rdCV_PLS_RF_ML . . . . .                                 | 88 |
| read_bruker_sample . . . . .                             | 89 |

|  |     |
|--|-----|
| regions_from_peak_table . . . . .                        | 90  |
| ROI_blood . . . . .                                      | 90  |
| ROI_cell . . . . .                                       | 91  |
| ROI_urine . . . . .                                      | 91  |
| save_files_to_rDolphin . . . . .                         | 91  |
| save_profiling_output . . . . .                          | 92  |
| SummarizedExperiment_to_nmr_dataset_peak_table . . . . . | 93  |
| SummarizedExperiment_to_nmr_data_1r . . . . .            | 94  |
| to_ChemoSpec . . . . .                                   | 94  |
| validate_nmr_dataset . . . . .                           | 95  |
| validate_nmr_dataset_family . . . . .                    | 96  |
| validate_nmr_dataset_peak_table . . . . .                | 97  |
| [.nmr_dataset . . . . .                                  | 102 |
| [.nmr_dataset_1D . . . . .                               | 103 |
| [.nmr_dataset_peak_table . . . . .                       | 104 |

|              |            |
|--------------|------------|
| <b>Index</b> | <b>105</b> |
|--------------|------------|

---

|                 |  |
|-----------------|--|
| AlpsNMR-package | <i>AlpsNMR: Automated spectral Processing System for NMR</i> |
|-----------------|--|

---

## Description

AlpsNMR allows you to import NMR spectra into R and provides automated and efficient signal processing for untargeted NMR metabolomics.

## Details

The following functions can be combined with the pipe. They create or modify the `nmr_dataset` object.

- `nmr_read_samples_dir()` or `nmr_read_samples()`
- `nmr_interpolate_1D()`
- `nmr_exclude_region()`
- `nmr_normalize()`
- `plot()`

There are also functions to extract the metadata and submit the samples to irods, see the example below.

The `nmr_dataset` object is essentially a list, so it is easy to access its components for further analysis.

## Author(s)

**Maintainer:** Luis Fernandez <lfernandez@ibecbarcelona.eu> ([ORCID](#))

Authors:

- Ivan Montoliu Roura <Ivan.MontoliuRoura@rd.nestle.com>
- Sergio Oller Moreno <soller@ibecbarcelona.eu> ([ORCID](#))
- Francisco Madrid Gambin <fmadrid@ibecbarcelona.eu> ([ORCID](#))
- Héctor Gracia Cabrera <hgracia@ibecbarcelona.eu>

Other contributors:

- Nestlé Institute of Health Sciences [copyright holder]
- Institute for Bioengineering of Catalonia [copyright holder]

### Examples

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
my_nmr_dataset <- dataset %>%
  nmr_interpolate_1D(axis = c(0.4, 10)) %>%
  nmr_exclude_region(exclude = list(water = c(4.6, 5))) %>%
  nmr_normalize(method = "pqn") %>%
  plot
```

---

AUC\_model

*Deprecated function Extracts AUC value*

---

### Description

The function extracts the AUC value from the middle MUVR model

### Usage

```
AUC_model(MVObj)
```

### Arguments

MVObj            a MUVR model

### Value

the AUC value of the middle model

### Examples

```
message("MUVR is not compatible with Bioconductor,
use bp_kfold_VIP_analysis method instead")
```

---

bp\_kfold\_VIP\_analysis    *K-fold bootstrap and permutation over PLS-VIP*

---

### Description

Bootstrap and permutation over PLS-VIP on AlpsNMR can be performed on both [nmr\\_dataset\\_1D](#) full spectra as well as [nmr\\_dataset\\_peak\\_table](#) peak tables.

### Usage

```
bp_kfold_VIP_analysis(dataset, y_column, k = 4, ncomp = 3, nbootstrap = 300)
```

**Arguments**

|            |   |
|------------|---|
| dataset    | An <code>nmr_dataset_family</code> object                                       |
| y_column   | A string with the name of the y column (present in the metadata of the dataset) |
| k          | Number of folds, recommended between 4 to 10                                    |
| ncomp      | number of components for the bootstrap models                                   |
| nbootstrap | number of bootstrap dataset   |

**Details**

Use of the bootstrap and permutation methods for a more robust variable importance in the projection metric for partial least squares regression, in a k-fold cross validation

**Value**

A list with the following elements:

- `important_vips`: A list with the important vips selected
- `relevant_vips`: List of vips with some relevance
- `wilcoxon_vips`: List of vips that pass a wilcoxon test
- `vip_means`: Means of the vips scores
- `vip_score_plot`: plot of the vips scores
- `kfold_results`: results of the k `bp_VIP_analysis`
- `kfold_index`: list of index of partitions of the folds

**Examples**

```
# Data analysis for a table of integrated peaks

## Generate an artificial nmr_dataset_peak_table:
### Generate artificial metadata:
num_samples <- 64 # use an even number in this example
num_peaks <- 20
metadata <- data.frame(
  NMRExperiment = as.character(1:num_samples),
  Condition = rep(c("A", "B"), times = num_samples/2),
  stringsAsFactors = FALSE
)

### The matrix with peaks
peak_means <- runif(n = num_peaks, min = 300, max = 600)
peak_sd <- runif(n = num_peaks, min = 30, max = 60)
peak_matrix <- mapply(function(mu, sd) rnorm(num_samples, mu, sd),
  mu = peak_means, sd = peak_sd)
colnames(peak_matrix) <- paste0("Peak", 1:num_peaks)
rownames(peak_matrix) <- paste0("Sample", 1:num_samples)

## Artificial differences depending on the condition:
peak_matrix[metadata$Condition == "A", "Peak2"] <-
  peak_matrix[metadata$Condition == "A", "Peak2"] + 70

peak_matrix[metadata$Condition == "A", "Peak6"] <-
  peak_matrix[metadata$Condition == "A", "Peak6"] - 60
```

```

### The nmr_dataset_peak_table
peak_table <- new_nmr_dataset_peak_table(
  peak_table = peak_matrix,
  metadata = list(external = metadata)
)

## We will use bootstrap and permutation method for VIPs selection
## in a k-fold cross validation
bp_results <- bp_kfold_VIP_analysis(peak_table, # Data to be analyzed
  y_column = "Condition", # Label
  k = 3,
  nbootstrap = 10)

message("Selected VIPs are: ", bp_results$importantn_vips)

```

bp\_VIP\_analysis

*Bootstrap and permutation over PLS-VIP***Description**

Bootstrap and permutation over PLS-VIP on AlpsNMR can be performed on both [nmr\\_dataset\\_1D](#) full spectra as well as [nmr\\_dataset\\_peak\\_table](#) peak tables.

**Usage**

```
bp_VIP_analysis(dataset, train_index, y_column, ncomp, nbootstrap = 300)
```

**Arguments**

|             |   |
|-------------|---|
| dataset     | An <a href="#">nmr_dataset_family</a> object                                    |
| train_index | set of index used to generate the bootstrap datasets                            |
| y_column    | A string with the name of the y column (present in the metadata of the dataset) |
| ncomp       | number of components used in the plsda models                                   |
| nbootstrap  | number of bootstrap dataset   |

**Details**

Use of the bootstrap and permutation methods for a more robust variable importance in the projection metric for partial least squares regression

**Value**

A list with the following elements:

- `important_vips`: A list with the important vips selected
- `relevant_vips`: List of vips with some relevance
- `pls_vip`: Pls-VIPs of every bootstrap
- `pls_vip_perm`: Pls-VIPs of every bootstrap with permuted variables
- `pls_vip_means`: Pls-VIPs normalized differences means

- `pls_vip_score_diff`: Differences of `pls_vip` and `pls_vip_perm`
- `pls_models`: pls models of the different bootstraps
- `pls_perm_models`: pls permuted models of the different bootstraps
- `classif_rate`: classification rate of the bootstrap models
- `general_model`: pls model trained with all train data
- `general_CR`: classification rate of the `general_model`
- `vips_model`: pls model trained with vips selection over all train data
- `vips_CR`: classification rate of the `vips_model`
- `error`: error expected in a t distribution
- `lower_bound`: lower bound of the confidence interval
- `upper_bound`: upper bound of the confidence interval

### Examples

```
# Data analysis for a table of integrated peaks

## Generate an artificial nmr_dataset_peak_table:
### Generate artificial metadata:
num_samples <- 32 # use an even number in this example
num_peaks <- 20
metadata <- data.frame(
  NMRExperiment = as.character(1:num_samples),
  Condition = rep(c("A", "B"), times = num_samples/2),
  stringsAsFactors = FALSE
)

### The matrix with peaks
peak_means <- runif(n = num_peaks, min = 300, max = 600)
peak_sd <- runif(n = num_peaks, min = 30, max = 60)
peak_matrix <- mapply(function(mu, sd) rnorm(num_samples, mu, sd),
                      mu = peak_means, sd = peak_sd)
colnames(peak_matrix) <- paste0("Peak", 1:num_peaks)

## Artificial differences depending on the condition:
peak_matrix[metadata$Condition == "A", "Peak2"] <-
  peak_matrix[metadata$Condition == "A", "Peak2"] + 70

peak_matrix[metadata$Condition == "A", "Peak6"] <-
  peak_matrix[metadata$Condition == "A", "Peak6"] - 60

### The nmr_dataset_peak_table
peak_table <- new_nmr_dataset_peak_table(
  peak_table = peak_matrix,
  metadata = list(external = metadata)
)

## We will use a double cross validation, splitting the samples with random
## subsampling both in the external and internal validation.
## The classification model will be a PLSDA, exploring at maximum 3 latent
## variables.
## The best model will be selected based on the area under the ROC curve
methodology <- plsda_auroc_vip_method(ncomp = 3)
model <- nmr_data_analysis(
```



```

    peak_table,
    y_column = "Condition",
    identity_column = NULL,
    external_val = list(iterations = 1, test_size = 0.25),
    internal_val = list(iterations = 3, test_size = 0.25),
    data_analysis_method = methodology
  )
## Area under ROC for each outer cross-validation iteration:
model$outer_cv_results_digested$auROC

## The number of components for the bootstrap models is selected
ncomps <- model$outer_cv_results`1`$model$ncomp
train_index <- model$train_test_partitions$outer`1`$outer_train

# Bootstrap and permutation for VIP selection
bp_VIPS <- bp_VIP_analysis(peak_table, # Data to be analyzed
                          train_index,
                          y_column = "Condition",
                          ncomp = ncomps,
                          nbootstrap = 10)

```

---

computes\_peak\_width\_ppm

*Peak width estimation for integration*

---

## Description

Estimates the peak width (ppm width) to perform peak integration using `nmr_integrate_peak_positions`. For this purpose, the full width at half maximum of a peak from alanine doublet is considered.

## Usage

```
computes_peak_width_ppm(nmr_dataset)
```

## Arguments

`nmr_dataset` An [nmr\\_dataset\\_1D](#).

## Value

Numerical. A peak width (ppm) that may be set in `nmr_integrate_peak_positions`

## See Also

Other peak integration functions: [Pipelines](#), [nmr\\_identify\\_regions\\_blood\(\)](#), [nmr\\_identify\\_regions\\_cell\(\)](#), [nmr\\_identify\\_regions\\_urine\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#)

Other `nmr_dataset_1D` functions: [\[.nmr\\_dataset\\_1D\(\)](#), [file\\_lister\(\)](#), [files\\_to\\_rDolphin\(\)](#), [format.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [nmr\\_align\\_find\\_ref\(\)](#), [nmr\\_baseline\\_removal\(\)](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_exclude\\_region\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_normalize\(\)](#), [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_outliers\(\)](#), [nmr\\_ppm\\_resolution\(\)](#),

`plot.nmr_dataset_1D()`, `plot_webgl()`, `print.nmr_dataset_1D()`, `rdCV_PLS_RF_ML()`, `rdCV_PLS_RF()`,  
`save_files_to_rDolphin()`, `to_ChemoSpec()`, `validate_nmr_dataset_peak_table()`, `validate_nmr_dataset()`

`confusion_matrix`      *Deprecated function Confusion matrix of the MUVR model*

### Description

After creating a model with the `rdCV_PLS_RF` function, you can run `confusion_matrix` on the model to make a confusion matrix from MUVR. This gives information about the model performance (e.g. classification rate).

### Usage

```
confusion_matrix(MVObj, model = "mid")
```

### Value

A confusion matrix of the model comparing actual vs predicted class

### Examples

```
message("MUVR is not compatible with Bioconductor,  
use bp_kfold_VIP_analysis method instead")
```

`files_to_rDolphin`      *Files to rDolphin*

### Description

The `rDolphin` family functions are introduced to perform automatic targeted metabolite profiling. Therefore, ensure that you interpolated from -0.1 ppm in order to consider the TSP/DSS signal at 0.0 ppm. The function generates a list with the files required by `to_rDolphin` function. Then, it is required to save them with the `save_files_to_rDolphin`. `to_rDolphin` function will read the generated "parameters.csv" file. function.

### Usage

```
files_to_rDolphin(nmr_dataset, biological_origin)
```

### Arguments

`nmr_dataset`      An `nmr_dataset` object  
`biological_origin`      String specify the type of sample (blood, urine, cell)

**Value**

a list containing:

- meta\_rDolphin: metadata in rDolphin format,
- NMR\_spectra: spectra matrix
- ROI: ROI template
- Parameters: parameters file

**See Also**

Other import/export functions: [Pipelines](#), [load\\_and\\_save\\_functions](#), [nmr\\_data\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_read\\_bruker\\_fid\(\)](#), [nmr\\_read\\_samples\(\)](#), [nmr\\_zip\\_bruker\\_samples\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [save\\_profiling\\_output\(\)](#), [to\\_ChemoSpec\(\)](#)

Other nmr\_dataset\_1D functions: [\[.nmr\\_dataset\\_1D\(\)\]](#), [computes\\_peak\\_width\\_ppm\(\)](#), [file\\_lister\(\)](#), [format.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [nmr\\_align\\_find\\_ref\(\)](#), [nmr\\_baseline\\_removal\(\)](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_exclude\\_region\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_normalize\(\)](#), [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_outliers\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [plot.nmr\\_dataset\\_1D\(\)](#), [plot\\_webgl\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [rdCV\\_PLS\\_RF\\_ML\(\)](#), [rdCV\\_PLS\\_RF\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [to\\_ChemoSpec\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

**Examples**

```
## Not run:
# Set the directory in which rDolphin files will be saved
output_dir_10_rDolphin <- file.path(your_path, "10-rDolphin")
fs::dir_create(output_dir_10_rDolphin)

# Generate the files (for plasma/serum)
files_rDolphin = files_to_rDolphin(nmr_dataset_0_10_ppm, blood)

# Save the files
save_files_to_rDolphin(files_rDolphin, output_dir_10_rDolphin)

# Build the rDolphin object. Do not forget to set the directory
setwd(output_dir_10_rDolphin)
rDolphin_object = to_rDolphin("Parameters.csv")

# Visualize your spectra
rDolphin_plot(rDolphin_object)

# Run the main profiling function (it takes a while)
targeted_profiling = Automatic_targeted_profiling(rDolphin_object)

# Save results
save_profiling_output(targeted_profiling, output_dir_10_rDolphin)

save_profiling_plots(output_dir_10_rDolphin, targeted_profiling$final_output,
targeted_profiling$reproducibility_data)

#Additionally, you can run some stats
intensities = targeted_profiling$final_output$intensity
group = as.factor(rDolphin_object$Metadata$type)
```

```
model_PLS <- rdCV_PLS_RF(X = intensities, Y = group)

## End(Not run)
```

---

|             |                        |
|-------------|------------------------|
| file_lister | <i>NMR file lister</i> |
|-------------|------------------------|

---

## Description

The function lists samples from the chosen folder required to import and create a `nmr_dataset_1D` object. The function is based on the `fs::dir_ls()` function.

## Usage

```
file_lister(dataset_path_nmr, glob)
```

## Arguments

|                  |  |
|------------------|--|
| dataset_path_nmr | A character vector of the path where samples are.  |
| glob             | A wildcard or globbing pattern common for the samples to be read, for example ending with *0 (spectra acquired by a NOESY sequence often end by 0: 10, 20, 30...) or *s (for example, samples from the tutorial in this package) passed on to <code>grep()</code> to filter paths. |

## Value

lists of samples from the chosen folder

## See Also

Other `nmr_dataset_1D` functions: `[.nmr_dataset_1D()`, `computes_peak_width_ppm()`, `files_to_rDolphin()`, `format.nmr_dataset_1D()`, `is.nmr_dataset_1D()`, `load_and_save_functions`, `new_nmr_dataset_1D()`, `nmr_align_find_ref()`, `nmr_baseline_removal()`, `nmr_baseline_threshold()`, `nmr_exclude_region()`, `nmr_integrate_regions()`, `nmr_interpolate_1D()`, `nmr_meta_add()`, `nmr_meta_export()`, `nmr_meta_get_column`, `nmr_meta_get()`, `nmr_normalize()`, `nmr_pca_build_model()`, `nmr_pca_outliers_filter()`, `nmr_pca_outliers_plot()`, `nmr_pca_outliers_robust()`, `nmr_pca_outliers()`, `nmr_ppm_resolution()`, `plot.nmr_dataset_1D()`, `plot_webgl()`, `print.nmr_dataset_1D()`, `rdCV_PLS_RF_ML()`, `rdCV_PLS_RF()`, `save_files_to_rDolphin()`, `to_ChemoSpec()`, `validate_nmr_dataset_peak_table()`, `validate_nmr_dataset()`

## Examples

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
lists_of_samples <- file_lister(dir_to_demo_dataset, "*0")
```

---

```
filter.nmr_dataset_family
```

*Keep samples based on metadata column criteria*

---

### Description

Keep samples based on metadata column criteria

### Usage

```
## S3 method for class 'nmr_dataset_family'  
filter(.data, ...)
```

### Arguments

```
.data      An nmr\_dataset\_family object  
...       conditions, as in dplyr
```

### Value

The same object, with the matching rows

### See Also

Other subsetting functions: [\[.nmr\\_dataset\\_1D\(\)](#), [\[.nmr\\_dataset\\_peak\\_table\(\)](#), [\[.nmr\\_dataset\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#)

### Examples

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")  
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)  
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))  
  
## example 1  
sample_10 <- filter(dataset_1D, NMRExperiment == "10")  
  
## example 2  
#test_samples <- dataset_1D %>% filter(nmr_peak_table$metadata$external$Group == "placebo")
```

---

```
format.nmr_dataset      Format for nmr_dataset
```

---

### Description

Format for nmr\_dataset

### Usage

```
## S3 method for class 'nmr_dataset'  
format(x, ...)
```

**Arguments**

x                    an `nmr_dataset` object  
...                  for future use

**Value**

Format for `nmr_dataset`

**See Also**

Other class helper functions: `format.nmr_dataset_1D()`, `format.nmr_dataset_peak_table()`, `is.nmr_dataset_1D()`, `is.nmr_dataset_peak_table()`, `new_nmr_dataset_1D()`, `new_nmr_dataset_peak_table()`, `new_nmr_dataset()`, `print.nmr_dataset_1D()`, `print.nmr_dataset_peak_table()`, `print.nmr_dataset()`, `validate_nmr_dataset_family()`, `validate_nmr_dataset_peak_table()`, `validate_nmr_dataset()`

Other `nmr_dataset` functions: `[.nmr_dataset()`, `load_and_save_functions`, `new_nmr_dataset()`, `nmr_interpolate_1D()`, `nmr_meta_add()`, `nmr_meta_export()`, `nmr_meta_get_column()`, `nmr_meta_get()`, `nmr_ppm_resolution()`, `nmr_read_samples()`, `print.nmr_dataset()`, `validate_nmr_dataset()`

**Examples**

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
format(dataset)
```

---

`format.nmr_dataset_1D` *format for nmr\_dataset\_1D*

---

**Description**

format for `nmr_dataset_1D`

**Usage**

```
## S3 method for class 'nmr_dataset_1D'
format(x, ...)
```

**Arguments**

x                    an `nmr_dataset_1D` object  
...                  for future use

**Value**

format for `nmr_dataset_1D`

**See Also**

Other class helper functions: `format.nmr_dataset_peak_table()`, `format.nmr_dataset()`, `is.nmr_dataset_1D()`, `is.nmr_dataset_peak_table()`, `new_nmr_dataset_1D()`, `new_nmr_dataset_peak_table()`, `new_nmr_dataset()`, `print.nmr_dataset_1D()`, `print.nmr_dataset_peak_table()`, `print.nmr_dataset()`, `validate_nmr_dataset_f`, `validate_nmr_dataset_peak_table()`, `validate_nmr_dataset()`

Other `nmr_dataset_1D` functions: `[.nmr_dataset_1D()`, `computes_peak_width_ppm()`, `file_lister()`, `files_to_rDolphin()`, `is.nmr_dataset_1D()`, `load_and_save_functions`, `new_nmr_dataset_1D()`, `nmr_align_find_ref()`, `nmr_baseline_removal()`, `nmr_baseline_threshold()`, `nmr_exclude_region()`, `nmr_integrate_regions()`, `nmr_interpolate_1D()`, `nmr_meta_add()`, `nmr_meta_export()`, `nmr_meta_get_column`, `nmr_meta_get()`, `nmr_normalize()`, `nmr_pca_build_model()`, `nmr_pca_outliers_filter()`, `nmr_pca_outliers_plot()`, `nmr_pca_outliers_robust()`, `nmr_pca_outliers()`, `nmr_ppm_resolution()`, `plot.nmr_dataset_1D()`, `plot_webgl()`, `print.nmr_dataset_1D()`, `rdCV_PLS_RF_ML()`, `rdCV_PLS_RF()`, `save_files_to_rDolphin()`, `to_ChemoSpec()`, `validate_nmr_dataset_peak_table()`, `validate_nmr_dataset()`

**Examples**

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
format(dataset_1D)
```

---

```
format.nmr_dataset_peak_table
```

*Format for nmr\_dataset\_peak\_table*

---

**Description**

Format for `nmr_dataset_peak_table`

**Usage**

```
## S3 method for class 'nmr_dataset_peak_table'
format(x, ...)
```

**Arguments**

|                  |   |
|------------------|---|
| <code>x</code>   | an <code>nmr_dataset_peak_table</code> object |
| <code>...</code> | for future use                                |

**Value**

Format for `nmr_dataset_peak_table`

**See Also**

Other `nmr_dataset_peak_table` functions: `[.nmr_dataset_peak_table()`, `is.nmr_dataset_peak_table()`, `load_and_save_functions`, `new_nmr_dataset_peak_table()`, `nmr_meta_add()`, `nmr_meta_export()`, `nmr_meta_get_column()`, `nmr_meta_get()`, `print.nmr_dataset_peak_table()`, `validate_nmr_dataset_peak_t`

Other class helper functions: `format.nmr_dataset_1D()`, `format.nmr_dataset()`, `is.nmr_dataset_1D()`, `is.nmr_dataset_peak_table()`, `new_nmr_dataset_1D()`, `new_nmr_dataset_peak_table()`, `new_nmr_dataset()`, `print.nmr_dataset_1D()`, `print.nmr_dataset_peak_table()`, `print.nmr_dataset()`, `validate_nmr_dataset_f`, `validate_nmr_dataset_peak_table()`, `validate_nmr_dataset()`

**Examples**

```

dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
meta <- file.path(dir_to_demo_dataset, "dummy_metadata.xlsx")
metadata <- readxl::read_excel(meta, sheet = 1)
dataset_1D <- nmr_meta_add(dataset_1D, metadata = metadata, by = "NMRExperiment")
metadata <- list(external = dataset_1D[["metadata"]][["external"]])
peak_table <- nmr_data(dataset_1D)
new <- new_nmr_dataset_peak_table(peak_table, metadata)
format(new)

```

hmdb

*The Human Metabolome DataBase multiplet table***Description**

The Human Metabolome DataBase multiplet table

**References**

<https://hmdb.ca/>

HMDB\_blood

*The Human Metabolome DataBase multiplet table: blood metabolites normally found in NMR-based metabolomics***Description**

The Human Metabolome DataBase multiplet table: blood metabolites normally found in NMR-based metabolomics

**References**

<https://hmdb.ca/>

HMDB\_cell

*The Human Metabolome DataBase multiplet table: cell metabolites normally found in NMR-based metabolomics***Description**

The Human Metabolome DataBase multiplet table: cell metabolites normally found in NMR-based metabolomics

**References**

<https://hmdb.ca/>



---

|            |  |
|------------|--|
| HMDB_urine | <i>The Human Metabolome DataBase multiplet table: urine metabolites normally found in NMR-based metabolomics</i> |
|------------|--|

---

### Description

The Human Metabolome DataBase multiplet table: urine metabolites normally found in NMR-based metabolomics

### References

<https://hmdb.ca/>

---

|                |   |
|----------------|---|
| is.nmr_dataset | <i>Object is of <a href="#">nmr_dataset</a> class</i> |
|----------------|---|

---

### Description

Object is of [nmr\\_dataset](#) class

### Usage

```
is.nmr_dataset(x)
```

### Arguments

x                    An object

### Value

TRUE if the object is an [nmr\\_dataset](#), FALSE otherwise

### Examples

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
is(dataset)
```

---

is.nmr\_dataset\_1D      *Object is of [nmr\\_dataset\\_1D](#) class*

---

### Description

Object is of [nmr\\_dataset\\_1D](#) class

### Usage

```
is.nmr_dataset_1D(x)
```

### Arguments

x                      an [nmr\\_dataset\\_1D](#) object

### Value

TRUE if the object is an [nmr\\_dataset\\_1D](#), FALSE otherwise

### See Also

Other class helper functions: [format.nmr\\_dataset\\_1D\(\)](#), [format.nmr\\_dataset\\_peak\\_table\(\)](#), [format.nmr\\_dataset\(\)](#), [is.nmr\\_dataset\\_peak\\_table\(\)](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [new\\_nmr\\_dataset\\_peak\\_table\(\)](#), [new\\_nmr\\_dataset\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [print.nmr\\_dataset\\_peak\\_table\(\)](#), [print.nmr\\_dataset\(\)](#), [validate\\_nmr\\_dataset\\_family\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

Other [nmr\\_dataset\\_1D](#) functions: [\[.nmr\\_dataset\\_1D\(\)](#), [computes\\_peak\\_width\\_ppm\(\)](#), [file\\_lister\(\)](#), [files\\_to\\_rDolphin\(\)](#), [format.nmr\\_dataset\\_1D\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [nmr\\_align\\_find\\_ref\(\)](#), [nmr\\_baseline\\_removal\(\)](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_exclude\\_region\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_normalize\(\)](#), [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_outliers\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [plot.nmr\\_dataset\\_1D\(\)](#), [plot\\_webgl\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [rdCV\\_PLS\\_RF\\_ML\(\)](#), [rdCV\\_PLS\\_RF\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [to\\_ChemoSpec\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

### Examples

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
result <- is(dataset_1D)
```

---

is.nmr\_dataset\_peak\_table  
                           *Object is of [nmr\\_dataset\\_peak\\_table](#) class*

---

### Description

Object is of [nmr\\_dataset\\_peak\\_table](#) class

**Usage**

```
is.nmr_dataset_peak_table(x)
```

**Arguments**

x                    an [nmr\\_dataset\\_peak\\_table](#) object

**Value**

TRUE if the object is an `nmr_dataset_peak_table`, FALSE otherwise

**See Also**

Other `nmr_dataset_peak_table` functions: [\[.nmr\\_dataset\\_peak\\_table\(\)](#), [format.nmr\\_dataset\\_peak\\_table\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_peak\\_table\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [print.nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#)

Other class helper functions: [format.nmr\\_dataset\\_1D\(\)](#), [format.nmr\\_dataset\\_peak\\_table\(\)](#), [format.nmr\\_dataset\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [new\\_nmr\\_dataset\\_peak\\_table\(\)](#), [new\\_nmr\\_dataset\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [print.nmr\\_dataset\\_peak\\_table\(\)](#), [print.nmr\\_dataset\(\)](#), [validate\\_nmr\\_dataset\\_family\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

**Examples**

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
meta <- file.path(dir_to_demo_dataset, "dummy_metadata.xlsx")
metadata <- readxl::read_excel(meta, sheet = 1)
dataset_1D <- nmr_meta_add(dataset_1D, metadata = metadata, by = "NMRExperiment")
metadata <- list(external = dataset_1D[["metadata"]][["external"]])
peak_table <- nmr_data(dataset_1D)
new <- new_nmr_dataset_peak_table(peak_table, metadata)
is(new)
```

---

load\_and\_save\_functions

*nmr\_dataset\_load*

---

**Description**

`nmr_dataset_load`

`nmr_dataset_save`

**Usage**

```
nmr_dataset_load(file_name)
```

```
nmr_dataset_save(nmr_dataset, file_name, ...)
```

**Arguments**

|             |  |
|-------------|--|
| file_name   | The file name to load or save to                         |
| nmr_dataset | An object from the <a href="#">nmr_dataset_family</a>    |
| ...         | Additional arguments passed to <a href="#">saveRDS</a> . |

**Value**

Functions to load and save nmr\_dataset objects

load nmr dataset

save nmr dataset

**See Also**

Other nmr\_dataset functions: [[.nmr\\_dataset\(\)](#), [format.nmr\\_dataset\(\)](#), [new\\_nmr\\_dataset\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [nmr\\_read\\_samples\(\)](#), [print.nmr\\_dataset\(\)](#), [validate\\_nmr\\_dataset\(\)](#)]

Other nmr\_dataset\_1D functions: [[.nmr\\_dataset\\_1D\(\)](#), [computes\\_peak\\_width\\_ppm\(\)](#), [file\\_lister\(\)](#), [files\\_to\\_rDolphin\(\)](#), [format.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [nmr\\_align\\_find\\_ref\(\)](#), [nmr\\_baseline\\_removal\(\)](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_exclude\\_region\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_normalize\(\)](#), [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_outliers\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [plot.nmr\\_dataset\\_1D\(\)](#), [plot\\_webgl\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [rdCV\\_PLS\\_RF\\_ML\(\)](#), [rdCV\\_PLS\\_RF\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [to\\_ChemoSpec\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)]

Other nmr\_dataset\_peak\_table functions: [[.nmr\\_dataset\\_peak\\_table\(\)](#), [format.nmr\\_dataset\\_peak\\_table\(\)](#), [is.nmr\\_dataset\\_peak\\_table\(\)](#), [new\\_nmr\\_dataset\\_peak\\_table\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [print.nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#)]

Other import/export functions: [Pipelines](#), [files\\_to\\_rDolphin\(\)](#), [nmr\\_data\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_read\\_bruker\\_fid\(\)](#), [nmr\\_read\\_samples\(\)](#), [nmr\\_zip\\_bruker\\_samples\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [save\\_profiling\\_output\(\)](#), [to\\_ChemoSpec\(\)](#)

**Examples**

```
# dataset <- nmr_dataset_load("test")
nmr_dataset <- nmr_dataset_load(system.file("extdata", "nmr_dataset.rds", package = "AlpsNMR"))

dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
# nmr_dataset_save(dataset, "test")
```

---

models\_stability\_plot\_bootstrap

*Models stability plot*

---

**Description**

Plot stability among models of the external cross validation

**Usage**

```
models_stability_plot_bootstrap(bp_results)
```

**Arguments**

```
bp_results      bp_kfold_VIP_analysis results
```

**Value**

A plot of models stability

**Examples**

```
# Data analysis for a table of integrated peaks

## Generate an artificial nmr_dataset_peak_table:
### Generate artificial metadata:
num_samples <- 64 # use an even number in this example
num_peaks <- 20
metadata <- data.frame(
  NMRExperiment = as.character(1:num_samples),
  Condition = rep(c("A", "B"), times = num_samples/2),
  stringsAsFactors = FALSE
)

### The matrix with peaks
peak_means <- runif(n = num_peaks, min = 300, max = 600)
peak_sd <- runif(n = num_peaks, min = 30, max = 60)
peak_matrix <- mapply(function(mu, sd) rnorm(num_samples, mu, sd),
                      mu = peak_means, sd = peak_sd)
colnames(peak_matrix) <- paste0("Peak", 1:num_peaks)

## Artificial differences depending on the condition:
peak_matrix[metadata$Condition == "A", "Peak2"] <-
  peak_matrix[metadata$Condition == "A", "Peak2"] + 70

peak_matrix[metadata$Condition == "A", "Peak6"] <-
  peak_matrix[metadata$Condition == "A", "Peak6"] - 60

### The nmr_dataset_peak_table
peak_table <- new_nmr_dataset_peak_table(
  peak_table = peak_matrix,
  metadata = list(external = metadata)
)

## We will use bootstrap and permutation method for VIPs selection
## in a a k-fold cross validation
#bp_results <- bp_kfold_VIP_analysis(peak_table, # Data to be analyzed
#                                y_column = "Condition", # Label
#                                k = 3,
#                                nbootstrap = 10)

#message("Selected VIPs are: ", bp_results$important_vips)

#models_stability_plot_bootstrap(bp_results)
```

---

```
models_stability_plot_plsda
      Models stability plot
```

---

**Description**

Plot stability among models of the external cross validation

**Usage**

```
models_stability_plot_plsda(model)
```

**Arguments**

```
model          A nmr_data_analysis_model
```

**Value**

A plot of models stability

**Examples**

```
# Data analysis for a table of integrated peaks

## Generate an artificial nmr_dataset_peak_table:
### Generate artificial metadata:
num_samples <- 32 # use an even number in this example
num_peaks <- 20
metadata <- data.frame(
  NMRExperiment = as.character(1:num_samples),
  Condition = rep(c("A", "B"), times = num_samples/2),
  stringsAsFactors = FALSE
)

### The matrix with peaks
peak_means <- runif(n = num_peaks, min = 300, max = 600)
peak_sd <- runif(n = num_peaks, min = 30, max = 60)
peak_matrix <- mapply(function(mu, sd) rnorm(num_samples, mu, sd),
                      mu = peak_means, sd = peak_sd)
colnames(peak_matrix) <- paste0("Peak", 1:num_peaks)

## Artificial differences depending on the condition:
peak_matrix[metadata$Condition == "A", "Peak2"] <-
  peak_matrix[metadata$Condition == "A", "Peak2"] + 70

peak_matrix[metadata$Condition == "A", "Peak6"] <-
  peak_matrix[metadata$Condition == "A", "Peak6"] - 60

### The nmr_dataset_peak_table
peak_table <- new_nmr_dataset_peak_table(
  peak_table = peak_matrix,
  metadata = list(external = metadata)
)
```

```

methodology <- plsda_auroc_vip_method(ncomp = 3)
model <- nmr_data_analysis(
  peak_table,
  y_column = "Condition",
  identity_column = NULL,
  external_val = list(iterations = 3, test_size = 0.25),
  internal_val = list(iterations = 3, test_size = 0.25),
  data_analysis_method = methodology
)

#models_stability_plot_plsda(model)

```

---

model\_VIP

*Deprecated function Model VIP values*


---

### Description

Once, the MVObj is created and validated, this function extracts autoselected ranked variables from the model (MUVR object). See rdCV\_PLS\_RF function.

### Usage

```
model_VIP(MVObj, model = "mid")
```

### Arguments

MVObj                    a MUVR model

### Value

a data frame with the order, name and average rank of selected variables

### Examples

```

message("MUVR is not compatible with Bioconductor,
use bp_kfold_VIP_analysis method instead")

```

---

MUVR\_model\_plot

*Deprecated function Model plot*


---

### Description

Plot the model (a MUVR object) obtained from rdCV\_PLS\_RF function. For more information about the multivariate modelling with minimally biased variable selection (MUVR) from the MUVR package, see Shi et al., 2018 (DOI: 10.1093/bioinformatics/bty710).

### Usage

```
MUVR_model_plot(MVObj, model = "mid", factCols, sampLabels, ylim = NULL)
```

**Value**

A plot with the model performance

**Examples**

```
message("MUVR is not compatible with Bioconductor,
use bp_kfold_VIP_analysis method instead")
```

---

|                 |                                     |
|-----------------|-------------------------------------|
| new_nmr_dataset | <i>Create an nmr_dataset object</i> |
|-----------------|-------------------------------------|

---

**Description**

Create an nmr\_dataset object

**Usage**

```
new_nmr_dataset(metadata, data_fields, axis)
```

**Arguments**

|             |                                  |
|-------------|----------------------------------|
| metadata    | A named list of data frames      |
| data_fields | A named list. Check the examples |
| axis        | A list. Check the examples       |

**Value**

Create an nmr\_dataset object

Create an nmr\_dataset object

**See Also**

Other class helper functions: [format.nmr\\_dataset\\_1D\(\)](#), [format.nmr\\_dataset\\_peak\\_table\(\)](#), [format.nmr\\_dataset\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_peak\\_table\(\)](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [new\\_nmr\\_dataset\\_peak\\_table\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [print.nmr\\_dataset\\_peak\\_table\(\)](#), [print.nmr\\_dataset\(\)](#), [validate\\_nmr\\_dataset\\_family\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

Other nmr\_dataset functions: [\[.nmr\\_dataset\(\)](#), [format.nmr\\_dataset\(\)](#), [load\\_and\\_save\\_functions](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [nmr\\_read\\_samples\(\)](#), [print.nmr\\_dataset\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

**Examples**

```
#
metadata_1D <- list(external = data.frame(NMRExperiment = c("10", "20")))
# Sample 10 and Sample 20 can have different lengths (due to different setups)
data_fields_1D <- list(data_1r = list(runif(16), runif(32)))
# Each sample has its own axis list, with one element (because this example is 1D)
axis_1D <- list(list(1:16), list(1:32))
my_1D_data <- new_nmr_dataset(metadata_1D, data_fields_1D, axis_1D)
```



```
# Example for 2D samples
metadata_2D <- list(external = data.frame(NMRExperiment = c("11", "21")))
data_fields_2D <- list(data_2rr = list(matrix(runif(16*3), nrow=16, ncol=3),
      runif(32*3), nrow=32, ncol=3))
# Each sample has its own axis list, with one element (because this example is 1D)
axis_2D <- list(list(1:16, 1:3), list(1:32, 1:3))
my_2D_data <- new_nmr_dataset(metadata_2D, data_fields_2D, axis_2D)
```

---

new\_nmr\_dataset\_1D      *Creates a new 1D nmr\_dataset object from scratch*

---

## Description

Creates a new 1D nmr\_dataset object from scratch

## Usage

```
new_nmr_dataset_1D(ppm_axis, data_1r, metadata)
```

## Arguments

|          |   |
|----------|---|
| ppm_axis | A numeric vector with the ppm values for the columns of data_1r |
| data_1r  | A numeric matrix with one NMR spectrum on each row              |
| metadata | A list of data frames with at least the NMRExperiment column    |

## Value

Creates a new 1D nmr\_dataset object from scratch

## See Also

Other class helper functions: [format.nmr\\_dataset\\_1D\(\)](#), [format.nmr\\_dataset\\_peak\\_table\(\)](#), [format.nmr\\_dataset\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_peak\\_table\(\)](#), [new\\_nmr\\_dataset\\_peak\\_table\(\)](#), [new\\_nmr\\_dataset\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [print.nmr\\_dataset\\_peak\\_table\(\)](#), [print.nmr\\_dataset\(\)](#), [validate\\_nmr\\_dataset\\_family\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

Other nmr\_dataset\_1D functions: [\[.nmr\\_dataset\\_1D\(\)](#), [computes\\_peak\\_width\\_ppm\(\)](#), [file\\_lister\(\)](#), [files\\_to\\_rDolphin\(\)](#), [format.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [load\\_and\\_save\\_functions](#), [nmr\\_align\\_find\\_ref\(\)](#), [nmr\\_baseline\\_removal\(\)](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_exclude\\_region\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_normalize\(\)](#), [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_outliers\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [plot.nmr\\_dataset\\_1D\(\)](#), [plot\\_webgl\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [rdCV\\_PLS\\_RF\\_ML\(\)](#), [rdCV\\_PLS\\_RF\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [to\\_ChemoSpec\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

**Examples**

```
# Create a random spectra matrix
nsamp <- 12
npoints <- 20
dummy_ppm_axis <- seq(from = 0.2, to = 10, length.out = npoints)
dummy_spectra_matrix <- matrix(runif(nsamp*npoints), nrow = nsamp, ncol = npoints)
metadata <- list(external = data.frame(NMRExperiment = paste0("Sample", 1:12),
                                     DummyClass = c("a", "b"),
                                     stringsAsFactors = FALSE))
dummy_nmr_dataset_1D <- new_nmr_dataset_1D(ppm_axis = dummy_ppm_axis,
                                           data_1r = dummy_spectra_matrix,
                                           metadata = metadata)
```

---

new\_nmr\_dataset\_peak\_table

*Creates a new nmr\_dataset\_peak\_table object from scratch*

---

**Description**

Creates a new nmr\_dataset\_peak\_table object from scratch

**Usage**

```
new_nmr_dataset_peak_table(peak_table, metadata)
```

**Arguments**

|            |  |
|------------|--|
| peak_table | A numeric matrix with one NMR spectrum on each row           |
| metadata   | A list of data frames with at least the NMRExperiment column |

**Value**

Creates a new nmr\_dataset\_peak\_table object from scratch

**See Also**

Other nmr\_dataset\_peak\_table functions: [\[.nmr\\_dataset\\_peak\\_table\(\)](#), [format.nmr\\_dataset\\_peak\\_table\(\)](#), [is.nmr\\_dataset\\_peak\\_table\(\)](#), [load\\_and\\_save\\_functions](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [print.nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#)

Other class helper functions: [format.nmr\\_dataset\\_1D\(\)](#), [format.nmr\\_dataset\\_peak\\_table\(\)](#), [format.nmr\\_dataset\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_peak\\_table\(\)](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [new\\_nmr\\_dataset\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [print.nmr\\_dataset\\_peak\\_table\(\)](#), [print.nmr\\_dataset\(\)](#), [validate\\_nmr\\_dataset\\_family\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

## Examples

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
meta <- file.path(dir_to_demo_dataset, "dummy_metadata.xlsx")
metadata <- readxl::read_excel(meta, sheet = 1)
dataset_1D <- nmr_meta_add(dataset_1D, metadata = metadata, by = "NMRExperiment")
metadata <- list(external = dataset_1D[["metadata"]][["external"]])
peak_table <- nmr_data(dataset_1D)
new <- new_nmr_dataset_peak_table(peak_table, metadata)
```

---

nmr\_align\_find\_ref      *Find alignment reference*

---

## Description

Find alignment reference

## Usage

```
nmr_align_find_ref(nmr_dataset, peak_data)
```

## Arguments

nmr\_dataset      An [nmr\\_dataset\\_1D](#)

peak\_data        The detected peak data given by [nmr\\_detect\\_peaks](#).

## Value

The NMRExperiment of the reference sample

## See Also

Other alignment functions: [Pipelines](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#)

Other peak alignment functions: [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#)

Other `nmr_dataset_1D` functions: [.nmr\\_dataset\\_1D\(\)](#), [computes\\_peak\\_width\\_ppm\(\)](#), [file\\_lister\(\)](#), [files\\_to\\_rDolphin\(\)](#), [format.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [nmr\\_baseline\\_removal\(\)](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_exclude\\_region\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_normalize\(\)](#), [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_outliers\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [plot.nmr\\_dataset\\_1D\(\)](#), [plot\\_webgl\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [rdCV\\_PLS\\_RF\\_ML\(\)](#), [rdCV\\_PLS\\_RF\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [to\\_ChemoSpec\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

---

nmr\_baseline\_removal *Baseline Removal NMR*

---

## Description

Removes the baseline on an `nmr_dataset_1D` object, using `baseline::baseline.als`.

## Usage

```
nmr_baseline_removal(nmr_dataset, lambda = 6, p = 0.05, maxit = 20)
```

## Arguments

|                          |                                  |
|--------------------------|----------------------------------|
| <code>nmr_dataset</code> | An <code>nmr_dataset_1D</code> . |
| <code>lambda</code>      | 2nd derivative constraint        |
| <code>p</code>           | Weighting of positive residuals  |
| <code>maxit</code>       | Maximum number of iterations     |

## Value

The same `nmr_dataset_1D` object after baseline removal.

## See Also

[baseline::baseline.als](#)

Other `nmr_dataset_1D` functions: [\[.nmr\\_dataset\\_1D\(\)](#), [computes\\_peak\\_width\\_ppm\(\)](#), [file\\_lister\(\)](#), [files\\_to\\_rDolphin\(\)](#), [format.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [nmr\\_align\\_find\\_ref\(\)](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_exclude\\_region\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_normalize\(\)](#), [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_outliers\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [plot.nmr\\_dataset\\_1D\(\)](#), [plot\\_webgl\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [rdCV\\_PLS\\_RF\\_ML\(\)](#), [rdCV\\_PLS\\_RF\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [to\\_ChemoSpec\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

## Examples

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
dataset_no_base_line <- nmr_baseline_removal(dataset_1D, lambda = 6, p = 0.01)
```

---

`nmr_baseline_threshold`*Threshold estimation for peak detection*

---

## Description

Estimates the threshold value for peak detection on an `nmr_dataset_1D` object. This is performed computing the mean and the standard deviation of each spectrum beyond 9.5 ppm. The threshold is then averaged of means and adding 3 times the mean of the standard deviations

## Usage

```
nmr_baseline_threshold(nmr_dataset)
```

## Arguments

`nmr_dataset` An `nmr_dataset_1D`.

## Value

Numerical. A threshold value in intensity below that no peak is detected.

## See Also

Other peak detection functions: [Pipelines](#), [nmr\\_identify\\_regions\\_blood\(\)](#), [nmr\\_identify\\_regions\\_cell\(\)](#), [nmr\\_identify\\_regions\\_urine\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [regions\\_from\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#)

Other `nmr_dataset_1D` functions: [\[.nmr\\_dataset\\_1D\(\)](#), [computes\\_peak\\_width\\_ppm\(\)](#), [file\\_lister\(\)](#), [files\\_to\\_rDolphin\(\)](#), [format.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [nmr\\_align\\_find\\_ref\(\)](#), [nmr\\_baseline\\_removal\(\)](#), [nmr\\_exclude\\_region\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_normalize\(\)](#), [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_outliers\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [plot.nmr\\_dataset\\_1D\(\)](#), [plot\\_webgl\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [rdCV\\_PLS\\_RF\\_ML\(\)](#), [rdCV\\_PLS\\_RF\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [to\\_ChemoSpec\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

## Examples

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
bl_threshold <- nmr_baseline_threshold(dataset_1D)
```

---

nmr\_batman

*Batman helpers*

---

## Description

Batman helpers

## Usage

```
nmr_batman_write_options(  
  bopts,  
  batman_dir = "BatmanInput",  
  filename = "batmanOptions.txt"  
)  
  
nmr_batman_export_dataset(  
  nmr_dataset,  
  batman_dir = "BatmanInput",  
  filename = "NMRdata.txt"  
)  
  
nmr_batman_multi_data_user_hmdb(  
  batman_dir = "BatmanInput",  
  filename = "multi_data_user.csv"  
)  
  
nmr_batman_multi_data_user(  
  multiplet_table,  
  batman_dir = "BatmanInput",  
  filename = "multi_data_user.csv"  
)  
  
nmr_batman_metabolites_list(  
  metabolite_names,  
  batman_dir = "BatmanInput",  
  filename = "metabolitesList.csv"  
)
```

## Arguments

|                  |  |
|------------------|--|
| bopts            | Batman options   |
| batman_dir       | Batman input directory                                 |
| filename         | Filename to use, inside batman_dir                     |
| nmr_dataset      | An <a href="#">nmr_dataset_1D</a> object               |
| multiplet_table  | A data frame, like the <a href="#">hmdb</a> dataset    |
| metabolite_names | A character vector of the metabolite names to consider |

**Value**

These are helper functions to make Batman tests easier

**See Also**

Other batman functions: [nmr\\_batman\\_options\(\)](#)

**Examples**

```
bopts <- nmr_batman_options()
# nmr_batman_write_options(bopts)

dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
#nmr_batman_export_dataset(dataset_1D)

message("Use of multi_data_user_hmdb")
#multi_data_user_hmdb <- nmr_batman_multi_data_user_hmdb()
hmdb <- NULL
#utils::data("hmdb", package = "AlpsNMR", envir = environment())
#hmdb <- nmr_batman_multi_data_user(hmdb)

metabolite_names <- c("alanine", "glucose")
#metabolite_names <- nmr_batman_metabolites_list(metabolite_names)
```

---

|                    |                              |
|--------------------|------------------------------|
| nmr_batman_options | <i>Batman Options helper</i> |
|--------------------|------------------------------|

---

**Description**

Batman Options helper

**Usage**

```
nmr_batman_options(
  ppmRange = matrix(c(3, 3.1, 3.6, 3.7, 3.9, 4, 4, 4.1, 6.95, 7.05, 7.6, 7.7, 7.8,
    7.9), ncol = 2, byrow = TRUE),
  specNo = "1",
  paraProc = 4L,
  negThresh = -0.5,
  scaleFac = 1e+06,
  downSamp = 1,
  hiresFlag = 1,
  randSeed = 100025L,
  nItBurnin = 200L,
  nItPostBurnin = 5000L,
  multFile = 2L,
  thinning = 50L,
  cfeFlag = 0,
  nItRerun = 5000L,
```

```

startTemp = 1000,
specFreq = 600,
a = 1e-05,
b = 1e-09,
muMean = 1.1,
muVar = 0.2,
muVar_prop = 0.002,
nuMVar = 0.0025,
nuMVarProp = 0.1,
tauMean = -0.05,
tauPrec = 2,
rdelta = 0.02,
csFlag = 0
)

```

### Arguments

|               |  |
|---------------|--|
| ppmRange      | Range of ppm to process  |
| specNo        | Index of spectra to process  |
| paraProc      | Number of cores to use   |
| negThresh     | Truncation threshold for negative intensities  |
| scaleFac      | Divide each spectrum by this number  |
| downSamp      | Decimate each spectrum by this factor  |
| hiresFlag     | Keep High Resolution deconvolved spectra   |
| randSeed      | A random seed  |
| nItBurnin     | Number of burn-in iterations   |
| nItPostBurnin | Number of iterations after burn-in   |
| multFile      | Multiplet file (integer)   |
| thinning      | Save MCMC state every thinning iterations  |
| cfeFlag       | Same concentration for all spectra (fixed effect)  |
| nItRerun      | Number of iterations for a batman rerun  |
| startTemp     | Start temperature  |
| specFreq      | NMR Spectrometer frequency   |
| a             | Shape parameter for the gamma distribution (used for lambda, the precision)                |
| b             | Rate distribution parameter for the gamma distribution (used for lambda, the precision)    |
| muMean        | Peak width mean in ln(Hz)  |
| muVar         | Peak width variance in ln(Hz)  |
| muVar_prop    | Peak width proposed variance in ln(Hz)   |
| nuMVar        | Peak width metabolite variance in ln(Hz)   |
| nuMVarProp    | Peak width metabolite proposed variance in ln(Hz)  |
| tauMean       | mean of the prior on tau   |
| tauPrec       | inverse of variance of prior on tau  |
| rdelta        | Truncation of the prior on peak shift (ppm)  |
| csFlag        | Specify chemical shift for each multiplet in each spectrum? (chemShiftperSpectra.csv file) |



**Value**

A batman\_options object with the Batman Options

**See Also**

Other batman functions: [nmr\\_batman](#)

**Examples**

```
bopts <- nmr_batman_options()
```

---

|          |   |
|----------|---|
| nmr_data | <i>Set/Return the full spectra matrix</i> |
|----------|---|

---

**Description**

Set/Return the full spectra matrix

**Usage**

```
nmr_data(nmr_dataset, ...)

nmr_data(nmr_dataset) <- value

## S3 replacement method for class 'nmr_dataset_1D'
nmr_data(nmr_dataset) <- value
```

**Arguments**

|             |  |
|-------------|--|
| nmr_dataset | An object from the <a href="#">nmr_dataset_family</a> to get the raw data from |
| ...         | Unused and left for future compatibility                                       |
| value       | A matrix   |

**Value**

a matrix  
The given nmr\_dataset

**See Also**

Other import/export functions: [Pipelines](#), [files\\_to\\_rDolphin\(\)](#), [load\\_and\\_save\\_functions](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_read\\_bruker\\_fid\(\)](#), [nmr\\_read\\_samples\(\)](#), [nmr\\_zip\\_bruker\\_samples\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [save\\_profiling\\_output\(\)](#), [to\\_ChemoSpec\(\)](#)

**Examples**

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
dataset_data <- nmr_data(dataset_1D)
```

---

|             |                               |
|-------------|-------------------------------|
| nmr_dataset | <i>nmr_dataset (S3 class)</i> |
|-------------|-------------------------------|

---

### Description

An `nmr_dataset` represents a set of NMR samples. It is defined as an S3 class, and it can be treated as a regular list.

### Details

It currently has the following elements:

- `metadata`: A list of data frames. Each data frame contains metadata of a given area (acquisition parameters, preprocessing parameters, general sample information...)
- `axis`: A list with length equal to the dimensionality of the data. For 1D spectra it is a list with a numeric vector
- `data_*`: Data arrays with the actual spectra. The first index represents the sample, the rest of the indices match the length of each axis. Typically `data_1r` is a matrix with one sample on each row and the chemical shifts in the columns.
- `num_samples`: The number of samples in the dataset

### See Also

[Functions to save and load these objects](#)

Other AlpsNMR dataset objects: [nmr\\_dataset\\_1D](#), [nmr\\_dataset\\_family](#)

---

|                |                                  |
|----------------|----------------------------------|
| nmr_dataset_1D | <i>nmr_dataset_1D (S3 class)</i> |
|----------------|----------------------------------|

---

### Description

An `nmr_dataset_1D` represents a set of 1D interpolated NMR samples. It is defined as an S3 class, and it can be treated as a regular list.

### Details

It currently has the following elements:

- `metadata`: A list of data frames. Each data frame contains metadata of a given area (acquisition parameters, preprocessing parameters, general sample information...)
- `axis`: A numeric vector with the chemical shift axis in ppm.
- `data_1r`: A matrix with one sample on each row and the chemical shifts in the columns.

### See Also

Other AlpsNMR dataset objects: [nmr\\_dataset\\_family](#), [nmr\\_dataset](#)

---

nmr\_dataset\_family     *nmr\_dataset like objects (S3 classes)*

---

### Description

The AlpsNMR package defines and uses several objects to manage NMR Data.

### Details

These objects share some structure and functions, so it makes sense to have an abstract class to ensure that the shared structures are compatible

### See Also

[Functions to save and load these objects](#)

Other AlpsNMR dataset objects: [nmr\\_dataset\\_1D](#), [nmr\\_dataset](#)

---

nmr\_dataset\_peak\_table  
*nmr\_dataset\_peak\_table (S3 class)*

---

### Description

An `nmr_dataset_peak_table` represents a peak table with metadata. It is defined as an S3 class, and it can be treated as a regular list.

### Details

- `metadata`: A list of data frames. Each data frame contains metadata. Usually the list only has one data frame named "external".
- `peak_table`: A matrix with one sample on each row and the peaks in the columns

---

nmr\_dataset\_peak\_table\_to\_SummarizedExperiment  
*Export nmr\_dataset\_peak\_table to SummarizedExperiment*

---

### Description

Export `nmr_dataset_peak_table` to `SummarizedExperiment`

### Usage

```
nmr_dataset_peak_table_to_SummarizedExperiment(nmr_peak_table)
```

### Arguments

`nmr_peak_table` An [nmr\\_dataset\\_peak\\_table](#) object

**Value**

SummarizedExperiment object (unmodified)

**Examples**

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
meta <- file.path(dir_to_demo_dataset, "dummy_metadata.xlsx")
metadata <- readxl::read_excel(meta, sheet = 1)
dataset_1D <- nmr_meta_add(dataset_1D, metadata = metadata, by = "NMRExperiment")
metadata <- list(external = dataset_1D[["metadata"]][["external"]])
peak_table <- nmr_data(dataset_1D)
nmr_peak_table <- new_nmr_dataset_peak_table(peak_table, metadata)
se <- nmr_dataset_peak_table_to_SummarizedExperiment(nmr_peak_table)
```

---

nmr\_data\_1r\_to\_SummarizedExperiment

*Export 1D NMR data to SummarizedExperiment*

---

**Description**

Export 1D NMR data to SummarizedExperiment

**Usage**

```
nmr_data_1r_to_SummarizedExperiment(nmr_dataset)
```

**Arguments**

nmr\_dataset     An [nmr\\_dataset\\_1D](#) object

**Value**

SummarizedExperiment An SummarizedExperiment object (unmodified)

**Examples**

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
se <- nmr_data_1r_to_SummarizedExperiment(dataset_1D)
```

---

nmr\_data\_analysis      *Data analysis*

---

### Description

Data analysis on AlpsNMR can be performed on both [nmr\\_dataset\\_ID](#) full spectra as well as [nmr\\_dataset\\_peak\\_table](#) peak tables.

### Usage

```
nmr_data_analysis(  
  dataset,  
  y_column,  
  identity_column,  
  external_val,  
  internal_val,  
  data_analysis_method  
)
```

### Arguments

|                            |  |
|----------------------------|--|
| dataset                    | An <a href="#">nmr_dataset_family</a> object   |
| y_column                   | A string with the name of the y column (present in the metadata of the dataset)                                |
| identity_column            | NULL or a string with the name of the identity column (present in the metadata of the dataset).                |
| external_val, internal_val | A list with two elements: iterations and test_size. See <a href="#">random_subsampling</a> for further details |
| data_analysis_method       | An <a href="#">nmr_data_analysis_method</a> object   |

### Details

The workflow consists of a double cross validation strategy using random subsampling for splitting into train and test sets. The classification model and the metric to choose the best model can be customized (see [new\\_nmr\\_data\\_analysis\\_method\(\)](#)), but for now only a PLSDA classification model with a best area under ROC curve metric is implemented (see the examples here and [plsda\\_auroc\\_vip\\_method](#))

### Value

A list with the following elements:

- `train_test_partitions`: A list with the indices used in train and test on each of the cross-validation iterations
- `inner_cv_results`: The output returned by `train_evaluate_model` on each inner cross-validation
- `inner_cv_results_digested`: The output returned by `choose_best_inner`.
- `outer_cv_results`: The output returned by `train_evaluate_model` on each outer cross-validation
- `outer_cv_results_digested`: The output returned by `train_evaluate_model_digest_outer`.

**Examples**

```

# Data analysis for a table of integrated peaks

## Generate an artificial nmr_dataset_peak_table:
### Generate artificial metadata:
num_samples <- 32 # use an even number in this example
num_peaks <- 20
metadata <- data.frame(
  NMRExperiment = as.character(1:num_samples),
  Condition = rep(c("A", "B"), times = num_samples/2),
  stringsAsFactors = FALSE
)

### The matrix with peaks
peak_means <- runif(n = num_peaks, min = 300, max = 600)
peak_sd <- runif(n = num_peaks, min = 30, max = 60)
peak_matrix <- mapply(function(mu, sd) rnorm(num_samples, mu, sd),
                      mu = peak_means, sd = peak_sd)

colnames(peak_matrix) <- paste0("Peak", 1:num_peaks)

## Artificial differences depending on the condition:
peak_matrix[metadata$Condition == "A", "Peak2"] <-
  peak_matrix[metadata$Condition == "A", "Peak2"] + 70

peak_matrix[metadata$Condition == "A", "Peak6"] <-
  peak_matrix[metadata$Condition == "A", "Peak6"] - 60

### The nmr_dataset_peak_table
peak_table <- new_nmr_dataset_peak_table(
  peak_table = peak_matrix,
  metadata = list(external = metadata)
)

## We will use a double cross validation, splitting the samples with random
## subsampling both in the external and internal validation.
## The classification model will be a PLSDA, exploring at maximum 3 latent
## variables.
## The best model will be selected based on the area under the ROC curve
methodology <- plsda_auroc_vip_method(ncomp = 3)
model <- nmr_data_analysis(
  peak_table,
  y_column = "Condition",
  identity_column = NULL,
  external_val = list(iterations = 3, test_size = 0.25),
  internal_val = list(iterations = 3, test_size = 0.25),
  data_analysis_method = methodology
)

## Area under ROC for each outer cross-validation iteration:
model$outter_cv_results_digested$auroc
## Rank Product of the Variable Importance in the Projection
## (Lower means more important)
sort(model$outter_cv_results_digested$vip_rankproducts)

```

---

nmr\_data\_analysis\_method

*Create method for NMR data analysis*


---

## Description

Create method for NMR data analysis

## Usage

```
new_nmr_data_analysis_method(
  train_evaluate_model,
  train_evaluate_model_params_inner,
  choose_best_inner,
  train_evaluate_model_params_outer,
  train_evaluate_model_digest_outer
)
```

## Arguments

`train_evaluate_model`

A function. The `train_evaluate_model` must have the following signature:

```
function(x_train, y_train, identity_train, x_test, y_test, identity_test, ...)
```

The `x_train` and `y_train` (and their test counterparts) are self-explanatory.

The `identity_` arguments are expected to be factors. They can be used for instance with a callback that uses `mixOmics::plsda` in a multilevel approach for longitudinal studies. In those studies the `identity` would be an identifier of the subject.

The `...` arguments are free to be defined for each `train_evaluate_model`.

`train_evaluate_model_params_inner`, `train_evaluate_model_params_outer`

A list with additional arguments to pass to `train_evaluate_model` either in the inner cv loop or in the outer cv loop.

`choose_best_inner`

A function with a single argument:

```
function(inner_cv_results)
```

The argument is a list of `train_evaluate_model` outputs. The return value of must be a list with at least an element named `train_evaluate_model_args`. `train_evaluate_model_args` must be a named list.

- Each element must be named as one of the `train_evaluate_model` arguments.
- Each element must be a vector as long as the number of outer cross-validations.
- The values of each vector must be the values that the `train_evaluate_model` argument must take on each outer cross-validation iteration. Additional list elements can be returned and will be given back to the user.

`train_evaluate_model_digest_outer`

A function with a single argument:

```
function(outer_cv_results)
```

The argument is a list of `train_evaluate_model` outputs in outer cross-validation. The return value is returned by `nmr_data_analysis`

**Value**

An object encapsulating the method dependent functions that can be used with [nmr\\_data\\_analysis](#)

**Examples**

```
help(new_nmr_data_analysis_method)
```

---

|                    |                                    |
|--------------------|------------------------------------|
| nmr_exclude_region | <i>Exclude region from samples</i> |
|--------------------|------------------------------------|

---

**Description**

Excludes a given region (for instance to remove the water peak)

**Usage**

```
nmr_exclude_region(samples, exclude = list(water = c(4.7, 5)))

## S3 method for class 'nmr_dataset_1D'
nmr_exclude_region(samples, exclude = list(water = c(4.7, 5)))
```

**Arguments**

|         |   |
|---------|---|
| samples | An object   |
| exclude | A list with regions to be removed Typically: <code>exclude = list(water = c(4.7, 5.0))</code> |

**Value**

The same object, with the regions excluded

**See Also**

Other `nmr_dataset_1D` functions: [\[.nmr\\_dataset\\_1D\(\)](#), [computes\\_peak\\_width\\_ppm\(\)](#), [file\\_lister\(\)](#), [files\\_to\\_rDolphin\(\)](#), [format.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [nmr\\_align\\_find\\_ref\(\)](#), [nmr\\_baseline\\_removal\(\)](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_normalize\(\)](#), [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_outliers\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [plot.nmr\\_dataset\\_1D\(\)](#), [plot\\_webgl\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [rdCV\\_PLS\\_RF\\_ML\(\)](#), [rdCV\\_PLS\\_RF\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [to\\_ChemoSpec\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

**Examples**

```
nmr_dataset <- nmr_dataset_load(system.file("extdata", "nmr_dataset.rds", package = "AlpsNMR"))
exclude_regions <- list(water = c(5.1, 4.5))
nmr_dataset <- nmr_exclude_region(nmr_dataset, exclude = exclude_regions)

nmr_dataset <- nmr_dataset_load(system.file("extdata", "nmr_dataset.rds", package = "AlpsNMR"))
exclude_regions <- list(water = c(5.1, 4.5))
nmr_dataset <- nmr_exclude_region(nmr_dataset, exclude = exclude_regions)
```



---

nmr\_export\_data\_1r      *Export 1D NMR data to a CSV file*

---

### Description

Export 1D NMR data to a CSV file

### Usage

```
nmr_export_data_1r(nmr_dataset, filename)
```

### Arguments

nmr\_dataset      An [nmr\\_dataset\\_1D](#) object  
filename          The csv filename

### Value

The nmr\_dataset object (unmodified)

### Examples

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")  
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)  
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))  
#nmr_export_data_1r(dataset_1D, "exported_nmr_dataset")
```

---

nmr\_identify\_regions\_blood  
                          *NMR peak identification (plasma/serum samples)*

---

### Description

Identify given regions and return a data frame with plausible assignments in human plasma/serum samples.

### Usage

```
nmr_identify_regions_blood(  
  ppm_to_assign,  
  num_proposed_compounds = 3,  
  verbose = FALSE  
)
```

### Arguments

ppm\_to\_assign    A vector with the ppm regions to assign  
num\_proposed\_compounds      set the number of proposed metabolites sorted by the number times reported in the HMDB: HMDB\_blood.  
verbose            Logical value. Set it to TRUE to print additional information

**Value**

a data frame with plausible assignments.

**See Also**

Other peak detection functions: [Pipelines](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_identify\\_regions\\_cell\(\)](#), [nmr\\_identify\\_regions\\_urine\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [regions\\_from\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#)

Other peak integration functions: [Pipelines](#), [computes\\_peak\\_width\\_ppm\(\)](#), [nmr\\_identify\\_regions\\_cell\(\)](#), [nmr\\_identify\\_regions\\_urine\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#)

**Examples**

```
# We identify regions from from the corresponding ppm stored in a vector.
ppm_to_assign <- c(4.060960203, 3.048970634, 2.405935596,
3.24146865, 0.990616851, 1.002075066, 0.955325548)
identification <- nmr_identify_regions_blood (ppm_to_assign)
```

---

nmr\_identify\_regions\_cell

*NMR peak identification (cell samples)*

---

**Description**

Identify given regions and return a data frame with plausible assignments in cell samples.

**Usage**

```
nmr_identify_regions_cell(  
  ppm_to_assign,  
  num_proposed_compounds = 3,  
  verbose = FALSE  
)
```

**Arguments**

ppm\_to\_assign A vector with the ppm regions to assign  
num\_proposed\_compounds set the number of proposed metabolites in HMDB\_cell.  
verbose Logical value. Set it to TRUE to print additional information

**Value**

a data frame with plausible assignments.

**See Also**

Other peak detection functions: [Pipelines](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_identify\\_regions\\_blood\(\)](#), [nmr\\_identify\\_regions\\_urine\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [regions\\_from\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#)

Other peak integration functions: [Pipelines](#), [computes\\_peak\\_width\\_ppm\(\)](#), [nmr\\_identify\\_regions\\_blood\(\)](#), [nmr\\_identify\\_regions\\_urine\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#)

## Examples

```
# We identify regions from from the corresponding ppm stored in a vector.
ppm_to_assign <- c(4.060960203, 3.048970634, 2.405935596,
3.24146865, 0.990616851, 1.002075066, 0.955325548)
identification <- nmr_identify_regions_cell (ppm_to_assign, num_proposed_compounds = 3)
```

---

nmr\_identify\_regions\_urine

*NMR peak identification (urine samples)*

---

## Description

Identify given regions and return a data frame with plausible assignments in human urine samples. The data frame contains the column "Bouatra\_2013" showing if the proposed metabolite was reported in this publication as regular urinary metabolite.

## Usage

```
nmr_identify_regions_urine(  
  ppm_to_assign,  
  num_proposed_compounds = 5,  
  verbose = FALSE  
)
```

## Arguments

**ppm\_to\_assign** A vector with the ppm regions to assign

**num\_proposed\_compounds** set the number of proposed metabolites sorted by the number times reported in the HMDB: HMDB\_urine.

**verbose** Logical value. Set it to TRUE to print additional information

## Value

a data frame with plausible assignments.

## See Also

Other peak detection functions: [Pipelines](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_identify\\_regions\\_blood\(\)](#), [nmr\\_identify\\_regions\\_cell\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [regions\\_from\\_peak\\_table\(\)](#), [validate\\_nmr\\_data](#)

Other peak integration functions: [Pipelines](#), [computes\\_peak\\_width\\_ppm\(\)](#), [nmr\\_identify\\_regions\\_blood\(\)](#), [nmr\\_identify\\_regions\\_cell\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#)

## Examples

```
# We identify regions from from the corresponding ppm stored in a vector.
ppm_to_assign <- c(4.060960203, 3.048970634, 2.405935596,
3.24146865, 0.990616851, 1.002075066, 0.955325548)
identification <- nmr_identify_regions_urine (ppm_to_assign, num_proposed_compounds = 5)
```

---

nmr\_integrate\_regions *Integrate regions*

---

## Description

Integrate given regions and return a data frame with them

## Usage

```
nmr_integrate_regions(samples, regions, ...)

## S3 method for class 'nmr_dataset_1D'
nmr_integrate_regions(
  samples,
  regions,
  fix_baseline = TRUE,
  excluded_regions_as_zero = FALSE,
  set_negative_areas_to_zero = FALSE,
  ...
)
```

## Arguments

|                            |  |
|----------------------------|--|
| samples                    | A <a href="#">nmr_dataset</a> object   |
| regions                    | A named list. Each element of the list is a region, given as a named numeric vector of length two with the range to integrate. The name of the region will be the name of the column   |
| ...                        | Keep for compatibility   |
| fix_baseline               | A logical. If TRUE it removes the baseline. See details below  |
| excluded_regions_as_zero   | A logical. It determines the behaviour of the integration when integrating regions that have been excluded. If TRUE, it will treat those regions as zero. If FALSE (the default) it will return NA values.<br>If <code>fix_baseline</code> is TRUE, then the region boundaries are used to estimate a baseline. The baseline is estimated "connecting the boundaries with a straight line". Only when the spectrum is above the baseline the area is integrated (negative contributions due to the baseline estimation are ignored). |
| set_negative_areas_to_zero | A logical. Ignored if <code>fix_baseline</code> is FALSE. When set to TRUE negative areas are set to zero.   |

## Value

An [nmr\\_dataset\\_peak\\_table](#) object

## See Also

Other peak detection functions: [Pipelines](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_identify\\_regions\\_blood\(\)](#), [nmr\\_identify\\_regions\\_cell\(\)](#), [nmr\\_identify\\_regions\\_urine\(\)](#), [regions\\_from\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#)

Other peak integration functions: `Pipelines`, `computes_peak_width_ppm()`, `nmr_identify_regions_blood()`, `nmr_identify_regions_cell()`, `nmr_identify_regions_urine()`, `validate_nmr_dataset_peak_table()`

Other `nmr_dataset_1D` functions: `[.nmr_dataset_1D()`, `computes_peak_width_ppm()`, `file_lister()`, `files_to_rDolphin()`, `format.nmr_dataset_1D()`, `is.nmr_dataset_1D()`, `load_and_save_functions`, `new_nmr_dataset_1D()`, `nmr_align_find_ref()`, `nmr_baseline_removal()`, `nmr_baseline_threshold()`, `nmr_exclude_region()`, `nmr_interpolate_1D()`, `nmr_meta_add()`, `nmr_meta_export()`, `nmr_meta_get_column()`, `nmr_meta_get()`, `nmr_normalize()`, `nmr_pca_build_model()`, `nmr_pca_outliers_filter()`, `nmr_pca_outliers_plot()`, `nmr_pca_outliers_robust()`, `nmr_pca_outliers()`, `nmr_ppm_resolution()`, `plot.nmr_dataset_1D()`, `plot_webgl()`, `print.nmr_dataset_1D()`, `rdCV_PLS_RF_ML()`, `rdCV_PLS_RF()`, `save_files_to_rDolphin()`, `to_ChemoSpec()`, `validate_nmr_dataset_peak_table()`, `validate_nmr_dataset()`

## Examples

```
#Creating a dataset
dataset <- new_nmr_dataset_1D(ppm_axis = 1:10,
  data_1r = matrix(sample(0:99,replace = TRUE), nrow = 10),
  metadata = list(external = data.frame(NMRExperiment = c("10",
    "20", "30", "40", "50", "60", "70", "80", "90", "100"))))

# Integrating selected regions
peak_table_integration = nmr_integrate_regions(
  samples = dataset,
  regions = list(ppm = c(2,5)),
  fix_baseline = TRUE)

#Creating a dataset
dataset <- new_nmr_dataset_1D(ppm_axis = 1:10,
  data_1r = matrix(sample(0:99,replace = TRUE), nrow = 10),
  metadata = list(external = data.frame(NMRExperiment = c("10",
    "20", "30", "40", "50", "60", "70", "80", "90", "100"))))

# Integrating selected regions
peak_table_integration = nmr_integrate_regions(
  samples = dataset,
  regions = list(ppm = c(2,5)),
  fix_baseline = TRUE)
```

---

`nmr_interpolate_1D`      *Interpolate a set of 1D NMR Spectra*

---

## Description

Interpolate a set of 1D NMR Spectra

## Usage

```
nmr_interpolate_1D(samples, axis = c(min = 0.2, max = 10, by = 8e-04))

## S3 method for class 'nmr_dataset'
nmr_interpolate_1D(samples, axis = c(min = 0.2, max = 10, by = 8e-04))
```

**Arguments**

|                      |  |
|----------------------|--|
| <code>samples</code> | An NMR dataset                                 |
| <code>axis</code>    | The ppm axis range and optionally the ppm step |

**Value**

Interpolate a set of 1D NMR Spectra

**See Also**

Other `nmr_dataset` functions: [[.nmr\\_dataset\(\)](#)], [format.nmr\\_dataset\(\)](#), [load\\_and\\_save\\_functions.new\\_nmr\\_dataset\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [nmr\\_read\\_samples\(\)](#), [print.nmr\\_dataset\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

Other `nmr_dataset_1D` functions: [[.nmr\\_dataset\\_1D\(\)](#)], [computes\\_peak\\_width\\_ppm\(\)](#), [file\\_lister\(\)](#), [files\\_to\\_rDolphin\(\)](#), [format.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [load\\_and\\_save\\_functions.new\\_nmr\\_dataset\\_1D\(\)](#), [nmr\\_align\\_find\\_ref\(\)](#), [nmr\\_baseline\\_removal\(\)](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_exclude\\_region\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_normalize\(\)](#), [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_outliers\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [plot.nmr\\_dataset\\_1D\(\)](#), [plot\\_webgl\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [rdCV\\_PLS\\_RF\\_ML\(\)](#), [rdCV\\_PLS\\_RF\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [to\\_ChemoSpec\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

**Examples**

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))

dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
```

---

`nmr_meta_add`

*Add metadata to an nmr\_dataset object*

---

**Description**

This is useful to add metadata to datasets that can be later used for plotting spectra or further analysis (PCA...).

**Usage**

```
nmr_meta_add(nmr_data, metadata, by = "NMRExperiment")

nmr_meta_add_tidy_excel(nmr_data, excel_file)
```

**Arguments**

|            |  |
|------------|--|
| nmr_data   | an <code>nmr_dataset_family</code> object  |
| metadata   | A data frame with metadata to add  |
| by         | A column name of both the <code>nmr_data\$metadata\$external</code> and the <code>metadata</code> data.frame. If you want to merge two columns with different headers you can use a named character vector <code>c("NMRExperiment" = "ExperimentNMR")</code> where the left side is the column name of the <code>nmr_data\$metadata\$external</code> and the right side is the column name of the <code>metadata</code> data frame.  |
| excel_file | Path to a tidy Excel file name. The Excel can consist of multiple sheets, that are added sequentially. The first column of the first sheet <b>MUST</b> be named as one of the metadata already present in the dataset, typically will be "NMRExperiment". The rest of the columns of the first sheet can be named at will. Similarly, the first column of the second sheet must be named as one of the metadata already present in the dataset, typically "NMRExperiment" or any of the columns of the first sheet. The rest of the columns of the second sheet can be named at will. See the package vignette for an example. |

**Value**

The `nmr_dataset_family` object with the added metadata

**See Also**

Other metadata functions: `Pipelines`, `nmr_meta_export()`, `nmr_meta_get_column()`, `nmr_meta_get()`

Other `nmr_dataset` functions: `[.nmr_dataset()`, `format.nmr_dataset()`, `load_and_save_functions`, `new_nmr_dataset()`, `nmr_interpolate_1D()`, `nmr_meta_export()`, `nmr_meta_get_column()`, `nmr_meta_get()`, `nmr_ppm_resolution()`, `nmr_read_samples()`, `print.nmr_dataset()`, `validate_nmr_dataset()`

Other `nmr_dataset_1D` functions: `[.nmr_dataset_1D()`, `computes_peak_width_ppm()`, `file_lister()`, `files_to_rDolphin()`, `format.nmr_dataset_1D()`, `is.nmr_dataset_1D()`, `load_and_save_functions`, `new_nmr_dataset_1D()`, `nmr_align_find_ref()`, `nmr_baseline_removal()`, `nmr_baseline_threshold()`, `nmr_exclude_region()`, `nmr_integrate_regions()`, `nmr_interpolate_1D()`, `nmr_meta_export()`, `nmr_meta_get_column()`, `nmr_meta_get()`, `nmr_normalize()`, `nmr_pca_build_model()`, `nmr_pca_outliers_filt`, `nmr_pca_outliers_plot()`, `nmr_pca_outliers_robust()`, `nmr_pca_outliers()`, `nmr_ppm_resolution()`, `plot.nmr_dataset_1D()`, `plot_webgl()`, `print.nmr_dataset_1D()`, `rdCV_PLS_RF_ML()`, `rdCV_PLS_RF()`, `save_files_to_rDolphin()`, `to_ChemoSpec()`, `validate_nmr_dataset_peak_table()`, `validate_nmr_dataset()`

Other `nmr_dataset_peak_table` functions: `[.nmr_dataset_peak_table()`, `format.nmr_dataset_peak_table()`, `is.nmr_dataset_peak_table()`, `load_and_save_functions`, `new_nmr_dataset_peak_table()`, `nmr_meta_export()`, `nmr_meta_get_column()`, `nmr_meta_get()`, `print.nmr_dataset_peak_table()`, `validate_nmr_dataset_peak_table()`

**Examples**

```
# Load a demo dataset with four samples:
dataset <- system.file("dataset-demo", package = "AlpsNMR")
nmr_dataset <- nmr_read_samples_dir(dataset)

# At first we just have the NMRExperiment column
nmr_meta_get(nmr_dataset, groups = "external")
# Get a table with NMRExperiment -> SubjectID
dummy_metadata <- system.file("dataset-demo", "dummy_metadata.xlsx", package = "AlpsNMR")
NMRExp_SubjID <- readxl::read_excel(dummy_metadata, sheet = 1)
```

```

NMRExp_SubjID
# We can link the SubjectID column of the first excel into the dataset
nmr_dataset <- nmr_meta_add(nmr_dataset, NMRExp_SubjID, by = "NMRExperiment")
nmr_meta_get(nmr_dataset, groups = "external")
# The second excel can use the SubjectID:
SubjID_Age <- readxl::read_excel(dummy_metadata, sheet = 2)
SubjID_Age
# Add the metadata by its SubjectID:
nmr_dataset <- nmr_meta_add(nmr_dataset, SubjID_Age, by = "SubjectID")
# The final loaded metadata:
nmr_meta_get(nmr_dataset, groups = "external")

# Read a tidy excel file:

dataset <- system.file("dataset-demo", package = "AlpsNMR")
nmr_dataset <- nmr_read_samples_dir(dataset)

# At first we just have the NMRExperiment column
nmr_meta_get(nmr_dataset, groups = "external")
# Get a table with NMRExperiment -> SubjectID
dummy_metadata <- system.file("dataset-demo", "dummy_metadata.xlsx", package = "AlpsNMR")

nmr_dataset <- nmr_meta_add_tidy_excel(nmr_dataset, dummy_metadata)
# Updated Metadata:
nmr_meta_get(nmr_dataset, groups = "external")

```

---

nmr\_meta\_export

*Export Metadata to an Excel file*


---

## Description

Export Metadata to an Excel file

## Usage

```

nmr_meta_export(
  nmr_dataset,
  xlsx_file,
  groups = c("info", "orig", "title", "external")
)

```

## Arguments

|             |   |
|-------------|---|
| nmr_dataset | An <a href="#">nmr_dataset_family</a> object  |
| xlsx_file   | "The .xlsx excel file"  |
| groups      | A character vector. Use "external" for the external metadata or the default for a more generic solution |

## Value

The Excel file name



**See Also**

Other metadata functions: [Pipelines](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#)

Other nmr\_dataset functions: [\[.nmr\\_dataset\(\)\]](#), [format.nmr\\_dataset\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [nmr\\_read\\_samples\(\)](#), [print.nmr\\_dataset\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

Other nmr\_dataset\_1D functions: [\[.nmr\\_dataset\\_1D\(\)\]](#), [computes\\_peak\\_width\\_ppm\(\)](#), [file\\_lister\(\)](#), [files\\_to\\_rDolphin\(\)](#), [format.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [nmr\\_align\\_find\\_ref\(\)](#), [nmr\\_baseline\\_removal\(\)](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_exclude\\_region\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_normalize\(\)](#), [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_filt](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_outliers\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [plot.nmr\\_dataset\\_1D\(\)](#), [plot\\_webgl\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [rdCV\\_PLS\\_RF\\_ML\(\)](#), [rdCV\\_PLS\\_RF\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [to\\_ChemoSpec\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

Other nmr\_dataset\_peak\_table functions: [\[.nmr\\_dataset\\_peak\\_table\(\)\]](#), [format.nmr\\_dataset\\_peak\\_table\(\)](#), [is.nmr\\_dataset\\_peak\\_table\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_peak\\_table\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [print.nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#)

Other import/export functions: [Pipelines](#), [files\\_to\\_rDolphin\(\)](#), [load\\_and\\_save\\_functions](#), [nmr\\_data\(\)](#), [nmr\\_read\\_bruker\\_fid\(\)](#), [nmr\\_read\\_samples\(\)](#), [nmr\\_zip\\_bruker\\_samples\(\)](#), [save\\_files\\_to\\_rDolp](#), [save\\_profiling\\_output\(\)](#), [to\\_ChemoSpec\(\)](#)

**Examples**

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
#nmr_meta_export(dataset, "metadata.xlsx")
```

nmr\_meta\_get

*Get metadata***Description**

Get metadata

**Usage**

```
nmr_meta_get(samples, columns = NULL, groups = NULL)
```

**Arguments**

|         |  |
|---------|--|
| samples | a <a href="#">nmr_dataset_family</a> object  |
| columns | Columns to get. By default gets all the columns.   |
| groups  | Groups to get. Groups are predefined of columns. Typically "external" for metadata added with <a href="#">nmr_meta_add</a> . |
|         | Both groups and columns can't be given simultaneously.   |

**Value**

a data frame with the injection metadata

**See Also**

Other metadata functions: [Pipelines](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#)

Other nmr\_dataset functions: [\[.nmr\\_dataset\(\)\]](#), [format.nmr\\_dataset\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [nmr\\_read\\_samples\(\)](#), [print.nmr\\_dataset\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

Other nmr\_dataset\_1D functions: [\[.nmr\\_dataset\\_1D\(\)\]](#), [computes\\_peak\\_width\\_ppm\(\)](#), [file\\_lister\(\)](#), [files\\_to\\_rDolphin\(\)](#), [format.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [nmr\\_align\\_find\\_ref\(\)](#), [nmr\\_baseline\\_removal\(\)](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_exclude\\_region\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_normalize\(\)](#), [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_f](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_outliers\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [plot.nmr\\_dataset\\_1D\(\)](#), [plot\\_webgl\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [rdCV\\_PLS\\_RF\\_ML\(\)](#), [rdCV\\_PLS\\_RF\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [to\\_ChemoSpec\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

Other nmr\_dataset\_peak\_table functions: [\[.nmr\\_dataset\\_peak\\_table\(\)\]](#), [format.nmr\\_dataset\\_peak\\_table\(\)](#), [is.nmr\\_dataset\\_peak\\_table\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_peak\\_table\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [print.nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#)

**Examples**

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
metadata <- nmr_meta_get(dataset)
```

---

nmr\_meta\_get\_column     *Get a single metadata column*

---

**Description**

Get a single metadata column

**Usage**

```
nmr_meta_get_column(samples, column = "NMRExperiment")
```

**Arguments**

`samples`            a [nmr\\_dataset\\_family](#) object  
`column`             A column to get

**Value**

A vector with the column

**See Also**

Other metadata functions: [Pipelines](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\(\)](#)

Other `nmr_dataset` functions: [\[.nmr\\_dataset\(\)\]](#), [format.nmr\\_dataset\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [nmr\\_read\\_samples\(\)](#), [print.nmr\\_dataset\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

Other `nmr_dataset_1D` functions: [\[.nmr\\_dataset\\_1D\(\)\]](#), [computes\\_peak\\_width\\_ppm\(\)](#), [file\\_lister\(\)](#), [files\\_to\\_rDolphin\(\)](#), [format.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [nmr\\_align\\_find\\_ref\(\)](#), [nmr\\_baseline\\_removal\(\)](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_exclude\\_region\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_normalize\(\)](#), [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_outliers\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [plot.nmr\\_dataset\\_1D\(\)](#), [plot\\_webgl\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [rdCV\\_PLS\\_RF\\_ML\(\)](#), [rdCV\\_PLS\\_RF\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [to\\_ChemoSpec\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

Other `nmr_dataset_peak_table` functions: [\[.nmr\\_dataset\\_peak\\_table\(\)\]](#), [format.nmr\\_dataset\\_peak\\_table\(\)](#), [is.nmr\\_dataset\\_peak\\_table\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_peak\\_table\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\(\)](#), [print.nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#)

**Examples**

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
metadata_column <- nmr_meta_get_column(dataset)
```

---

|               |   |
|---------------|---|
| nmr_normalize | <i>Normalize nmr_dataset_1D samples</i> |
|---------------|---|

---

**Description**

The `nmr_normalize` function is used to normalize all the samples according to a given criteria.

**Usage**

```
nmr_normalize(
  samples,
  method = c("area", "max", "value", "region", "pqn", "none"),
  ...
)

nmr_normalize_extra_info(samples)
```

**Arguments**

|         |  |
|---------|--|
| samples | A <a href="#">nmr_dataset_1D</a> object  |
| method  | The criteria to be used for normalization - area: Normalize to the total area - max: Normalize to the maximum intensity - value: Normalize each sample to a user defined value - region: Integrate a region and normalize each sample to that region - pqn: Use Probabalistic Quotient Normalization for normalization - none: Do not normalize at all |

... Method dependent arguments: - method == "value": - value: A numeric vector with the normalization values to use - method == "region": - ppm\_range: A chemical shift region to integrate - ...: Other arguments passed on to [nmr\\_integrate\\_regions](#)

### Details

The `nmr_normalize_extra_info` function is used to extract additional information after the normalization. Typically, we want to know what was the actual normalization factor applied to each sample. The extra information includes a plot, representing the dispersion of the normalization factor for each sample.

### Value

The `nmr_dataset_1D` object, with the samples normalized. Further information for diagnostic of the normalization process is also saved and can be extracted by calling `nmr_normalize_extra_info()` afterwards.

### See Also

Other `nmr_dataset_1D` functions: [[.nmr\\_dataset\\_1D\(\)](#), [computes\\_peak\\_width\\_ppm\(\)](#), [file\\_lister\(\)](#), [files\\_to\\_rDolphin\(\)](#), [format.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [nmr\\_align\\_find\\_ref\(\)](#), [nmr\\_baseline\\_removal\(\)](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_exclude\\_region\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_fit](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_outliers\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [plot.nmr\\_dataset\\_1D\(\)](#), [plot\\_webgl\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [rdCV\\_PLS\\_RF\\_ML\(\)](#), [rdCV\\_PLS\\_RF\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [to\\_ChemoSpec\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)]

### Examples

```
nmr_dataset <- nmr_dataset_load(system.file("extdata", "nmr_dataset.rds", package = "AlpsNMR"))
nmr_dataset <- nmr_normalize(nmr_dataset, method = "area")
norm_dataset <- nmr_normalize(nmr_dataset)
norm_dataset$plot
nmr_dataset <- nmr_dataset_load(system.file("extdata", "nmr_dataset.rds", package = "AlpsNMR"))
nmr_dataset <- nmr_normalize(nmr_dataset, method = "area")
norm_extra_info <- nmr_normalize_extra_info(nmr_dataset)
norm_extra_info$plot
```

---

`nmr_pca_build_model`     *Build a PCA on for an nmr\_dataset*

---

### Description

This function builds a PCA model with all the NMR spectra. Regions with zero values (excluded regions) or near-zero variance regions are automatically excluded from the analysis.

**Usage**

```

nmr_pca_build_model(
  nmr_dataset,
  ncomp = NULL,
  center = TRUE,
  scale = FALSE,
  ...
)

## S3 method for class 'nmr_dataset_1D'
nmr_pca_build_model(
  nmr_dataset,
  ncomp = NULL,
  center = TRUE,
  scale = FALSE,
  ...
)

```

**Arguments**

|                          |   |
|--------------------------|---|
| <code>nmr_dataset</code> | a <a href="#">nmr_dataset_1D</a> object   |
| <code>ncomp</code>       | Integer, if data is complete <code>ncomp</code> decides the number of components and associated eigenvalues to display from the <code>pcasvd</code> algorithm and if the data has missing values, <code>ncomp</code> gives the number of components to keep to perform the reconstitution of the data using the NIPALS algorithm. If <code>NULL</code> , function sets <code>ncomp = min(nrow(X), ncol(X))</code>         |
| <code>center</code>      | (Default= <code>TRUE</code> ) Logical, whether the variables should be shifted to be zero centered. Alternatively, a vector of length equal the number of columns of <code>X</code> can be supplied. The value is passed to <a href="#">scale</a> .   |
| <code>scale</code>       | (Default= <code>FALSE</code> ) Logical indicating whether the variables should be scaled to have unit variance before the analysis takes place. The default is <code>FALSE</code> for consistency with <code>prcomp</code> function, but in general scaling is advisable. Alternatively, a vector of length equal the number of columns of <code>X</code> can be supplied. The value is passed to <a href="#">scale</a> . |
| <code>...</code>         | Additional arguments passed on to <a href="#">mixOmics::pca</a>   |

**Value**

A PCA model as given by [mixOmics::pca](#) with two additional attributes:

- `nmr_data_axis` containing the full ppm axis
- `nmr_included` with the data points included in the model These attributes are used internally by `AlpsNMR` to create loading plots

**See Also**

Other PCA related functions: [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_outliers\(\)](#), [nmr\\_pca\\_plots](#)

Other `nmr_dataset_1D` functions: [\[.nmr\\_dataset\\_1D\(\)](#), [computes\\_peak\\_width\\_ppm\(\)](#), [file\\_lister\(\)](#), [files\\_to\\_rDolphin\(\)](#), [format.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [nmr\\_align\\_find\\_ref\(\)](#), [nmr\\_baseline\\_removal\(\)](#), [nmr\\_baseline\\_threshold\(\)](#),

```
nmr_exclude_region(), nmr_integrate_regions(), nmr_interpolate_1D(), nmr_meta_add(),
nmr_meta_export(), nmr_meta_get_column(), nmr_meta_get(), nmr_normalize(), nmr_pca_outliers_filter(),
nmr_pca_outliers_plot(), nmr_pca_outliers_robust(), nmr_pca_outliers(), nmr_ppm_resolution(),
plot.nmr_dataset_1D(), plot_webgl(), print.nmr_dataset_1D(), rdCV_PLS_RF_ML(), rdCV_PLS_RF(),
save_files_to_rDolphin(), to_ChemoSpec(), validate_nmr_dataset_peak_table(), validate_nmr_dataset()
```

## Examples

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
model <- nmr_pca_build_model(dataset_1D)
```

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
model <- nmr_pca_build_model(dataset_1D)
```

---

|                  |   |
|------------------|---|
| nmr_pca_outliers | <i>Compute PCA residuals and score distance for each sample</i> |
|------------------|---|

---

## Description

Compute PCA residuals and score distance for each sample

## Usage

```
nmr_pca_outliers(
  nmr_dataset,
  pca_model,
  ncomp = NULL,
  quantile_critical = 0.975
)
```

## Arguments

|                   |   |
|-------------------|---|
| nmr_dataset       | An <a href="#">nmr_dataset_1D</a> object                                    |
| pca_model         | A <a href="#">pca</a> model returned by <a href="#">nmr_pca_build_model</a> |
| ncomp             | Number of components to use. Use NULL for 90% of the variance               |
| quantile_critical | critical quantile   |

## Value

A list with:

- outlier\_info: A data frame with the NMRExperiment, the Q residuals and T scores
- ncomp: Number of components used to compute Q and T
- Tscore\_critical, QResidual\_critical: Critical values, given a quantile, for both Q and T.

## See Also

Other PCA related functions: [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_plots](#)

Other outlier detection functions: [Pipelines](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#)

Other `nmr_dataset_1D` functions: [\[.nmr\\_dataset\\_1D\(\)\]](#), [computes\\_peak\\_width\\_ppm\(\)](#), [file\\_lister\(\)](#), [files\\_to\\_rDolphin\(\)](#), [format.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [nmr\\_align\\_find\\_ref\(\)](#), [nmr\\_baseline\\_removal\(\)](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_exclude\\_region\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_normalize\(\)](#), [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [plot.nmr\\_dataset\\_1D\(\)](#), [plot\\_webgl\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [rdCV\\_PLS\\_RF\\_ML\(\)](#), [rdCV\\_PLS\\_RF\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [to\\_ChemoSpec\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

## Examples

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
model <- nmr_pca_build_model(dataset_1D)
outliers_info <- nmr_pca_outliers(dataset_1D, model)
```

---

`nmr_pca_outliers_filter`

*Exclude outliers*

---

## Description

Exclude outliers

## Usage

```
nmr_pca_outliers_filter(nmr_dataset, pca_outliers)
```

## Arguments

`nmr_dataset` An [nmr\\_dataset\\_1D](#) object

`pca_outliers` The output from [nmr\\_pca\\_outliers\(\)](#)

## Value

An [nmr\\_dataset\\_1D](#) without the detected outliers

**See Also**

Other PCA related functions: `nmr_pca_build_model()`, `nmr_pca_outliers_plot()`, `nmr_pca_outliers_robust()`, `nmr_pca_outliers()`, `nmr_pca_plots`

Other outlier detection functions: `Pipelines`, `nmr_pca_outliers_plot()`, `nmr_pca_outliers_robust()`, `nmr_pca_outliers()`

Other `nmr_dataset_1D` functions: `[.nmr_dataset_1D()]`, `computes_peak_width_ppm()`, `file_lister()`, `files_to_rDolphin()`, `format.nmr_dataset_1D()`, `is.nmr_dataset_1D()`, `load_and_save_functions`, `new_nmr_dataset_1D()`, `nmr_align_find_ref()`, `nmr_baseline_removal()`, `nmr_baseline_threshold()`, `nmr_exclude_region()`, `nmr_integrate_regions()`, `nmr_interpolate_1D()`, `nmr_meta_add()`, `nmr_meta_export()`, `nmr_meta_get_column()`, `nmr_meta_get()`, `nmr_normalize()`, `nmr_pca_build_model()`, `nmr_pca_outliers_plot()`, `nmr_pca_outliers_robust()`, `nmr_pca_outliers()`, `nmr_ppm_resolution()`, `plot.nmr_dataset_1D()`, `plot_webgl()`, `print.nmr_dataset_1D()`, `rdCV_PLS_RF_ML()`, `rdCV_PLS_RF()`, `save_files_to_rDolphin()`, `to_ChemoSpec()`, `validate_nmr_dataset_peak_table()`, `validate_nmr_dataset()`

Other subsetting functions: `[.nmr_dataset_1D()]`, `[.nmr_dataset_peak_table()]`, `[.nmr_dataset()]`, `filter.nmr_dataset_family()`

**Examples**

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
model <- nmr_pca_build_model(dataset_1D)
outliers_info <- nmr_pca_outliers(dataset_1D, model)
dataset_whitout_outliers <- nmr_pca_outliers_filter(dataset_1D, outliers_info)
```

---

`nmr_pca_outliers_plot` *Plot for outlier detection diagnostic*

---

**Description**

Plot for outlier detection diagnostic

**Usage**

```
nmr_pca_outliers_plot(nmr_dataset, pca_outliers, ...)
```

**Arguments**

|                           |   |
|---------------------------|---|
| <code>nmr_dataset</code>  | An <code>nmr_dataset_1D</code> object                                 |
| <code>pca_outliers</code> | The output from <code>nmr_pca_outliers()</code>                       |
| <code>...</code>          | Additional parameters passed on to <code>ggplot2::aes_string()</code> |

**Value**

A plot for the outlier detection



**See Also**

Other PCA related functions: [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_outliers\(\)](#), [nmr\\_pca\\_plots](#)

Other outlier detection functions: [Pipelines](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_outliers\(\)](#)

Other `nmr_dataset_1D` functions: [\[.nmr\\_dataset\\_1D\(\)](#), [computes\\_peak\\_width\\_ppm\(\)](#), [file\\_lister\(\)](#), [files\\_to\\_rDolphin\(\)](#), [format.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [nmr\\_align\\_find\\_ref\(\)](#), [nmr\\_baseline\\_removal\(\)](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_exclude\\_region\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_normalize\(\)](#), [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_outliers\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [plot.nmr\\_dataset\\_1D\(\)](#), [plot\\_webgl\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [rdCV\\_PLS\\_RF\\_ML\(\)](#), [rdCV\\_PLS\\_RF\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [to\\_ChemoSpec\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

**Examples**

```
#dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
#dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
#dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
#model <- nmr_pca_build_model(dataset_1D)
#outliers_info <- nmr_pca_outliers(dataset_1D, model)
#nmr_pca_outliers_plot(dataset_1D, outliers_info)
```

---

nmr\_pca\_outliers\_robust

*Outlier detection through robust PCA*

---

**Description**

Outlier detection through robust PCA

**Usage**

```
nmr_pca_outliers_robust(nmr_dataset, ncomp = 5)
```

**Arguments**

`nmr_dataset` An `nmr_dataset_1D` object

`ncomp` Number of rPCA components to use

We have observed that the statistical test used as a threshold for outlier detection usually flags as outliers too many samples, due possibly to a lack of gaussianity. As a workaround, a heuristic method has been implemented: We know that in the Q residuals vs T scores plot from [nmr\\_pca\\_outliers\\_plot\(\)](#) outliers are on the right or on the top of the plot, and quite separated from non-outlier samples.

To determine the critical value, both for Q and T, we find the biggest gap between samples in the plot and use as critical value the center of the gap.

This approach seems to work well when there are outliers, but it fails when there isn't any outlier. For that case, the gap would be placed anywhere and that is not

desirable as many samples would be incorrectly flagged. The second assumption that we use is that no more than 10% the samples may pass our critical value. If more than 10% pass the critical value, then we assume that our heuristics are not reasonable and we don't set any critical limit.

### Value

A list similar to [nmr\\_pca\\_outliers](#)

### See Also

Other PCA related functions: [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\(\)](#), [nmr\\_pca\\_plots](#)

Other outlier detection functions: [Pipelines](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\(\)](#)

Other `nmr_dataset_1D` functions: [\[.nmr\\_dataset\\_1D\(\)\]](#), [computes\\_peak\\_width\\_ppm\(\)](#), [file\\_lister\(\)](#), [files\\_to\\_rDolphin\(\)](#), [format.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [nmr\\_align\\_find\\_ref\(\)](#), [nmr\\_baseline\\_removal\(\)](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_exclude\\_region\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_normalize\(\)](#), [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [plot.nmr\\_dataset\\_1D\(\)](#), [plot\\_webgl\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [rdCV\\_PLS\\_RF\\_ML\(\)](#), [rdCV\\_PLS\\_RF\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [to\\_ChemoSpec\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

### Examples

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
outliers_info <- nmr_pca_outliers_robust(dataset_1D)
```

---

nmr\_pca\_plots

*Plotting functions for PCA*

---

### Description

Plotting functions for PCA

### Usage

```
nmr_pca_plot_variance(pca_model)
```

```
nmr_pca_scoreplot(nmr_dataset, pca_model, comp = seq_len(2), ...)
```

```
nmr_pca_loadingplot(pca_model, comp)
```

### Arguments

|                          |   |
|--------------------------|---|
| <code>pca_model</code>   | A PCA model trained with <a href="#">nmr_pca_build_model</a>                              |
| <code>nmr_dataset</code> | an <a href="#">nmr_dataset_1D</a> object  |
| <code>comp</code>        | Components to represent   |
| <code>...</code>         | Additional aesthetics passed on to <a href="#">ggplot2::aes</a> (use bare unquoted names) |

**Value**

Plot of PCA

**See Also**

Other PCA related functions: [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_outliers\(\)](#)

**Examples**

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
model <- nmr_pca_build_model(dataset_1D)
nmr_pca_plot_variance(model)
```

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
model <- nmr_pca_build_model(dataset_1D)
nmr_pca_scoreplot(dataset_1D, model)
```

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
model <- nmr_pca_build_model(dataset_1D)
nmr_pca_loadingplot(model, 1)
```

---

nmr\_ppm\_resolution      *PPM resolution of the spectra*

---

**Description**

The function gets the ppm resolution of the dataset using the median of the difference of data points.

**Usage**

```
nmr_ppm_resolution(nmr_dataset)

## S3 method for class 'nmr_dataset'
nmr_ppm_resolution(nmr_dataset)

## S3 method for class 'nmr_dataset_1D'
nmr_ppm_resolution(nmr_dataset)
```

**Arguments**

nmr\_dataset      An object containing NMR samples

**Value**

Numeric (the ppm resolution, measured in ppms)

**See Also**

Other `nmr_dataset` functions: [[.nmr\\_dataset\(\)](#), [format.nmr\\_dataset\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_read\\_samples\(\)](#), [print.nmr\\_dataset\(\)](#), [validate\\_nmr\\_dataset\(\)](#)]

Other `nmr_dataset_1D` functions: [[.nmr\\_dataset\\_1D\(\)](#), [computes\\_peak\\_width\\_ppm\(\)](#), [file\\_lister\(\)](#), [files\\_to\\_rDolphin\(\)](#), [format.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [nmr\\_align\\_find\\_ref\(\)](#), [nmr\\_baseline\\_removal\(\)](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_exclude\\_region\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_normalize\(\)](#), [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_outliers\(\)](#), [plot.nmr\\_dataset\\_1D\(\)](#), [plot\\_webgl\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [rdCV\\_PLS\\_RF\\_ML\(\)](#), [rdCV\\_PLS\\_RF\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [to\\_ChemoSpec\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)]

**Examples**

```
nmr_dataset <- nmr_dataset_load(system.file("extdata", "nmr_dataset.rds", package = "AlpsNMR"))
nmr_ppm_resolution(nmr_dataset)
message("the ppm resolution of this dataset is ", nmr_ppm_resolution(nmr_dataset), " ppm")

nmr_dataset <- nmr_dataset_load(system.file("extdata", "nmr_dataset.rds", package = "AlpsNMR"))
nmr_ppm_resolution(nmr_dataset)
message("the ppm resolution of this dataset is ", nmr_ppm_resolution(nmr_dataset), " ppm")

nmr_dataset <- nmr_dataset_load(system.file("extdata", "nmr_dataset.rds", package = "AlpsNMR"))
nmr_ppm_resolution(nmr_dataset)
message("the ppm resolution of this dataset is ", nmr_ppm_resolution(nmr_dataset), " ppm")
```

---

`nmr_read_bruker_fid`     *Read Free Induction Decay file*

---

**Description**

Reads an FID file. This is a very simple function.

**Usage**

```
nmr_read_bruker_fid(sample_name, endian = "little")
```

**Arguments**

|                          |  |
|--------------------------|--|
| <code>sample_name</code> | A single sample name   |
| <code>endian</code>      | Endianness of the fid file ("little" by default, use "big" if <code>acqus\$BYTORDA == 1</code> ) |

**Value**

A numeric vector with the free induction decay values

**See Also**

Other import/export functions: [Pipelines](#), [files\\_to\\_rDolphin\(\)](#), [load\\_and\\_save\\_functions](#), [nmr\\_data\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_read\\_samples\(\)](#), [nmr\\_zip\\_bruker\\_samples\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [save\\_profiling\\_output\(\)](#), [to\\_ChemoSpec\(\)](#)

**Examples**

```
fid <- nmr_read_bruker_fid("sample.fid")
```

---

|                  |                         |
|------------------|-------------------------|
| nmr_read_samples | <i>Read NMR samples</i> |
|------------------|-------------------------|

---

**Description**

These functions load samples from files and return a [nmr\\_dataset](#).

**Usage**

```
nmr_read_samples_dir(
  samples_dir,
  format = "bruker",
  pulse_sequence = NULL,
  metadata_only = FALSE,
  ...
)

nmr_read_samples(
  sample_names,
  format = "bruker",
  pulse_sequence = NULL,
  metadata_only = FALSE,
  ...
)
```

**Arguments**

|                             |   |
|-----------------------------|---|
| <code>samples_dir</code>    | A directory that contains multiple samples  |
| <code>format</code>         | Either "bruker" or "jdx"  |
| <code>pulse_sequence</code> | If it is set to a pulse sequence ("NOESY", "JRES", "CPMG"...) it will only load the samples that match that pulse sequence. |
| <code>metadata_only</code>  | A logical, to load only metadata (default: FALSE)   |
| <code>...</code>            | Arguments passed to <a href="#">read_bruker_sample()</a> for data loading   |
| <code>sample_names</code>   | A character vector with file or directory names.  |

**Value**

a [nmr\\_dataset](#) object

**See Also**

Other [nmr\\_dataset](#) functions: [.nmr\\_dataset\(\)](#), [format.nmr\\_dataset\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [print.nmr\\_dataset\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

Other import/export functions: [Pipelines](#), [files\\_to\\_rDolphin\(\)](#), [load\\_and\\_save\\_functions](#), [nmr\\_data\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_read\\_bruker\\_fid\(\)](#), [nmr\\_zip\\_bruker\\_samples\(\)](#), [save\\_files\\_to\\_rDolphin](#), [save\\_profiling\\_output\(\)](#), [to\\_ChemoSpec\(\)](#)

## Examples

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
```

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
zip_files <- fs::dir_ls(dir_to_demo_dataset, glob = "*.zip")
dataset <- nmr_read_samples(sample_names = zip_files)
```

---

nmr\_zip\_bruker\_samples

*Create one zip file for each brucker sample path*

---

## Description

Create one zip file for each brucker sample path

## Usage

```
nmr_zip_bruker_samples(path, workdir, overwrite = FALSE, ...)
```

## Arguments

|           |   |
|-----------|---|
| path      | Character vector with sample directories  |
| workdir   | Directory to store zip files              |
| overwrite | Should existing zip files be overwritten? |
| ...       | Passed to <a href="#">utils::zip</a>      |

## Value

A character vector of the same length as path, with the zip file names

## See Also

Other import/export functions: [Pipelines](#), [files\\_to\\_rDolphin\(\)](#), [load\\_and\\_save\\_functions](#), [nmr\\_data\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_read\\_bruker\\_fid\(\)](#), [nmr\\_read\\_samples\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [save\\_profiling\\_output\(\)](#), [to\\_ChemoSpec\(\)](#)

## Examples

```
outpaths <- nmr_zip_bruker_samples(".", getwd())
```

---

Parameters\_blood      *to rDolphin*

---

**Description**

Parameters for blood (plasma/serum) samples profiling

**Details**

The template Parameters\_blood contains the chosen normalization approach (by default, PQN), the Spectrometer Frequency (by default, 600.04MHz), alignment (by default, TSP 0.00 ppm), bucket resolution (by default, 0.00023)

**References**

[github.com/danielcanueto/rDolphin](https://github.com/danielcanueto/rDolphin)

---

Parameters\_cell      *Parameters for cell samples profiling*

---

**Description**

The template Parameters\_cell contains the chosen normalization approach (by default, PQN), the Spectrometer Frequency (by default, 600.04MHz), alignment (by default, TSP 0.00 ppm), bucket resolution (by default, 0.00023)

**References**

[github.com/danielcanueto/rDolphin](https://github.com/danielcanueto/rDolphin)

---

Parameters\_urine      *Parameters for urine samples profiling*

---

**Description**

The template Parameters\_urine contains the chosen normalization approach (by default, PQN), the Spectrometer Frequency (by default, 600.04MHz), alignment (by default, TSP 0.00 ppm), bucket resolution (by default, 0.00023)

**References**

[github.com/danielcanueto/rDolphin](https://github.com/danielcanueto/rDolphin)

---

`permutation_test_model`*Permutation test*

---

### Description

Make permutations with data and default settings from an `nmr_data_analysis_method`

### Usage

```
permutation_test_model(  
  dataset,  
  y_column,  
  identity_column,  
  external_val,  
  internal_val,  
  data_analysis_method,  
  nPerm = 50  
)
```

### Arguments

|                                   |   |
|-----------------------------------|---|
| <code>dataset</code>              | An <a href="#">nmr_dataset_family</a> object  |
| <code>y_column</code>             | A string with the name of the y column (present in the metadata of the dataset)   |
| <code>identity_column</code>      | NULL or a string with the name of the identity column (present in the metadata of the dataset).   |
| <code>external_val</code>         | A list with two elements: <code>iterations</code> and <code>test_size</code> . See <a href="#">random_subsampling</a> for further details |
| <code>internal_val</code>         | A list with two elements: <code>iterations</code> and <code>test_size</code> . See <a href="#">random_subsampling</a> for further details |
| <code>data_analysis_method</code> | An <a href="#">nmr_data_analysis_method</a> object  |
| <code>nPerm</code>                | number of permutations  |

### Value

A permutation matrix with permuted values

### Examples

```
# Data analysis for a table of integrated peaks  
  
## Generate an artificial nmr_dataset_peak_table:  
### Generate artificial metadata:  
num_samples <- 32 # use an even number in this example  
num_peaks <- 20  
metadata <- data.frame(  
  NMRExperiment = as.character(1:num_samples),  
  Condition = rep(c("A", "B"), times = num_samples/2),  
  stringsAsFactors = FALSE
```



```

)

### The matrix with peaks
peak_means <- runif(n = num_peaks, min = 300, max = 600)
peak_sd <- runif(n = num_peaks, min = 30, max = 60)
peak_matrix <- mapply(function(mu, sd) rnorm(num_samples, mu, sd),
                      mu = peak_means, sd = peak_sd)

colnames(peak_matrix) <- paste0("Peak", 1:num_peaks)

## Artificial differences depending on the condition:
peak_matrix[metadata$Condition == "A", "Peak2"] <-
  peak_matrix[metadata$Condition == "A", "Peak2"] + 70

peak_matrix[metadata$Condition == "A", "Peak6"] <-
  peak_matrix[metadata$Condition == "A", "Peak6"] - 60

### The nmr_dataset_peak_table
peak_table <- new_nmr_dataset_peak_table(
  peak_table = peak_matrix,
  metadata = list(external = metadata)
)

methodology <- plsda_auroc_vip_method(ncomp = 3)
model <- nmr_data_analysis(
  peak_table,
  y_column = "Condition",
  identity_column = NULL,
  external_val = list(iterations = 3, test_size = 0.25),
  internal_val = list(iterations = 3, test_size = 0.25),
  data_analysis_method = methodology
)

p = permutation_test_model(peak_table,
  y_column = "Condition",
  identity_column = NULL,
  external_val = list(iterations = 3, test_size = 0.25),
  internal_val = list(iterations = 3, test_size = 0.25),
  data_analysis_method = methodology,
  nPerm = 10)

```

---

permutation\_test\_plot *Permutation test plot*

---

## Description

Plot permutation test using actual model and permuted models

## Usage

```

permutation_test_plot(
  nmr_data_analysis_model,
  permMatrix,
  xlab = "AUCs",

```

```

    xlim,
    ylim = NULL,
    breaks = "Sturges",
    main = "Permutation test"
  )

```

### Arguments

|                                      |  |
|--------------------------------------|--|
| <code>nmr_data_analysis_model</code> | A <code>nmr_data_analysis_model</code>                                 |
| <code>permMatrix</code>              | A permutation fitness outcome from <code>permutation_test_model</code> |
| <code>xlab</code>                    | optional xlabel  |
| <code>xlim</code>                    | optional x-range   |
| <code>ylim</code>                    | optional y-range   |
| <code>breaks</code>                  | optional custom histogram breaks (defaults to 'sturges')               |
| <code>main</code>                    | optional plot title (or TRUE for autoname)                             |

### Value

A plot with the comparison between the actual model versus the permuted models

### Examples

```

# Data analysis for a table of integrated peaks

## Generate an artificial nmr_dataset_peak_table:
### Generate artificial metadata:
num_samples <- 32 # use an even number in this example
num_peaks <- 20
metadata <- data.frame(
  NMRExperiment = as.character(1:num_samples),
  Condition = rep(c("A", "B"), times = num_samples/2),
  stringsAsFactors = FALSE
)

### The matrix with peaks
peak_means <- runif(n = num_peaks, min = 300, max = 600)
peak_sd <- runif(n = num_peaks, min = 30, max = 60)
peak_matrix <- mapply(function(mu, sd) rnorm(num_samples, mu, sd),
                      mu = peak_means, sd = peak_sd)
colnames(peak_matrix) <- paste0("Peak", 1:num_peaks)

## Artificial differences depending on the condition:
peak_matrix[metadata$Condition == "A", "Peak2"] <-
  peak_matrix[metadata$Condition == "A", "Peak2"] + 70

peak_matrix[metadata$Condition == "A", "Peak6"] <-
  peak_matrix[metadata$Condition == "A", "Peak6"] - 60

### The nmr_dataset_peak_table
peak_table <- new_nmr_dataset_peak_table(
  peak_table = peak_matrix,
  metadata = list(external = metadata)
)

```

```

methodology <- plsda_auroc_vip_method(ncomp = 3)
model <- nmr_data_analysis(
  peak_table,
  y_column = "Condition",
  identity_column = NULL,
  external_val = list(iterations = 3, test_size = 0.25),
  internal_val = list(iterations = 3, test_size = 0.25),
  data_analysis_method = methodology
)

p = permutation_test_model(peak_table,
  y_column = "Condition",
  identity_column = NULL,
  external_val = list(iterations = 3, test_size = 0.25),
  internal_val = list(iterations = 3, test_size = 0.25),
  data_analysis_method = methodology,
  nPerm = 10)

permutation_test_plot(model, p)

```

---

Pipelines

*Pipelines*


---

### Description

Uses [nmr\\_pca\\_outliers\\_robust](#) to perform the detection of outliers

Normalize the full spectra to the internal calibrant region, then exclude that region and finally perform PQN normalization.

### Usage

```
pipe_load_samples(samples_dir, glob = "*0", output_dir = NULL)
```

```
pipe_add_metadata(nmr_dataset_rds, excel_file, output_dir)
```

```
pipe_interpolate_1D(nmr_dataset_rds, axis, output_dir)
```

```
pipe_exclude_regions(nmr_dataset_rds, exclude, output_dir)
```

```
pipe_outlier_detection(nmr_dataset_rds, output_dir)
```

```
pipe_filter_samples(nmr_dataset_rds, conditions, output_dir)
```

```

pipe_peakdet_align(
  nmr_dataset_rds,
  nDivRange_ppm = 0.1,
  scales = seq(1, 16, 2),
  baselineThresh = 0.01,
  SNR.Th = -1,
  maxShift_ppm = 0.0015,

```

```

    acceptLostPeak = FALSE,
    output_dir = NULL
  )

  pipe_peak_integration(
    nmr_dataset_rds,
    peak_det_align_dir,
    peak_width_ppm,
    output_dir
  )

  pipe_normalization(
    nmr_dataset_rds,
    internal_calibrant = NULL,
    output_dir = NULL
  )

```

### Arguments

|                              |   |
|------------------------------|---|
| <code>samples_dir</code>     | The directory where the samples are   |
| <code>glob</code>            | A wildcard aka globbing pattern (e.g. *.csv) passed on to <code>grep()</code> to filter paths.  |
| <code>output_dir</code>      | The output directory for this pipe element  |
| <code>nmr_dataset_rds</code> | The <code>nmr_dataset.rds</code> file name coming from previous nodes   |
| <code>excel_file</code>      | <p>An excel file name. See details for the requirements</p> <p>The excel file can have one or more sheets. The excel sheets need to be as simple as possible: One header column on the first row and values below.</p> <p>Each of the sheets contain metadata that has to be integrated. The merge (technically a left join) is done using the first column of each sheet as key.</p> <p>In practical terms this means that the first sheet of the excel file MUST start with an "NMRExperiment" column, and as many additional columns to add (e.g. FluidXBarcode, SampleCollectionDate, TimePoint and SubjectID).</p> <p>The second sheet can have as the first column any of the already added columns, for instance the "SubjectID", and any additional columns (e.g. Gender, Age).</p> <p>The first column on each sheet, named the key column, MUST have unique values. For instance, a sheet starting with "SubjectID" MUST specify each subject ID only once (without repetitions).</p> |
| <code>axis</code>            | The ppm axis range and optionally the ppm step  |
| <code>exclude</code>         | A list with regions to be removed Typically: <code>exclude = list(water = c(4.7, 5.0))</code>   |
| <code>conditions</code>      | <p>A character vector with conditions to filter metadata. The <code>conditions</code> parameter should be a character vector of valid R logical conditions. Some examples:</p> <ul style="list-style-type: none"> <li><code>conditions &lt;- 'Gender == "Female"'</code></li> <li><code>conditions &lt;- 'Cohort == "Chuv"'</code></li> <li><code>conditions &lt;- 'TimePoint %in% c("T0", "T31")'</code></li> <li><code>conditions &lt;- c(Cohort == "Chuv", 'TimePoint %in% c("T0", "T31")')</code></li> </ul> <p>Only samples fulfilling all the given conditions are kept in further analysis.</p>  |
| <code>nDivRange_ppm</code>   | Segment size, in ppms, to divide the spectra and search for peaks.  |
| <code>scales</code>          | The parameter of <code>peakDetectionCWT</code> function of <code>MassSpecWavelet</code> package, look it up in the original function.   |

|                    |  |
|--------------------|--|
| baselineThresh     | It will remove all peaks under an intensity set by baselineThresh. If you set it to 'NULL', nmr_detect_peaks will automatically compute an approximate value considering baseline between 9.5 and 10.0 ppm (automatically calculation using baselineThresh = NULL will not work if spectra were not interpolated up to 10.0 ppm) |
| SNR.Th             | The parameter of peakDetectionCWT function of MassSpecWavelet package, look it up in the original function. If you set -1, the function will itself recompute this value.  |
| maxShift_ppm       | The maximum shift allowed, in ppm  |
| acceptLostPeak     | This is an option for users, TRUE is the default value. If the users believe that all the peaks in the peak list are true positive, change it to FALSE.  |
| peak_det_align_dir | Output directory from <a href="#">pipe_peakdet_align</a>   |
| peak_width_ppm     | A peak width in ppm  |
| internal_calibrant | A ppm range where the internal calibrant is, or NULL.  |

### Details

If there is no internal calibrant, only the PQN normalization is done.

### Value

This function saves the result to the output directory

This function saves the result to the output directory

This function saves the result to the output directory

This function saves the result to the output directory

This function saves the result to the output directory

Pipeline: Filter samples according to metadata conditions

Pipeline: Peak detection and Alignment

Pipeline: Peak integration

Pipe: Full spectra normalization

### See Also

Other import/export functions: [files\\_to\\_rDolphin\(\)](#), [load\\_and\\_save\\_functions](#), [nmr\\_data\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_read\\_bruker\\_fid\(\)](#), [nmr\\_read\\_samples\(\)](#), [nmr\\_zip\\_bruker\\_samples\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [save\\_profiling\\_output\(\)](#), [to\\_ChemoSpec\(\)](#)

Other metadata functions: [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#)

Other outlier detection functions: [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_outliers\(\)](#)

Other peak detection functions: [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_identify\\_regions\\_blood\(\)](#), [nmr\\_identify\\_regions\\_cell\(\)](#), [nmr\\_identify\\_regions\\_urine\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [regions\\_from\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#)

Other alignment functions: [nmr\\_align\\_find\\_ref\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#)

Other peak integration functions: [computes\\_peak\\_width\\_ppm\(\)](#), [nmr\\_identify\\_regions\\_blood\(\)](#), [nmr\\_identify\\_regions\\_cell\(\)](#), [nmr\\_identify\\_regions\\_urine\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#)

## Examples

```
## Example of pipeline usage
## There are differet ways of load the dataset
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
#excel_file <- system.file("dataset-demo",
#                           "dummy_metadata.xlsx",
#                           package = "AlpsNMR")
#output_dir <- tempdir()

## Load samples with pipes
#pipe_load_samples(dir_to_demo_dataset,
#                  glob = "*.zip",
#                  output_dir = "../pipe_output")

## Another way to load it
#nmr_dataset <- nmr_read_samples_dir(dir_to_demo_dataset)

## Saving the dataset in a .rds file
#nmr_dataset_rds <- tempfile(fileext = ".rds")
#nmr_dataset_save(nmr_dataset, nmr_dataset_rds)

## Interpolation
#pipe_interpolate_1D(nmr_dataset_rds,
#                   axis = c(min = -0.5, max = 10, by = 2.3E-4),
#                   output_dir)

## Get the new path, based in output_dir
#nmr_dataset_rds <- paste(output_dir, "\", "nmr_dataset.rds", sep = "", collapse = NULL)

## Adding metadata to samples
#pipe_add_metadata(nmr_dataset_rds = nmr_dataset_rds, output_dir = output_dir,
#                  excel_file = excel_file)

## Filtering samples
#conditions <- 'SubjectID == "Ana"'
#pipe_filter_samples(nmr_dataset_rds, conditions, output_dir)

## Outlier detection
#pipe_outlier_detection(nmr_dataset_rds, output_dir)

## Exclude regions
#exclude_regions <- list(water = c(5.1, 4.5))
#pipe_exclude_regions(nmr_dataset_rds, exclude_regions, output_dir)

## peak aling
#pipe_peakdet_align(nmr_dataset_rds, output_dir = output_dir)

## peak integration
#pipe_peak_integration(nmr_dataset_rds,
#                      peak_det_align_dir = output_dir,
#                      peak_width_ppm = 0.006, output_dir)

## Normalization
#pipe_normalization(nmr_dataset_rds, output_dir = output_dir)
```

---

```
plot.nmr_dataset_1D Plot an nmr_dataset_1D
```

---

## Description

Plot an nmr\_dataset\_1D

## Usage

```
## S3 method for class 'nmr_dataset_1D'
plot(
  x,
  NMRExperiment = NULL,
  chemshift_range = NULL,
  interactive = FALSE,
  quantile_plot = NULL,
  quantile_colors = NULL,
  ...
)
```

## Arguments

|                 |  |
|-----------------|--|
| x               | a <a href="#">nmr_dataset_1D</a> object  |
| NMRExperiment   | A character vector with the NMRExperiments to include. Use "all" to include all experiments.   |
| chemshift_range | range of the chemical shifts to be included. Can be of length 3 to include the resolution in the third element (e.g. <code>c(0.2, 0.8, 0.005)</code> ) |
| interactive     | if TRUE return an interactive plotly plot, otherwise return a ggplot one.  |
| quantile_plot   | If TRUE plot the 10\ If two numbers between 0 and 1 are given then a custom percentile can be plotted  |
| quantile_colors | A vector with the colors for each of the quantiles   |
| ...             | arguments passed to <a href="#">ggplot2::aes_string</a> .  |

## Value

The plot

## See Also

Other plotting functions: [plot\\_interactive\(\)](#)

Other nmr\_dataset\_1D functions: [\[.nmr\\_dataset\\_1D\(\)](#), [computes\\_peak\\_width\\_ppm\(\)](#), [file\\_lister\(\)](#), [files\\_to\\_rDolphin\(\)](#), [format.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [nmr\\_align\\_find\\_ref\(\)](#), [nmr\\_baseline\\_removal\(\)](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_exclude\\_region\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_normalize\(\)](#), [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_outliers\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [plot\\_webgl\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [rdCV\\_PLS\\_RF\\_ML\(\)](#), [rdCV\\_PLS\\_RF\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [to\\_ChemoSpec\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

**Examples**

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
#dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
#dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
#plot(dataset_1D)
```

---

```
plot_bootstrap_multimodel
      Bootstrap plot predictions
```

---

**Description**

Bootstrap plot predictions

**Usage**

```
plot_bootstrap_multimodel(bp_results, dataset, y_column, plot = TRUE)
```

**Arguments**

|            |   |
|------------|---|
| bp_results | bp_kfold_VIP_analysis results   |
| dataset    | An <a href="#">nmr_dataset_family</a> object                                    |
| y_column   | A string with the name of the y column (present in the metadata of the dataset) |
| plot       | A boolean that indicate if results are plotted or not                           |

**Value**

A plot of the results or a ggplot object

**Examples**

```
# Data analysis for a table of integrated peaks

## Generate an artificial nmr_dataset_peak_table:
### Generate artificial metadata:
num_samples <- 64 # use an even number in this example
num_peaks <- 20
metadata <- data.frame(
  NMRExperiment = as.character(1:num_samples),
  Condition = rep(c("A", "B"), times = num_samples/2),
  stringsAsFactors = FALSE
)

### The matrix with peaks
peak_means <- runif(n = num_peaks, min = 300, max = 600)
peak_sd <- runif(n = num_peaks, min = 30, max = 60)
peak_matrix <- mapply(function(mu, sd) rnorm(num_samples, mu, sd),
                      mu = peak_means, sd = peak_sd)
colnames(peak_matrix) <- paste0("Peak", 1:num_peaks)

## Artificial differences depending on the condition:
```



```
peak_matrix[metadata$Condition == "A", "Peak2"] <-  
  peak_matrix[metadata$Condition == "A", "Peak2"] + 70  
  
peak_matrix[metadata$Condition == "A", "Peak6"] <-  
  peak_matrix[metadata$Condition == "A", "Peak6"] - 60  
  
### The nmr_dataset_peak_table  
peak_table <- new_nmr_dataset_peak_table(  
  peak_table = peak_matrix,  
  metadata = list(external = metadata)  
)  
  
## We will use bootstrap and permutation method for VIPs selection  
## in a a k-fold cross validation  
#bp_results <- bp_kfold_VIP_analysis(peak_table, # Data to be analyzed  
#                                     y_column = "Condition", # Label  
#                                     k = 3,  
#                                     nbootstrap = 10)  
  
#message("Selected VIPs are: ", bp_results$important_vips)  
  
#plot_bootstrap_multimodel(bp_results, peak_table, "Condition")
```

---

plot\_interactive      *Plots in WebGL*

---

## Description

Plots in WebGL

## Usage

```
plot_interactive(plt, html_filename)
```

## Arguments

plt                    A plot created with plotly or ggplot2  
html\_filename        The file name where the plot will be saved

## Value

The html\_filename

## See Also

Other plotting functions: [plot.nmr\\_dataset\\_1D\(\)](#)

**Examples**

```

dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
# plot <- plot(dataset_1D)
# html_plot_interactive <- plot_interactive(plot, "html_plot_interactive.html")

```

---

plot\_plsda\_multimodel *Multi PLDSA model plot predictions*

---

**Description**

Multi PLDSA model plot predictions

**Usage**

```
plot_plsda_multimodel(model, plot = TRUE)
```

**Arguments**

|       |   |
|-------|---|
| model | A nmr_data_analysis_model                             |
| plot  | A boolean that indicate if results are plotted or not |

**Value**

A plot of the results or a ggplot object

**Examples**

```

#' # Data analysis for a table of integrated peaks

## Generate an artificial nmr_dataset_peak_table:
### Generate artificial metadata:
num_samples <- 32 # use an even number in this example
num_peaks <- 20
metadata <- data.frame(
  NMRExperiment = as.character(1:num_samples),
  Condition = rep(c("A", "B"), times = num_samples/2),
  stringsAsFactors = FALSE
)

### The matrix with peaks
peak_means <- runif(n = num_peaks, min = 300, max = 600)
peak_sd <- runif(n = num_peaks, min = 30, max = 60)
peak_matrix <- mapply(function(mu, sd) rnorm(num_samples, mu, sd),
                      mu = peak_means, sd = peak_sd)
colnames(peak_matrix) <- paste0("Peak", 1:num_peaks)

## Artificial differences depending on the condition:
peak_matrix[metadata$Condition == "A", "Peak2"] <-
  peak_matrix[metadata$Condition == "A", "Peak2"] + 70

```

```
peak_matrix[metadata$Condition == "A", "Peak6"] <-  
  peak_matrix[metadata$Condition == "A", "Peak6"] - 60  
  
### The nmr_dataset_peak_table  
peak_table <- new_nmr_dataset_peak_table(  
  peak_table = peak_matrix,  
  metadata = list(external = metadata)  
)  
  
## We will use a double cross validation, splitting the samples with random  
## subsampling both in the external and internal validation.  
## The classification model will be a PLSDA, exploring at maximum 3 latent  
## variables.  
## The best model will be selected based on the area under the ROC curve  
methodology <- plsda_auroc_vip_method(ncomp = 1)  
model <- nmr_data_analysis(  
  peak_table,  
  y_column = "Condition",  
  identity_column = NULL,  
  external_val = list(iterations = 2, test_size = 0.25),  
  internal_val = list(iterations = 2, test_size = 0.25),  
  data_analysis_method = methodology  
)  
  
#plot_plsda_multimodel(model)
```

---

plot\_plsda\_samples      *Plot PLSDA predictions*

---

## Description

Plot PLSDA predictions

## Usage

```
plot_plsda_samples(model, newdata = NULL, plot = TRUE)
```

## Arguments

|         |  |
|---------|--|
| model   | A plsda model  |
| newdata | newdata to predict, if not included model\$X_test will be used |
| plot    | A boolean that indicate if results are plotted or not          |

## Value

A plot of the samples or a ggplot object

**Examples**

```

#' # Data analysis for a table of integrated peaks

## Generate an artificial nmr_dataset_peak_table:
### Generate artificial metadata:
num_samples <- 32 # use an even number in this example
num_peaks <- 20
metadata <- data.frame(
  NMRExperiment = as.character(1:num_samples),
  Condition = rep(c("A", "B"), times = num_samples/2),
  stringsAsFactors = FALSE
)

### The matrix with peaks
peak_means <- runif(n = num_peaks, min = 300, max = 600)
peak_sd <- runif(n = num_peaks, min = 30, max = 60)
peak_matrix <- mapply(function(mu, sd) rnorm(num_samples, mu, sd),
  mu = peak_means, sd = peak_sd)
colnames(peak_matrix) <- paste0("Peak", 1:num_peaks)

## Artificial differences depending on the condition:
peak_matrix[metadata$Condition == "A", "Peak2"] <-
  peak_matrix[metadata$Condition == "A", "Peak2"] + 70

peak_matrix[metadata$Condition == "A", "Peak6"] <-
  peak_matrix[metadata$Condition == "A", "Peak6"] - 60

### The nmr_dataset_peak_table
peak_table <- new_nmr_dataset_peak_table(
  peak_table = peak_matrix,
  metadata = list(external = metadata)
)

## We will use a double cross validation, splitting the samples with random
## subsampling both in the external and internal validation.
## The classification model will be a PLSDA, exploring at maximum 3 latent
## variables.
## The best model will be selected based on the area under the ROC curve
methodology <- plsda_auroc_vip_method(ncomp = 1)
model <- nmr_data_analysis(
  peak_table,
  y_column = "Condition",
  identity_column = NULL,
  external_val = list(iterations = 1, test_size = 0.25),
  internal_val = list(iterations = 1, test_size = 0.25),
  data_analysis_method = methodology
)

#plot_plsda_samples(model$outer_cv_results[[1]]$model)

```

**Description**

Plot vip scores of bootstrap

**Usage**

```
plot_vip_scores(vip_means, error, nbootstrap, plot = TRUE)
```

**Arguments**

|            |   |
|------------|---|
| vip_means  | vips means values of bootstraps                       |
| error      | error tolerated, calculated in the bootstrap          |
| nbootstrap | number of bootstraps realized                         |
| plot       | A boolean that indicate if results are plotted or not |

**Value**

A plot of the results or a ggplot object

**Examples**

```
# Data analysis for a table of integrated peaks

## Generate an artificial nmr_dataset_peak_table:
### Generate artificial metadata:
num_samples <- 64 # use an even number in this example
num_peaks <- 20
metadata <- data.frame(
  NMRExperiment = as.character(1:num_samples),
  Condition = rep(c("A", "B"), times = num_samples/2),
  stringsAsFactors = FALSE
)

### The matrix with peaks
peak_means <- runif(n = num_peaks, min = 300, max = 600)
peak_sd <- runif(n = num_peaks, min = 30, max = 60)
peak_matrix <- mapply(function(mu, sd) rnorm(num_samples, mu, sd),
  mu = peak_means, sd = peak_sd)
colnames(peak_matrix) <- paste0("Peak", 1:num_peaks)

## Artificial differences depending on the condition:
peak_matrix[metadata$Condition == "A", "Peak2"] <-
  peak_matrix[metadata$Condition == "A", "Peak2"] + 70

peak_matrix[metadata$Condition == "A", "Peak6"] <-
  peak_matrix[metadata$Condition == "A", "Peak6"] - 60

### The nmr_dataset_peak_table
peak_table <- new_nmr_dataset_peak_table(
  peak_table = peak_matrix,
  metadata = list(external = metadata)
)

## We will use bootstrap and permutation method for VIPs selection
## in a a k-fold cross validation
#bp_results <- bp_kfold_VIP_analysis(peak_table, # Data to be analyzed
```

```

#           y_column = "Condition", # Label
#           k = 3,
#           nbootstrap = 10)

#message("Selected VIPs are: ", bp_results$important_vips)

#plot_vip_scores(bp_results$fold_results[[1]]$vip_means,
#               bp_results$fold_results[[1]]$error[1],
#               nbootstrap = 10)

```

---

plot\_webgl

*Plot a dataset into a HTML file*


---

## Description

Uses WebGL for performance

## Usage

```
plot_webgl(nmr_dataset, html_filename, ...)
```

## Arguments

|                 |  |
|-----------------|--|
| nmr_dataset     | An <a href="#">nmr_dataset_1D</a>  |
| html_filename   | The output HTML filename to be created   |
| ...             | Arguments passed on to <a href="#">plot.nmr_dataset_1D</a>   |
|                 | x a <a href="#">nmr_dataset_1D</a> object  |
| chemshift_range | range of the chemical shifts to be included. Can be of length 3 to include the resolution in the third element (e.g. <code>c(0.2, 0.8, 0.005)</code> ) |
| NMRExperiment   | A character vector with the NMRExperiments to include. Use "all" to include all experiments.   |
| quantile_plot   | If TRUE plot the 10\ If two numbers between 0 and 1 are given then a custom percentile can be plotted  |
| quantile_colors | A vector with the colors for each of the quantiles   |
| interactive     | if TRUE return an interactive plotly plot, otherwise return a ggplot one.  |

## Value

the html filename created

## See Also

Other [nmr\\_dataset\\_1D](#) functions: [\[.nmr\\_dataset\\_1D\(\)](#), [computes\\_peak\\_width\\_ppm\(\)](#), [file\\_lister\(\)](#), [files\\_to\\_rDolphin\(\)](#), [format.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [nmr\\_align\\_find\\_ref\(\)](#), [nmr\\_baseline\\_removal\(\)](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_exclude\\_region\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_normalize\(\)](#), [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_outliers\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [plot.nmr\\_dataset\\_1D\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [rdCV\\_PLS\\_RF\\_ML\(\)](#), [rdCV\\_PLS\\_RF\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [to\\_ChemoSpec\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

**Examples**

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
#dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
#dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
#html_plot <- plot_webgl(dataset_1D, "html_plot.html")
```

---

```
plsda_auroc_vip_compare
```

*Compare PLSDA auroc VIP results*

---

**Description**

Compare PLSDA auroc VIP results

**Usage**

```
plsda_auroc_vip_compare(...)
```

**Arguments**

... Results of [nmr\\_data\\_analysis](#) to be combined. Give each result a name.

**Value**

A plot of the AUC for each method

**Examples**

```
# Data analysis for a table of integrated peaks

## Generate an artificial nmr_dataset_peak_table:
### Generate artificial metadata:
num_samples <- 32 # use an even number in this example
num_peaks <- 20
metadata <- data.frame(
  NMRExperiment = as.character(1:num_samples),
  Condition = rep(c("A", "B"), times = num_samples/2),
  stringsAsFactors = FALSE
)

### The matrix with peaks
peak_means <- runif(n = num_peaks, min = 300, max = 600)
peak_sd <- runif(n = num_peaks, min = 30, max = 60)
peak_matrix <- mapply(function(mu, sd) rnorm(num_samples, mu, sd),
  mu = peak_means, sd = peak_sd)
colnames(peak_matrix) <- paste0("Peak", 1:num_peaks)

## Artificial differences depending on the condition:
peak_matrix[metadata$Condition == "A", "Peak2"] <-
  peak_matrix[metadata$Condition == "A", "Peak2"] + 70

peak_matrix[metadata$Condition == "A", "Peak6"] <-
```

```

peak_matrix[metadata$Condition == "A", "Peak6"] - 60

### The nmr_dataset_peak_table
peak_table <- new_nmr_dataset_peak_table(
  peak_table = peak_matrix,
  metadata = list(external = metadata)
)

## We will use a double cross validation, splitting the samples with random
## subsampling both in the external and internal validation.
## The classification model will be a PLSDA, exploring at maximum 3 latent
## variables.
## The best model will be selected based on the area under the ROC curve
methodology <- plsda_auroc_vip_method(ncomp = 1)
model1 <- nmr_data_analysis(
  peak_table,
  y_column = "Condition",
  identity_column = NULL,
  external_val = list(iterations = 1, test_size = 0.25),
  internal_val = list(iterations = 1, test_size = 0.25),
  data_analysis_method = methodology
)

methodology2 <- plsda_auroc_vip_method(ncomp = 2)
model2 <- nmr_data_analysis(
  peak_table,
  y_column = "Condition",
  identity_column = NULL,
  external_val = list(iterations = 1, test_size = 0.25),
  internal_val = list(iterations = 1, test_size = 0.25),
  data_analysis_method = methodology2
)

plsda_auroc_vip_compare(model1 = model1, model2 = model2)

```

---

plsda\_auroc\_vip\_method

*Method for nmr\_data\_analysis (PLSDA model with AUROC and VIP outputs)*

---

## Description

Method for nmr\_data\_analysis (PLSDA model with AUROC and VIP outputs)

## Usage

```
plsda_auroc_vip_method(ncomp, auc_increment_threshold = 0.05)
```

## Arguments

ncomp                    Max. number of latent variables to explore in the PLSDA analysis

auc\_increment\_threshold                    Choose the number of latent variables when the AUC does not increment more than this threshold.



**Value**

Returns an object to be used with `nmr_data_analysis` to perform a (optionally multilevel) PLS-DA model, using the area under the ROC curve as figure of merit to determine the optimum number of latent variables.

**Examples**

```
method <- plsda_auroc_vip_method(3)
```

---

|                |                                |
|----------------|--------------------------------|
| ppm_resolution | <i>Unlisted PPM resolution</i> |
|----------------|--------------------------------|

---

**Description**

A wrapper to unlist the output from the function `nmr_ppm_resolution(nmr_dataset)` when no interpolation has been applied.

**Usage**

```
ppm_resolution(nmr_dataset)
```

**Arguments**

`nmr_dataset` An object containing NMR samples

**Value**

A number (the ppm resolution, measured in ppms)

Numeric (the ppm resolution, measured in ppms)

**Examples**

```
nmr_dataset <- nmr_dataset_load(system.file("extdata", "nmr_dataset.rds", package = "AlpsNMR"))
nmr_ppm_resolution(nmr_dataset)
```

---

|                |  |
|----------------|--|
| ppm_VIP_vector | <i>Feature selection and validation in multivariate analysis</i> |
|----------------|--|

---

**Description**

Numeric VIPs vector

**Usage**

```
ppm_VIP_vector(VIPs)
```

**Arguments**

`VIPs` a dataframe from the `model_VIP(MVObj)` function. It requires a "ppms" variable

**Details**

The function extracts the VIPs vector (numeric) from the `model_VIP(MVObj)` function. It is not necessary if you have the ppm values in a numeric vector. This is needed in case that an automated pipeline is applied, connecting the output from `model_VIP(MVObj)` to `nmr_identify_regions` family functions.

**Value**

a numeric ppm vector ready to be identified with `nmr_identify_regions_blood`, `nmr_identify_cell` or `nmr_identify_regions_urine`

**Examples**

```
message("MUVR is not compatible with Bioconductor,
use bp_kfold_VIP_analysis method instead")

## Example of MUVR usage
# 1.Build a model with the X data from your nmr object and your class:
#MVObj <- rdCV_PLS_RF(nmr_data(nmr_peak_table),
#Y = nmr_peak_table_completed$Timepoint)

# 2.Model performance
#confusion_matrix(MVObj)

# 3.Plotting the model
#MUVR_model_plot(MVObj)

# 4.Permutation test
#permutations <- permutation_test_model(MVObj, nPerm = 50)

# 5.Plotting permutation test results
#permutation_test_plot(MVObj, permutations, model = "Mid", type = "t")

# 6.p-Value
#p.value <- p_value_perm(MVObj$miss[["mid"]], permutations[, "Mid"])

# 7.Significant variables
#VIPs <- model_VIP(MVObj)

# 8.Identification
#results <- nmr_identify_regions_blood(ppm_VIP_vector(VIPs))
```

---

```
print.nmr_dataset      Print for nmr_dataset
```

---

**Description**

Print for `nmr_dataset`

**Usage**

```
## S3 method for class 'nmr_dataset'
print(x, ...)
```

**Arguments**

x                    an `nmr_dataset` object  
...                    for future use

**Value**

Print for `nmr_dataset`

**See Also**

Other class helper functions: `format.nmr_dataset_1D()`, `format.nmr_dataset_peak_table()`, `format.nmr_dataset()`, `is.nmr_dataset_1D()`, `is.nmr_dataset_peak_table()`, `new_nmr_dataset_1D()`, `new_nmr_dataset_peak_table()`, `new_nmr_dataset()`, `print.nmr_dataset_1D()`, `print.nmr_dataset_peak_table()`, `validate_nmr_dataset_family()`, `validate_nmr_dataset_peak_table()`, `validate_nmr_dataset()`

Other `nmr_dataset` functions: `[.nmr_dataset()`, `format.nmr_dataset()`, `load_and_save_functions`, `new_nmr_dataset()`, `nmr_interpolate_1D()`, `nmr_meta_add()`, `nmr_meta_export()`, `nmr_meta_get_column()`, `nmr_meta_get()`, `nmr_ppm_resolution()`, `nmr_read_samples()`, `validate_nmr_dataset()`

**Examples**

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
print(dataset)
```

---

`print.nmr_dataset_1D`    *print for nmr\_dataset\_1D*

---

**Description**

print for `nmr_dataset_1D`

**Usage**

```
## S3 method for class 'nmr_dataset_1D'
print(x, ...)
```

**Arguments**

x                    an `nmr_dataset_1D` object  
...                    for future use

**Value**

print for `nmr_dataset_1D`

**See Also**

Other class helper functions: `format.nmr_dataset_1D()`, `format.nmr_dataset_peak_table()`, `format.nmr_dataset()`, `is.nmr_dataset_1D()`, `is.nmr_dataset_peak_table()`, `new_nmr_dataset_1D()`, `new_nmr_dataset_peak_table()`, `new_nmr_dataset()`, `print.nmr_dataset_peak_table()`, `print.nmr_dataset()`, `validate_nmr_dataset_family()`, `validate_nmr_dataset_peak_table()`, `validate_nmr_dataset()`

Other `nmr_dataset_1D` functions: `[.nmr_dataset_1D()`, `computes_peak_width_ppm()`, `file_lister()`, `files_to_rDolphin()`, `format.nmr_dataset_1D()`, `is.nmr_dataset_1D()`, `load_and_save_functions`, `new_nmr_dataset_1D()`, `nmr_align_find_ref()`, `nmr_baseline_removal()`, `nmr_baseline_threshold()`, `nmr_exclude_region()`, `nmr_integrate_regions()`, `nmr_interpolate_1D()`, `nmr_meta_add()`, `nmr_meta_export()`, `nmr_meta_get_column()`, `nmr_meta_get()`, `nmr_normalize()`, `nmr_pca_build_model()`, `nmr_pca_outliers_filter()`, `nmr_pca_outliers_plot()`, `nmr_pca_outliers_robust()`, `nmr_pca_outliers()`, `nmr_ppm_resolution()`, `plot.nmr_dataset_1D()`, `plot_webgl()`, `rdCV_PLS_RF_ML()`, `rdCV_PLS_RF()`, `save_files_to_rDolphin()`, `to_ChemoSpec()`, `validate_nmr_dataset_peak_table()`, `validate_nmr_dataset()`

**Examples**

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
print(dataset_1D)
```

---

```
print.nmr_dataset_peak_table
      print for nmr_dataset_peak_table
```

---

**Description**

print for `nmr_dataset_peak_table`

**Usage**

```
## S3 method for class 'nmr_dataset_peak_table'
print(x, ...)
```

**Arguments**

```
x          an nmr_dataset_peak_table object
...        for future use
```

**Value**

print for `nmr_dataset_peak_table`

**See Also**

Other `nmr_dataset_peak_table` functions: `[.nmr_dataset_peak_table()`, `format.nmr_dataset_peak_table()`, `is.nmr_dataset_peak_table()`, `load_and_save_functions`, `new_nmr_dataset_peak_table()`, `nmr_meta_add()`, `nmr_meta_export()`, `nmr_meta_get_column()`, `nmr_meta_get()`, `validate_nmr_dataset_peak_table()`

Other class helper functions: `format.nmr_dataset_1D()`, `format.nmr_dataset_peak_table()`, `format.nmr_dataset()`, `is.nmr_dataset_1D()`, `is.nmr_dataset_peak_table()`, `new_nmr_dataset_1D()`, `new_nmr_dataset_peak_table()`, `new_nmr_dataset()`, `print.nmr_dataset_1D()`, `print.nmr_dataset()`, `validate_nmr_dataset_family()`, `validate_nmr_dataset_peak_table()`, `validate_nmr_dataset()`

**Examples**

```

dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
meta <- file.path(dir_to_demo_dataset, "dummy_metadata.xlsx")
metadata <- readxl::read_excel(meta, sheet = 1)
dataset_1D <- nmr_meta_add(dataset_1D, metadata = metadata, by = "NMRExperiment")
metadata <- list(external = dataset_1D[["metadata"]][["external"]])
peak_table <- nmr_data(dataset_1D)
new <- new_nmr_dataset_peak_table(peak_table, metadata)
new

```

---

p\_value\_perm

*Deprecated function p-Value from permutation test*


---

**Description**

The function calculates the cumulative (1-tailed) probability of 'actual' belonging to 'h0' (permutation\_object from the permutation\_test\_model function).

**Usage**

```
p_value_perm(model_actual, permutation_object)
```

**Arguments**

model\_actual    The actual model performance (e.g. misclassifications or Q2)

permutation\_object

Null hypothesis distribution from permutation test from permutation\_test\_model function

**Value**

The p-value indicating if there is significant differences between the model performance and the null hypothesis distribution from permutation test test

**Examples**

```

message("MUVR is not compatible with Bioconductor,
use bp_kfold_VIP_analysis method instead")

```

---

|                    |                           |
|--------------------|---------------------------|
| random_subsampling | <i>Random subsampling</i> |
|--------------------|---------------------------|

---

### Description

Random subsampling

### Usage

```
random_subsampling(  
  sample_idx,  
  iterations = 10L,  
  test_size = 0.25,  
  keep_together = NULL,  
  balance_in_train = NULL  
)
```

### Arguments

|                  |   |
|------------------|---|
| sample_idx       | Typically a numeric vector with sample index to be separated. A character vector with sample IDs could also be used   |
| iterations       | An integer, the number of iterations in the random subsampling  |
| test_size        | A number between 0 and 1. The samples to be included in the test set on each iteration.   |
| keep_together    | Either NULL or a factor with the same length as sample_idx. keep_together can be used to ensure that groups of samples are kept in together in all iterations (either on training or on test, but never split). A typical use case for this is when you have sample replicates and you want to keep all replicates together to prevent overoptimistic results (having one sample on the train subset and its replicate on the test subset would make the prediction easier to guess). Another use case for this is when you have a longitudinal study and you want to keep some subjects in the same train or test group, because you want to use some information in a longitudinal way (e.g. a multilevel plsda model). |
| balance_in_train | Either NULL or a factor with the same length as sample_idx. balance_in_train can be used to force that on each iteration, the train partition contains the same number of samples of the given factor levels. For instance, if we have a dataset with 40 samples of class "A" and 20 samples of class "B", using a test_size = 0.25, we can force to always have 16 samples of class "A" and 16 samples of class "B" in the training subset. This is beneficial to those algorithms that require that the training groups are balanced.   |

### Value

A list of length equal to iterations. Each element of the list is a list with two entries (training and test) containing the sample\_idx values that will belong to each subset.

**Examples**

```
random_subsampling(1:100, iterations = 4, test_size = 0.25)

subject_id <- c("Alice", "Bob", "Charlie", "Eve")
random_subsampling(1:4, iterations = 2, test_size = 0.25, keep_together = subject_id)
```

---

rdCV\_PLS\_RF

*Deprecated function*

---

**Description**

Deprecated function

**Usage**

```
rdCV_PLS_RF(  
  X,  
  Y,  
  ID,  
  scale = TRUE,  
  nRep = 10,  
  nOuter = 5,  
  nInner,  
  varRatio = 0.75,  
  DA = FALSE,  
  fitness = "MISS",  
  method = "PLS",  
  nCompMax,  
  methParam,  
  ML = FALSE,  
  modReturn = FALSE,  
  logg = FALSE,  
  parallel = TRUE  
)
```

**Value**

a MUVR model containing selection parameters, validation and fitness

**References**

Shi,L. et al. (2018) Variable selection and validation in multivariate modelling. *Bioinformatics*.

**See Also**

`nmr_data_analysis` Feature selection and validation in multivariate analysis

Statistical analysis and feature selection in a repeated double cross-validation frame based on the partial least squares (PLS) or random forest (RF) analyses using an algorithm for multivariate modelling with minimally biased variable selection (MUVR) from the MUVR package. If your work with a `nmr_peak_table` object from `AlpsNMR`, first you need to extract the X data from the main

nmr\_dataset object (e.g. your peak table) with the nmr\_data function, otherwise you would try to set a list on the X. You also need to set the class from this object, or just set it from another Y vector.

Other nmr\_dataset\_1D functions: [.nmr\_dataset\_1D(), computes\_peak\_width\_ppm(), file\_lister(), files\_to\_rDolphin(), format.nmr\_dataset\_1D(), is.nmr\_dataset\_1D(), load\_and\_save\_functions, new\_nmr\_dataset\_1D(), nmr\_align\_find\_ref(), nmr\_baseline\_removal(), nmr\_baseline\_threshold(), nmr\_exclude\_region(), nmr\_integrate\_regions(), nmr\_interpolate\_1D(), nmr\_meta\_add(), nmr\_meta\_export(), nmr\_meta\_get\_column(), nmr\_meta\_get(), nmr\_normalize(), nmr\_pca\_build\_model(), nmr\_pca\_outliers\_filter(), nmr\_pca\_outliers\_plot(), nmr\_pca\_outliers\_robust(), nmr\_pca\_outliers(), nmr\_ppm\_resolution(), plot.nmr\_dataset\_1D(), plot\_webgl(), print.nmr\_dataset\_1D(), rdCV\_PLS\_RF\_ML(), save\_files\_to\_rDolphin(), to\_ChemoSpec(), validate\_nmr\_dataset\_peak\_table(), validate\_nmr\_dataset()]

### Examples

```
message("MUVR is not compatible with Bioconductor,
use bp_kfold_VIP_analysis method instead")
```

---

|                |   |
|----------------|---|
| rdCV_PLS_RF_ML | <i>Deprecated function Feature selection and validation in MULTI-LEVEL analysis</i> |
|----------------|---|

---

### Description

Statistical analysis and feature selection in a repeated double cross-validation frame based on the partial least squares (PLS) or random forest (RF) analyses using an algorithm for multivariate modelling with minimally biased variable selection (MUVR) from the MUVR package. The function rdCV\_PLS\_RF\_ML allows the multilevel comparison, especially useful in crossover or longitudinal studies (2 timepoints) considering the same individual (it requires 2 samples of the same observation).

### Usage

```
rdCV_PLS_RF_ML(
  nmr_peak_table,
  label,
  scale = TRUE,
  nRep = 10,
  nOuter = 5,
  nInner,
  varRatio = 0.75,
  DA = FALSE,
  fitness = "MISS",
  method = "PLS",
  ML = TRUE,
  modReturn = FALSE,
  logg = FALSE,
  parallel = TRUE
)
```



**Arguments**

nmr\_peak\_table an AlpsNMR integration object (2 classes)  
 label the name of the variable to test (e.g. "Timepoint")

**Value**

a MUVR model containing selection parameters, validation and fitness

**References**

Shi,L. et al. (2018) Variable selection and validation in multivariate modelling. *Bioinformatics*.

**See Also**

Other nmr\_dataset\_1D functions: [[.nmr\\_dataset\\_1D\(\)](#), [computes\\_peak\\_width\\_ppm\(\)](#), [file\\_lister\(\)](#), [files\\_to\\_rDolphin\(\)](#), [format.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [nmr\\_align\\_find\\_ref\(\)](#), [nmr\\_baseline\\_removal\(\)](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_exclude\\_region\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_normalize\(\)](#), [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_outliers\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [plot.nmr\\_dataset\\_1D\(\)](#), [plot\\_webgl\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [rdCV\\_PLS\\_RF\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [to\\_ChemoSpec\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)]

**Examples**

```
message("MUVR is not compatible with Bioconductor,
use bp_kfold_VIP_analysis method instead")
```

---

|                    |                                       |
|--------------------|---------------------------------------|
| read Bruker sample | <i>Read a Bruker sample directory</i> |
|--------------------|---------------------------------------|

---

**Description**

Read a Bruker sample directory

**Usage**

```
read Bruker sample(
  sample_path,
  pdata_file = NULL,
  pdata_path = "pdata/1",
  all_components = FALSE
)
```

**Arguments**

sample\_path A character path of the sample directory  
 pdata\_file File name of the binary NMR data to load. Usually "1r". If it is null it is autodetected and all files are loaded.  
 pdata\_path Path from sample\_path to the preprocessed data  
 all\_components If FALSE load only the real component. Otherwise load all of them

**Value**

a list with all the bruker sample information

---

regions\_from\_peak\_table

*Build list of regions for peak integration*

---

**Description**

Build list of regions for peak integration

**Usage**

```
regions_from_peak_table(peak_pos_ppm, peak_width_ppm)
```

**Arguments**

peak\_pos\_ppm The peak positions, in ppm

peak\_width\_ppm The peak widths (or a single peak width for all peaks)

**Value**

A list of regions suitable for [nmr\\_integrate\\_regions](#)

**See Also**

Other peak detection functions: [Pipelines](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_identify\\_regions\\_blood\(\)](#), [nmr\\_identify\\_regions\\_cell\(\)](#), [nmr\\_identify\\_regions\\_urine\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#)

---

ROI\_blood

*ROIs for blood (plasma/serum) samples*

---

**Description**

The template ROI\_blood contains the targeted list of metabolites to be quantified (blood samples)

**References**

[github.com/danielcanueto/rDolphin](https://github.com/danielcanueto/rDolphin)

---

|          |                              |
|----------|------------------------------|
| ROI_cell | <i>ROIs for cell samples</i> |
|----------|------------------------------|

---

**Description**

The template ROI\_cell contains the targeted list of metabolites to be quantified (cell samples)

**References**

[github.com/danielcanueto/rDolphin](https://github.com/danielcanueto/rDolphin)

---

|           |                               |
|-----------|-------------------------------|
| ROI_urine | <i>ROIs for urine samples</i> |
|-----------|-------------------------------|

---

**Description**

The template ROI\_urine contains the targeted list of metabolites to be quantified (urine samples)

**References**

[github.com/danielcanueto/rDolphin](https://github.com/danielcanueto/rDolphin)

---

|                        |                               |
|------------------------|-------------------------------|
| save_files_to_rDolphin | <i>Save files to rDolphin</i> |
|------------------------|-------------------------------|

---

**Description**

The function saves the CSV files required by to\_rDolphin and Automatic\_targeted\_profiling functions for metabolite profiling.

**Usage**

```
save_files_to_rDolphin(files_rDolphin, output_directory)
```

**Arguments**

files\_rDolphin a list containing 4 elements from files\_to\_rDolphin

- meta\_rDolphin: metadata in rDolphin format,
- NMR\_spectra: spectra matrix
- ROI: ROI template
- Parameters\_blood: parameters file

output\_directory

a directory in which the CSV files are saved

**Value**

CSV files containing:

**See Also**

Other import/export functions: [Pipelines](#), [files\\_to\\_rDolphin\(\)](#), [load\\_and\\_save\\_functions](#), [nmr\\_data\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_read\\_bruker\\_fid\(\)](#), [nmr\\_read\\_samples\(\)](#), [nmr\\_zip\\_bruker\\_samples\(\)](#), [save\\_profiling\\_output\(\)](#), [to\\_ChemoSpec\(\)](#)

Other `nmr_dataset_1D` functions: [\[.nmr\\_dataset\\_1D\(\)](#), [computes\\_peak\\_width\\_ppm\(\)](#), [file\\_lister\(\)](#), [files\\_to\\_rDolphin\(\)](#), [format.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [nmr\\_align\\_find\\_ref\(\)](#), [nmr\\_baseline\\_removal\(\)](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_exclude\\_region\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_normalize\(\)](#), [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_outliers\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [plot.nmr\\_dataset\\_1D\(\)](#), [plot\\_webgl\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [rdCV\\_PLS\\_RF\\_ML\(\)](#), [rdCV\\_PLS\\_RF\(\)](#), [to\\_ChemoSpec\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

Other `to_rDolphin_blood` functions: [save\\_profiling\\_output\(\)](#)

**Examples**

```
## Not run:
dataset <- system.file("dataset-demo", package = "AlpsNMR")
excel_file <- system.file("dataset-demo", "dummy_metadata.xlsx", package = "AlpsNMR")
nmr_dataset <- nmr_read_samples_dir(dataset)
files_rDolphin = files_to_rDolphin_blood(nmr_dataset)
output_directory = "."
save_files_to_rDolphin(files_rDolphin, output_directory)

## End(Not run)
```

---

save\_profiling\_output *Save rDolphin output*

---

**Description**

The function saves the output from `Automatic_targeted_profiling` function in CSV format.

**Usage**

```
save_profiling_output(targeted_profiling, output_directory)
```

**Arguments**

`targeted_profiling`  
A list from `Automatic_targeted_profiling` function

`output_directory`  
a directory in which the CSV files are saved

**Value**

rDolphin output from `Automatic_targeted_profiling` function:

- `metabolites_intensity`
- `metabolites_quantification`

- ROI\_profiles\_used
- chemical\_shift
- fitting\_error
- half\_bandwidth
- signal\_area\_ratio

### See Also

Other import/export functions: [Pipelines](#), [files\\_to\\_rDolphin\(\)](#), [load\\_and\\_save\\_functions](#), [nmr\\_data\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_read\\_bruker\\_fid\(\)](#), [nmr\\_read\\_samples\(\)](#), [nmr\\_zip\\_bruker\\_samples\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [to\\_ChemoSpec\(\)](#)

Other to\_rDolphin\_blood functions: [save\\_files\\_to\\_rDolphin\(\)](#)

### Examples

```
## Not run:
rDolphin_object = to_rDolphin(parameters)
targeted_profiling = Automatic_targeted_profiling(rDolphin)
save_profiling_output(targeted_profiling, output_directory)

## End(Not run)
```

---

SummarizedExperiment\_to\_nmr\_dataset\_peak\_table

*Import SummarizedExperiment as mr\_dataset\_peak\_table*

---

### Description

Import SummarizedExperiment as mr\_dataset\_peak\_table

### Usage

```
SummarizedExperiment_to_nmr_dataset_peak_table(se)
```

### Arguments

se                    An SummarizedExperiment object

### Value

nmr\_dataset\_peak\_table An [nmr\\_dataset\\_peak\\_table](#) object (unmodified)

### Examples

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
meta <- file.path(dir_to_demo_dataset, "dummy_metadata.xlsx")
metadata <- readxl::read_excel(meta, sheet = 1)
dataset_1D <- nmr_meta_add(dataset_1D, metadata = metadata, by = "NMRExperiment")
metadata <- list(external = dataset_1D[["metadata"]][["external"]])
peak_table <- nmr_data(dataset_1D)
```

```
nmr_peak_table <- new_nmr_dataset_peak_table(peak_table, metadata)
se <- nmr_dataset_peak_table_to_SummarizedExperiment(nmr_peak_table)
nmr_peak_table <- SummarizedExperiment_to_nmr_dataset_peak_table(se)
```

---

SummarizedExperiment\_to\_nmr\_data\_1r

*Import SummarizedExperiment as 1D NMR data*

---

### Description

Import SummarizedExperiment as 1D NMR data

### Usage

```
SummarizedExperiment_to_nmr_data_1r(se)
```

### Arguments

se                    An SummarizedExperiment object

### Value

nmr\_dataset An [nmr\\_dataset\\_1D](#) object (unmodified)

### Examples

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
se <- nmr_data_1r_to_SummarizedExperiment(dataset_1D)
dataset_1D <- SummarizedExperiment_to_nmr_data_1r(se)
```

---

to\_ChemoSpec

*Convert to ChemoSpec Spectra class*

---

### Description

Convert to ChemoSpec Spectra class

### Usage

```
to_ChemoSpec(nmr_dataset, desc = "A nmr_dataset")
```

### Arguments

nmr\_dataset        An [nmr\\_dataset\\_1D](#) object  
 desc                a description for the dataset

### Value

A Spectra object from the ChemoSpec package

**See Also**

Other import/export functions: [Pipelines](#), [files\\_to\\_rDolphin\(\)](#), [load\\_and\\_save\\_functions](#), [nmr\\_data\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_read\\_bruker\\_fid\(\)](#), [nmr\\_read\\_samples\(\)](#), [nmr\\_zip\\_bruker\\_samples\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [save\\_profiling\\_output\(\)](#)

Other `nmr_dataset_1D` functions: [\[.nmr\\_dataset\\_1D\(\)](#), [computes\\_peak\\_width\\_ppm\(\)](#), [file\\_lister\(\)](#), [files\\_to\\_rDolphin\(\)](#), [format.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [nmr\\_align\\_find\\_ref\(\)](#), [nmr\\_baseline\\_removal\(\)](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_exclude\\_region\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_normalize\(\)](#), [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_outliers\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [plot.nmr\\_dataset\\_1D\(\)](#), [plot\\_webgl\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [rdCV\\_PLS\\_RF\\_ML\(\)](#), [rdCV\\_PLS\\_RF\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

**Examples**

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
chemo_spectra <- to_ChemoSpec(dataset_1D)
```

---

`validate_nmr_dataset`    *Validate nmr\_dataset objects*

---

**Description**

Validate `nmr_dataset` objects

Validate 1D `nmr` datasets

**Usage**

```
validate_nmr_dataset(samples)
```

```
validate_nmr_dataset_1D(nmr_dataset_1D)
```

**Arguments**

`samples`            An `nmr_dataset` object

`nmr_dataset_1D`    An [nmr\\_dataset\\_1D](#) object

**Value**

Validate `nmr_dataset` objects

The [nmr\\_dataset\\_1D](#) unchanged

This function is useful for its side-effects. Stopping in case of error

**See Also**

Other class helper functions: `format.nmr_dataset_1D()`, `format.nmr_dataset_peak_table()`, `format.nmr_dataset()`, `is.nmr_dataset_1D()`, `is.nmr_dataset_peak_table()`, `new_nmr_dataset_1D()`, `new_nmr_dataset_peak_table()`, `new_nmr_dataset()`, `print.nmr_dataset_1D()`, `print.nmr_dataset_peak_table()`, `print.nmr_dataset()`, `validate_nmr_dataset_family()`, `validate_nmr_dataset_peak_table()`

Other `nmr_dataset` functions: `[.nmr_dataset()`, `format.nmr_dataset()`, `load_and_save_functions`, `new_nmr_dataset()`, `nmr_interpolate_1D()`, `nmr_meta_add()`, `nmr_meta_export()`, `nmr_meta_get_column()`, `nmr_meta_get()`, `nmr_ppm_resolution()`, `nmr_read_samples()`, `print.nmr_dataset()`

Other class helper functions: `format.nmr_dataset_1D()`, `format.nmr_dataset_peak_table()`, `format.nmr_dataset()`, `is.nmr_dataset_1D()`, `is.nmr_dataset_peak_table()`, `new_nmr_dataset_1D()`, `new_nmr_dataset_peak_table()`, `new_nmr_dataset()`, `print.nmr_dataset_1D()`, `print.nmr_dataset_peak_table()`, `print.nmr_dataset()`, `validate_nmr_dataset_family()`, `validate_nmr_dataset_peak_table()`

Other `nmr_dataset_1D` functions: `[.nmr_dataset_1D()`, `computes_peak_width_ppm()`, `file_lister()`, `files_to_rDolphin()`, `format.nmr_dataset_1D()`, `is.nmr_dataset_1D()`, `load_and_save_functions`, `new_nmr_dataset_1D()`, `nmr_align_find_ref()`, `nmr_baseline_removal()`, `nmr_baseline_threshold()`, `nmr_exclude_region()`, `nmr_integrate_regions()`, `nmr_interpolate_1D()`, `nmr_meta_add()`, `nmr_meta_export()`, `nmr_meta_get_column()`, `nmr_meta_get()`, `nmr_normalize()`, `nmr_pca_build_model()`, `nmr_pca_outliers_filter()`, `nmr_pca_outliers_plot()`, `nmr_pca_outliers_robust()`, `nmr_pca_outliers()`, `nmr_ppm_resolution()`, `plot.nmr_dataset_1D()`, `plot_webgl()`, `print.nmr_dataset_1D()`, `rdCV_PLS_RF_ML()`, `rdCV_PLS_RF()`, `save_files_to_rDolphin()`, `to_ChemoSpec()`, `validate_nmr_dataset_peak`

**Examples**

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
validate_nmr_dataset(dataset)
```

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
dataset_1D_validated <- validate_nmr_dataset_1D(dataset_1D)
```

---

validate\_nmr\_dataset\_family

*Validate nmr\_dataset\_family objects*

---

**Description**

Validate `nmr_dataset_family` objects

**Usage**

```
validate_nmr_dataset_family(nmr_dataset_family)
```

**Arguments**

`nmr_dataset_family`

An `nmr_dataset_family` object



**Value**

The `nmr_dataset_family` unchanged

This function is useful for its side-effects: Stopping in case of error

**See Also**

Other class helper functions: `format.nmr_dataset_1D()`, `format.nmr_dataset_peak_table()`, `format.nmr_dataset()`, `is.nmr_dataset_1D()`, `is.nmr_dataset_peak_table()`, `new_nmr_dataset_1D()`, `new_nmr_dataset_peak_table()`, `new_nmr_dataset()`, `print.nmr_dataset_1D()`, `print.nmr_dataset_peak_table()`, `print.nmr_dataset()`, `validate_nmr_dataset_peak_table()`, `validate_nmr_dataset()`

**Examples**

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))
validate_nmr_dataset_family(dataset_1D)
```

---

```
validate_nmr_dataset_peak_table
```

*Validate nmr\_dataset\_peak\_table objects*

---

**Description**

The function detects peaks on an `nmr_dataset_1D` object, using `speaq::detectSpecPeaks`. `detectSpecPeaks` divides the whole spectra into smaller segments and uses `MassSpecWavelet::peakDetectionCWT` for peak detection.

This function is based on `speaq::dohCluster`.

The function allows the integration of a given ppm vector with a specific width.

**Usage**

```
validate_nmr_dataset_peak_table(nmr_dataset_peak_table)
```

```
nmr_detect_peaks(
  nmr_dataset,
  nDivRange_ppm = 0.1,
  scales = seq(1, 16, 2),
  baselineThresh = NULL,
  SNR.Th = 3
)
```

```
nmr_detect_peaks_plot(nmr_dataset, peak_data, NMRExperiment, ...)
```

```
nmr_detect_peaks_tune_snr(
  ds,
  NMRExperiment = NULL,
  SNR_thresholds = seq(from = 2, to = 6, by = 0.1)
)
```

```

nmr_align(
  nmr_dataset,
  peak_data,
  NMRExp_ref = NULL,
  maxShift_ppm = 0.0015,
  acceptLostPeak = FALSE
)

nmr_integrate_peak_positions(
  samples,
  peak_pos_ppm,
  peak_width_ppm = 0.006,
  ...
)

get_integration_with_metadata(integration_object)

```

### Arguments

|                        |  |
|------------------------|--|
| nmr_dataset_peak_table | An <a href="#">nmr_dataset_peak_table</a> object   |
| nmr_dataset            | An <a href="#">nmr_dataset_1D</a>  |
| nDivRange_ppm          | Segment size, in ppms, to divide the spectra and search for peaks.   |
| scales                 | The parameter of peakDetectionCWT function of MassSpecWavelet package, look it up in the original function.  |
| baselineThresh         | It will remove all peaks under an intensity set by baselineThresh. If you set it to 'NULL', nmr_detect_peaks will automatically compute an approximate value considering baseline between 9.5 and 10.0 ppm (automatically calculation using baselineThresh = NULL will not work if spectra were not interpolated up to 10.0 ppm) |
| SNR.Th                 | The parameter of peakDetectionCWT function of MassSpecWavelet package, look it up in the original function. If you set -1, the function will itself re-compute this value.   |
| peak_data              | The detected peak data given by <a href="#">nmr_detect_peaks</a> .   |
| NMRExperiment          | A string with the single NMRExperiment used explore the SNR thresholds. If not given, use the first one.   |
| ...                    | Arguments passed on to <a href="#">nmr_integrate_regions</a>   |
| regions                | A named list. Each element of the list is a region, given as a named numeric vector of length two with the range to integrate. The name of the region will be the name of the column   |
| ds                     | An <a href="#">nmr_dataset_1D</a> dataset  |
| SNR_thresholds         | A numeric vector with the SNR thresholds to explore  |
| NMRExp_ref             | NMRExperiment of the reference to use for alignment  |
| maxShift_ppm           | The maximum shift allowed, in ppm  |
| acceptLostPeak         | This is an option for users, TRUE is the default value. If the users believe that all the peaks in the peak list are true positive, change it to FALSE.  |
| samples                | A <a href="#">nmr_dataset</a> object   |
| peak_pos_ppm           | The peak positions, in ppm   |

peak\_width\_ppm The peak widths (or a single peak width for all peaks)  
 integration\_object  
     A [nmr\\_dataset](#) object

## Value

The [nmr\\_dataset\\_peak\\_table](#) unchanged

This function is useful for its side-effects: Stopping in case of error

A data frame with the NMRExperiment, the sample index, the position in ppm and index and the peak intensity

Plot peak detection results

A list with the following elements:

- peaks\_detected: A data frame with the columns from the [nmr\\_detect\\_peaks](#) output and an additional column SNR\_threshold with the threshold used on each row.
- num\_peaks\_per\_region: A summary of the peaks\_detected table, with the number of peaks detected on each chemical shift region
- plot\_num\_peaks\_per\_region: A visual representation of num\_peaks\_per\_region
- plot\_spectrum\_and\_detections: A visual representation of the spectrum and the peaks detected with each SNR threshold. Use [plotly::ggplotly](#) or [plot\\_interactive](#) on this to zoom and explore the results.

An [nmr\\_dataset\\_1D](#), with the spectra aligned

Integrate peak positions

Get integrals with metadata from integrate peak positions

integration dataframe

## Parallelization

This function accepts parallelization with future strategies. You can use `plan(multiprocess)` or `plan(sequential)` before calling this function to determine if it should be parallelized or not.

## See Also

[nmr\\_align](#) for peak alignment with the detected peak table

Other [nmr\\_dataset\\_peak\\_table](#) functions: [\[.nmr\\_dataset\\_peak\\_table\(\)](#), [format.nmr\\_dataset\\_peak\\_table\(\)](#), [is.nmr\\_dataset\\_peak\\_table\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_peak\\_table\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [print.nmr\\_dataset\\_peak\\_table\(\)](#)

Other class helper functions: [format.nmr\\_dataset\\_1D\(\)](#), [format.nmr\\_dataset\\_peak\\_table\(\)](#), [format.nmr\\_dataset\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_peak\\_table\(\)](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [new\\_nmr\\_dataset\\_peak\\_table\(\)](#), [new\\_nmr\\_dataset\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [print.nmr\\_dataset\\_peak\\_table\(\)](#), [print.nmr\\_dataset\(\)](#), [validate\\_nmr\\_dataset\\_family\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

Other peak detection functions: [Pipelines](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_identify\\_regions\\_blood\(\)](#), [nmr\\_identify\\_regions\\_cell\(\)](#), [nmr\\_identify\\_regions\\_urine\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [regions\\_from\\_peak\\_table\(\)](#)

Other [nmr\\_dataset\\_1D](#) functions: [\[.nmr\\_dataset\\_1D\(\)](#), [computes\\_peak\\_width\\_ppm\(\)](#), [file\\_lister\(\)](#), [files\\_to\\_rDolphin\(\)](#), [format.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [nmr\\_align\\_find\\_ref\(\)](#), [nmr\\_baseline\\_removal\(\)](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_exclude\\_region\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#),

nmr\_meta\_export(), nmr\_meta\_get\_column(), nmr\_meta\_get(), nmr\_normalize(), nmr\_pca\_build\_model(), nmr\_pca\_outliers\_filter(), nmr\_pca\_outliers\_plot(), nmr\_pca\_outliers\_robust(), nmr\_pca\_outliers(), nmr\_ppm\_resolution(), plot.nmr\_dataset\_1D(), plot\_webgl(), print.nmr\_dataset\_1D(), rdCV\_PLS\_RF\_ML(), rdCV\_PLS\_RF(), save\_files\_to\_rDolphin(), to\_ChemoSpec(), validate\_nmr\_dataset()

Other peak detection functions: Pipelines, nmr\_baseline\_threshold(), nmr\_identify\_regions\_blood(), nmr\_identify\_regions\_cell(), nmr\_identify\_regions\_urine(), nmr\_integrate\_regions(), regions\_from\_peak\_table()

Other peak detection functions: Pipelines, nmr\_baseline\_threshold(), nmr\_identify\_regions\_blood(), nmr\_identify\_regions\_cell(), nmr\_identify\_regions\_urine(), nmr\_integrate\_regions(), regions\_from\_peak\_table()

Other nmr\_dataset\_1D functions: [.nmr\_dataset\_1D(), computes\_peak\_width\_ppm(), file\_lister(), files\_to\_rDolphin(), format.nmr\_dataset\_1D(), is.nmr\_dataset\_1D(), load\_and\_save\_functions, new\_nmr\_dataset\_1D(), nmr\_align\_find\_ref(), nmr\_baseline\_removal(), nmr\_baseline\_threshold(), nmr\_exclude\_region(), nmr\_integrate\_regions(), nmr\_interpolate\_1D(), nmr\_meta\_add(), nmr\_meta\_export(), nmr\_meta\_get\_column(), nmr\_meta\_get(), nmr\_normalize(), nmr\_pca\_build\_model(), nmr\_pca\_outliers\_filter(), nmr\_pca\_outliers\_plot(), nmr\_pca\_outliers\_robust(), nmr\_pca\_outliers(), nmr\_ppm\_resolution(), plot.nmr\_dataset\_1D(), plot\_webgl(), print.nmr\_dataset\_1D(), rdCV\_PLS\_RF\_ML(), rdCV\_PLS\_RF(), save\_files\_to\_rDolphin(), to\_ChemoSpec(), validate\_nmr\_dataset()

Other peak detection functions: Pipelines, nmr\_baseline\_threshold(), nmr\_identify\_regions\_blood(), nmr\_identify\_regions\_cell(), nmr\_identify\_regions\_urine(), nmr\_integrate\_regions(), regions\_from\_peak\_table()

Other nmr\_dataset\_1D functions: [.nmr\_dataset\_1D(), computes\_peak\_width\_ppm(), file\_lister(), files\_to\_rDolphin(), format.nmr\_dataset\_1D(), is.nmr\_dataset\_1D(), load\_and\_save\_functions, new\_nmr\_dataset\_1D(), nmr\_align\_find\_ref(), nmr\_baseline\_removal(), nmr\_baseline\_threshold(), nmr\_exclude\_region(), nmr\_integrate\_regions(), nmr\_interpolate\_1D(), nmr\_meta\_add(), nmr\_meta\_export(), nmr\_meta\_get\_column(), nmr\_meta\_get(), nmr\_normalize(), nmr\_pca\_build\_model(), nmr\_pca\_outliers\_filter(), nmr\_pca\_outliers\_plot(), nmr\_pca\_outliers\_robust(), nmr\_pca\_outliers(), nmr\_ppm\_resolution(), plot.nmr\_dataset\_1D(), plot\_webgl(), print.nmr\_dataset\_1D(), rdCV\_PLS\_RF\_ML(), rdCV\_PLS\_RF(), save\_files\_to\_rDolphin(), to\_ChemoSpec(), validate\_nmr\_dataset()

Other alignment functions: Pipelines, nmr\_align\_find\_ref()

Other peak alignment functions: nmr\_align\_find\_ref()

Other nmr\_dataset\_1D functions: [.nmr\_dataset\_1D(), computes\_peak\_width\_ppm(), file\_lister(), files\_to\_rDolphin(), format.nmr\_dataset\_1D(), is.nmr\_dataset\_1D(), load\_and\_save\_functions, new\_nmr\_dataset\_1D(), nmr\_align\_find\_ref(), nmr\_baseline\_removal(), nmr\_baseline\_threshold(), nmr\_exclude\_region(), nmr\_integrate\_regions(), nmr\_interpolate\_1D(), nmr\_meta\_add(), nmr\_meta\_export(), nmr\_meta\_get\_column(), nmr\_meta\_get(), nmr\_normalize(), nmr\_pca\_build\_model(), nmr\_pca\_outliers\_filter(), nmr\_pca\_outliers\_plot(), nmr\_pca\_outliers\_robust(), nmr\_pca\_outliers(), nmr\_ppm\_resolution(), plot.nmr\_dataset\_1D(), plot\_webgl(), print.nmr\_dataset\_1D(), rdCV\_PLS\_RF\_ML(), rdCV\_PLS\_RF(), save\_files\_to\_rDolphin(), to\_ChemoSpec(), validate\_nmr\_dataset()

Other peak integration functions: Pipelines, computes\_peak\_width\_ppm(), nmr\_identify\_regions\_blood(), nmr\_identify\_regions\_cell(), nmr\_identify\_regions\_urine(), nmr\_integrate\_regions()

Other nmr\_dataset\_1D functions: [.nmr\_dataset\_1D(), computes\_peak\_width\_ppm(), file\_lister(), files\_to\_rDolphin(), format.nmr\_dataset\_1D(), is.nmr\_dataset\_1D(), load\_and\_save\_functions, new\_nmr\_dataset\_1D(), nmr\_align\_find\_ref(), nmr\_baseline\_removal(), nmr\_baseline\_threshold(), nmr\_exclude\_region(), nmr\_integrate\_regions(), nmr\_interpolate\_1D(), nmr\_meta\_add(), nmr\_meta\_export(), nmr\_meta\_get\_column(), nmr\_meta\_get(), nmr\_normalize(), nmr\_pca\_build\_model(), nmr\_pca\_outliers\_filter(), nmr\_pca\_outliers\_plot(), nmr\_pca\_outliers\_robust(), nmr\_pca\_outliers(), nmr\_ppm\_resolution(), plot.nmr\_dataset\_1D(), plot\_webgl(), print.nmr\_dataset\_1D(), rdCV\_PLS\_RF\_ML(), rdCV\_PLS\_RF(), save\_files\_to\_rDolphin(), to\_ChemoSpec(), validate\_nmr\_dataset()

Other peak integration functions: `Pipelines`, `computes_peak_width_ppm()`, `nmr_identify_regions_blood()`, `nmr_identify_regions_cell()`, `nmr_identify_regions_urine()`, `nmr_integrate_regions()`

Other `nmr_dataset_1D` functions: `[.nmr_dataset_1D()`, `computes_peak_width_ppm()`, `file_lister()`, `files_to_rDolphin()`, `format.nmr_dataset_1D()`, `is.nmr_dataset_1D()`, `load_and_save_functions`, `new_nmr_dataset_1D()`, `nmr_align_find_ref()`, `nmr_baseline_removal()`, `nmr_baseline_threshold()`, `nmr_exclude_region()`, `nmr_integrate_regions()`, `nmr_interpolate_1D()`, `nmr_meta_add()`, `nmr_meta_export()`, `nmr_meta_get_column()`, `nmr_meta_get()`, `nmr_normalize()`, `nmr_pca_build_model()`, `nmr_pca_outliers_filter()`, `nmr_pca_outliers_plot()`, `nmr_pca_outliers_robust()`, `nmr_pca_outliers()`, `nmr_ppm_resolution()`, `plot.nmr_dataset_1D()`, `plot_webgl()`, `print.nmr_dataset_1D()`, `rdCV_PLS_RF_ML()`, `rdCV_PLS_RF()`, `save_files_to_rDolphin()`, `to_ChemoSpec()`, `validate_nmr_dataset()`

## Examples

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")
nmr_dataset <- nmr_read_samples_dir(dir_to_demo_dataset)
dataset_1D <- nmr_interpolate_1D(nmr_dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))

sample_10 <- filter(dataset_1D, NMRExperiment == "10")

# 1. Peak detection in the dataset.
peak_data <- nmr_detect_peaks(dataset_1D,
                             nDivRange_ppm = 0.1, # Size of detection segments
                             scales = seq(1, 16, 2),
                             baselineThresh = 0, # Minimum peak intensity
                             SNR.Th = 4) # Signal to noise ratio

#nmr_detect_peaks_plot(sample_10, peak_data, "NMRExp_ref")

peaks_detected <- nmr_detect_peaks_tune_snr(sample_10,
                                             SNR_thresholds = seq(from = 2,
                                                                    to = 3, by = 0.5))

# 2. Find the reference spectrum to align with.
NMRExp_ref <- nmr_align_find_ref(dataset_1D, peak_data)

# 3. Spectra alignment using the ref spectrum and a maximum alignment shift
nmr_dataset <- nmr_align(dataset_1D, # the dataset
                        peak_data, # detected peaks
                        NMRExp_ref = NMRExp_ref, # ref spectrum
                        maxShift_ppm = 0.0015, # max alignment shift
                        acceptLostPeak = FALSE) # lost peaks

# 4. PEAK INTEGRATION (please, consider previous normalization step).
# First we take the peak table from the reference spectrum
peak_data_ref <- filter(peak_data, NMRExperiment == NMRExp_ref)

# Then we integrate spectra considering the peaks from the ref spectrum
nmr_peak_table <- nmr_integrate_peak_positions(
  samples = nmr_dataset,
  peak_pos_ppm = peak_data_ref$ppm,
  peak_width_ppm = NULL)

validate_nmr_dataset_peak_table(nmr_peak_table)

#If you wanted the final peak table before machine learning you can run
```

```
nmr_peak_table_completed <- get_integration_with_metadata(nmr_peak_table)
```

---

[.nmr\_dataset]      *Extract parts of an nmr\_dataset*

---

## Description

Extract parts of an nmr\_dataset

## Usage

```
## S3 method for class 'nmr_dataset'  
x[i]
```

## Arguments

x                    an [nmr\\_dataset](#) object  
i                    indices of the samples to keep

## Value

an nmr\_dataset with the extracted samples

## See Also

Other subsetting functions: [\[.nmr\\_dataset\\_1D\(\)](#), [\[.nmr\\_dataset\\_peak\\_table\(\)](#), [filter.nmr\\_dataset\\_family\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#)

Other nmr\_dataset functions: [format.nmr\\_dataset\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column\(\)](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [nmr\\_read\\_samples\(\)](#), [print.nmr\\_dataset\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

## Examples

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")  
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)  
dataset2 <- dataset[1:3] # get the first 3 samples
```

---

[.nmr\_dataset\_1D      *Extract parts of an nmr\_dataset\_1D*

---

## Description

Extract parts of an nmr\_dataset\_1D

## Usage

```
## S3 method for class 'nmr_dataset_1D'  
x[i]
```

## Arguments

x                    an [nmr\\_dataset\\_1D](#) object  
i                    indices of the samples to keep

## Value

an nmr\_dataset\_1D with the extracted samples

## See Also

Other subsetting functions: [\[.nmr\\_dataset\\_peak\\_table\(\)](#), [\[.nmr\\_dataset\(\)](#), [filter.nmr\\_dataset\\_family\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#)

Other nmr\_dataset\_1D functions: [computes\\_peak\\_width\\_ppm\(\)](#), [file\\_lister\(\)](#), [files\\_to\\_rDolphin\(\)](#), [format.nmr\\_dataset\\_1D\(\)](#), [is.nmr\\_dataset\\_1D\(\)](#), [load\\_and\\_save\\_functions](#), [new\\_nmr\\_dataset\\_1D\(\)](#), [nmr\\_align\\_find\\_ref\(\)](#), [nmr\\_baseline\\_removal\(\)](#), [nmr\\_baseline\\_threshold\(\)](#), [nmr\\_exclude\\_region\(\)](#), [nmr\\_integrate\\_regions\(\)](#), [nmr\\_interpolate\\_1D\(\)](#), [nmr\\_meta\\_add\(\)](#), [nmr\\_meta\\_export\(\)](#), [nmr\\_meta\\_get\\_column](#), [nmr\\_meta\\_get\(\)](#), [nmr\\_normalize\(\)](#), [nmr\\_pca\\_build\\_model\(\)](#), [nmr\\_pca\\_outliers\\_filter\(\)](#), [nmr\\_pca\\_outliers\\_plot\(\)](#), [nmr\\_pca\\_outliers\\_robust\(\)](#), [nmr\\_pca\\_outliers\(\)](#), [nmr\\_ppm\\_resolution\(\)](#), [plot.nmr\\_dataset\\_1D\(\)](#), [plot\\_webgl\(\)](#), [print.nmr\\_dataset\\_1D\(\)](#), [rdCV\\_PLS\\_RF\\_ML\(\)](#), [rdCV\\_PLS\\_RF\(\)](#), [save\\_files\\_to\\_rDolphin\(\)](#), [to\\_ChemoSpec\(\)](#), [validate\\_nmr\\_dataset\\_peak\\_table\(\)](#), [validate\\_nmr\\_dataset\(\)](#)

## Examples

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")  
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)  
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))  
dataset_1D[0]
```

---

`[.nmr_dataset_peak_table`*Extract parts of an nmr\_dataset\_peak\_table*

---

### Description

Extract parts of an `nmr_dataset_peak_table`

### Usage

```
## S3 method for class 'nmr_dataset_peak_table'  
x[i]
```

### Arguments

`x` an `nmr_dataset_peak_table` object  
`i` indices of the samples to keep

### Value

an `nmr_dataset_peak_table` with the extracted samples

### See Also

Other subsetting functions: `[.nmr_dataset_1D()`, `[.nmr_dataset()`, `filter.nmr_dataset_family()`, `nmr_pca_outliers_filter()`

Other `nmr_dataset_peak_table` functions: `format.nmr_dataset_peak_table()`, `is.nmr_dataset_peak_table()`, `load_and_save_functions`, `new_nmr_dataset_peak_table()`, `nmr_meta_add()`, `nmr_meta_export()`, `nmr_meta_get_column()`, `nmr_meta_get()`, `print.nmr_dataset_peak_table()`, `validate_nmr_dataset_peak_t`

### Examples

```
dir_to_demo_dataset <- system.file("dataset-demo", package = "AlpsNMR")  
dataset <- nmr_read_samples_dir(dir_to_demo_dataset)  
dataset_1D <- nmr_interpolate_1D(dataset, axis = c(min = -0.5, max = 10, by = 2.3E-4))  
meta <- file.path(dir_to_demo_dataset, "dummy_metadata.xlsx")  
metadata <- readxl::read_excel(meta, sheet = 1)  
dataset_1D <- nmr_meta_add(dataset_1D, metadata = metadata, by = "NMRExperiment")  
metadata <- list(external = dataset_1D[["metadata"]][["external"]])  
peak_table <- nmr_data(dataset_1D)  
new <- new_nmr_dataset_peak_table(peak_table, metadata)  
new[0]
```



# Index

- \* **AlpsNMR dataset objects**
  - [nmr\\_dataset](#), [34](#)
  - [nmr\\_dataset\\_1D](#), [34](#)
  - [nmr\\_dataset\\_family](#), [35](#)
- \* **PCA related functions**
  - [nmr\\_pca\\_build\\_model](#), [52](#)
  - [nmr\\_pca\\_outliers](#), [54](#)
  - [nmr\\_pca\\_outliers\\_filter](#), [55](#)
  - [nmr\\_pca\\_outliers\\_plot](#), [56](#)
  - [nmr\\_pca\\_outliers\\_robust](#), [57](#)
  - [nmr\\_pca\\_plots](#), [58](#)
- \* **alignment functions**
  - [nmr\\_align\\_find\\_ref](#), [27](#)
  - [Pipelines](#), [67](#)
  - [validate\\_nmr\\_dataset\\_peak\\_table](#), [97](#)
- \* **baseline removal functions**
  - [nmr\\_baseline\\_removal](#), [28](#)
- \* **batman functions**
  - [nmr\\_batman](#), [30](#)
  - [nmr\\_batman\\_options](#), [31](#)
- \* **class helper functions**
  - [format.nmr\\_dataset](#), [13](#)
  - [format.nmr\\_dataset\\_1D](#), [14](#)
  - [format.nmr\\_dataset\\_peak\\_table](#), [15](#)
  - [is.nmr\\_dataset\\_1D](#), [18](#)
  - [is.nmr\\_dataset\\_peak\\_table](#), [18](#)
  - [new\\_nmr\\_dataset](#), [24](#)
  - [new\\_nmr\\_dataset\\_1D](#), [25](#)
  - [new\\_nmr\\_dataset\\_peak\\_table](#), [26](#)
  - [print.nmr\\_dataset](#), [82](#)
  - [print.nmr\\_dataset\\_1D](#), [83](#)
  - [print.nmr\\_dataset\\_peak\\_table](#), [84](#)
  - [validate\\_nmr\\_dataset](#), [95](#)
  - [validate\\_nmr\\_dataset\\_family](#), [96](#)
  - [validate\\_nmr\\_dataset\\_peak\\_table](#), [97](#)
- \* **data**
  - [hmdb](#), [16](#)
  - [HMDB\\_blood](#), [16](#)
  - [HMDB\\_cell](#), [16](#)
  - [HMDB\\_urine](#), [17](#)
  - [Parameters\\_blood](#), [63](#)
  - [Parameters\\_cell](#), [63](#)
  - [Parameters\\_urine](#), [63](#)
  - [ROI\\_blood](#), [90](#)
  - [ROI\\_cell](#), [91](#)
  - [ROI\\_urine](#), [91](#)
- \* **import/export functions**
  - [files\\_to\\_rDolphin](#), [10](#)
  - [load\\_and\\_save\\_functions](#), [19](#)
  - [nmr\\_data](#), [33](#)
  - [nmr\\_meta\\_export](#), [48](#)
  - [nmr\\_read\\_bruker\\_fid](#), [60](#)
  - [nmr\\_read\\_samples](#), [61](#)
  - [nmr\\_zip\\_bruker\\_samples](#), [62](#)
  - [Pipelines](#), [67](#)
  - [save\\_files\\_to\\_rDolphin](#), [91](#)
  - [save\\_profiling\\_output](#), [92](#)
  - [to\\_ChemoSpec](#), [94](#)
- \* **metadata functions**
  - [nmr\\_meta\\_add](#), [46](#)
  - [nmr\\_meta\\_export](#), [48](#)
  - [nmr\\_meta\\_get](#), [49](#)
  - [nmr\\_meta\\_get\\_column](#), [50](#)
  - [Pipelines](#), [67](#)
- \* **nmr\_dataset functions**
  - [\[.nmr\\_dataset](#), [102](#)
  - [format.nmr\\_dataset](#), [13](#)
  - [load\\_and\\_save\\_functions](#), [19](#)
  - [new\\_nmr\\_dataset](#), [24](#)
  - [nmr\\_interpolate\\_1D](#), [45](#)
  - [nmr\\_meta\\_add](#), [46](#)
  - [nmr\\_meta\\_export](#), [48](#)
  - [nmr\\_meta\\_get](#), [49](#)
  - [nmr\\_meta\\_get\\_column](#), [50](#)
  - [nmr\\_ppm\\_resolution](#), [59](#)
  - [nmr\\_read\\_samples](#), [61](#)
  - [print.nmr\\_dataset](#), [82](#)
  - [validate\\_nmr\\_dataset](#), [95](#)
- \* **nmr\_dataset manipulation functions**
  - [is.nmr\\_dataset](#), [17](#)
- \* **nmr\_dataset\_1D functions**
  - [\[.nmr\\_dataset\\_1D](#), [103](#)
  - [computes\\_peak\\_width\\_ppm](#), [9](#)
  - [file\\_lister](#), [12](#)

- files\_to\_rDolphin, 10
- format.nmr\_dataset\_1D, 14
- is.nmr\_dataset\_1D, 18
- load\_and\_save\_functions, 19
- new\_nmr\_dataset\_1D, 25
- nmr\_align\_find\_ref, 27
- nmr\_baseline\_removal, 28
- nmr\_baseline\_threshold, 29
- nmr\_exclude\_region, 40
- nmr\_integrate\_regions, 44
- nmr\_interpolate\_1D, 45
- nmr\_meta\_add, 46
- nmr\_meta\_export, 48
- nmr\_meta\_get, 49
- nmr\_meta\_get\_column, 50
- nmr\_normalize, 51
- nmr\_pca\_build\_model, 52
- nmr\_pca\_outliers, 54
- nmr\_pca\_outliers\_filter, 55
- nmr\_pca\_outliers\_plot, 56
- nmr\_pca\_outliers\_robust, 57
- nmr\_ppm\_resolution, 59
- plot.nmr\_dataset\_1D, 71
- plot\_webgl, 78
- print.nmr\_dataset\_1D, 83
- rdCV\_PLS\_RF, 87
- rdCV\_PLS\_RF\_ML, 88
- save\_files\_to\_rDolphin, 91
- to\_ChemoSpec, 94
- validate\_nmr\_dataset, 95
- validate\_nmr\_dataset\_peak\_table, 97
- \* **nmr\_dataset\_family functions**
  - validate\_nmr\_dataset\_family, 96
- \* **nmr\_dataset\_family manipulation functions**
  - filter.nmr\_dataset\_family, 13
- \* **nmr\_dataset\_peak\_table functions**
  - [.nmr\_dataset\_peak\_table, 104
  - format.nmr\_dataset\_peak\_table, 15
  - is.nmr\_dataset\_peak\_table, 18
  - load\_and\_save\_functions, 19
  - new\_nmr\_dataset\_peak\_table, 26
  - nmr\_meta\_add, 46
  - nmr\_meta\_export, 48
  - nmr\_meta\_get, 49
  - nmr\_meta\_get\_column, 50
  - print.nmr\_dataset\_peak\_table, 84
  - validate\_nmr\_dataset\_peak\_table, 97
- \* **outlier detection functions**
  - nmr\_pca\_outliers, 54
  - nmr\_pca\_outliers\_filter, 55
  - nmr\_pca\_outliers\_plot, 56
  - nmr\_pca\_outliers\_robust, 57
  - Pipelines, 67
- \* **peak alignment functions**
  - nmr\_align\_find\_ref, 27
  - validate\_nmr\_dataset\_peak\_table, 97
- \* **peak detection functions**
  - nmr\_baseline\_threshold, 29
  - nmr\_identify\_regions\_blood, 41
  - nmr\_identify\_regions\_cell, 42
  - nmr\_identify\_regions\_urine, 43
  - nmr\_integrate\_regions, 44
  - Pipelines, 67
  - regions\_from\_peak\_table, 90
  - validate\_nmr\_dataset\_peak\_table, 97
- \* **peak integration functions**
  - computes\_peak\_width\_ppm, 9
  - nmr\_identify\_regions\_blood, 41
  - nmr\_identify\_regions\_cell, 42
  - nmr\_identify\_regions\_urine, 43
  - nmr\_integrate\_regions, 44
  - Pipelines, 67
  - validate\_nmr\_dataset\_peak\_table, 97
- \* **pipeline functions**
  - Pipelines, 67
- \* **plotting functions**
  - plot.nmr\_dataset\_1D, 71
  - plot\_interactive, 73
- \* **plotting nmr datasets**
  - plot\_webgl, 78
- \* **subsetting functions**
  - [.nmr\_dataset, 102
  - [.nmr\_dataset\_1D, 103
  - [.nmr\_dataset\_peak\_table, 104
  - filter.nmr\_dataset\_family, 13
  - nmr\_pca\_outliers\_filter, 55
- \* **to\_rDolphin functions**
  - files\_to\_rDolphin, 10
- \* **to\_rDolphin\_blood functions**
  - save\_files\_to\_rDolphin, 91
  - save\_profiling\_output, 92
- [.nmr\_dataset, 13, 14, 20, 24, 46, 47, 49–51, 56, 60, 61, 83, 96, 102, 103, 104
- [.nmr\_dataset\_1D, 9, 11–13, 15, 18, 20, 25, 27–29, 40, 45–47, 49–53, 55–58, 60, 71, 78, 84, 88, 89, 92, 95, 96, 99–102, 103, 104
- [.nmr\_dataset\_peak\_table, 13, 15, 19, 20,

- 26, 47, 49–51, 56, 84, 99, 102, 103, 104
- AlpsNMR (AlpsNMR-package), 4
- AlpsNMR-package, 4
- AUC\_model, 5
- baseline::baseline.als, 28
- bp\_kfold\_VIP\_analysis, 5
- bp\_VIP\_analysis, 6, 7
- computes\_peak\_width\_ppm, 9, 11, 12, 15, 18, 20, 25, 27–29, 40, 42, 43, 45–47, 49–53, 55–58, 60, 69, 71, 78, 84, 88, 89, 92, 95, 96, 99–101, 103
- confusion\_matrix, 10
- dplyr, 13
- file\_list, 9, 11, 12, 15, 18, 20, 25, 27–29, 40, 45–47, 49–53, 55–58, 60, 71, 78, 84, 88, 89, 92, 95, 96, 99–101, 103
- files\_to\_rDolphin, 9, 10, 12, 15, 18, 20, 25, 27–29, 33, 40, 45–47, 49–53, 55–58, 60–62, 69, 71, 78, 84, 88, 89, 92, 93, 95, 96, 99–101, 103
- filter.nmr\_dataset\_family, 13, 56, 102–104
- format.nmr\_dataset, 13, 15, 18–20, 24–26, 46, 47, 49–51, 60, 61, 83, 84, 96, 97, 99, 102
- format.nmr\_dataset\_1D, 9, 11, 12, 14, 14, 15, 18–20, 24–29, 40, 45–47, 49–53, 55–58, 60, 71, 78, 83, 84, 88, 89, 92, 95–97, 99–101, 103
- format.nmr\_dataset\_peak\_table, 14, 15, 15, 18–20, 24–26, 47, 49–51, 83, 84, 96, 97, 99, 104
- fs::dir\_ls(), 12
- Functions to save and load these objects, 34, 35
- get\_integration\_with\_metada  
(validate\_nmr\_dataset\_peak\_table), 97
- get\_integration\_with\_metadadata  
(validate\_nmr\_dataset\_peak\_table), 97
- ggplot2::aes, 58
- ggplot2::aes\_string, 71
- ggplot2::aes\_string(), 56
- grep(), 68
- hmdb, 16, 30
- HMDB\_blood, 16
- HMDB\_cell, 16
- HMDB\_urine, 17
- is.nmr\_dataset, 17
- is.nmr\_dataset\_1D, 9, 11, 12, 14, 15, 18, 19, 20, 24–29, 40, 45–47, 49–53, 55–58, 60, 71, 78, 83, 84, 88, 89, 92, 95–97, 99–101, 103
- is.nmr\_dataset\_peak\_table, 14, 15, 18, 18, 20, 24–26, 47, 49–51, 83, 84, 96, 97, 99, 104
- load\_and\_save\_functions, 9, 11, 12, 14, 15, 18, 19, 19, 24–29, 33, 40, 45–47, 49–53, 55–58, 60–62, 69, 71, 78, 83, 84, 88, 89, 92, 93, 95, 96, 99–104
- MassSpecWavelet::peakDetectionCWT, 97
- mixOmics::pca, 53
- mixOmics::plsda, 39
- model\_VIP, 23
- models\_stability\_plot\_bootstrap, 20
- models\_stability\_plot\_plsda, 22
- MUVR\_model\_plot, 23
- new\_nmr\_data\_analysis\_method  
(nmr\_data\_analysis\_method), 39
- new\_nmr\_data\_analysis\_method(), 37
- new\_nmr\_dataset, 14, 15, 18–20, 24, 25, 26, 46, 47, 49–51, 60, 61, 83, 84, 96, 97, 99, 102
- new\_nmr\_dataset\_1D, 9, 11, 12, 14, 15, 18–20, 24, 25, 26–29, 40, 45–47, 49–53, 55–58, 60, 71, 78, 83, 84, 88, 89, 92, 95–97, 99–101, 103
- new\_nmr\_dataset\_peak\_table, 14, 15, 18–20, 24, 25, 26, 47, 49–51, 83, 84, 96, 97, 99, 104
- nmr\_align, 99
- nmr\_align  
(validate\_nmr\_dataset\_peak\_table), 97
- nmr\_align\_find\_ref, 9, 11, 12, 15, 18, 20, 25, 27, 28, 29, 40, 45–47, 49–53, 55–58, 60, 69, 71, 78, 84, 88, 89, 92, 95, 96, 99–101, 103
- nmr\_baseline\_removal, 9, 11, 12, 15, 18, 20, 25, 27, 28, 29, 40, 45–47, 49–53, 55–58, 60, 71, 78, 84, 88, 89, 92, 95, 96, 99–101, 103
- nmr\_baseline\_threshold, 9, 11, 12, 15, 18, 20, 25, 27, 28, 29, 40, 42–47, 49–53,

- 55–58, 60, 69, 71, 78, 84, 88–90, 92, 95, 96, 99–101, 103
- `nmr_batman`, 30, 33
- `nmr_batman_export_dataset` (`nmr_batman`), 30
- `nmr_batman_metabolites_list` (`nmr_batman`), 30
- `nmr_batman_multi_data_user` (`nmr_batman`), 30
- `nmr_batman_multi_data_user_hmndb` (`nmr_batman`), 30
- `nmr_batman_options`, 31, 31
- `nmr_batman_write_options` (`nmr_batman`), 30
- `nmr_data`, 11, 20, 33, 49, 60–62, 69, 92, 93, 95
- `nmr_data<-` (`nmr_data`), 33
- `nmr_data_1r_to_SummarizedExperiment`, 36
- `nmr_data_analysis`, 37, 40, 79, 81
- `nmr_data_analysis_method`, 37, 39, 64
- `nmr_dataset`, 4, 10, 14, 17, 34, 34, 35, 44, 61, 83, 98, 99, 102
- `nmr_dataset_1D`, 5, 7, 9, 12, 14, 18, 27–30, 34, 34, 35–37, 41, 51–56, 58, 71, 78, 83, 94, 95, 97–99, 103
- `nmr_dataset_family`, 6, 7, 13, 20, 33, 34, 35, 37, 47–50, 64, 72, 96, 97
- `nmr_dataset_load` (`load_and_save_functions`), 19
- `nmr_dataset_peak_table`, 5, 7, 15, 18, 19, 35, 35, 37, 44, 84, 93, 98, 99, 104
- `nmr_dataset_peak_table_to_SummarizedExperiment`, 35
- `nmr_dataset_save` (`load_and_save_functions`), 19
- `nmr_detect_peaks`, 27, 98, 99
- `nmr_detect_peaks` (`validate_nmr_dataset_peak_table`), 97
- `nmr_detect_peaks_plot` (`validate_nmr_dataset_peak_table`), 97
- `nmr_detect_peaks_tune_snr` (`validate_nmr_dataset_peak_table`), 97
- `nmr_exclude_region`, 9, 11, 12, 15, 18, 20, 25, 27–29, 40, 45–47, 49–52, 54–58, 60, 71, 78, 84, 88, 89, 92, 95, 96, 99–101, 103
- `nmr_exclude_region()`, 4
- `nmr_export_data_1r`, 41
- `nmr_identify_regions_blood`, 9, 29, 41, 42–45, 69, 90, 99–101
- `nmr_identify_regions_cell`, 9, 29, 42, 42, 43–45, 69, 90, 99–101
- `nmr_identify_regions_urine`, 9, 29, 42, 43, 44, 45, 69, 90, 99–101
- `nmr_integrate_peak_positions` (`validate_nmr_dataset_peak_table`), 97
- `nmr_integrate_regions`, 9, 11, 12, 15, 18, 20, 25, 27–29, 40, 42, 43, 44, 46, 47, 49–52, 54–58, 60, 69, 71, 78, 84, 88–90, 92, 95, 96, 98–101, 103
- `nmr_interpolate_1D`, 9, 11, 12, 14, 15, 18, 20, 24, 25, 27–29, 40, 45, 45, 47, 49–52, 54–58, 60, 61, 71, 78, 83, 84, 88, 89, 92, 95, 96, 99–103
- `nmr_interpolate_1D()`, 4
- `nmr_meta_add`, 9, 11, 12, 14, 15, 18–20, 24–29, 40, 45, 46, 46, 49–52, 54–58, 60, 61, 69, 71, 78, 83, 84, 88, 89, 92, 95, 96, 99–104
- `nmr_meta_add_tidy_excel` (`nmr_meta_add`), 46
- `nmr_meta_export`, 9, 11, 12, 14, 15, 18–20, 24–29, 33, 40, 45–47, 48, 50–52, 54–58, 60–62, 69, 71, 78, 83, 84, 88, 89, 92, 93, 95, 96, 99–104
- `nmr_meta_get`, 9, 11, 12, 14, 15, 18–20, 24–29, 40, 45–47, 49, 49, 51, 52, 54–58, 60, 61, 69, 71, 78, 83, 84, 88, 89, 92, 95, 96, 99–104
- `nmr_meta_get_column`, 9, 11, 12, 14, 15, 18–20, 24–29, 40, 45–47, 49, 50, 50, 52, 54–58, 60, 61, 69, 71, 78, 83, 84, 88, 89, 92, 95, 96, 99–104
- `nmr_normalize`, 9, 11, 12, 15, 18, 20, 25, 27–29, 40, 45–47, 49–51, 51, 54–58, 60, 71, 78, 84, 88, 89, 92, 95, 96, 100, 101, 103
- `nmr_normalize()`, 4
- `nmr_normalize_extra_info` (`nmr_normalize`), 51
- `nmr_pca_build_model`, 9, 11, 12, 15, 18, 20, 25, 27–29, 40, 45–47, 49–52, 52, 54–60, 71, 78, 84, 88, 89, 92, 95, 96, 100, 101, 103
- `nmr_pca_loadingplot` (`nmr_pca_plots`), 58
- `nmr_pca_outliers`, 9, 11, 12, 15, 18, 20, 25, 27–29, 40, 45–47, 49–54, 54, 56–60, 69, 71, 78, 84, 88, 89, 92, 95, 96, 100, 101, 103
- `nmr_pca_outliers()`, 55, 56

- nmr\_pca\_outliers\_filter, 9, 11–13, 15, 18, 20, 25, 27–29, 40, 45–47, 49–55, 55, 57–60, 69, 71, 78, 84, 88, 89, 92, 95, 96, 100–104  
 nmr\_pca\_outliers\_plot, 9, 11, 12, 15, 18, 20, 25, 27–29, 40, 45–47, 49–56, 56, 58–60, 69, 71, 78, 84, 88, 89, 92, 95, 96, 100, 101, 103  
 nmr\_pca\_outliers\_plot(), 57  
 nmr\_pca\_outliers\_robust, 9, 11, 12, 15, 18, 20, 25, 27–29, 40, 45–47, 49–57, 57, 59, 60, 67, 69, 71, 78, 84, 88, 89, 92, 95, 96, 100, 101, 103  
 nmr\_pca\_plot\_variance (nmr\_pca\_plots), 58  
 nmr\_pca\_plots, 53, 55–58, 58  
 nmr\_pca\_scoreplot (nmr\_pca\_plots), 58  
 nmr\_ppm\_resolution, 9, 11, 12, 14, 15, 18, 20, 24, 25, 27–29, 40, 45–47, 49–52, 54–58, 59, 61, 71, 78, 83, 84, 88, 89, 92, 95, 96, 100–103  
 nmr\_read\_bruker\_fid, 11, 20, 33, 49, 60, 61, 62, 69, 92, 93, 95  
 nmr\_read\_samples, 11, 14, 20, 24, 33, 46, 47, 49–51, 60, 61, 62, 69, 83, 92, 93, 95, 96, 102  
 nmr\_read\_samples(), 4  
 nmr\_read\_samples\_dir (nmr\_read\_samples), 61  
 nmr\_read\_samples\_dir(), 4  
 nmr\_zip\_bruker\_samples, 11, 20, 33, 49, 60, 61, 62, 69, 92, 93, 95  
  
 p\_value\_perm, 85  
 Parameters\_blood, 63  
 Parameters\_cell, 63  
 Parameters\_urine, 63  
 Peak\_detection (validate\_nmr\_dataset\_peak\_table), 97  
 permutation\_test\_model, 64  
 permutation\_test\_plot, 65  
 pipe\_add\_metadata (Pipelines), 67  
 pipe\_exclude\_regions (Pipelines), 67  
 pipe\_filter\_samples (Pipelines), 67  
 pipe\_interpolate\_1D (Pipelines), 67  
 pipe\_load\_samples (Pipelines), 67  
 Pipe\_normalization (Pipelines), 67  
 pipe\_normalization (Pipelines), 67  
 pipe\_outlier\_detection (Pipelines), 67  
 pipe\_pakdet\_align (Pipelines), 67  
 pipe\_peak\_integration (Pipelines), 67  
 pipe\_peakdet\_align, 69  
 pipe\_peakdet\_align (Pipelines), 67  
 Pipelines, 9, 11, 20, 27, 29, 33, 42–45, 47, 49–51, 55–58, 60–62, 67, 90, 92, 93, 95, 99–101  
 plot(), 4  
 plot\_nmr\_dataset\_1D, 10–12, 15, 18, 20, 25, 27–29, 40, 45–47, 49–52, 54–58, 60, 71, 73, 78, 84, 88, 89, 92, 95, 96, 100, 101, 103  
 plot\_bootstrap\_multimodel, 72  
 plot\_interactive, 71, 73, 99  
 plot\_plsda\_multimodel, 74  
 plot\_plsda\_samples, 75  
 plot\_vip\_scores, 76  
 plot\_webgl, 10–12, 15, 18, 20, 25, 27–29, 40, 45–47, 49–52, 54–58, 60, 71, 78, 84, 88, 89, 92, 95, 96, 100, 101, 103  
 plotly::ggplotly, 99  
 plsda\_aucroc\_vip\_compare, 79  
 plsda\_aucroc\_vip\_method, 37, 80  
 ppm\_resolution, 81  
 ppm\_VIP\_vector, 81  
 print.nmr\_dataset, 14, 15, 18–20, 24–26, 46, 47, 49–51, 60, 61, 82, 84, 96, 97, 99, 102  
 print.nmr\_dataset\_1D, 10–12, 14, 15, 18–20, 24–29, 40, 45–47, 49–52, 54–58, 60, 71, 78, 83, 83, 84, 88, 89, 92, 95–97, 99–101, 103  
 print.nmr\_dataset\_peak\_table, 14, 15, 18–20, 24–26, 47, 49–51, 83, 84, 84, 96, 97, 99, 104  
  
 random\_subsampling, 37, 64, 86  
 rdCV\_PLS\_RF, 10–12, 15, 18, 20, 25, 27–29, 40, 45–47, 49–52, 54–58, 60, 71, 78, 84, 87, 89, 92, 95, 96, 100, 101, 103  
 rdCV\_PLS\_RF\_ML, 10–12, 15, 18, 20, 25, 27–29, 40, 45–47, 49–52, 54–58, 60, 71, 78, 84, 88, 88, 92, 95, 96, 100, 101, 103  
 read\_bruker\_sample, 89  
 read\_bruker\_sample(), 61  
 regions\_from\_peak\_table, 29, 42–44, 69, 90, 99, 100  
 ROI\_blood, 90  
 ROI\_cell, 91  
 ROI\_urine, 91  
  
 save\_files\_to\_rDolphin, 10–12, 15, 18, 20, 25, 27–29, 33, 40, 45–47, 49–52, 54–58, 60–62, 69, 71, 78, 84, 88, 89, 91, 93, 95, 96, 100, 101, 103

save\_profiling\_output, [11](#), [20](#), [33](#), [49](#),  
[60–62](#), [69](#), [92](#), [92](#), [95](#)  
saveRDS, [20](#)  
scale, [53](#)  
speaq::detectSpecPeaks, [97](#)  
speaq::dohCluster, [97](#)  
SummarizedExperiment\_to\_nmr\_data\_1r,  
[94](#)  
SummarizedExperiment\_to\_nmr\_dataset\_peak\_table,  
[93](#)

to\_ChemoSpec, [10–12](#), [15](#), [18](#), [20](#), [25](#), [27–29](#),  
[33](#), [40](#), [45–47](#), [49–52](#), [54–58](#), [60–62](#),  
[69](#), [71](#), [78](#), [84](#), [88](#), [89](#), [92](#), [93](#), [94](#), [96](#),  
[100](#), [101](#), [103](#)

utils::zip, [62](#)

validate\_nmr\_dataset, [10–12](#), [14](#), [15](#),  
[18–20](#), [24–29](#), [40](#), [45–47](#), [49–52](#),  
[54–58](#), [60](#), [61](#), [71](#), [78](#), [83](#), [84](#), [88](#), [89](#),  
[92](#), [95](#), [95](#), [97](#), [99–103](#)  
validate\_nmr\_dataset\_1D  
(validate\_nmr\_dataset), [95](#)  
validate\_nmr\_dataset\_family, [14](#), [15](#), [18](#),  
[19](#), [24–26](#), [83](#), [84](#), [96](#), [96](#), [99](#)  
validate\_nmr\_dataset\_peak\_table, [9–12](#),  
[14](#), [15](#), [18–20](#), [24–29](#), [40](#), [42–47](#),  
[49–52](#), [54–58](#), [60](#), [69](#), [71](#), [78](#), [83](#), [84](#),  
[88–90](#), [92](#), [95–97](#), [97](#), [103](#), [104](#)