# affy: Custom Processing Methods (HowTo)

#### Laurent

#### April 13, 2011

#### Contents

1	Introduction	1
2	How-to	1
3	Examples	3

## 1 Introduction

This document describes briefly how to customize the affy package by adding one's own processing methods. The types of processing methods are background correction, normalization, perfect match correction and summary expression value computation. We tried our best to make this as easy as we could, but we are aware that it is far from being perfect. We are still working on things to improve them. Hopefully this document should let you extend the package with supplementary processing methods easily.

As usual, loading the package in your R session is required.

R> library(affy) ##load the affy package

#### 2 How-to

For each processing step, labels for the methods known to the package are stored in variables.

> normalize.AffyBatch.methods()

[1]	"constant"	"contrasts"	"invariantset"	"loess"
[5]	"methods"	"qspline"	"quantiles"	"quantiles.robust"

> bgcorrect.methods()

variable for labels	naming convention
bgcorrect.methods	bg.correct. <label></label>
normalize.AffyBatch.methods	normalize.AffyBatch. <label></label>
pmcorrect.methods	pmcorrect. < label >
express.summary.stat.methods	generateExprset.methods. <label></label>

Table 1: Summary table for the processing methods.

step	$\operatorname{argument}(s)$	returns
background correction	AffyBatch	AffyBatch
normalization	AffyBatch	AffyBatch
pm correction	ProbeSet	a matrix of corrected PM
		values (one probe per row,
		one column per chip).
expression value	a matrix of PM values	a list of two elements exprs
		and se.exprs

Table 2: Summary table for the processing methods.

- [1] "bg.correct" "mas" "none" "rma"
- > pmcorrect.methods()
- [1] "mas" "methods" "pmonly" "subtractmm"
- > express.summary.stat.methods()
- [1] "avgdiff" "liwong" "mas" "medianpolish" "playerout"

We would recommend the use of the method normalize.methods to access the list of available normalization methods (as a scheme for normalization methods that would go beyond 'affy' is thought).

- > library(affydata)
- > data(Dilution)
- > normalize.methods(Dilution)
- [1] "constant" "contrasts" "invariantset" "loess"
- [5] "methods" "qspline" "quantiles" "quantiles.robust"

For each processing step, a naming convention exists between the method label and the function name in R (see table 1). Each processing methods should be passed objects (and return objects) corresponding to the processing step (see table 2).

Practically, this means that to add a method one has to

- 1. create an appropriate method with a name satisfying the convention.
- 2. register the method by adding the label name to the corresponding variable.

### 3 Examples

As an example we show how to add two simple new methods. The first one does pm correction. The method subtract mm values to the pm values, except when the result is negative. In this case the pm value remains unchanged.

We create a function using the label name subtractmmsometimes.

Once the method defined, we just register the *label name* in the corresponding variable.

```
> upDate.pmcorrect.methods(c(pmcorrect.methods(), "subtractmmsometimes"))
```

The second new method intends to be an robust alternative to the summary expression value computation <code>avgdiff</code>. The idea is to use the function <code>huber</code> of the package <code>MASS</code>.

```
> huber <- function(y, k = 1.5, tol = 1e-06) {
      y <- y[!is.na(y)]</pre>
      n <- length(y)
      mu <- median(y)</pre>
      s \leftarrow mad(y)
+
      if (s == 0)
           stop("cannot estimate scale: MAD is zero for this sample")
+
           yy \leftarrow pmin(pmax(mu - k * s, y), mu + k * s)
           mu1 <- sum(yy)/n
           if (abs(mu - mu1) < tol * s)
               break
           mu <- mu1
      }
      list(mu = mu, s = s)
+ }
```

This method returns P.J. Huber's location estimate with MAD scale. You think this is what you want to compute the summary expression value from the probe intensities. What is needed to have as a processing method is a simple wrapper:

```
> computeExprVal.huber <- function(probes) {
+ res <- apply(probes, 2, huber)</pre>
```

```
+ mu <- unlist(lapply(res, function(x) x$mu))
+ s <- unlist(lapply(res, function(x) x$s))
+ return(list(exprs = mu, se.exprs = s))
+ }
> upDate.generateExprSet.methods(c(generateExprSet.methods(), "huber"))
```

From now the package is aware of the two new methods...in this session.

The code for the methods included in the package can be informative if you plan to develop methods. An example that demonstrates how a normalization method can be added is given by the function normalize.AffyBatch.vsn in the package vsn; see its help file.