# Package 'missMethyl'

October 17, 2020

Type Package

Title Analysing Illumina HumanMethylation BeadChip Data

**Version** 1.22.0 **Date** 2019-10-04

Author Belinda Phipson and Jovana Maksimovic

Maintainer Belinda Phipson <br/>
belinda.phipson@mcri.edu.au>, Jovana Maksimovic <jovana.maksimovic@petermac.org>, Andrew Lonsdale <andrew.lonsdale@petermac.org>

**Depends** R (>= 3.6.0), IlluminaHumanMethylation450kanno.ilmn12.hg19, IlluminaHumanMethylationEPICanno.ilm10b4.hg19

Imports AnnotationDbi, BiasedUrn, Biobase, BiocGenerics,
GenomicRanges, GO.db, IlluminaHumanMethylation450kmanifest,
IlluminaHumanMethylation450kanno.ilmn12.hg19,
IlluminaHumanMethylationEPICmanifest,
IlluminaHumanMethylationEPICanno.ilm10b4.hg19, IRanges, limma,
methods, methylumi, minfi, org.Hs.eg.db, ruv, S4Vectors,
statmod, stringr, SummarizedExperiment

### VignetteBuilder knitr

**Suggests** BiocStyle, edgeR, knitr, minfiData, rmarkdown, tweeDEseqCountData

Description Normalisation, testing for differential variability and differential methylation and gene set testing for data from Illumina's Infinium HumanMethylation arrays. The normalisation procedure is subset-quantile within-array normalisation (SWAN), which allows Infinium I and II type probes on a single array to be normalised together. The test for differential variability is based on an empirical Bayes version of Levene's test. Differential methylation testing is performed using RUV, which can adjust for systematic errors of unknown origin in high-dimensional data by using negative control probes. Gene ontology analysis is performed by taking into account the number of probes per gene on the array.

License GPL-2

**biocViews** Normalization, DNAMethylation, MethylationArray, GenomicVariation, GeneticVariability, DifferentialMethylation, GeneSetEnrichment 2 missMethyl-package

RoxygenNote 7.0.2
Encoding UTF-8
${\bf git\_url}\ \ https://git.bioconductor.org/packages/missMethyl$
git_branch RELEASE_3_11
git_last_commit f805189
git_last_commit_date 2020-04-27
Date/Publication 2020-10-16

# **R** topics documented:

	missMethyl-package	2
	contrasts.varFit	3
	densityByProbeType	4
	getAdj	6
	getINCs	7
	getLeveneResiduals	8
	getMappedEntrezIDs	9
	gometh	10
	goregion	13
	gsameth	16
	gsaregion	19
	RUVadj	21
	RUVfit	24
	SWAN	26
	topGSA	28
	topRUV	29
	topVar	31
	varFit	32
Index		36
missM	ethyl-package Introduction to the missMethyl package	

# Description

missMethyl is a library for the analysis of Illumina's 450K human methylation BeadChip. Specifically, functions for SWAN normalisation and differential variability analysis are provided. SWAN normalisation uses probe specific information, and the differential variability procedure uses linear models which can handle any designed experiment.

# **Details**

Package: missMethyl Type: Package Version: 0.99.0 Date: 2014-06-30 License: GPL-2 contrasts.varFit 3

Normalisation of the 450K arrays can be performed using the function SWAN.

Differential variability analysis can be performed by calling varFit followed by topVar for a list of the top ranked differentially variable CpGs between conditions.

More detailed help documentation is provided in each function's help page.

### Author(s)

Belinda Phipson and Jovana Maksimovic

Maintainer: Belinda Phipson <br/> belinda.phipson@mcri.edu.au>, Jovana Maksimovic <jovana.maksimovic@petermac.org

#### References

Maksimovic, J., Gordon, L., Oshlack, A. (2012). SWAN: Subset-quantile within array normalization for illumina infinium HumanMethylation450 BeadChips. *Genome Biology*, 13:R44.

Phipson, B., and Oshlack, A. (2014). DiffVar: A new method for detecting differential variability with application to methylation in cancer and aging. *Genome Biology*, **15**:465.

contrasts.varFit

Compute contrasts for a varFit object.

# **Description**

Compute estimated coefficients, standard errors and LogVarRatios for a given set of contrasts.

# Usage

```
contrasts.varFit(fit, contrasts = NULL)
```

### **Arguments**

fit List containing a linear model fit produced by varFit. The fit object should be

of class MArrayLM.

containing contrasts.

### **Details**

This function calls the contrasts.fit function in limma to compute coefficients and standard errors for the specified contrasts corresponding to a linear model fit obtained from the varFit function. LogVarRatios are also computed in terms of the contrasts. A contrasts matrix can be computed using the makeContrasts function.

# Value

A list object of the same class as fit.

# Author(s)

Belinda Phipson

### See Also

```
varFit, contrasts.fit, makeContrasts
```

# **Examples**

```
# Randomly generate data for a 3 group problem with 100 CpG sites and 4
# arrays in each group.
library(limma)
y<-matrix(rnorm(1200),ncol=12)
group<-factor(rep(c(1,2,3),each=4))</pre>
design<-model.matrix(~0+group)</pre>
colnames(design)<-c("grp1", "grp2", "grp3")</pre>
# Fit linear model for differential variability
vfit<-varFit(y,design)</pre>
# Specify contrasts
contr<-makeContrasts(grp2-grp1,grp3-grp1,grp3-grp2,levels=colnames(design))</pre>
# Compute contrasts from fit object
vfit.contr<-contrasts.varFit(vfit,contrasts=contr)</pre>
summary(decideTests(vfit.contr))
# Look at top table of results for first contrast
topVar(vfit.contr,coef=1)
```

 ${\tt densityByProbeType}$ 

Plot the beta value distributions of the Infinium I and II probe types relative to the overall beta value distribution.

# Description

Plot the overall density distribution of beta values and the density distributions of the Infinium I and II probe types.

# Usage

```
densityByProbeType(
  data,
  legendPos = "top",
  colors = c("black", "red", "blue"),
  main = "",
  lwd = 3,
  cex.legend = 1
)
```

densityByProbeType 5

### **Arguments**

A MethylSet or a matrix or a vector. We either use the getBeta function to get Beta values (in the first case) or we assume the matrix or vector contains Beta values.

The x and y co-ordinates to be used to position the legend. They can be specified by keyword or in any way which is accepted by xy.coords. See legend for details.

Colors

Colors to be used for the different beta value density distributions. Must be a vector of length 3.

main Plot title.

The line width to be used for the different beta value density distributions.

cex.legend The character expansion factor for the legend text.

### **Details**

The density distribution of the beta values for a single sample is plotted. The density distributions of the Infinium I and II probes are then plotted individually, showing how they contribute to the overall distribution. This is useful for visualising how using SWAN affects the data.

### Value

No return value. Plot is produced as a side-effect.

# Author(s)

Jovana Maksimovic

### See Also

```
densityPlot, densityBeanPlot, par, legend
```

```
if (require(minfi) & require(minfiData)) {
  dat <- preprocessRaw(RGsetEx)
  datSwan <- SWAN(dat)
  par(mfrow=c(1,2))
  densityByProbeType(dat[,1], main="Raw")
  densityByProbeType(datSwan[,1], main="SWAN")
}</pre>
```

6 getAdj

getAdj

Extract values adjusted for unwanted variation by RUVm

# **Description**

Extract values adjusted for unwanted variation by RUVm.

### Usage

```
getAdj(Y, fit)
```

# **Arguments**

Y A matrix of M-values.

fit The list list object produced by RUVfit.

### **Details**

This function extracts values adjusted for unwanted variation by RUVm. These values are ONLY intendeded to be used for visualisation purposes. It is NOT recommended that they are used for any further analysis.

### Value

An matrix of M-values.

# Author(s)

Jovana Maksimovic

# See Also

MArrayLM

```
if(require(minfi) & require(minfiData) & require(limma)) {

# Get methylation data for a 2 group comparison
meth <- getMeth(MsetEx)
unmeth <- getUnmeth(MsetEx)
Mval <- log2((meth + 100)/(unmeth + 100))

group <- factor(pData(MsetEx)$Sample_Group, labels=c(0,1))
design <- model.matrix(~group)

# Perform initial analysis to empirically identify negative control features
# when not known a priori
lFit <- lmFit(Mval,design)
lFit2 <- eBayes(lFit)
lTop <- topTable(lFit2,coef=2,num=Inf)</pre>
```

getINCs 7

```
# The negative control features should *not* be associated with factor of
# interest but *should* be affected by unwanted variation
ctl <- rownames(Mval) %in% rownames(lTop[lTop$adj.P.Val > 0.5,])
# Perform RUV adjustment and fit
fit <- RUVfit(Y=Mval, X=group, ctl=ctl)
fit2 <- RUVadj(Y=Mval, fit=fit)
# get adjusted values
Madj <- getAdj(Y=Mval,fit=fit)
}</pre>
```

getINCs

Extract intensity data for Illumina negative controls found on 450k or EPIC arrays.

# **Description**

Extracts the intensity data for the Illumina negative controls found on 450k or EPIC arrays and returns a matrix of M-values (log2 ratio of the green to red intensities).

### Usage

```
getINCs(rgSet)
```

# **Arguments**

rgSet

An object of class RGChannelSet.

# **Details**

The getINCs function extracts the intensity data for the INCs from an RGChannelSet object. The function retrieves both the green and red channel intensity values and returns the data as the log2 ratio of the green to red intensities. Essentially, the INCs are being treated like Type II probes for which the M-values are also given as the log2 ratio of the green to red intensities.

### Value

An matrix of M-values.

# Author(s)

Jovana Maksimovic

# See Also

**RGChannelSet** 

8 getLeveneResiduals

### **Examples**

```
if (require(minfi) & require(minfiData)) {
   INCs <- getINCs(RGsetEx)
   head(INCs)
   dim(INCs)
}</pre>
```

getLeveneResiduals

Obtain Levene residuals

# **Description**

Obtain absolute or squared Levene residuals for each CpG given a series of methylation arrays

# Usage

```
getLeveneResiduals(data, design = NULL, type = NULL)
```

# **Arguments**

data	Object of class matrix of M values, with rows corresponding to features of interest such as CpG sites and columns corresponding to samples or arrays.
design	The design matrix of the experiment, with rows corresponding to arrays/samples and columns to coefficients to be estimated. Defaults to the unit vector.
type	Character string, "AD" for absolute residuals or "SQ" for squared residuals. Default is "AD".

# **Details**

This function will return absolute or squared Levene residuals given a matrix of M values and a design matrix. This can be used for graphing purposes or for downstream analysis such a gene set testing based on differential variability rather than differential methylation. If no design matrix is given, the residuals are determined by treating all samples as coming from one group.

# Value

Returns a list with three components. data contains a matrix of absolute or squared residuals, AvgVar is a vector of sample variances and LogVarRatio corresponds to the columns of the design matrix and is usually the ratios of the log of the group variances.

# Author(s)

Belinda Phipson

# References

Phipson, B., and Oshlack, A. (2014). A method for detecting differential variability in methylation data shows CpG islands are highly variably methylated in cancers. *Genome Biology*, **15**:465.

getMappedEntrezIDs 9

### See Also

varFit

### **Examples**

```
# Randomly generate data for a 2 group problem with 100 CpG sites and 5
# arrays in each group
y <- matrix(rnorm(1000),ncol=10)

group <- factor(rep(c(1,2),each=5))
design <- model.matrix(~group)

# Get absolute Levene Residuals
resid <- getLeveneResiduals(y,design)

# Plot the first CpG
barplot(resid$data[1,],col=rep(c(2,4),each=5),
ylab="Absolute Levene Residuals",names=group)</pre>
```

getMappedEntrezIDs

Get mapped Entrez Gene IDs from CpG probe names

### **Description**

Given a set of CpG probe names and optionally all the CpG sites tested, this function outputs a list containing the mapped Entrez Gene IDs as well as the numbers of probes per gene, and a vector indicating significance.

# Usage

```
getMappedEntrezIDs(sig.cpg, all.cpg = NULL, array.type, anno = NULL)
```

# **Arguments**

sig.cpg	Character vector of significant CpG sites used for testing gene set enrichment.
all.cpg	Character vector of all CpG sites tested. Defaults to all CpG sites on the array.
array.type	The Illumina methylation array used. Options are "450K" or "EPIC".
anno	Optional. A DataFrame object containing the complete array annotation as generated by the minfi getAnnotation function. Speeds up execution, if provided.

### **Details**

This function is used by the gene set testing functions gometh and gsameth. It maps the significant CpG probe names to Entrez Gene IDs, as well as all the CpG sites tested. It also calculates the numbers of probes for gene.

Genes associated with each CpG site are obtained from the annotation package IlluminaHumanMethylation450kanno.: if the array type is "450K". For the EPIC array, the annotation package IlluminaHumanMethylationEPICanno.ilm10b4 is used. To use a different annotation package, please supply it using the anno argument.

10 gometh

### Value

A list with the following elements

sig.eg mapped Entrez Gene IDs for the significant probes

universe mapped Entrez Gene IDs for all probes on the array, or for all the CpG probes

tested.

freq table output with numbers of probes associated with each gene

de a vector of ones and zeroes of the same length of universe indicating which

genes in the universe are significantly differentially methylated.

fract.counts a dataframe with 2 columns corresponding to the Entrez Gene IDS for the sig-

nificant probes and the associated weight to account for multi-mapping probes.

# Author(s)

Belinda Phipson

### See Also

```
gometh, gsameth
```

# **Examples**

```
## Not run: # to avoid timeout on Bioconductor build
library(IlluminaHumanMethylation450kanno.ilmn12.hg19)
library(org.Hs.eg.db)
library(limma)
ann <- getAnnotation(IlluminaHumanMethylation450kanno.ilmn12.hg19)

# Randomly select 1000 CpGs to be significantly differentially methylated
sigcpgs <- sample(rownames(ann),1000,replace=FALSE)

# All CpG sites tested
allcpgs <- rownames(ann)

mappedEz <- getMappedEntrezIDs(sigcpgs,allcpgs,array.type="450K")
mappedEz$sig.eg[1:10]
mappedEz$freq[1:10]
mappedEz$freq[1:10]
mappedEz$de[1:10]

## End(Not run)</pre>
```

gometh

Gene ontology testing for Ilumina methylation array data

# **Description**

Tests gene ontology enrichment for significant CpGs from Illumina's Infinium HumanMethylation450 or MethylationEPIC array, taking into account the differing number of probes per gene present on the array.

gometh 11

### Usage

```
gometh(
    sig.cpg,
    all.cpg = NULL,
    collection = c("GO", "KEGG"),
    array.type = c("450K", "EPIC"),
    plot.bias = FALSE,
    prior.prob = TRUE,
    anno = NULL,
    equiv.cpg = TRUE,
    fract.counts = TRUE
)
```

## **Arguments**

sig.cpg	Character vector of significant CpG sites to test for GO term enrichment.
all.cpg	Character vector of all CpG sites tested. Defaults to all CpG sites on the array.
collection	The collection of pathways to test. Options are "GO" and "KEGG". Defaults to "GO".
array.type	The Illumina methylation array used. Options are "450K" or "EPIC". Defaults to "450K".
plot.bias	Logical, if true a plot showing the bias due to the differing numbers of probes per gene will be displayed.
prior.prob	Logical, if true will take into account the probability of significant differentially methylation due to numbers of probes per gene. If false, a hypergeometric test is performed ignoring any bias in the data.
anno	Optional. A DataFrame object containing the complete array annotation as generated by the minfi. getAnnotation function. Speeds up execution, if provided.
equiv.cpg	Logical, if true then equivalent numbers of cpgs are used for odds calculation rather than total number cpgs. Only used if prior.prob=TRUE.
fract.counts	Logical, if true then fractional counting of cpgs is used to account for cgps that map to multiple genes. Only used if prior.prob=TRUE.

# **Details**

This function takes a character vector of significant CpG sites, maps the CpG sites to Entrez Gene IDs, and tests for GO term or KEGG pathway enrichment using a hypergeometric test, taking into account the number of CpG sites per gene on the 450K/EPIC array. Geeleher et al. (2013) showed that a severe bias exists when performing gene set analysis for genome-wide methylation data that occurs due to the differing numbers of CpG sites profiled for each gene. gometh is based on the goseq method (Young et al., 2010) and calls the goana function for GO testing, or the kegga function for KEGG testing, both of which are from the limma package (Ritchie et al. 2015). If prior.prob is set to FALSE, then prior probabilities are not used and it is assumed that each gene is equally likely to have a significant CpG site associated with it.

The testing now also takes into account that some CpGs map to multiple genes. For a small number of gene families, this previously caused their associated GO categories/gene sets to be erroneously overrepresented and thus highly significant. If fract.counts=FALSE then CpGs are allowed to map to multiple genes (this is NOT recommended).

12 gometh

Genes associated with each CpG site are obtained from the annotation package IlluminaHumanMethylation450kanno.: if the array type is "450K". For the EPIC array, the annotation package IlluminaHumanMethylationEPICanno.ilm10b4 is used. To use a different annotation package, please supply it using the anno argument.

In order to get a list which contains the mapped Entrez gene IDs, please use the getMappedEntrezIDs function. gometh tests all GO or KEGG terms, and false discovery rates are calculated using the method of Benjamini and Hochberg (1995). The limma functions topGO and topKEGG can be used to display the top 20 most enriched pathways.

For more generalised gene set testing where the user can specify the gene set/s of interest to be tested, please use the gsameth function.

#### Value

A data frame with a row for each GO or KEGG term and the following columns:

Term GO term if testing GO pathways

Ont ontology that the GO term belongs to if testing GO pathways. "BP" - biological

process, "CC" - cellular component, "MF" - molecular function.

Pathway the KEGG pathway being tested if testing KEGG terms.

N number of genes in the GO or KEGG term

DE number of genes that are differentially methylated

P.DE p-value for over-representation of the GO or KEGG term term

FDR False discovery rate

# Author(s)

Belinda Phipson

### References

Phipson, B., Maksimovic, J., and Oshlack, A. (2016). missMethyl: an R package for analysing methylation data from Illuminas HumanMethylation450 platform. *Bioinformatics*, **15**;32(2), 286–8.

Geeleher, P., Hartnett, L., Egan, L. J., Golden, A., Ali, R. A. R., and Seoighe, C. (2013). Gene-set analysis is severely biased when applied to genome-wide methylation data. *Bioinformatics*, **29**(15), 1851–1857.

Young, M. D., Wakefield, M. J., Smyth, G. K., and Oshlack, A. (2010). Gene ontology analysis for RNA-seq: accounting for selection bias. *Genome Biology*, 11, R14.

Ritchie, M. E., Phipson, B., Wu, D., Hu, Y., Law, C. W., Shi, W., and Smyth, G. K. (2015). limma powers differential expression analyses for RNA-sequencing and microarray studies. *Nucleic Acids Research*, gkv007.

Benjamini, Y., and Hochberg, Y. (1995). Controlling the false discovery rate: a practical and powerful approach to multiple testing. *Journal of the Royal Statistical Society Series*, B, **57**, 289-300.

# See Also

goana, kegga, gsameth

goregion 13

# **Examples**

```
## Not run: # to avoid timeout on Bioconductor build
library(IlluminaHumanMethylation450kanno.ilmn12.hg19)
library(limma)
ann <- getAnnotation(IlluminaHumanMethylation450kanno.ilmn12.hg19)</pre>
# Randomly select 1000 CpGs to be significantly differentially methylated
sigcpgs <- sample(rownames(ann),1000,replace=FALSE)</pre>
# All CpG sites tested
allcpgs <- rownames(ann)</pre>
# GO testing with prior probabilities taken into account
# Plot of bias due to differing numbers of CpG sites per gene
gst <- gometh(sig.cpg = sigcpgs, all.cpg = allcpgs, collection = "GO",</pre>
              plot.bias = TRUE, prior.prob = TRUE, anno = ann)
# Total number of GO categories significant at 5% FDR
table(gst$FDR<0.05)
# Table of top GO results
topGSA(gst)
# GO testing ignoring bias
gst.bias <- gometh(sig.cpg = sigcpgs, all.cpg = allcpgs, collection = "GO",</pre>
                     prior.prob=FALSE, anno = ann)
\mbox{\#} Total number of GO categories significant at 5% FDR ignoring bias
table(gst.bias$FDR<0.05)</pre>
# Table of top GO results ignoring bias
topGSA(gst.bias)
# KEGG testing
kegg <- gometh(sig.cpg = sigcpgs, all.cpg = allcpgs, collection = "KEGG",</pre>
                prior.prob=TRUE, anno = ann)
# Table of top KEGG results
topGSA(kegg)
## End(Not run)
```

goregion

Gene ontology testing for Ilumina methylation array data

### **Description**

Tests gene ontology or KEGG pathway enrichment for differentially methylated regions (DMRs) identified from Illumina's Infinium HumanMethylation450 or MethylationEPIC array, taking into account the differing number of probes per gene present on the array.

# Usage

```
goregion(
  regions,
  all.cpg = NULL,
  collection = c("GO", "KEGG"),
  array.type = c("450K", "EPIC"),
  plot.bias = FALSE,
  prior.prob = TRUE,
  anno = NULL,
```

14 goregion

```
equiv.cpg = TRUE,
fract.counts = TRUE
)
```

### **Arguments**

regions	GRanges object of DMR coordinates to test for GO term enrichment.
all.cpg	Character vector of all CpG sites tested. Defaults to all CpG sites on the array.
collection	The collection of pathways to test. Options are "GO" and "KEGG". Defaults to "GO".
array.type	The Illumina methylation array used. Options are "450K" or "EPIC". Defaults to "450K".
plot.bias	Logical, if true a plot showing the bias due to the differing numbers of probes per gene will be displayed.
prior.prob	Logical, if true will take into account the probability of significant differentially methylation due to numbers of probes per gene. If false, a hypergeometric test is performed ignoring any bias in the data.
anno	Optional. A DataFrame object containing the complete array annotation as generated by the minfi. getAnnotation function. Speeds up execution, if provided.
equiv.cpg	Logical, if true then equivalent numbers of cpgs are used for odds calculation rather than total number cpgs. Only used if prior.prob=TRUE.
fract.counts	Logical, if true then fractional counting of cpgs is used to account for cgps that map to multiple genes. Only used if prior.prob=TRUE.

# **Details**

This function takes a GRanges object of DMR coordinates, maps them to CpG sites on the array and then to Entrez Gene IDs, and tests for GO term or KEGG pathway enrichment using a hypergeometric test, taking into account the number of CpG sites per gene on the 450K/EPIC array. Geeleher et al. (2013) showed that a severe bias exists when performing gene set analysis for genome-wide methylation data that occurs due to the differing numbers of CpG sites profiled for each gene. gometh is based on the goseq method (Young et al., 2010) and calls the goana function for GO testing, or the kegga function for KEGG testing, both of which are from the limma package (Ritchie et al. 2015). If prior.prob is set to FALSE, then prior probabilities are not used and it is assumed that each gene is equally likely to have a significant CpG site associated with it.

The testing now also takes into account that some CpGs map to multiple genes. For a small number of gene families, this previously caused their associated GO categories/gene sets to be erroneously overrepresented and thus highly significant. If fract.counts=FALSE then CpGs are allowed to map to multiple genes (this is NOT recommended).

Genes associated with each CpG site are obtained from the annotation package IlluminaHumanMethylation450kanno.: if the array type is "450K". For the EPIC array, the annotation package IlluminaHumanMethylationEPICanno.ilm10b4 is used. To use a different annotation package, please supply it using the anno argument.

In order to get a list which contains the mapped Entrez gene IDS, please use the getMappedEntrezIDs function. gometh tests all GO or KEGG terms, and false discovery rates are calculated using the method of Benjamini and Hochberg (1995). The limma functions topGO and topKEGG can be used to display the top 20 most enriched pathways.

For more generalised gene set testing where the user can specify the gene set/s of interest to be tested, please use the gsameth function.

goregion 15

#### Value

A data frame with a row for each GO or KEGG term and the following columns:

Term GO term if testing GO pathways

Ont ontology that the GO term belongs to if testing GO pathways. "BP" - biological

process, "CC" - cellular component, "MF" - molecular function.

Pathway the KEGG pathway being tested if testing KEGG terms.

N number of genes in the GO or KEGG term

DE number of genes that are differentially methylated

P.DE p-value for over-representation of the GO or KEGG term term

FDR False discovery rate

### Author(s)

Jovana Maksimovic

### References

Phipson, B., Maksimovic, J., and Oshlack, A. (2016). missMethyl: an R package for analysing methylation data from Illuminas HumanMethylation450 platform. *Bioinformatics*, **15**;32(2), 286–8.

Geeleher, P., Hartnett, L., Egan, L. J., Golden, A., Ali, R. A. R., and Seoighe, C. (2013). Gene-set analysis is severely biased when applied to genome-wide methylation data. *Bioinformatics*, **29**(15), 1851–1857.

Young, M. D., Wakefield, M. J., Smyth, G. K., and Oshlack, A. (2010). Gene ontology analysis for RNA-seq: accounting for selection bias. *Genome Biology*, 11, R14.

Ritchie, M. E., Phipson, B., Wu, D., Hu, Y., Law, C. W., Shi, W., and Smyth, G. K. (2015). limma powers differential expression analyses for RNA-sequencing and microarray studies. *Nucleic Acids Research*, gkv007.

Benjamini, Y., and Hochberg, Y. (1995). Controlling the false discovery rate: a practical and powerful approach to multiple testing. *Journal of the Royal Statistical Society Series*, B, **57**, 289-300.

# See Also

```
goana, kegga, gsameth
```

16 gsameth

```
dmrcoutput <- dmrcate(myannotation, lambda=1000)</pre>
regions <- extractRanges(dmrcoutput)</pre>
ann <- getAnnotation(IlluminaHumanMethylation450kanno.ilmn12.hg19)</pre>
# All CpG sites tested
allcpgs <- rownames(ann)</pre>
# GO testing with prior probabilities taken into account
# Plot of bias due to differing numbers of CpG sites per gene
gst <- goregion(regions = regions, all.cpg = allcpgs, collection = "GO",</pre>
              plot.bias = TRUE, prior.prob = TRUE, anno = ann)
# Total number of GO categories significant at 5% FDR
table(gst$FDR<0.05)
# Table of top GO results
topGSA(gst)
# GO testing ignoring bias
gst.bias <- goregion(regions = regions, all.cpg = allcpgs, collection = "GO",</pre>
                    prior.prob=FALSE, anno = ann)
\mbox{\#} Total number of GO categories significant at 5% FDR ignoring bias
table(gst.bias$FDR<0.05)
# Table of top GO results ignoring bias
topGSA(gst.bias)
# KEGG testing
kegg <- goregion(regions = regions, all.cpg = allcpgs, collection = "KEGG",</pre>
                prior.prob=TRUE, anno = ann)
# Table of top KEGG results
topGSA(kegg)
## End(Not run)
```

gsameth

Generalised gene set testing for Illumina's methylation array data

# **Description**

Given a user specified list of gene sets to test, grameth tests whether significantly differentially methylated CpG sites are enriched in these gene sets.

### Usage

```
gsameth(
    sig.cpg,
    all.cpg = NULL,
    collection,
    array.type = c("450K", "EPIC"),
    plot.bias = FALSE,
    prior.prob = TRUE,
    anno = NULL,
    equiv.cpg = TRUE,
    fract.counts = TRUE
)
```

gsameth 17

# **Arguments**

sig.cpg	Character vector of significant CpG sites to test for gene set enrichment.
all.cpg	Character vector of all CpG sites tested. Defaults to all CpG sites on the array.
collection	A list of user specified gene sets to test. Can also be a single character vector gene set. Gene identifiers must be Entrez Gene IDs.
array.type	The Illumina methylation array used. Options are "450K" or "EPIC". Defaults to "450K".
plot.bias	Logical, if true a plot showing the bias due to the differing numbers of probes per gene will be displayed
prior.prob	Logical, if true will take into account the probability of significant differentially methylation due to numbers of probes per gene. If false, a hypergeometric test is performed ignoring any bias in the data.
anno	Optional. A DataFrame object containing the complete array annotation as generated by the minfi getAnnotation function. Speeds up execution, if provided.
equiv.cpg	Logical, if true then equivalent numbers of cpgs are used for odds calculation rather than total number cpgs. Only used if $prior.prob=TRUE$ .
fract.counts	Logical, if true then fractional counting of cpgs is used to account for cgps that map to multiple genes. Only used if prior.prob=TRUE.

### **Details**

This function extends gometh, which only tests GO and KEGG pathways. gsameth can take a list of user specified gene sets and test whether the significant CpG sites are enriched in these pathways. gsameth maps the CpG sites to Entrez Gene IDs and tests for pathway enrichment using a hypergeometric test, taking into account the number of CpG sites per gene on the 450K/EPIC arrays. Please note the gene ids for the collection of gene sets must be Entrez Gene IDs. Geeleher et al. (2013) showed that a severe bias exists when performing gene set analysis for genome-wide methylation data that occurs due to the differing numbers of CpG sites profiled for each gene. gsameth and gometh is based on the goseq method (Young et al., 2010). If prior.prob is set to FALSE, then prior probabilities are not used and it is assumed that each gene is equally likely to have a significant CpG site associated with it.

The testing now also takes into account that some CpGs map to multiple genes. For a small number of gene families, this previously caused their associated GO categories/gene sets to be erroneously overrepresented and thus highly significant. If fract.counts=FALSE then CpGs are allowed to map to multiple genes (this is NOT recommended).

Genes associated with each CpG site are obtained from the annotation package IlluminaHumanMethylation450kanno.: if the array type is "450K". For the EPIC array, the annotation package IlluminaHumanMethylationEPICanno.ilm10b4 is used. To use a different annotation package, please supply it using the anno argument.

In order to get a list which contains the mapped Entrez gene IDS, please use the getMappedEntrezIDs function.

### Value

A data frame with a row for each gene set and the following columns:

N number of genes in the gene set

DE number of genes that are differentially methylated P.DE p-value for over-representation of the gene set

FDR False discovery rate, calculated using the method of Benjamini and Hochberg

(1995).

18 gsameth

### Author(s)

Belinda Phipson

#### References

Phipson, B., Maksimovic, J., and Oshlack, A. (2016). missMethyl: an R package for analysing methylation data from Illuminas HumanMethylation450 platform. *Bioinformatics*, **15**;32(2), 286–8.

Geeleher, P., Hartnett, L., Egan, L. J., Golden, A., Ali, R. A. R., and Seoighe, C. (2013). Gene-set analysis is severely biased when applied to genome-wide methylation data. *Bioinformatics*, **29**(15), 1851–1857.

Young, M. D., Wakefield, M. J., Smyth, G. K., and Oshlack, A. (2010). Gene ontology analysis for RNA-seq: accounting for selection bias. *Genome Biology*, 11, R14.

Ritchie, M. E., Phipson, B., Wu, D., Hu, Y., Law, C. W., Shi, W., and Smyth, G. K. (2015). limma powers differential expression analyses for RNA-sequencing and microarray studies. *Nucleic Acids Research*, gkv007.

Benjamini, Y., and Hochberg, Y. (1995). Controlling the false discovery rate: a practical and powerful approach to multiple testing. *Journal of the Royal Statistical Society Series*, B, **57**, 289-300.

### See Also

gometh,getMappedEntrezIDs

```
## Not run: # to avoid timeout on Bioconductor build
library(IlluminaHumanMethylation450kanno.ilmn12.hg19)
library(org.Hs.eg.db)
library(limma)
ann <- getAnnotation(IlluminaHumanMethylation450kanno.ilmn12.hg19)</pre>
# Randomly select 1000 CpGs to be significantly differentially methylated
sigcpgs <- sample(rownames(ann),1000,replace=FALSE)</pre>
# All CpG sites tested
allcpgs <- rownames(ann)</pre>
# Use org.Hs.eg.db to extract a GO term
GOtoID <- suppressMessages(select(org.Hs.eg.db, keys=keys(org.Hs.eg.db),</pre>
                                    columns=c("ENTREZID", "GO"),
                                    keytype="ENTREZID"))
setname1 <- GOtoID$GO[1]</pre>
setname1
keep.set1 <- GOtoID$GO %in% setname1
set1 <- GOtoID$ENTREZID[keep.set1]</pre>
setname2 <- GOtoID$GO[2]</pre>
setname2
keep.set2 <- GOtoID$GO %in% setname2
set2 <- GOtoID$ENTREZID[keep.set2]</pre>
# Make the gene sets into a list
sets <- list(set1, set2)</pre>
names(sets) <- c(setname1, setname2)</pre>
# Testing with prior probabilities taken into account
# Plot of bias due to differing numbers of CpG sites per gene
gst <- gsameth(sig.cpg = sigcpgs, all.cpg = allcpgs, collection = sets,</pre>
                 plot.bias = TRUE, prior.prob = TRUE)
```

gsaregion 19

gsaregion

Generalised gene set testing for Illumina's methylation array data

# **Description**

Given a user specified list of gene sets to test, gsameth tests whether differentially methylated regions (DMRs) identified from Illumina's Infinium HumanMethylation450 or MethylationEPIC array are enriched, taking into account the differing number of probes per gene present on the array.

# Usage

```
gsaregion(
  regions,
  all.cpg = NULL,
  collection,
  array.type = c("450K", "EPIC"),
  plot.bias = FALSE,
  prior.prob = TRUE,
  anno = NULL,
  equiv.cpg = TRUE,
  fract.counts = TRUE
)
```

# **Arguments**

regions	GRanges Object of DMR coordinates to test for GO term enrichment.
all.cpg	Character vector of all CpG sites tested. Defaults to all CpG sites on the array.
collection	A list of user specified gene sets to test. Can also be a single character vector gene set. Gene identifiers must be Entrez Gene IDs.
array.type	The Illumina methylation array used. Options are " $450K$ " or "EPIC". Defaults to " $450K$ ".
plot.bias	Logical, if true a plot showing the bias due to the differing numbers of probes per gene will be displayed.
prior.prob	Logical, if true will take into account the probability of significant differentially methylation due to numbers of probes per gene. If false, a hypergeometric test is performed ignoring any bias in the data.
anno	Optional. A DataFrame object containing the complete array annotation as generated by the minfi getAnnotation function. Speeds up execution, if provided.
equiv.cpg	Logical, if true then equivalent numbers of cpgs are used for odds calculation rather than total number cpgs. Only used if prior.prob=TRUE.
fract.counts	Logical, if true then fractional counting of cpgs is used to account for cgps that map to multiple genes. Only used if prior.prob=TRUE.

20 gsaregion

#### **Details**

This function extends gometh, which only tests GO and KEGG pathways. gsameth can take a list of user specified gene sets and test whether the significant DMRs are enriched in these pathways. This function takes a GRanges object of DMR coordinates, maps them to CpG sites on the array and then to Entrez Gene IDs, and tests for enrichment using a hypergeometric test, taking into account the number of CpG sites per gene on the 450K/EPIC array. Geeleher et al. (2013) showed that a severe bias exists when performing gene set analysis for genome-wide methylation data that occurs due to the differing numbers of CpG sites profiled for each gene. gsameth and gometh is based on the goseq method (Young et al., 2010). If prior.prob is set to FALSE, then prior probabilities are not used and it is assumed that each gene is equally likely to have a significant CpG site associated with it.

The testing now also takes into account that some CpGs map to multiple genes. For a small number of gene families, this previously caused their associated GO categories/gene sets to be erroneously overrepresented and thus highly significant. If fract.counts=FALSE then CpGs are allowed to map to multiple genes (this is NOT recommended).

Genes associated with each CpG site are obtained from the annotation package IlluminaHumanMethylation450kanno.: if the array type is "450K". For the EPIC array, the annotation package IlluminaHumanMethylationEPICanno.ilm10b4 is used. To use a different annotation package, please supply it using the anno argument.

In order to get a list which contains the mapped Entrez gene IDS, please use the getMappedEntrezIDs function.

### Value

A data frame with a row for each gene set and the following columns:

N number of genes in the gene set

DE number of genes that are differentially methylated
P.DE p-value for over-representation of the gene set

FDR False discovery rate, calculated using the method of Benjamini and Hochberg

(1995).

# Author(s)

Belinda Phipson

# References

Phipson, B., Maksimovic, J., and Oshlack, A. (2016). missMethyl: an R package for analysing methylation data from Illuminas HumanMethylation450 platform. *Bioinformatics*, **15**;32(2), 286–8.

Geeleher, P., Hartnett, L., Egan, L. J., Golden, A., Ali, R. A. R., and Seoighe, C. (2013). Gene-set analysis is severely biased when applied to genome-wide methylation data. *Bioinformatics*, **29**(15), 1851–1857.

Young, M. D., Wakefield, M. J., Smyth, G. K., and Oshlack, A. (2010). Gene ontology analysis for RNA-seq: accounting for selection bias. *Genome Biology*, 11, R14.

Ritchie, M. E., Phipson, B., Wu, D., Hu, Y., Law, C. W., Shi, W., and Smyth, G. K. (2015). limma powers differential expression analyses for RNA-sequencing and microarray studies. *Nucleic Acids Research*, gkv007.

Benjamini, Y., and Hochberg, Y. (1995). Controlling the false discovery rate: a practical and powerful approach to multiple testing. *Journal of the Royal Statistical Society Series*, B, **57**, 289-300.

RUVadj 21

### See Also

```
gometh,getMappedEntrezIDs
```

### **Examples**

```
## Not run: # to avoid timeout on Bioconductor build
library(IlluminaHumanMethylation450kanno.ilmn12.hg19)
library(org.Hs.eg.db)
library(limma)
library(DMRcate)
data(dmrcatedata)
myMs <- logit2(myBetas)</pre>
myMs.noSNPs <- rmSNPandCH(myMs, dist=2, mafcut=0.05)</pre>
patient <- factor(sub("-.*", "", colnames(myMs)))</pre>
type <- factor(sub(".*-", "", colnames(myMs)))</pre>
design <- model.matrix(~patient + type)</pre>
myannotation <- cpg.annotate("array", myMs.noSNPs, what="M", arraytype = "450K",</pre>
                               analysis.type="differential", design=design, coef=39)
dmrcoutput <- dmrcate(myannotation, lambda=1000)</pre>
regions <- extractRanges(dmrcoutput)</pre>
ann <- getAnnotation(IlluminaHumanMethylation450kanno.ilmn12.hg19)</pre>
# All CpG sites tested
allcpgs <- rownames(ann)</pre>
# Use org.Hs.eg.db to extract a GO term
GOtoID <- suppressMessages(select(org.Hs.eg.db, keys=keys(org.Hs.eg.db),</pre>
                                    columns=c("ENTREZID","GO"), keytype="ENTREZID"))
keep.set1 <- GOtoID$GO %in% "GO:0010951"
set1 <- GOtoID$ENTREZID[keep.set1]</pre>
keep.set2 <- GOtoID$GO %in% "GO:0042742"
set2 <- GOtoID$ENTREZID[keep.set2]</pre>
# Make the gene sets into a list
sets <- list(set1, set2)</pre>
names(sets) <- c("GO:0010951", "GO:0042742")
# Testing with prior probabilities taken into account
# Plot of bias due to differing numbers of CpG sites per gene
gst <- gsaregion(regions = regions, all.cpg = allcpgs, collection = sets,</pre>
                 plot.bias = TRUE, prior.prob = TRUE)
topGSA(gst)
# Testing ignoring bias
gst.bias <- gsaregion(regions = regions, all.cpg = allcpgs, collection = sets,</pre>
                     prior.prob = FALSE)
topGSA(gst.bias)
## End(Not run)
```

RUVadi

RUV adjust

# **Description**

Post-process and summarize the results of call to RUVfit.

22 RUVadj

# Usage

```
RUVadj(
   Y,
   fit,
   var.type = c("ebayes", "standard", "pooled"),
   p.type = c("standard", "rsvar", "evar"),
   cpginfo = NULL,
   ...
)
```

# **Arguments**

Υ	The original data matrix used in the call to RUVfit.
fit	A RUV model fit (a list) as returned by RUVfit.
var.type	Which type of estimate for sigma2 should be used from the call to variance_adjust? The options are "ebayes", "standard", or "pooled." See variance_adjust for details.
p.type	Which type of p-values should be used from the call to variance_adjust? The options are "standard", "rsvar", or "evar".
cpginfo	A matrix or dataframe containing information about the CpGs. This information is included in the summary that is returned.
	Other parameters that can be passed to ruv function ruv_summary.

### **Details**

This function post-processes the results of a call to RUVfit and then summarizes the output. The post-processing step primarily consists of a call to ruv\_summary and variance\_adjust, which computes various adjustments to variances, t-statistics, and and p-values. See variance\_adjust for details. The var.type and p.type options determine which of these adjustments are used.

After post-processing, the results are summarized into a list containing 4 objects: 1) the data matrix Y; 2) a dataframe R containing information about the rows (samples); 3) a dataframe C containing information about the columns (features, e.g. genes), and 4) a list misc of other information returned by RUVfit.

# Value

An list containing:

Υ	The original data matrix
R	A dataframe of sample-wise information, including $X$ , $Z$ , and any other data passed in with rowinfo.
С	A dataframe of cpg-wise information, including p-values, estimated regression coefficients, estimated variances, column means, an index of the negative controls, and any other data passed in with cpginfo.
misc	A list of additional information returned by RUVfit.

# Author(s)

Jovana Maksimovic < jovana.maksimovic@mcri.edu.au>

RUVadj 23

#### References

Benjamini, Y., and Hochberg, Y. (1995). Controlling the false discovery rate: a practical and powerful approach to multiple testing. *Journal of the Royal Statistical Society Series*, B, **57**, 289-300.

Gagnon-Bartsch JA, Speed TP. (2012). Using control genes to correct for unwanted variation in microarray data. *Biostatistics*. **13**(3), 539-52. Available at: http://biostatistics.oxfordjournals.org/content/13/3/539.full.

Gagnon-Bartsch, Jacob, and Speed. 2013. Removing Unwanted Variation from High Dimensional Data with Negative Controls. Available at: http://statistics.berkeley.edu/tech-reports/820.

Smyth, G. K. (2004). Linear models and empirical Bayes methods for assessing differential expression in microarray experiments. *Statistical Applications in Genetics and Molecular Biology*, Volume 3, Article 3. http://www.statsci.org/smyth/pubs/ebayes.pdf.

### See Also

MArrayLM, RUV2, RUV4, RUVinv, RUVrinv, p.adjust, get\_empirical\_variances, sigmashrink

```
if(require(minfi) & require(minfiData) & require(limma)) {
# Get methylation data for a 2 group comparison
meth <- getMeth(MsetEx)</pre>
unmeth <- getUnmeth(MsetEx)</pre>
Mval <- log2((meth + 100)/(unmeth + 100))
group<-factor(pData(MsetEx)$Sample_Group)</pre>
design<-model.matrix(~group)</pre>
# Perform initial analysis to empirically identify negative control features
# when not known a priori
lFit <- lmFit(Mval,design)</pre>
lFit2 <- eBayes(lFit)</pre>
lTop <- topTable(lFit2,coef=2,num=Inf)</pre>
# The negative control features should *not* be associated with factor of
# interest but *should* be affected by unwanted variation
ctl <- rownames(Mval) %in% rownames(lTop[lTop$adj.P.Val > 0.5,])
# Perform RUV adjustment and fit
fit <- RUVfit(Y=Mval, X=group, ctl=ctl)</pre>
fit2 <- RUVadj(Y=Mval, fit=fit)</pre>
# Look at table of top results
top <- topRUV(fit2)</pre>
}
```

24 RUVfit

RUVfit

Remove unwanted variation when testing for differential methylation

# **Description**

Provides an interface similar to lmFit from limma to the RUV2, RUV4, RUVinv and RUVrinv functions from the ruv package, which facilitates the removal of unwanted variation in a differential methylation analysis. A set of negative control variables, as described in the references, must be specified.

# Usage

```
RUVfit(
    Y,
    X,
    ctl,
    Z = 1,
    k = NULL,
    method = c("inv", "rinv", "ruv4", "ruv2"),
    ...
)
```

### **Arguments**

Υ	numeric matrix with rows corresponding to the features of interest such as CpG sites and columns corresponding to samples or arrays.
X	The factor(s) of interest. A m by p matrix, where m is the number of samples and p is the number of factors of interest. Very often $p = 1$ . Factors and dataframes are also permissible, and converted to a matrix by design.matrix.
ctl	logical vector, length == nrow(Y). Features that are to be used as negative control variables are indicated as TRUE, all other features are FALSE.
Z	Any additional covariates to include in the model, typically a m by q matrix. Factors and dataframes are also permissible, and converted to a matrix by design.matrix. Alternatively, may simply be 1 (the default) for an intercept term. May also be NULL.
k	integer, required if method is "ruv2" or "ruv4". Indicates the number of unwanted factors to use. Can be 0.
method	character string, indicates which ruv method should be used.
	additional arguments that can be passed to RUV2, RUV4, RUVinv and RUVrinv. See linked function documentation for details.

### **Details**

This function depends on the ruv package and is used to estimate and adjust for unwanted variation in a differential methylation analysis. Briefly, the unwanted factors W are estimated using negative control variables. Y is then regressed on the variables X, Z, and W. For methylation data, the analysis is performed on the M-values, defined as the log base 2 ratio of the methylated signal to the unmethylated signal.

RUVfit 25

#### Value

### A list containing:

betahat The estimated coefficients of the factor(s) of interest. A p by n matrix.

t t statistics for the factor(s) of interest. A p by n matrix.

P-values for the factor(s) of interest. A p by n matrix.

Fstats F statistics for testing all of the factors in X simultaneously.

Fpvals P-values for testing all of the factors in X simultaneously.

multiplier The constant by which sigma2 must be multiplied in order to get an estimate of

the variance of betahat.

df The number of residual degrees of freedom.

W The estimated unwanted factors.

alpha The estimated coefficients of W.

byx The coefficients in a regression of Y on X (after both Y and X have been "ad-

justed" for Z). Useful for projection plots.

bwx The coefficients in a regression of W on X (after X has been "adjusted" for Z).

Useful for projection plots.

X X. Included for reference.
k k. Included for reference.
ctl ctl. Included for reference.
Z Z. Included for reference.

fullW0 Can be used to speed up future calls of RUVfit.

include.intercept

include.intercept. Included for reference.

method Character variable with value indicating which RUV method was used. Included

for reference.

# Author(s)

Jovana Maksimovic

### References

Gagnon-Bartsch JA, Speed TP. (2012). Using control genes to correct for unwanted variation in microarray data. *Biostatistics*. **13**(3), 539-52. Available at: http://biostatistics.oxfordjournals.org/content/13/3/539.full.

Gagnon-Bartsch, Jacob, and Speed. 2013. Removing Unwanted Variation from High Dimensional Data with Negative Controls. Available at: http://statistics.berkeley.edu/tech-reports/820.

# See Also

RUV2, RUV4, RUVinv, RUVrinv, topRUV

26 SWAN

### **Examples**

```
if(require(minfi) & require(minfiData) & require(limma)) {
# Get methylation data for a 2 group comparison
meth <- getMeth(MsetEx)</pre>
unmeth <- getUnmeth(MsetEx)</pre>
Mval <- log2((meth + 100)/(unmeth + 100))
group <- factor(pData(MsetEx)$Sample_Group)</pre>
design <- model.matrix(~group)</pre>
# Perform initial analysis to empirically identify negative control features
# when not known a priori
lFit <- lmFit(Mval,design)</pre>
lFit2 <- eBayes(lFit)</pre>
lTop <- topTable(lFit2,coef=2,num=Inf)</pre>
# The negative control features should *not* be associated with factor of
# interest but *should* be affected by unwanted variation
ctl <- rownames(Mval) %in% rownames(lTop[lTop$adj.P.Val > 0.5,])
# Perform RUV adjustment and fit
fit <- RUVfit(Y=Mval, X=group, ctl=ctl)</pre>
fit2 <- RUVadj(Y=Mval, fit=fit)</pre>
# Look at table of top results
top <- topRUV(fit2)</pre>
}
```

**SWAN** 

Subset-quantile Within Array Normalisation for Illumina Infinium HumanMethylation450 BeadChips

# **Description**

Subset-quantile Within Array Normalisation (SWAN) is a within array normalisation method for the Illumina Infinium HumanMethylation450 platform. It allows Infinium I and II type probes on a single array to be normalized together.

# Usage

```
SWAN(data, verbose = FALSE)
## S3 method for class 'MethyLumiSet'
SWAN(data, verbose = FALSE)
## S3 method for class 'RGChannelSet'
SWAN(data, verbose = FALSE)
## Default S3 method:
SWAN(data, verbose = FALSE)
```

### **Arguments**

data An object of class either MethylSet, RGChannelSet or MethyLumiSet.

verbose Should the function be verbose?

SWAN 27

#### **Details**

The SWAN method has two parts. First, an average quantile distribution is created using a subset of probes defined to be biologically similar based on the number of CpGs underlying the probe body. This is achieved by randomly selecting N Infinium I and II probes that have 1, 2 and 3 underlying CpGs, where N is the minimum number of probes in the 6 sets of Infinium I and II probes with 1, 2 or 3 probe body CpGs. If no probes have previously been filtered out e.g. sex chromosome probes, etc. N=11,303. This results in a pool of 3N Infinium I and 3N Infinium II probes. The subset for each probe type is then sorted by increasing intensity. The value of each of the 3N pairs of observations is subsequently assigned to be the mean intensity of the two probe types for that row or 'quantile'. This is the standard quantile procedure. The intensities of the remaining probes are then separately adjusted for each probe type using linear interpolation between the subset probes.

### Value

An object of class MethylSet.

NULL

NULL

NULL

### Note

SWAN uses a random subset of probes to perform the within-array normalization. In order to achive reproducible results, the seed needs to be set using set.seed.

### Author(s)

Jovana Maksimovic

### References

J Maksimovic, L Gordon and A Oshlack (2012). SWAN: Subset quantile Within-Array Normalization for Illumina Infinium HumanMethylation450 BeadChips. Genome Biology 13, R44.

# See Also

RGChannelSet and MethylSet as well as MethyLumiSet and IlluminaMethylationManifest.

```
if (require(minfi) & require(minfiData)) {
    set.seed(100)
    datSwan1 <- SWAN(RGsetEx)

    dat <- preprocessRaw(RGsetEx)
    set.seed(100)
    datSwan2 <- SWAN(dat)

    head(getMeth(datSwan2)) == head(getMeth(datSwan1))
}</pre>
```

28 topGSA

topGSA	Get table of top 20 enriched pathways

# **Description**

After using gsameth, calling topGSA will output the top 20 most significantly enriched pathways.

# Usage

```
topGSA(gsa, number = 20, sort = TRUE)
```

# **Arguments**

gsa Matrix, from output of gsameth.

number Scalar, number of pathway results to output. Default is 20.

sort Logical, should the table be ordered by p-value. Default is TRUE.

### **Details**

This function will output the top 20 most significant pathways from a pathway analysis using the gsameth function. The output is ordered by p-value.

### Value

A matrix ordered by P.DE, with a row for each gene set and the following columns:

N number of genes in the gene set

DE number of genes that are differentially methylated
P.DE p-value for over-representation of the gene set

FDR False discovery rate, calculated using the method of Benjamini and Hochberg

(1995).

# Author(s)

Belinda Phipson

# See Also

gsameth

```
library(IlluminaHumanMethylation450kanno.ilmn12.hg19)
library(org.Hs.eg.db)
library(limma)
ann <- getAnnotation(IlluminaHumanMethylation450kanno.ilmn12.hg19)

# Randomly select 1000 CpGs to be significantly differentially methylated sigcpgs <- sample(rownames(ann),1000,replace=FALSE)

# All CpG sites tested</pre>
```

topRUV 29

```
allcpgs <- rownames(ann)</pre>
# Use org.Hs.eg.db to extract a GO term
GOtoID <- toTable(org.Hs.egG02EG)</pre>
setname1 <- GOtoID$go_id[1]</pre>
setname1
keep.set1 <- GOtoID\$go\_id \%in\% \ setname1
set1 <- GOtoID$gene_id[keep.set1]</pre>
setname2 <- GOtoID$go_id[2]</pre>
setname2
keep.set2 <- GOtoID$go_id %in% setname2</pre>
set2 <- GOtoID$gene_id[keep.set2]</pre>
# Make the gene sets into a list
sets <- list(set1, set2)</pre>
names(sets) <- c(setname1,setname2)</pre>
# Testing with prior probabilities taken into account
# Plot of bias due to differing numbers of CpG sites per gene
gst <- gsameth(sig.cpg = sigcpgs, all.cpg = allcpgs, collection = sets,</pre>
plot.bias = TRUE, prior.prob = TRUE)
topGSA(gst)
# Testing ignoring bias
gst.bias <- gsameth(sig.cpg = sigcpgs, all.cpg = allcpgs, collection = sets,</pre>
                      prior.prob = FALSE)
topGSA(gst.bias)
```

topRUV

Table of top-ranked differentially methylated CpGs obatained from a differential methylation analysis using RUV

# **Description**

Extract a table of the top-ranked CpGs from a linear model fit after performing a differential methylation analysis using RUVfit and RUVadj.

### Usage

```
topRUV(fitsum, number = 10, sort.by = c("p", "F.p"), p.BH = 1)
```

# Arguments

fitsum	An object containing the summary fit object produced by RUVadj. The object should be a list.
number	integer, maximum number of genes to list. Default is 10.
sort.by	character string, what the results should be sorted by. Default is unadjusted p-value.
p.BH	numeric, cutoff value for Benjamini-Hochberg adjusted p-values. Only features with lower p-values are listed. Must be between 0 and 1. Default is 1.

30 topRUV

#### **Details**

This function summarises the results of a differential methylation analysis performed using RUVfit, followed by RUVadj. The top ranked CpGs are sorted by p-value.

### Value

Produces a dataframe with rows corresponding to the top number CpGs and the following columns: F.p F.p.BH p\_X1 p.BH\_X1 b\_X1 sigma2 var.b\_X1 fit.ctl mean

F.p	P-values for testing all of the factors of interest simultaneously.
F.p.BH	Benjamini-Hochberg adjusted p-values for testing all of the factors of interest simultaneously.
p_X1	p-values for the factor of interest.
p.BH_X1	Benjamini-Hochberg adjusted p-values for the factor of interest.
b_X1	The estimated coefficients of the factor of interest.
sigma2	Estimate of the methylation variance.
var.b_X1	Variance estimate of betahat.
fit.ctl	logical, indicating whether CpG was designated as a negative control.
mean	The mean methylation (M-value).

### Author(s)

Jovana Maksimovic < jovana.maksimovic@mcri.edu.au>

# References

Benjamini, Y., and Hochberg, Y. (1995). Controlling the false discovery rate: a practical and powerful approach to multiple testing. *Journal of the Royal Statistical Society Series*, B, **57**, 289-300.

Smyth, G. K. (2004). Linear models and empirical Bayes methods for assessing differential expression in microarray experiments. *Statistical Applications in Genetics and Molecular Biology*, Volume 3, Article 3. http://www.statsci.org/smyth/pubs/ebayes.pdf.

### See Also

```
RUVfit, RUVadj, MArrayLM
```

```
if(require(minfi) & require(minfiData) & require(limma)){

# Get methylation data for a 2 group comparison
meth <- getMeth(MsetEx)
unmeth <- getUnmeth(MsetEx)
Mval <- log2((meth + 100)/(unmeth + 100))

group <- factor(pData(MsetEx)$Sample_Group)
design <- model.matrix(~group)

# Perform initial analysis to empirically identify negative control features
# when *not* known a priori
lFit <- lmFit(Mval,design)</pre>
```

topVar 31

```
lFit2 <- eBayes(lFit)
lTop <- topTable(lFit2,coef=2,num=Inf)

# The negative control features should *not* be associated with factor of
# interest but *should* be affected by unwanted variation
ctl <- rownames(Mval) %in% rownames(lTop[lTop$adj.P.Val > 0.5,])

# Perform RUV adjustment and fit
fit <- RUVfit(Y=Mval, X=group, ctl=ctl)
fit2 <- RUVadj(Y=Mval, fit=fit)

# Look at table of top results
top <- topRUV(fit2)
}</pre>
```

topVar

Table of top-ranked differentially variable CpGs

### **Description**

Extract a table of the top-ranked CpGs from a linear model fit after a differential variability analysis.

### Usage

```
topVar(fit, coef = NULL, number = 10, sort = TRUE)
```

### **Arguments**

fit	List containing a linear model fit produced by varFit. The fit object should be of class MArrayLM.
coef	Column number or column name specifying which coefficient of the linear model fit is of interest. It should be the same coefficient that the differential variability testing was performed on. Default is last column of fit object.
number	Maximum number of genes to list. Default is 10.
sort	Logical, default is TRUE. Sorts output according the P-value. FALSE will return results in same order as fit object.

# Details

This function summarises the results of a differential variability analysis performed with varFit. The p-values from the comparison of interest are adjusted using Benjamini and Hochberg's false discovery rate with the function p.adjust. The top ranked CpGs are selected by first ranking the adjusted p-values, then ranking the raw p-values. At this time no other sorting option is catered for.

# Value

Produces a dataframe with rows corresponding to the top CpGs and the following columns:

genelist one or more columns of annotation for each CpG, if the gene information is

available in fit

AvgVar average of the absolute or squared Levene residuals across all samples

DiffVar estimate of the difference in the Levene residuals corresponding to the compar-

ison of interest

t moderated t-statistic

P.Value raw p-value Adj.P.Value adjusted p-value

### Author(s)

Belinda Phipson

### References

Phipson, B., and Oshlack, A. (2014). A method for detecting differential variability in methylation data shows CpG islands are highly variably methylated in cancers. *Genome Biology*, **15**:465.

Benjamini, Y., and Hochberg, Y. (1995). Controlling the false discovery rate: a practical and powerful approach to multiple testing. *Journal of the Royal Statistical Society Series*, B, **57**, 289-300.

### See Also

```
varFit, p.adjust
```

# **Examples**

```
# Randomly generate data for a 2 group problem with 100 CpG sites and 5
# arrays in each group.

y<-matrix(rnorm(1000),ncol=10)

group<-factor(rep(c(1,2),each=5))
design<-model.matrix(~group)

# Fit linear model for differential variability
vfit<-varFit(y,design)

# Look at top table of results
topVar(vfit,coef=2)</pre>
```

varFit

Testing for differential variability

# **Description**

Fit linear model on mean absolute or squared deviations for each CpG given a series of methylation arrays

### Usage

```
varFit(
  data,
  design = NULL,
  coef = NULL,
  type = NULL,
  trend = TRUE
  robust = TRUE,
  weights = NULL
## S3 method for class 'MethylSet'
varFit(
  data,
  design = NULL,
  coef = NULL,
  type = NULL,
  trend = TRUE,
  robust = TRUE,
  weights = NULL
## S3 method for class 'DGEList'
varFit(
  data,
  design = NULL,
  coef = NULL,
  type = NULL,
  trend = TRUE,
  robust = TRUE,
  weights = NULL
## Default S3 method:
varFit(
  data,
  design = NULL,
  coef = NULL,
  type = NULL,
  trend = TRUE,
  robust = TRUE,
  weights = NULL
)
```

# **Arguments**

coef

Object of class MethylSet or matrix with rows corresponding to the features data of interest such as CpG sites and columns corresponding to samples or arrays.

The design matrix of the experiment, with rows corresponding to arrays/samples design

and columns to coefficients to be estimated. Defaults to the unit vector.

The columns of the design matrix containing the comparisons to test for differential variability.

type Character string, "AD" for absolute residuals or "SQ" for squared residuals. De-

fault is absolute.

trend Logical, if true fits a mean variance trend on the absolute or squared deviations.

robust Logical, if true performs robust empirical Bayes shrinkage of the variances for

the moderated t statistics.

weights Non-negative observation weights. Can be a numeric matrix of individual weights,

of same size as the object matrix, or a numeric vector of array weights, or a nu-

meric vector of gene/feature weights.

### **Details**

This function depends on the limma package and is used to rank features such as CpG sites or genes in order of evidence of differential variability between different comparisons corresponding to the columns of the design matrix. A measure of variability is calculated for each CpG in each sample by subtracting out the group mean and taking the absolute or squared deviation. A linear model is then fitted to the absolute or squared deviations. The residuals of the linear model fit are subjected to empirical Bayes shrinkage and moderated t statistics (Smyth, 2004) calculated. False discovery rates are calculated using the method of Benjamini and Hochberg (1995).

If coef is not specified, then group means are estimated based on all the columns of the design matrix and subtracted out before testing for differential variability. If the design matrix contains nuisance parameters, then subsetting the design matrix columns by coef should remove these columns from the design matrix. If the design matrix includes an intercept term, this should be included in coef. The nuisance parameters are included in the linear model fit to the absolute or squared deviations, but should not be considered when subtracting group means to obtain the deviations. Note that design matrices without an intercept term are permitted, and specific contrasts tested using the function contrasts.varFit.

For methylation data, the analysis is performed on the M-values, defined as the log base 2 ratio of the methylated signal to the unmethylated signal. If a MethylSet object is supplied, M-values are extracted with an offset of 100 added to the numerator and denominator.

For testing differential variability on RNA-Seq data, a DGEList object can be supplied directly to the function. A voom transformation is applied before testing for differential variability. The weights calculated in voom are used in the linear model fit.

Since the output is of class MArrayLM, any functions that can be applied to fit objects from lmFit and eBayes can be applied, for example, topTable and decideTests.

### Value

Produces an object of class MArrayLM (see MArrayLM-class) containing everything found in a fitted model object produced by 1mFit and eBayes as well as a vector containing the sample CpG-wise variances and a matrix of LogVarRatios corresponding to the differential variability analysis.

NULL

NULL

NULL

# Author(s)

Belinda Phipson

#### References

Phipson, B., and Oshlack, A. (2014). A method for detecting differential variability in methylation data shows CpG islands are highly variably methylated in cancers. *Genome Biology*, **15**:465.

Smyth, G.K. (2004). Linear models and empirical Bayes methods for assessing differential expression in microarray experiments. *Statistical Applications in Genetics and Molecular Biology*, Volume **3**, Article 3.

Smyth, G. K. (2005). Limma: linear models for microarray data. In: *Bioinformatics and Computational Biology Solutions using R and Bioconductor*. R. Gentleman, V. Carey, S. Dudoit, R. Irizarry, W. Huber (eds), Springer, New York, 2005.

Benjamini, Y., and Hochberg, Y. (1995). Controlling the false discovery rate: a practical and powerful approach to multiple testing. *Journal of the Royal Statistical Society Series*, B, **57**, 289-300.

### See Also

 $\verb|contrasts.varFit|, topVar|, getLeveneResiduals|, lmFit|, eBayes|, topTable|, decideTests|, voom|$ 

```
# Randomly generate data for a 2 group problem with 100 CpG sites and 5
# arrays in each # group.

y<-matrix(rnorm(1000),ncol=10)

group<-factor(rep(c(1,2),each=5))
design<-model.matrix(~group)

# Fit linear model for differential variability
vfit<-varFit(y,design,coef=c(1,2))

# Look at top table of results
topVar(vfit,coef=2)</pre>
```

# Index

<pre>contrasts.varFit, 3, 35  decideTests, 35 densityBeanPlot, 5 densityByProbeType, 4 densityPlot, 5 design.matrix, 24</pre>	RUVadj, 21, 30 RUVfit, 21, 22, 24, 30 RUVinv, 23-25 RUVrinv, 23-25 sigmashrink, 23 SWAN, 5, 26
eBayes, 35  get_empirical_variances, 23  getAdj, 6  getAnnotation, 9, 11, 14, 17, 19	topGSA, 28 topRUV, 25, 29 topTable, 35 topVar, 31, 35
getINCs, 7 getLeveneResiduals, 8, 35 getMappedEntrezIDs, 9, 18, 21	varFit, 9, 32 variance_adjust, 22 voom, 35
goana, 12, 15 gometh, 10, 10, 18, 21 goregion, 13 gsameth, 10, 12, 15, 16, 28 gsaregion, 19	xy.coords,5
IlluminaMethylationManifest, 27	
kegga, <i>12</i> , <i>15</i>	
legend, 5 limma, 24 lmFit, 24, 35	
MArrayLM, 6, 23, 30 MethylSet, 27 MethyLumiSet, 27 minfi, 9, 11, 14, 17, 19 missMethyl (missMethyl-package), 2 missMethyl-package, 2	
p.adjust, 23 par, 5	
RGChannelSet, 7, 27 ruv, 22, 24 RUV2, 23-25 RUV4, 23-25 ruv_summary, 22	