Package 'snm'

April 15, 2017

Type Package

Version 1.22.0

Title Supervised Normalization of Microarrays

Author	Brig Mecham and John D. Storey <jstorey@princeton.edu></jstorey@princeton.edu>
noi pre we eith sta fac dat mo	ion SNM is a modeling strategy especially designed for malizing high-throughput genomic data. The underlying smise of our approach is that your data is a function of what refer to as study-specific variables. These variables are ner biological variables that represent the target of the tistical analysis, or adjustment variables that represent tors arising from the experimental or biological setting the a is drawn from. The SNM approach aims to simultaneously del all study-specific variables in order to more accurately tracterize the biological or clinical variables of interest.
Maintaiı	ner John D. Storey <jstorey@princeton.edu></jstorey@princeton.edu>
Depends	R (>= 2.12.0)
Imports	corpcor, lme4 (>= 1.0), splines
License	LGPL
Di: Tra	s Microarray, OneChannel, TwoChannel, MultiChannel, ferentialExpression, ExonArray, GeneExpression, unscription, MultipleComparison, Preprocessing, alityControl
NeedsCo	mpilation no
R topi	cs documented:
	sim.doubleChannel sim.preProcessed sim.refDesign sim.singleChannel snm snm.fitted snm.plot snm.summary 1
Index	1

2 sim.doubleChannel

sim.doubleChannel	Simulated data	for a double channe	l microarray	ovnoriment
SIII. UUUDIECHAIIIEI	simulalea aala	јот и идидне спитне	і тистойттаў	ехрентет.

Description

Simulates data used in snm function examples.

Usage

sim.doubleChannel(seed)

Arguments

seed Numeric value used to seed random number generator.

Details

Simulated data set influenced by a probe-specific biological variable, a probe-specific adjustment variables, and intensity-dependent array and dye effects. Data were simulated for a total of 25,000 probes and 20 arrays. The biological variable is a dichotomous variable specifying two groups (Group 1 and Group 2), with 10 channels per group. The continuous probe-specific adjustment variable is sampled from a Normal(1,0.1) and mimics an age effect. The baseline probe intensities were sampled from a chi(1,2) distribution. Any baseline intensities greater than 15 were set to a random value from the interval [15,16]. The random variation terms were sampled from a Normal(0,0.25) and the array and dye functions were defined by randomly sampling coefficients for a two-dimensional B-spline basis functions from a Normal(0,1).

Randomly selected subsets of 30% and 20% of the probes were defined as influenced by the biological groups and age variables, respectfully. The magnitude of the biological effects were sampled from a Normal(1,0.3) distribution, and the probe-specific age effects from a Normal(1,0.1). An instance of this simulated data can be generated using the code in the examples section below.

Value

raw.data	a 25,000 by 50 matrix of simulated data generated according to the description above.
true.nulls	a vector of indices corresponding to the rows in raw.data of the probes unaffected by the biological variable of interest
bio.var	a model matrix of the biological variable of interest.
adj.var	a model matrix of the adjustment variables
int.var	a data frame of the intensity-dependent adjustment variables

Author(s)

Brig Mecham

See Also

```
snm, sim.singleChannel, sim.preProcessed, sim.refDesign
```

sim.preProcessed 3

Examples

```
## Not run:
doubleChannel <- sim.doubleChannel(12346)
snm.obj <- snm(doubleChannel$raw.data, doubleChannel$bio.var,
doubleChannel$adj.var, doubleChannel$int.var)
ks.test(snm.obj$pval[doubleChannel$true.nulls],"punif")
## End(Not run)</pre>
```

sim.preProcessed

Simulate data from a microarray experiment without any intensity-dependent effects.

Description

Simulated data set influenced by a single probe-specific biological and two probe-specific adjustment variables. This parameters for this data are identical to single channel simulated data available as sim.singleChannel(seed) with the difference that this example does not include the intensity-dependent effects. Consult the corresponding help file for details on this simulation.

Usage

```
sim.preProcessed(seed)
```

Arguments

seed Numeric value used to seed random number generator.

Value

raw.data a 25,000 by 50 matrix of simulated data generated according to the description

above.

true.nulls a vector of indices corresponding to the rows in raw.data of the probes unaffected

by the biological variable of interest

bio.var a model matrix of the biological variable of interest.

adj.var a model matrix of the adjustment variables

int.var set to NULL

Author(s)

Brig Mecham

See Also

```
snm, sim.doubleChannel, sim.singleChannel, sim.refDesign
```

Examples

4 sim.refDesign

sim.refDesign	Simulates data from a two-color microarray experiment using a reference design.
	ence design.

Description

Simulates a two channel experiment with a reference design. Used as an example for snm function call.

Usage

```
sim.refDesign(seed)
```

Arguments

seed

Numeric value used to seed random number generator.

Details

Simulated data set influenced by a probe-specific biological variable, and intensity-dependent array and dye effects. Data were simulated assuming a uniform reference design for a total of 25,000 probes and 20 arrays, each consisting of two channels. The reference channel consists of a single reference RNA population. The experimental channel measures a dichotomous biological variable specifying two groups (Group 1 and Group 2), with 10 samples per group. The baseline probe intensities were sampled from a chi(1,2) distribution. Any baseline intensities greater than 15 were set to a random value from the interval [15,16]. The random variation terms were sampled from a Normal(0,0.25) and the array and dye functions were defined by randomly sampling coefficients for a two-dimensional B-spline basis functions from a Normal(0,1).

A randomly selected subset of 30% of the probes was defined as influenced by the biological group variable. The magnitude of the biological effects were sampled from a Normal(1,0.3) distribution. An instance of this simulated data can be generated using the code in the examples section below.

Value

raw.data	a 25,000 by 50 matrix of simulated data generated according to the description above.
true.nulls	a vector of indices corresponding to the rows in raw.data of the probes unaffected by the biological variable of interest
bio.var	a model matrix of the biological variable of interest.
adj.var	a model matrix of the probe-specific adjustment variables. In this case set to NULL.
int.var	a data frame of the intensity-dependent adjustment variables

Author(s)

Brig Mecham and John D. Storey <jstorey@princeton.edu>

See Also

```
snm, sim.singleChannel, sim.doubleChannel, sim.preProcessed
```

sim.singleChannel 5

Examples

```
# reference design on channel level data
# reference channel is included in bio.var
refChannel <- sim.refDesign(12347)</pre>
snm.obj <- snm(refChannel$raw.data,refChannel$bio.var, refChannel$int.var)</pre>
ks.test(snm.obj$pval[refChannel$true.nulls], "punif")
# this is a different model formulation
# where reference channel is not included
# in bio.var
bio.var2 = matrix(c(rep(-1,10), rep(1, 10), rep(0,20)), ncol=1, dimnames=list(NULL, "treatment"))
adj.var2 = matrix(c(rep(1,40), \, rep(1,20), \, rep(0,20)), \, ncol=2, \, byrow=FALSE, \, dimnames=list(NULL, \, c("intercept", \, color=1, \, byrow=FALSE))))
# compare bio.var2 to refChannel$bio.var and adj.var2 to refChannel$adj.var
snm.obj <- snm(refChannel$raw.data, bio.var2, adj.var2, refChannel$int.var)</pre>
ks.test(snm.obj$pval[refChannel$true.nulls],"punif")
# reference design on log ratio data
# that is, log(target/ref) data
refChannel <- sim.refDesign(12347)</pre>
refChannel$raw.data = refChannel$raw.data[,1:20]-refChannel$raw.data[,21:40]
# removing reference channel variables
refChannel$bio.var = refChannel$bio.var[1:20,-3]
refChannel$adj.var = matrix(refChannel$adj.var[1:20,], ncol=1)
refChannel$int.var = data.frame(refChannel$int.var[1:20,1])
snm.obj <- snm(refChannel$raw.data,refChannel$bio.var, refChannel$adj.var, refChannel$int.var)</pre>
ks.test(snm.obj$pval[refChannel$true.nulls],"punif")
```

sim.singleChannel

Simulate data from a single channel microarray experiment

Description

Simulates single channel data used as an example for snm function call.

Usage

```
sim.singleChannel(seed)
```

Arguments

seed

Numeric value used to seed random number generator.

Details

Simulated data set influenced by a single probe-specific biological, two probe-specific adjustment variables, and intensity-dependent array effects. Data were simulated for a total of 25,000 probes and 50 arrays. The biological variable is a dichotomous variable specifying two groups (Group 1 and Group 2), with 25 arrays sampled from each group. The dichtomous probe-specific adjustment variables has 5 different levels and mimics a batch effect. The 5 batches each contain 10 samples, and are balanced with respect to the biological grouping factor. The continuous probe-specific adjustment variable is sampled from a Normal(1,0.1) and mimics an age effect. The baseline probe intensities were sampled from a chi(1,2) distribution. Any baseline intensities greater than 15 were set to a random value from the interval [15,16]. The random variation terms were sampled from

6 snm

a Normal(0,0.25) and the array functions were defined by randomly sampling coefficients for a two-dimensional B-spline basis functions from a Normal(0,1).

Randomly selected subsets of 30%, 10%, and 20% of the probes were defined as influenced by the biological groups, batch, and age variables, respectfully. The magnitude of the biological effects were sampled from a Normal(1,0.3) distribution, the probe-specific batch effects from a Normal(0,0.3) and the probe-specific age effects from a Normal(1,0.1). An instance of this simulated data can be generated using the code in the examples section below.

Value

raw.data	a 25,000 by 50 matrix of simulated data generated according to the description above.
true.nulls	a vector of indices corresponding to the rows in raw.data of the probes unaffected by the biological variable of interest
bio.var	a model matrix of the biological variable of interest.
adj.var	a model matrix of the adjustment variables
int.var	a data frame of the intensity-dependent adjustment variables

Author(s)

Brig Mecham

See Also

```
snm, sim.doubleChannel, sim.preProcessed, sim.refDesign
```

Examples

snm

Perform a supervised normalization of microarray data

Description

Implement Supervised Normalization of Microarrays on a gene expression matrix. Requires a set of biological covariates of interest and at least one probe-specific or intensity-dependent adjustment variable.

Usage

```
snm(raw.dat, bio.var=NULL, adj.var=NULL, int.var=NULL,
    weights=NULL, spline.dim = 4, num.iter = 10, lmer.max.iter=1000,
    nbins=20, rm.adj=FALSE, verbose=TRUE, diagnose=TRUE)
```

snm 7

Arguments

An m probes by n arrays matrix of expression data. If the user wishes to remove raw.dat intensity-dependent effects, then we request the matrix corresponds to the raw, log transformed data. A model matrix (see model.matrix) or data frame with n rows of the biological bio.var variables. If NULL, then all probes are treated as "null" in the algorithm. adj.var A model matrix (see model.matrix) or data frame with n rows of the probespecific adjustment variables. If NULL, a model with an intercept term is used. int.var A data frame with n rows of type factor with the unique levels of intensitydependent effects. Each column parametrizes a unique source of intensitydependent effect (e.g., array effects for column 1 and dye effects for column weights A vector of length m. Values unchanged by algorithm, used to control the influence of each probe on the intensity-dependent array effects. spline.dim Dimension of basis spline used for array effects. Number of snm model fit iterations to run. num.iter Number of lmer iterations that are permitted. Set lmer.max.iter=NULL if no lmer.max.iter maximum is desired. Number of bins used by binning strategy. Array effects are calculated from a nbins $nbins \times n$ data matrix, where the (i, j) value is equal to that bin i's average intensity on array j. rm.adj If set to FALSE, then only the intensity dependent effects have been removed from the normalized data, implying the effects from the adjustment variables are still present. If TRUE, then the adjustment variables effects and the intensity dependent effects are both removed from the returned normalized data. verbose A flag telling the software whether or not to display a report after each iteration. TRUE produces the output. diagnose A flag telling the software whether or not to produce diagnostic output in the form of consecutive plots. TRUE produces the plot.

Details

This function implements the supervised normalization of microarrays algorithm described in Mecham, Nelson, and Storey (2010).

Value

norm.dat The matrix of normalized data. The default setting is rm.adj=FALSE, which

means that only the intensity-dependent effects have been subtracted from the data. If the user wants the adjustment variable effects removed as well, then set

rm.adj=TRUE when calling the snm function.

pvalues A vector of p-values testing the association of the biological variables with each

probe. These p-values are obtained from an ANOVA comparing models where the full model contains both the probe-specific biological and adjustment variables versus a reduced model that just contains the probe-specific adjustment variables. The data used for this comparison has the intensity-dependent variables removed. These returned p-values are calculated after the final iteration of

the algorithm.

8 snm

pi0	The estimated proportion of true null probes pi_0 , calculated after the final iteration of the algorithm.
iter.pi0s	A vector of length equal to num.iter containing the estimated pi_0 values at each iteration of the snm algorithm. These values should converge and any non-convergence suggests a problem with the data, the assumed model, or both
nulls	A vector indexing the probes utilized in estimating the intensity-dependent effects on the final iteration.
М	A matrix containing the estimated probe intensities for each array utilized in estimating the intensity-dependent effects on the final iteration. For memory parsimony, only a subset of values spanning the range is returned, currently nbins*100 values.
array.fx	A matrix of the final estimated intensity-dependent array effects. For memory parsimony, only a subset of values spanning the range is returned, currently nbins*100 values.
bio.var	The processed version of the same input variable.
adj.var	The processed version of the same input variable.
int.var	The processed version of the same input variable.
df0	Degrees of freedom of the adjustment variables.
df1	Degrees of freedom of the full model matrix, which includes the biological variables and the adjustment variables.
raw.dat	The input data.
rm.var	Same as the input (useful for later analyses).
call	Function call.

Note

It is necessary for adj.var and adj.var+bio.var to be valid model matrices (e.g., the models cannot be over-determined).

We suggest that the probe level data be analyzed on the log-transformed scale, particularly if the user wishes to remove intensity-dependent effects. It is recommended that the normalized data (and resulting inference) be inspected for latent structure using Surrogate Variable Analysis (Leek and Storey 2007, PLoS Genetics).

Author(s)

Brig Mecham and John D. Storey <jstorey@princeton.edu>

References

Mecham BH, Nelson PS, Storey JD (2010) Supervised normalization of microarrays. Bioinformatics, 26: 1308-1315.

See Also

model.matrix,plot.snm, fitted.snm, summary.snm, sim.singleChannel, sim.doubleChannel,
sim.preProcessed, sim.refDesign

snm.fitted 9

Examples

snm.fitted

Extract fitted values from an snm object

Description

Computes fitted values under models used in snm normalization.

Usage

```
snm.fitted(object, ...)
## S3 method for class 'snm'
fitted(object, ...)
```

Arguments

object Output from the snm function.
... Not used.

Details

Returns the fitted values under the "null model" (adjustment variables only) and the "full model" (adjustment variables + biological variables).

Value

fit0	Linear model fits when regressing each probe's normalized data on the null model, \sim adj.var.
fit1	Linear model fits when regressing each probe's normalized data on the full model, ~adi,var+bio,var.

Note

These fits are useful for investigating the quality of the study-specific model used in the normalization. For example, the residuals can be obtained from the full model fit and examined for latent structure.

Author(s)

John D. Storey <jstorey@princeton.edu>

10 snm.plot

References

Mecham BH, Nelson PS, Storey JD (2010) Supervised normalization of microarrays. Bioinformatics, 26: 1308-1315.

See Also

```
snm, sim.singleChannel
```

Examples

snm.plot

Display plots for an snm object

Description

Creates a diagnostic plot of the snm fit.

Usage

```
snm.plot(x, col.by=NULL, ...)
## S3 method for class 'snm'
plot(x, ...)
```

Arguments

X Output from the snm function.

col.by

A factor vector of length equal to the number of arrays providing a grouping by which to color the intensity-dependent effects. Instead of a factor vector, the input may also be a model matrix composed only of 0's and 1's with the number of rows equal to the number of arrays and number of columns less than or equal to the number of arrays.

... Arguments passed to the plot functions. Not recommended.

snm.summary 11

Details

A four panel plot composed of the following:

1. Convergence of π_0 estimates over model fitting iterations. The pi_0 estimates for each iteration are compared to the pi_0 estimate calculated during the final model fit.

- 2. A scree plot of the principal components analysis of the full model residual matrix.
- 3. A plot of the estimated intensity-dependent effects.
- 4. A histogram of the p-values testing each probe for an association with the biological variables (bio.var). All probes to the right of the vertical red line are the least $\hat{\pi}_0$ significant probes (i.e., those used in estimating intensity-dependent effects). The dashed horizontal line is the pi_0 estimate from the final model fit.

Value

Nothing of interest.

Author(s)

John D. Storey <jstorey@princeton.edu>

References

Mecham BH, Nelson PS, Storey JD (2010) Supervised normalization of microarrays. Bioinformatics, 26: 1308-1315.

See Also

```
snm, sim.singleChannel
```

Examples

snm.summary

Display summary information for an snm object

Description

Provides a summary of the snm fit.

12 snm.summary

Usage

```
snm.summary(object, cuts=c(0.0001, 0.001, 0.01, 0.025, 0.05, 0.10, 1), ...)
## S3 method for class 'snm'
summary(object, ...)
```

Arguments

object An object of class "snm", i.e., output from the snm function.

Cut points at which to calculate the number of significant probes.

Optional arguments to send to the qvalue function in calculating statistical sig-

nificance.

Details

A summary of the snm fit.

Value

Nothing of interest.

Author(s)

John D. Storey <jstorey@princeton.edu>

References

Mecham BH, Nelson PS, Storey JD (2010) Supervised normalization of microarrays. Bioinformatics, 26: 1308-1315.

See Also

```
snm, sim.singleChannel
```

Examples

Index

```
*Topic datagen
    sim.doubleChannel, 2
    sim.preProcessed, 3
    sim.refDesign,4
    sim.singleChannel,5
*Topic htest
    snm, 6
*Topic misc
    snm.fitted, 9
    snm.plot, 10
    snm.summary, 11
*Topic models
    snm, 6
fitted.snm, 8
fitted.snm(snm.fitted), 9
1mer, 7
model.matrix, 7, 8
plot.snm, 8
plot.snm(snm.plot), 10
qvalue, 12
sim.doubleChannel, 2, 3, 4, 6, 8
sim.preProcessed, 2, 3, 4, 6, 8
sim.refDesign, 2, 3, 4, 6, 8
sim.singleChannel, 2-4, 5, 8, 10-12
snm, 2-4, 6, 6, 10-12
snm.fitted, 9
snm.plot, 10
snm.summary, 11
summary.snm, 8
summary.snm(snm.summary), 11
```