# Package 'methylKit'

April 12, 2018

Type Package

**Title** DNA methylation analysis from high-throughput bisulfite sequencing results

Version 1.4.1

**Author** Altuna Akalin [aut, cre], Matthias Kormaksson [aut], Sheng Li [aut], Arsene Wabo [ctb], Adrian Bierling [aut], Alexander Gosdschan [aut]

**Maintainer** Altuna Akalin <aakalin@gmail.com>, Alexander Gosdschan <alex.gos90@gmail.com>

**Description** methylKit is an R package for DNA methylation analysis and annotation from high-throughput bisulfite sequencing. The package is designed to deal with sequencing data from RRBS and its variants, but also target-capture methods and whole genome bisulfite sequencing. It also has functions to analyze base-pair resolution 5hmC data from experimental protocols such as oxBS-Seq and TAB-Seq. Perl is needed to read SAM files only.

License Artistic-2.0

URL http://code.google.com/p/methylkit/

LazyLoad yes

**NeedsCompilation** yes

LinkingTo Rcpp, Rhtslib, zlibbioc

biocViews DNAMethylation, Sequencing, MethylSeq

**Depends** R (>= 3.3.0), GenomicRanges (>= 1.18.1), methods

Imports IRanges, data.table (>= 1.9.6), parallel, S4Vectors (>= 0.13.13), GenomeInfoDb, KernSmooth, qvalue, emdbook, Rsamtools, gtools, fastseg, rtracklayer, mclust, Rcpp, R.utils, limma, grDevices, graphics, stats, utils

Suggests testthat, knitr, rmarkdown, genomation

VignetteBuilder knitr

Collate 'methylKit.R' 'backbone.R' 'diffMeth.R' 'clusterSamples.R' 'regionalize.R' 'processBismarkAln.R' 'RcppExports.R' 'document\_data.R' 'bedgraph.R' 'reorganize.R' 'percMethylation.R' 'normalizeCoverage.R' 'pool.R' 'adjust.methylC.R' 'updateMethObject.R' 'batchControl.R'

2 R topics documented:

 $\label{local-control} 'dataSim.R' 'methylDBClasses.R' 'methylDBFunctions.R' 'tabix.functions.R' 'methSeg.R' 'diffMethDSS.R' 'deprecated_defunct.R' 'onUnload.R'$ 

RoxygenNote 6.0.1

# **R** topics documented:

adjustMethylC	3
assocComp	5
bedgraph	5
calculateDiffMeth	7
calculateDiffMethDSS	10
clusterSamples	11
dataSim	12
diffMethPerChr	14
extract	15
filterByCoverage	16
getAssembly	
getContext	19
getCorrelation	20
getCoverageStats	22
getData	23
	24
<del>-</del>	25
	26
· ·	
	30
	31
	32
	35
	36
	37
	38
	38
·	39
	40
	40
methylRawDB-class	
methylRawList-class	44
methylRawList.obj	
methylRawListDB-class	
normalizeCoverage	46
PCASamples	47
<u>.</u>	49
pool	50
processBismarkAln	51
reconstruct	53
regionCounts	55

adjustMethylC 3

adju	stMethylC	Adju	st 1	neo	ası	ıre	d £	5m	C	leı	el	s u	ısi	ng	51	hn	ıC	le	ve	ls										
Index																														69
	updateMethObject.			٠							•			•		•					•	•	•		•	•	•	•	•	68
	unite																													66
	tileMethylCounts .																													64
	show,methylBase-m	ethod																												63
	selectByOverlap																													62
	select																													60
	reorganize																													59
	removeComp																													57

# Description

Measured 5mC levels via bisulfite sequencing might be a combination of 5hmC and 5mC levels since bisulfite sequencing can not distinguish between the two. This function can adjust 5mC levels of a bisulfite sequencing experiment if the user supplies corresponding 5hmC levels from the same sample.

# Usage

```
adjustMethylC(mc,hmc,save.db,...,chunk.size)
## S4 method for signature 'methylRaw,methylRaw'
adjustMethylC(mc, hmc, save.db = FALSE, ...,
    chunk.size = 1e+06)
## S4 method for signature 'methylRawList,methylRawList'
adjustMethylC(mc, hmc,
    save.db = FALSE, ..., chunk.size = 1e+06)
## S4 method for signature 'methylRawDB,methylRawDB'
adjustMethylC(mc, hmc, save.db = TRUE,
    ..., chunk.size = 1e+06)
## S4 method for signature 'methylRawListDB,methylRawListDB'
adjustMethylC(mc, hmc,
    save.db = TRUE, ..., chunk.size = 1e+06)
```

## Arguments

mc	a methylRawList, methylRaw, methylRawDB or methylRawListDB containing 5mC levels of a sample or set of samples
hmc	a methylRawList, methylRaw, methylRawDB or methylRawListDB containing 5hmC levels of a sample or set of samples. If a methylRawList or methylRawListDB given the sample order should be same as "mc" methylRawList or methylRawListDB object.
save.db	A Logical to decide whether the resulting object should be saved as flat file database or not, default: explained in Details sections

4 adjustMethylC

... optional Arguments used when save.db is TRUE

suffix A character string to append to the name of the output flat file database, only used if save.db is true, default actions: append "\_filtered" to current file-name if database already exists or generate new file with filename "sampleID\_filtered" dbdir The directory where flat file database(s) should be stored, defaults to getwd(), working directory for newly stored databases and to same directory for already existing database

dbtype The type of the flat file database, currently only option is "tabix" (only used for newly stored databases)

used for newly stored databases)

Number of rows to be taken as a chunk for processing the methylRawDB or

methylRawListDB objects (default: 1e6)

#### Value

chunk.size

returns adjusted 5-methyl cytosine levels in the form of methylRawList, methylRaw, methylRawDB or methylRawListDB object depending on the input object

#### **Details**

The parameter chunk.size is only used when working with methylRawDB or methylRawListDB objects, as they are read in chunk by chunk to enable processing large-sized objects which are stored as flat file database. Per default the chunk.size is set to 1M rows, which should work for most systems. If you encounter memory problems or have a high amount of memory available feel free to adjust the chunk.size.

The parameter save. db is per default TRUE for methylDB objects as methylRawDB and methylRawListDB, while being per default FALSE for methylRaw and methylRawList. If you wish to save the result of an in-memory-calculation as flat file database or if the size of the database allows the calculation in-memory, then you might change the value of this parameter.

### References

- 1. Booth, Branco, et al. (2012). Quantitative Sequencing of 5-Methylcytosine and 5-Hydroxymethylcytosine at Single-Base Resolution. Science, 934
- 2. Yu, Hon, et al. (2012). Base-resolution analysis of 5-hydroxymethylcytosine in the Mammalian genome. Cell, 149(6), 1368-80.

```
# read 5hmC and 5mC files
hmc.file=system.file("extdata", "test1.myCpG.txt", package = "methylKit")
mc.file =system.file("extdata", "test2.myCpG.txt", package = "methylKit")
my5hmC=methRead( hmc.file,sample.id="hmc",assembly="hg18")
my5mC =methRead( mc.file,sample.id="mc",assembly="hg18")
# adjusting the 5mC levels using 5hmC levels
adjusted.5mC=adjustMethylC(my5mC,my5hmC)
```

assocComp 5

assocComp

Associate principal components with sample annotations

# **Description**

This function associates principal components with sample annotations such as age, gender, batch\_id. Can be used to detect which batch effects are associated with the variation in the methylation values.

### Usage

```
assocComp(mBase, sampleAnnotation)
```

### **Arguments**

mBase methylBase or methylBaseDB object with no NA values in the data part. sampleAnnotation

a data frame where columns are different annotations and rows are the samples, in the same order as in the methylBase object.

# Value

a named list of principal component matrix (named 'pcs'),

by principal components (named 'vars') and a p-value matrix showing association p-values between sample annotations and principal components (named 'association').

### Author(s)

Altuna Akalin

### **Examples**

```
data(methylKit)
sampleAnnotation=data.frame(batch_id=c("a","a","b","b"),age=c(19,34,23,40))
as=assocComp(mBase=methylBase.obj,sampleAnnotation)
```

bedgraph

Get bedgraph from methylRaw, methylRawList and methylDiff objects

# Description

The function converts methylRaw, methylRawDB, methylRawList, methylRawListDB, methylDiff or methylDiffDB object into a bedgraph format. It either writes as a file or returns a data.frame

6 bedgraph

#### Usage

```
bedgraph(methylObj, file.name = NULL, col.name, unmeth = FALSE,
  log.transform = FALSE, negative = FALSE, add.on = "",
  chunk.size = 1e+06)
## S4 method for signature 'methylDiff'
bedgraph(methylObj, file.name, col.name, unmeth,
  log.transform, negative, add.on)
## S4 method for signature 'methylRaw'
bedgraph(methylObj, file.name, col.name, unmeth,
  log.transform, negative, add.on)
## S4 method for signature 'methylRawList'
bedgraph(methylObj, file.name, col.name, unmeth,
 log.transform, negative, add.on)
## S4 method for signature 'methylRawDB'
bedgraph(methylObj, file.name = NULL, col.name,
  unmeth = FALSE, log.transform = FALSE, negative = FALSE, add.on = "",
  chunk.size = 1e+06)
## S4 method for signature 'methylRawListDB'
bedgraph(methylObj, file.name = NULL, col.name,
  unmeth = FALSE, log.transform = FALSE, negative = FALSE, add.on = "",
  chunk.size = 1e+06)
## S4 method for signature 'methylDiffDB'
bedgraph(methylObj, file.name, col.name, log.transform,
  negative, add.on, chunk.size)
```

or methylDiffDB object

# Arguments

methylObj

log.transform

negative

	or meetily 1511 155 coject
file.name	Default: NULL. if a string is given a bedgraph file will be written, if NULL a data.frame or a list of data frames will be returned
col.name	name of the column in methylRaw, methylRawDB, methylRawList, methylRawListDB, methylDiff or methylDiffDB objects to be used as a score for the bedgraph. For methylDiff or methylDiffDB, col.name must be one of the following 'pvalue', 'qvalue', 'meth.diff'. For methylRaw, methylRawDB, methylRawList and methylRawListDB it must be one of the following 'coverage', 'numCs', 'numTs', 'perc.meth'
unmeth	when working with methylRaw, methylRawDB, methylRawList and methylRawListDB objects should you output unmethylated C percentage this makes it easier to see the unmethylated bases because their methylation percentage values will be zero. Only invoked when file.name is not NULL.

scale. Ignored when col.name='perc.meth'

by -1. Ignored when col.name='perc.meth'

a methylRaw, methylRawDB, methylRawList, methylRawListDB, methylDiff

Default FALSE, If TRUE the score column of the bedgraph wil be in log10

Default FALSE, If TRUE, the score column of the bedgraph will be multiplied

calculateDiffMeth 7

add.on additional string to be add on the track line of bedgraph. can be viewlimits,priority etc. Check bedgraph track line options at UCSC browser

Chunk.size Number of rows to be taken as a chunk for processing the methylRawDB, methylRawListDB or methylDiffDB objects, default: 1e6

•

#### Value

Returns a data.frame or list of data.frames if file.name=NULL, if a file.name given appropriate bed file will be written to that file

#### **Details**

The parameter chunk.size is only used when working with methylRawDB, methylRawListDB or methylDiffDB objects, as they are read in chunk by chunk to enable processing large-sized objects which are stored as flat file database. Per default the chunk.size is set to 1M rows, which should work for most systems. If you encounter memory problems or have a high amount of memory available feel free to adjust the chunk.size.

### **Examples**

calculateDiffMeth

Calculate differential methylation statistics

### **Description**

The function calculates differential methylation statistics between two groups of samples. The function uses either logistic regression test or Fisher's Exact test to calculate differential methylation. See the rest of the help page and references for detailed explanation on statistics.

# Usage

```
calculateDiffMeth(.Object, covariates = NULL, overdispersion = c("none",
   "MN", "shrinkMN"), adjust = c("SLIM", "holm", "hochberg", "hommel",
   "bonferroni", "BH", "BY", "fdr", "none", "qvalue"), effect = c("wmean",
   "mean", "predicted"), parShrinkMN = list(), test = c("F", "Chisq"),
```

8 calculateDiffMeth

```
mc.cores = 1, slim = TRUE, weighted.mean = TRUE, chunk.size = 1e+06,
    save.db = FALSE, ...)

## S4 method for signature 'methylBaseDB'
calculateDiffMeth(.Object, covariates = NULL,
    overdispersion = c("none", "MN", "shrinkMN"), adjust = c("SLIM", "holm",
    "hochberg", "hommel", "bonferroni", "BH", "BY", "fdr", "none", "qvalue"),
    effect = c("wmean", "mean", "predicted"), parShrinkMN = list(),
    test = c("F", "Chisq"), mc.cores = 1, slim = TRUE,
    weighted.mean = TRUE, chunk.size = 1e+06, save.db = TRUE, ...)
```

#### **Arguments**

 $. Object \hspace{1.5cm} a \hspace{1mm} methyl Base \hspace{1mm} Ob \hspace{1mm} ject \hspace{1mm} to \hspace{1mm} calculate \hspace{1mm} differential \hspace{1mm} methyl ation$ 

covariates a data.frame containing covariates, which should be included in the test.

overdispersion If set to "none" (default), no overdispersion correction will be attempted. If set

to "MN", basic overdispersion correction, proposed by McCullagh and Nelder (1989) will be applied. This correction applies a scaling parameter to variance estimated by the model. EXPERIMENTAL: If set to "shrinkMN", scaling parameter will be shrunk towards a common value (not thoroughly tested as of

yet).

adjust different methods to correct the p-values for multiple testing. Default is "SLIM"

from methylKit. For "qvalue" please see qvalue and for all other methods see

p.adjust.

effect method to calculate the mean methylation different between groups using read

coverage as weights (default). When set to "mean", the generic mean is applied and when set to "predicted", predicted means from the logistic regression model

is used for calculating the effect.

parShrinkMN a list for squeezeVar(). (NOT IMPLEMENTED)

test the statistical test used to determine the methylation differences. The Chisq-test

is used by default, while the F-test can be chosen if overdispersion control ist

applied.

mc.cores integer denoting how many cores should be used for parallel differential methy-

lation calculations (can only be used in machines with multiple cores).

slim If set to FALSE, adjust will be set to "BH" (default behaviour of earlier ver-

sions)

weighted.mean If set to FALSE, effect will be set to "mean" (default behaviour of earlier

versions)

chunk.size Number of rows to be taken as a chunk for processing the methylBaseDB objects

(default: 1e6)

save.db A Logical to decide whether the resulting object should be saved as flat file

database or not, default: explained in Details section.

... optional Arguments used when save.db is TRUE

suffix A character string to append to the name of the output flat file database, only used if save.db is true, default actions: The default suffix is a 13-character

random string appended to the fixed prefix "methylDiff", e.g. "methylDiff\_16d3047c1a254.txt.bgz".

 $\label{thm:condition} \begin{tabular}{ll} dbdir\ The\ directory\ where\ flat\ file\ database(s)\ should\ be\ stored,\ defaults\ to\ getwd(),\ working\ directory\ for\ newly\ stored\ databases\ and\ to\ same\ directory\ for\ directory\ for\ defaults\ to\ getwd(),\ working\ directory\ for\ newly\ stored\ databases\ and\ to\ same\ directory\ for\ defaults\ to\ getwd().$ 

already existing database

calculateDiffMeth 9

dbtype The type of the flat file database, currently only option is "tabix" (only used for newly stored databases)

#### Value

a methylDiff object containing the differential methylation statistics and locations for regions or bases

### **Details**

Covariates can be included in the analysis. The function will then try to separate the influence of the covariates from the treatment effect via a linear model.

The Chisq-test is used per default only when no overdispersion correction is applied. If overdispersion correction is applied, the function automatically switches to the F-test. The Chisq-test can be manually chosen in this case as well, but the F-test only works with overdispersion correction switched on.

The parameter chunk.size is only used when working with methylBaseDB objects, as they are read in chunk by chunk to enable processing large-sized objects which are stored as flat file database. Per default the chunk.size is set to 1M rows, which should work for most systems. If you encounter memory problems or have a high amount of memory available feel free to adjust the chunk.size.

The parameter save.db is per default TRUE for methylDB objects as methylBaseDB, while being per default FALSE for methylBase. If you wish to save the result of an in-memory-calculation as flat file database or if the size of the database allows the calculation in-memory, then you might want to change the value of this parameter.

#### References

Altuna Akalin, Matthias Kormaksson, Sheng Li, Francine E. Garrett-Bakelman, Maria E. Figueroa, Ari Melnick, Christopher E. Mason. (2012). "methylKit: A comprehensive R package for the analysis of genome-wide DNA methylation profiles." Genome Biology.

McCullagh and Nelder. (1989). Generalized Linear Models. Chapman and Hall. London New York.

# See Also

```
pool, reorganize dataSim
```

10 calculateDiffMethDSS

calculateDiffMethDSS calculate Differential Methylation with DSS

# **Description**

This function provides an interface to the beta-binomial model from DSS package by Hao Wu. It calculates the differential methylation statistics using a beta-binomial model with parameter shrinkage. See the reference for details.

### Usage

```
calculateDiffMethDSS(meth, adjust = c("SLIM", "holm", "hochberg", "hommel",
   "bonferroni", "BH", "BY", "fdr", "none", "qvalue"), mc.cores = 1)

## S4 method for signature 'methylBase'
calculateDiffMethDSS(meth, adjust = c("SLIM", "holm",
   "hochberg", "hommel", "bonferroni", "BH", "BY", "fdr", "none", "qvalue"),
   mc.cores = 1)
```

### **Arguments**

meth a methylBase object

adjust different methods to correct the p-values for multiple testing. Default is "SLIM"

from methylKit. For "qvalue" please see qvalue and for all other methods see

p.adjust.

mc.cores integer denoting how many cores should be used for parallel differential methy-

lation calculations (can only be used in machines with multiple cores).

### Value

a methylDiff object

### References

Feng H, Conneely K and Wu H (2014). A bayesian hierarchical model to detect differentially methylated loci from single nucleotide resolution sequencing data. Nucleic acids research

clusterSamples 11

# **Examples**

```
data(methylKit)
dssDiffay <- calculateDiffMethDSS(methylBase.obj, adjust="SLIM", mc.cores=1)</pre>
```

clusterSamples

Hierarchical Clustering using methylation data The function clusters samples using hclust function and various distance metrics derived from percent methylation per base or per region for each sample.

# **Description**

Hierarchical Clustering using methylation data

The function clusters samples using hclust function and various distance metrics derived from percent methylation per base or per region for each sample.

## Usage

# Arguments

.Object	a methylBase or methylBaseDB object
dist	the distance measure to be used. This must be one of "correlation", "euclidean", "maximum", "manhattan", "canberra", "binary" or "minkowski". Any unambiguous abbreviation can be given. (default:"correlation")
method	the agglomeration method to be used. This should be (an unambiguous abbreviation of) one of "ward", "single", "complete", "average", "mcquitty", "median" or "centroid". (default:"ward")
sd.filter	If TRUE, the bases/regions with low variation will be discarded prior to clustering (default:TRUE)
sd.threshold	A numeric value. If filterByQuantile is TRUE, features whose standard deviations is less than the quantile denoted by sd.threshold will be removed. If filterByQuantile is FALSE, then features whose standard deviations is less than the value of sd.threshold will be removed.(default:0.5)

12 dataSim

filterByQuantile

A logical determining if sd.threshold is to be interpreted as a quantile of all Standard Deviation values from bases/regions (the default), or as an absolute

value

plot a logical value indicating whether to plot hierarchical clustering. (default:TRUE)

chunk.size Number of rows to be taken as a chunk for processing the methylBaseDB ob-

jects, default: 1e6

#### Value

a tree object of a hierarchical cluster analysis using a set of dissimilarities for the n objects being clustered.

#### Details

The parameter chunk.size is only used when working with methylBaseDB objects, as they are read in chunk by chunk to enable processing large-sized objects which are stored as flat file database. Per default the chunk.size is set to 1M rows, which should work for most systems. If you encounter memory problems or have a high amount of memory available feel free to adjust the chunk.size.

### **Examples**

```
data(methylKit)
clusterSamples(methylBase.obj, dist="correlation", method="ward", plot=TRUE)
```

dataSim

Simulate DNA methylation data

# Description

The function simulates DNA methylation data from multiple samples. See references for detailed explanation on statistics.

#### Usage

```
dataSim(replicates, sites, treatment, percentage = 10, effect = 25,
   alpha = 0.4, beta = 0.5, theta = 10, covariates = NULL,
   sample.ids = NULL, assembly = "hg18", context = "CpG",
   add.info = FALSE)
```

### **Arguments**

replicates the number of samples that should be simulated.

sites the number of CpG sites per sample.
treatment a vector containing treatment information.

percentage the proportion of sites which should be affected by the treatment.

dataSim 13

effect	a number between 0 and 100 specifying the effect size of the treatment. This is essentially describing the average percent methylation difference between differentially methylated bases. See 'Examples' and 'Details'.
alpha	shape1 parameter for beta distribution (used for initial sampling of methylation proportions)
beta	shape2 parameter for beta distribution (used for initial sampling of methylation proportions)
theta	dispersion parameter for beta distribution (initial sampling of methylation proportions)
covariates	a data.frame containing covariates (optional)
sample.ids	will be generated automatically from treatment, but can be overwritten by a character vector containing sample names.
assembly	the assembly description (e.g. "hg18"). Only needed for book keeping.
context	the experimanteal context of the data (e.g. "CpG"). Only needed for book keeping.
add.info	if set to TRUE, the output will be a list with the first element being the methyl- base object and a vector of indices that indicate which CpGs should be differ- entially methylated. This vector can be used to subset simulated methylBase or methylDiff object with differentially methylated bases.

#### Value

a methylBase object containing simulated methylation data, or if add.info=TRUE a list containing the methylbase object and the indices of all treated sites (differentially methylated bases or regions) as the second element.

# **Details**

The function uses a Beta distribution to simulate the methylation proportion background across all samples. The parameters alpha, beta used in a beta distribution to draw methylation proportions,  $\mu$ , from a typical bimodal distribution. For each initial methylation proportion drawn using the parameters above, a range of methylation proportions is distributed around the original  $\mu$  with overdispersion parameter  $\theta$ , this is using an alternative parameterization of Beta distribution:  $Beta(\mu,\theta)$ . The parameters percentage and effect determine the proportion of sites that are affected by the treatment (meaning differential sites) and the strength of this influence, respectively. effect is added on top of  $\mu$  for the CpGs that are affected by the treament. The affected group of samples for that particular CpG will now be distributed by  $Beta(\mu + effect, \theta)$ . The coverage is modeled with a negative binomial distribution, using rnbinom function with size=1 and prob=0.01. The additional information needed for a valid methylBase object, such as CpG start, end and strand, is generated as "dummy values", but can be overwritten as needed.

```
data(methylKit)
# Simulate data for 4 samples with 20000 sites each.
# The methylation in 10% of the sites are elevated by 25%.
my.methylBase=dataSim(replicates=4,sites=2000,treatment=c(1,1,0,0),
percentage=10,effect=25)
```

14 diffMethPerChr

diffMethPerChr	Get and plot the number of hyper/hypo methylated regions/bases per
	chromosome

## Description

This function gets number of hyper/hypo methylated regions/bases from methylDiff object. It can also plot percentages of differentially methylated bases per chromosome.

# Usage

```
diffMethPerChr(x, plot = TRUE, qvalue.cutoff = 0.01, meth.cutoff = 25,
   exclude = NULL, ...)

## S4 method for signature 'methylDiff'
diffMethPerChr(x, plot = TRUE, qvalue.cutoff = 0.01,
   meth.cutoff = 25, exclude = NULL, ...)

## S4 method for signature 'methylDiffDB'
diffMethPerChr(x, plot = TRUE,
   qvalue.cutoff = 0.01, meth.cutoff = 25, exclude = NULL, ...)
```

# Arguments

x	a methylDiff object
plot	TRUEIFALSE. If TRUE horizontal barplots for proportion of hypo/hyper methylated bases/regions
qvalue.cutoff	cutoff for q-value
meth.cutoff	cutoff for percent methylation difference
exclude	names of chromosomes to be excluded
	extra graphical parameters to be passed to barplot function

### Value

plots a piechart or a barplot for percentage of the target features overlapping with annotation

extract 15

extract	$extract\ parts\ of\ methyl Raw, methyl Raw DB, methyl Base, methyl Base DB$
	and methylDiff data

# Description

The function extracts part of the data and returns a new object.

# Usage

```
## S4 method for signature 'methylRaw,ANY,ANY'
x[i, j]

## S4 method for signature 'methylBase,ANY,ANY,ANY'
x[i, j]

## S4 method for signature 'methylDiff,ANY,ANY,ANY'
x[i, j]

## S4 method for signature 'methylRawDB,ANY,ANY,ANY'
x[i, j]

## S4 method for signature 'methylBaseDB,ANY,ANY,ANY'
x[i, j]

## S4 method for signature 'methylDiffDB,ANY,ANY,ANY'
x[i, j]
```

### **Arguments**

X	an $methylBase, methylBaseDB, methylRaw, methylRawDB$ or $methylDiff$ object
i	a numeric or logical vector. This vector corresponds to bases or regions contained in methylKit objects. The vector is used to subset the data.
j	This argument can not be used for the extraction of columns. As unintentional extraction of the columns will cause an error in the downstream analysis. Using this argument will cause an error. Use getData to access the data part of the objects.

```
data(methylKit)
# selects first hundred rows, returns a methylRaw object
subset1=methylRawList.obj[[1]][1:100]
# selects first hundred rows, returns a methylBase object
subset2=methylBase.obj[1:100,]
# selects first hundred rows, returns a methylDiff object
subset3=methylDiff.obj[1:100,]
```

16 filterByCoverage

```
# This will get chromomsomes, will return a factor
# That means the resulting object will ceases to be a methylKit object
chrs=methylDiff.obj[[2]]
```

filterByCoverage

Filter methylRaw, methylRawDB, methylRawList and methyl-RawListDB object based on read coverage

### **Description**

This function filters methylRaw, methylRawDB, methylRawList and methylRawListDB objects. You can filter based on lower read cutoff or high read cutoff. Higher read cutoff is usefull to eliminate PCR effects Lower read cutoff is usefull for doing better statistical tests.

### Usage

```
filterByCoverage(methylObj, lo.count = NULL, lo.perc = NULL,
 hi.count = NULL, hi.perc = NULL, chunk.size = 1e+06, save.db = FALSE,
  ...)
## S4 method for signature 'methylRaw'
filterByCoverage(methylObj, lo.count = NULL,
  lo.perc = NULL, hi.count = NULL, hi.perc = NULL, chunk.size = 1e+06,
  save.db = FALSE, ...)
## S4 method for signature 'methylRawList'
filterByCoverage(methylObj, lo.count = NULL,
  lo.perc = NULL, hi.count = NULL, hi.perc = NULL, chunk.size = 1e+06,
  save.db = FALSE, ...)
## S4 method for signature 'methylRawDB'
filterByCoverage(methylObj, lo.count = NULL,
 lo.perc = NULL, hi.count = NULL, hi.perc = NULL, chunk.size = 1e+06,
 save.db = TRUE, ...)
## S4 method for signature 'methylRawListDB'
filterByCoverage(methylObj, lo.count = NULL,
 lo.perc = NULL, hi.count = NULL, hi.perc = NULL, chunk.size = 1e+06,
 save.db = TRUE, ...)
```

# Arguments

methylObj	a methylRaw, methylRawDB, methylRawList or methylRawListDB object
lo.count	An integer for read counts.Bases/regions having lower coverage than this count is discarded
lo.perc	A double [0-100] for percentile of read counts. Bases/regions having lower coverage than this percentile is discarded

filterByCoverage 17

hi.count	An integer for read counts. Bases/regions having higher coverage than this is count discarded
hi.perc	A double [0-100] for percentile of read counts. Bases/regions having higher coverage than this percentile is discarded
chunk.size	Number of rows to be taken as a chunk for processing the methylRawDB or methylRawListDB objects, default: 1e6
save.db	A Logical to decide whether the resulting object should be saved as flat file database or not, default: explained in Details sections
	optional Arguments used when save.db is TRUE suffix A character string to append to the name of the output flat file database, only used if save.db is true, default actions: append "_filtered" to current file-name if database already exists or generate new file with filename "sampleID_filtered" dbdir The directory where flat file database(s) should be stored, defaults to getwd(), working directory for newly stored databases and to original directory for already existing database
	dbtype The type of the flat file database, currently only option is "tabix" (only used for newly stored databases)

#### Value

methylRaw, methylRawDB, methylRawList or methylRawListDB object depending on input object

### **Details**

The parameter chunk.size is only used when working with methylRawDB or methylRawListDB objects, as they are read in chunk by chunk to enable processing large-sized objects which are stored as flat file database. Per default the chunk.size is set to 1M rows, which should work for most systems. If you encounter memory problems or have a high amount of memory available feel free to adjust the chunk.size.

The parameter save. db is per default TRUE for methylDB objects as methylRawDB and methylRawListDB, while being per default FALSE for methylRaw and methylRawList. If you wish to save the result of an in-memory-calculation as flat file database or if the size of the database allows the calculation in-memory, then you might change the value of this parameter.

18 getAssembly

```
# tidy up
rm(filtered3)
unlink("methylDB",recursive=TRUE)
```

```
{\tt genomation-deprecated}\ \ \textit{Deprecated/Defunct functions}
```

# Description

These are deprecated or defunct functions. Most of them are replaced by genomation functions. See the vignette for examples on how to use genomation functions for annotation purposes.

# Usage

```
annotate.WithFeature()
annotate.WithFeature.Flank()
annotate.WithGenicParts()
read.bed()
read.feature.flank()
read.transcript.features()
getFeatsWithTargetsStats()
getFlanks()
getMembers()
getTargetAnnotationStats()
plotTargetAnnotation()
```

getAssembly

get assembly of the genome

# Description

The function returns the genome assembly stored in any of the methylBase,methylBaseDB,methylRaw,methylDiff objects

getContext 19

### Usage

```
getAssembly(x)
## S4 method for signature 'methylBase'
getAssembly(x)
## S4 method for signature 'methylRaw'
getAssembly(x)
## S4 method for signature 'methylDiff'
getAssembly(x)
## S4 method for signature 'methylRawDB'
getAssembly(x)
## S4 method for signature 'methylBaseDB'
getAssembly(x)
## S4 method for signature 'methylBaseDB'
getAssembly(x)
```

# **Arguments**

x an methylBase,methylBaseDB, methylRaw,methylRawDB or methylDiff object

# Value

the assembly string for the object

# **Examples**

```
data(methylKit)
getAssembly(methylBase.obj)
getAssembly(methylDiff.obj)
getAssembly(methylRawList.obj[[1]])
```

getContext

get the context of methylation

# **Description**

The function returns the context of methylation. For example: "CpG", "CHH" or "CHG"

20 getCorrelation

#### Usage

```
getContext(x)

## S4 method for signature 'methylBase'
getContext(x)

## S4 method for signature 'methylRaw'
getContext(x)

## S4 method for signature 'methylDiff'
getContext(x)

## S4 method for signature 'methylRawDB'
getContext(x)

## S4 method for signature 'methylBaseDB'
getContext(x)

## S4 method for signature 'methylBaseDB'
getContext(x)
```

# **Arguments**

Х

an methylBase,methylBaseDB, methylRaw,methylRawDB or an methylDiff object

### Value

a string for the context methylation

# **Examples**

```
data(methylKit)
getContext(methylBase.obj)
getContext(methylDiff.obj)
getContext(methylRawList.obj[[1]])
```

 ${\tt getCorrelation}$ 

get correlation between samples in methylBase or methylBaseDB object

# Description

The functions returns a matrix of correlation coefficients and/or a set of scatterplots showing the relationship between samples. The scatterplots will contain also fitted lines using lm() for linear regression and lowess for polynomial regression.

getCorrelation 21

#### Usage

```
getCorrelation(object,method="pearson",plot=FALSE,nrow)
## S4 method for signature 'methylBase'
getCorrelation(object, method, plot)

## S4 method for signature 'methylBaseDB'
getCorrelation(object, method = "pearson",
    plot = FALSE, nrow = 2e+06)
```

## **Arguments**

object a methylBase or methylBaseDB object

method a character string indicating which correlation coefficient (or covariance) is to

be computed (default: "pearson", other options are "kendall" and "spearman")

plot scatterPlot if TRUE (default:FALSE)

nrow a numeric giving the number of lines to read in of methylBaseDB object, de-

faults to 2e6

### Value

a correlation matrix object and plot scatterPlot

#### **Details**

The argument 'nrow' is only evaluated if the input is a methylBaseDB object. If 'nrow' is not specified getCorrelation will read the first 2M records of the given object, but if you want to read all records 'nrow' has to be NULL. You should change 'nrow' if using getCorrelation with all records of the methylBaseDB object would take too long.

If the scatter plot is plotted, the red line in the plot is from linear regression fit and the green line is from polynomial regression fit with stats::lowess.

```
data(methylKit)
getCorrelation(methylBase.obj,method="pearson",plot=FALSE)

# create methylBaseDB
methylBaseDB.obj <- unite(methylRawList.obj,save.db=TRUE,dbdir="methylDB")
getCorrelation(methylBaseDB.obj,method="pearson",plot=FALSE,nrow=10000)

# remove Database again
rm(methylBaseDB.obj)
unlink("methylDB",recursive=TRUE)</pre>
```

22 getCoverageStats

# Description

The function returns basic statistics about read coverage per base. It can also plot a histogram of read coverage values.

# Usage

```
getCoverageStats(object, plot = FALSE, both.strands = FALSE,
    labels = TRUE, ..., chunk.size = 1e+06)

## S4 method for signature 'methylRaw'
getCoverageStats(object, plot = FALSE,
    both.strands = FALSE, labels = TRUE, ..., chunk.size = 1e+06)

## S4 method for signature 'methylRawDB'
getCoverageStats(object, plot = FALSE,
    both.strands = FALSE, labels = TRUE, ..., chunk.size = 1e+06)
```

### **Arguments**

object a methylRaw or methylRawDB object

plot plot a histogram of coverage if TRUE (default:FALSE) both.strands do stats and plot for both strands if TRUE (default:FALSE)

labels should the bars of the histrogram have labels showing the percentage of values

in each bin (default:TRUE)

... options to be passed to hist function

chunk.size Number of rows to be taken as a chunk for processing the methylRawDB objects

(default: 1e6)

#### Value

a summary of coverage statistics or plot a histogram of coverage

### **Details**

The parameter chunk.size is only used when working with methylRawDB or methylRawListDB objects, as they are read in chunk by chunk to enable processing large-sized objects which are stored as flat file database. Per default the chunk.size is set to 1M rows, which should work for most systems. If you encounter memory problems or have a high amount of memory available feel free to adjust the chunk.size.

```
data(methylKit)
# gets coverage stats for the first sample in methylRawList.obj object
getCoverageStats(methylRawList.obj[[1]],plot=TRUE,
both.strands=FALSE,labels=TRUE)
```

getData 23

getData

get the data slot from the methylKit objects

### **Description**

The functions retrieves the table containing methylation information from methylKit Objects. The data retrived from this function is of a data. frame. This is basically containing all relevant methylation information per genomic region or base.

# Usage

```
getData(x)
## S4 method for signature 'methylBase'
getData(x)
## S4 method for signature 'methylRaw'
getData(x)
## S4 method for signature 'methylDiff'
getData(x)
## S4 method for signature 'methylRawDB'
getData(x)
## S4 method for signature 'methylRawDB'
getData(x)
## S4 method for signature 'methylBaseDB'
getData(x)
## S4 method for signature 'methylDiffDB'
getData(x)
```

# **Arguments**

x an methylBase,methylBaseDB, methylRaw,methylRawDB or methylDiff object

### Value

data frame for methylation events

```
data(methylKit)
# following commands show first lines of returned
# data.frames from getData() function
head(
getData(methylBase.obj)
)
head( getData(methylDiff.obj))
```

24 getDBPath

```
head(getData(methylRawList.obj[[1]]))
```

getDBPath

Get path to database of the methylDB objects

# Description

The function returns the path to the flat file database that stores the data of the methylRawDB, methylRawListDB, methylBaseDB or methylDiffDB objects.

# Usage

```
getDBPath(x)
## S4 method for signature 'methylRawListDB'
getDBPath(x)
## S4 method for signature 'methylBaseDB'
getDBPath(x)
## S4 method for signature 'methylRawDB'
getDBPath(x)
## S4 method for signature 'methylDiffDB'
getDBPath(x)
```

### **Arguments**

Х

an methylBaseDB,methylRawDB, methylRawListDB or methylDiffDB object

```
data(methylKit)
methylBaseDB.obj <- unite(methylRawList.obj,save.db=TRUE,dbdir="methylDB")
#The path to the database is returned
getDBPath(methylBaseDB.obj)
# remove Database again
rm(methylBaseDB.obj)
unlink("methylDB",recursive=TRUE)</pre>
```

getMethylationStats 25

getMethylationStats get Methylation stats from methylRaw or methylRawDB object

## **Description**

The function returns basic statistics about It can also plot a histogram of

# Usage

```
getMethylationStats(object, plot = FALSE, both.strands = FALSE,
  labels = TRUE, ..., chunk.size = 1e+06)

## S4 method for signature 'methylRaw'
getMethylationStats(object, plot = FALSE,
  both.strands = FALSE, labels = TRUE, ..., chunk.size = 1e+06)

## S4 method for signature 'methylRawDB'
getMethylationStats(object, plot = FALSE,
  both.strands = FALSE, labels = TRUE, ..., chunk.size = 1e+06)
```

# **Arguments**

object a methylRaw or methylRawDB object

plot plot a histogram of Methylation if TRUE (deafult:FALSE)

both.strands do plots and stats for both strands seperately if TRUE (deafult:FALSE)

labels should the bars of the histrogram have labels showing the percentage of values

in each bin (default:TRUE)

... options to be passed to hist function.

chunk.size Number of rows to be taken as a chunk for processing the methylRawDB objects

(default: 1e6)

# Value

a summary of Methylation statistics or plot a histogram of coverage

#### **Details**

The parameter chunk.size is only used when working with methylRawDB or methylRawListDB objects, as they are read in chunk by chunk to enable processing large-sized objects which are stored as flat file database. Per default the chunk.size is set to 1M rows, which should work for most systems. If you encounter memory problems or have a high amount of memory available feel free to adjust the chunk.size.

```
data(methylKit)
# gets Methylation stats for the first sample in methylRawList.obj object
getMethylationStats(methylRawList.obj[[1]],plot=TRUE,
both.strands=FALSE,labels=TRUE)
```

26 getMethylDiff

getMethylDiff

get differentially methylated regions/bases based on cutoffs

### **Description**

The function subsets a methylDiff or methylDiffDB object in order to get differentially methylated bases/regions satisfying thresholds.

### Usage

```
getMethylDiff(.Object, difference = 25, qvalue = 0.01, type = "all",
    chunk.size = 1e+06, save.db = FALSE, ...)

## S4 method for signature 'methylDiff'
getMethylDiff(.Object, difference = 25,
    qvalue = 0.01, type = "all", chunk.size = 1e+06, save.db = FALSE, ...)

## S4 method for signature 'methylDiffDB'
getMethylDiff(.Object, difference = 25,
    qvalue = 0.01, type = "all", chunk.size = 1e+06, save.db = TRUE, ...)
```

# **Arguments**

a methylDiff or methylDiffDB object .Object difference cutoff for absolute value of methylation percentage change between test and control (default:25) qvalue cutoff for qvalue of differential methylation statistic (default:0.01) one of the "hyper", "hypo" or "all" strings. Specifies what type of differentially type menthylated bases/regions should be returned. For retrieving Hyper-methylated regions/bases type="hyper", for hypo-methylated type="hypo" (default: "all") chunk.size Number of rows to be taken as a chunk for processing the methylDiffDB objects (default: 1e6) A Logical to decide whether the resulting object should be saved as flat file save.db database or not, default: see Details optional Arguments used when save.db is TRUE

suffix A character string to append to the name of the output flat file database, only used if save.db is true, default actions: append "\_type" to current filename if database already exists or generate new file with filename "methylDiff\_type" dbdir The directory where flat file database(s) should be stored, defaults to getwd(), working directory for newly stored databases and to same directory for already existing database

dbtype The type of the flat file database, currently only option is "tabix" (only

#### Value

a methylDiff or methylDiffDB object containing the differential methylated locations satisfying the criteria

used for newly stored databases)

getSampleID 27

#### **Details**

The parameter chunk.size is only used when working with methylDiffDB objects, as they are read in chunk by chunk to enable processing large-sized objects which are stored as flat file database. Per default the chunk.size is set to 1M rows, which should work for most systems. If you encounter memory problems or have a high amount of memory available feel free to adjust the chunk.size.

The parameter save.db is per default TRUE for methylDB objects as methylDiffDB, while being per default FALSE for methylDiff. If you wish to save the result of an in-memory-calculation as flat file database or if the size of the database allows the calculation in-memory, then you might want to change the value of this parameter.

### **Examples**

```
data(methylKit)
# get differentially methylated bases/regions with specific cutoffs
all.diff=getMethylDiff(methylDiff.obj,difference=25,qvalue=0.01,type="all")
# get hyper-methylated
hyper=getMethylDiff(methylDiff.obj,difference=25,qvalue=0.01,type="hyper")
# get hypo-methylated
hypo=getMethylDiff(methylDiff.obj,difference=25,qvalue=0.01,type="hypo")
```

getSampleID

Get or Set Sample-IDs of the methylKit objects

### **Description**

The function returns or replaces the sample-ids stored in any of the following methylKit objects: methylRaw, methylRawDB, methylBase, methylBaseDB, methylRawList, methylRawListDB, methylDiff, methylDiffDB.

### Usage

```
getSampleID(x)
getSampleID(x) <- value

getSampleID(x) <- value

## S4 method for signature 'methylRawList'
getSampleID(x)

## S4 replacement method for signature 'methylRawList'
getSampleID(x) <- value

## S4 method for signature 'methylBase'
getSampleID(x)

## S4 replacement method for signature 'methylBase'</pre>
```

28 getSampleID

```
getSampleID(x) \leftarrow value
## S4 method for signature 'methylRaw'
getSampleID(x)
## S4 replacement method for signature 'methylRaw'
getSampleID(x) \leftarrow value
## S4 method for signature 'methylDiff'
getSampleID(x)
## S4 replacement method for signature 'methylDiff'
getSampleID(x) \leftarrow value
## S4 method for signature 'methylRawListDB'
getSampleID(x)
## S4 replacement method for signature 'methylRawListDB'
getSampleID(x) \leftarrow value
## S4 method for signature 'methylBaseDB'
getSampleID(x)
## S4 replacement method for signature 'methylBaseDB'
getSampleID(x) \leftarrow value
## S4 method for signature 'methylRawDB'
getSampleID(x)
## S4 replacement method for signature 'methylRawDB'
getSampleID(x) \leftarrow value
## S4 method for signature 'methylDiffDB'
getSampleID(x)
## S4 replacement method for signature 'methylDiffDB'
getSampleID(x) \leftarrow value
```

### **Arguments**

x an methylBaseDB,methylRawListDB or methylDiffDB object value a valid replacement vector for the sample-ids of the object

```
data(methylKit)
#The Sample-Ids can be printed ..
getSampleID(methylBase.obj)
# .. or replaced.
newObj <- methylBase.obj</pre>
```

getTreatment 29

```
getSampleID(newObj) <- c("sample1","sample2","sample3","sample4")
getSampleID(newObj)</pre>
```

getTreatment

Get or Set treatment vector of methylKit object

## **Description**

The function returns or replaces the treatment vector stored in any of the following methylKit objects: methylBase,methylRawList,methylBaseDB, methylRawListDB,methylDiff,methylDiffDB.

### Usage

```
getTreatment(x)
getTreatment(x) <- value</pre>
getTreatment(x) <- value</pre>
## S4 method for signature 'methylRawList'
getTreatment(x)
## S4 replacement method for signature 'methylRawList'
getTreatment(x) <- value</pre>
## S4 method for signature 'methylBase'
getTreatment(x)
## S4 replacement method for signature 'methylBase'
getTreatment(x) <- value</pre>
## S4 method for signature 'methylDiff'
getTreatment(x)
## S4 replacement method for signature 'methylDiff'
getTreatment(x) <- value</pre>
## S4 method for signature 'methylRawListDB'
getTreatment(x)
## S4 replacement method for signature 'methylRawListDB'
getTreatment(x) <- value</pre>
## S4 method for signature 'methylBaseDB'
getTreatment(x)
## S4 replacement method for signature 'methylBaseDB'
getTreatment(x) <- value</pre>
## S4 method for signature 'methylDiffDB'
```

```
getTreatment(x)
## S4 replacement method for signature 'methylDiffDB'
getTreatment(x) <- value</pre>
```

# **Arguments**

x a methylKit object

value a valid replacement for the treatment vector of the object

# **Examples**

```
data(methylKit)
# The treatment vector can be printed ..
getTreatment(methylBase.obj)
# .. or replaced with a new one
newObj <- methylBase.obj
getTreatment(newObj) <- c(1,2,3,4)
getTreatment(newObj)</pre>
```

joinSegmentNeighbours Join directly neighbouring segments produced by methSeg

# **Description**

Segmentation and clustering are two distinct steps in methSeg(), leading to adjacent segments of the same class. This leads to a bias segment length distributions, which is removed by joining those neighbours.

# Usage

```
joinSegmentNeighbours(res)
```

## Arguments

res A GRanges object with segment classification and information prudoced by

methSeg

# Value

A GRanges object with segment classification and information.

### See Also

methSeg

makeMethylDB 31

makeMethylDB

coerce methylKit objects from memory to flat file database objects

# Description

The function converts in-memory methylKit objects to methylDB objects

# Usage

```
makeMethylDB(obj, dbdir = getwd())
## S4 method for signature 'methylBase'
makeMethylDB(obj, dbdir = getwd())
## S4 method for signature 'methylRaw'
makeMethylDB(obj, dbdir = getwd())
## S4 method for signature 'methylDiff'
makeMethylDB(obj, dbdir = getwd())
## S4 method for signature 'methylRawList'
makeMethylDB(obj, dbdir = getwd())
```

# **Arguments**

obj an methylBase, methylRaw, methylRawList or an methylDiff object dbdir directory where flat file database(s) should be stored, defaults to getwd(), work-

ing directory.

# Value

an methylBaseDB,methylRawDB, methylRawListDB or an methylDiffDB object

```
## Not run:
data(methylKit)

makeMethylDB(methylBase.obj,"my/path")
## End(Not run)
```

32 methRead

methRead

read file(s) to methylRaw or methylRawList objects

## **Description**

The function reads a list of files or single files with methylation information for bases/region in the genome and creates a methylrawList or methylraw object. The information can be stored as flat file database by creating a methylrawlistDB or methylrawDB object.

### Usage

```
methRead(location, sample.id, assembly, dbtype = NA, pipeline = "amp",
  header = TRUE, skip = 0, sep = "\t", context = "CpG",
  resolution = "base", treatment, dbdir = getwd(), mincov = 10)

## S4 method for signature 'character, character, character'
methRead(location, sample.id,
  assembly, dbtype, pipeline, header, skip, sep, context, resolution, dbdir,
  mincov)

## S4 method for signature 'list,list, character'
methRead(location, sample.id, assembly,
  dbtype = NA, pipeline = "amp", header = TRUE, skip = 0, sep = "\t",
  context = "CpG", resolution = "base", treatment, dbdir = getwd(),
  mincov = 10)
```

# **Arguments**

sep

location

location	file location(s), either a list of locations (each a character string) or one location string
sample.id	sample.id(s)
assembly	a string that defines the genome assembly such as hg18, mm9. this is just a string for book keeping. It can be any string. Although, when using multiple files from the same assembly, this string should be consistent in each object.
dbtype	type of the flat file database, currently only option other than NA is "tabix". When "tabix" is given the objects are stored in tabix files, which are compressed and indexed. The default value is NA, in which case the objects are stored in memory.
pipeline	name of the alignment pipeline, it can be either "amp", "bismark", "bismarkCoverage", "bismarkCytosineReport" or a list (default:'amp'). The methylation text files generated from other pipelines can be read as generic methylation text files by supplying a named list argument as "pipeline" argument. The named list should containt column numbers which denotes which column of the text file corresponds to values and genomic location of the methylation events. See Details for more on possible values for this argument.
header	if the input file has a header or not (default: TRUE)
skip	number of lines to skip when reading. Can be set to 1 for bed files with track

seperator between fields, same as read. table argument (default: "\t")

line (default: 0)

file location(s) either a list of locations (each a character string) or one location

methRead 33

context methylation context string, ex: CpG,CpH,CHH, etc. (default:CpG)

resolution designates whether methylation information is base-pair resolution or regional

resolution. allowed values 'base' or 'region'. Default 'base'

treatment a vector containing 0 and 1 denoting which samples are control which samples

are test

dbdir directory where flat file database(s) should be stored, defaults to getwd(), work-

ing directory.

mincov minimum read coverage to be included in the methylKit objects. defaults to

10. Any methylated base/region in the text files below the mincov value will be

ignored.

#### Value

returns methylRaw, methylRawList, methylRawDB, methylRawListDB object

#### **Details**

The output of methRead is determined by specific input arguments, as there are location, sample.id, assembly and dbtype. The first three are obligatory, while if the last argument is given database features are enabled. If then location refers to an uncompressed file the function will create a flat file database and the associated methylRawDB object will link to this database. If then location refers to an earlier created database file then the object will directly link to this database, skipping the preprocessing steps.

When pipeline argument is a list, it is exptected to provide a named list with following names. 'fraction' is a logical value, denoting if the column frequency of Cs has a range from [0-1] or [0-100]. If true it assumes range is [0-1]. 'chr.col" is the number of the column that has chrosome string. 'start.col' is the number of the column that has start coordinate of the base/region of the methylation event. 'end.col' is the number of the column that has end coordinate of the base/region of the methylation event. 'coverage.col' is the number of the column that has read coverage values. 'strand.col' is the number of the column that has strand information in the file has to be in the form of '+' or '-', 'freqC.col' is the number of the column that has the frequency of Cs. See examples to see how to read a generic methylation text file.

Other possible values for pipeline argument are 'amp', 'bismark', 'bismarkCoverage' and 'bismarkCytosineReport'. For 'amp' and 'bismark' the function expects a tabular format shown in the webpage (http://github.com/al2na/methylKit). "amp" and "bismark" expect identical input and are kept for historical reasons. 'amp' was a pipeline used in Akalin et al. 2012 Plos Genetics paper, publicly available in googlecode.

Bismark aligner can output methylation information per base in multiple formats. With pipeline='bismarkCoverage', the function reads bismark coverage files, which have chr,start,end, number of cytosines (methylated bases) and number of thymines (unmethylated bases) format. If bismark coverage files are used the function will not have the strand information,so beware of that fact. With pipeline='bismarkCytosineReport', the function expects cytosine report files from Bismark, which have chr,start, strand, number of cytosines (methylated bases), number of thymines (unmethylated bases),context and trinucletide context format.

The function can also read gzipped files. On unix systems, this is achieved by using zcat filename and feeding that into data.table::fread. On Windows, the file is first uncompressed then read into R using data.table::fread.

34 methRead

```
# this is a list of example files, ships with the package
# for your own analysis you will just need to provide set of paths to files
# you will not need the "system.file(..." part
file.list=list(
         system.file("extdata", "test1.myCpG.txt", package = "methylKit"),
         system.file("extdata", "test2.myCpG.txt", package = "methylKit"),
         system.file("extdata", "control1.myCpG.txt", package = "methylKit"),
         system.file("extdata", "control2.myCpG.txt", package = "methylKit")
# read the files to a methylRawList object: myobj
myobj=methRead(file.list,
            sample.id=list("test1","test2","ctrl1","ctrl2"),
            assembly="hg18", treatment=c(1,1,0,0))
# read one file as methylRaw object
myobj=methRead( file.list[[1]],
            sample.id="test1",assembly="hg18")
# read a generic text file containing CpG methylation values
# let's first look at the content of the file
generic.file=system.file("extdata", "generic1.CpG.txt",package = "methylKit")
read.table(generic.file,header=TRUE)
# And this is how you can read that generic file as a methylKit object
myobj=methRead( generic.file,
             pipeline=list(fraction=FALSE,chr.col=1,start.col=2,end.col=2,
                           coverage.col=4,strand.col=3,freqC.col=5),
             sample.id="test1",assembly="hg18")
# This creates tabix files that save methylation data
# Without specified dbdir first creates a folder named the following
# in working directory:
# paste("methylDB",Sys.Date(),paste(sample(c(0:9, letters, LETTERS),3,
# replace=TRUE),collapse=""))
# Then, saves tabix files from methylKit objects there
myobj=methRead( file.list,
               sample.id=list("test1","test2","ctrl1","ctrl2"),
               assembly="hg18", treatment=c(1,1,0,0),
               dbtype="tabix")
# This creates a single tabix files that saves methylation data
# first creates a "methylDB_objects" directory
# Then, saves tabix file from methylKit objects there
myobj=methRead(file.list[[1]],
               sample.id="test1",
               assembly="hg18",
               dbtype="tabix",dbdir="methylDB_objects")
 # tidy up
 rm(myobj)
 unlink(list.files(pattern = "methylDB",full.names = TRUE),recursive = TRUE)
```

methSeg 35

methSeg

Segment methylation or differential methylation profile

# Description

The function uses a segmentation algorithm (fastseg) to segment the methylation profiles. Following that, it uses gaussian mixture modelling to cluster the segments into k components. This process uses mean methylation value of each segment in the modeling phase. Each component ideally indicates quantitative classification of segments, such as high or low methylated regions.

### Usage

```
methSeg(obj, diagnostic.plot = TRUE, join.neighbours = FALSE, ...)
```

### **Arguments**

obj

GRanges, methylDiff, methylDiffDB, methylRaw or methylRawDB . If the object is a GRanges it should have one meta column with methylation scores and has to be sorted by position, i.e. ignoring the strand information.

diagnostic.plot

if TRUE a diagnostic plot is plotted. The plot shows methylation and length statistics per segment group. In addition, it shows diagnostics from mixture modeling: the density function estimated and BIC criterion used to decide the optimum number of components in mixture modeling.

join.neighbours

if TRUE neighbouring segments that cluster to the same seg.group are joined by extending the ranges, summing up num.marks and averaging over seg.means.

. . .

arguments to fastseg function in fastseg package, or to densityMclust in Mclust package, could be used to fine tune the segmentation algorithm. E.g. Increasing "alpha" will give more segments. Increasing "cyberWeight" will give also more segments. "maxInt" controls the segment extension around a breakpoint. "minSeg" controls the minimum segment length. "G" argument denotes number of components used in BIC selection in mixture modeling. For more details see fastseg and Mclust documentation.

#### **Details**

To be sure that the algorithm will work on your data, the object should have at least 5000 records

### Value

A GRanges object with segment classification and information. 'seg.mean' column shows the mean methylation per segment. 'seg.group' column shows the segment groups obtained by mixture modeling

### Author(s)

Altuna Akalin, contributions by Arsene Wabo

36 methSeg2bed

#### See Also

```
methSeg2bed, joinSegmentNeighbours
```

## **Examples**

methSeg2bed

Export segments to BED files

# Description

The segments are color coded based on their score (methylation or differential methylation value). They are named by segment group (components in mixture modeling) and the score in the BED file is obtained from 'seg.mean' column of segments object.

### Usage

```
methSeg2bed(segments, filename,
  trackLine = "track name='meth segments' description='meth segments' itemRgb=On",
  colramp = colorRamp(c("gray", "green", "darkgreen")))
```

# Arguments

segments	GRanges object with segment classification and information. This should be the result of methSeg function
filename	name of the output data
trackLine	UCSC browser trackline
colramp	color scale to be used in the BED display defaults to gray, green, darkgreen scale.

### Value

A BED files with the segmented data which can be visualized in the UCSC browser

methylBase-class 37

#### See Also

methSeg

methylBase-class

An S4 class for methylation events sampled in multiple experiments

# **Description**

This class is designed to contain methylation information such as coverage, number of methylated bases, etc.. The methylation events contained in the class must be sampled in multiple experiments (ex: only CpG bases covered in multiple experiments are stored in the object of this class). The class extends data.frame and creates an object that holds methylation information and genomic location. The object belonging to this class is produced by unite function.

#### **Slots**

sample.ids: character vector for ids of samples in the object

assembly: name of the genome assembly

context: context of methylation. Ex: CpG,CpH,CHH, etc

treatment: treatment vector denoting which samples are test and control

coverage.index: vector denoting which columns in the data correspons to coverage values

numCs.index: vector denoting which columns in the data correspons to number of methylatedCs values

numTs.index: vector denoting which columns in the data correspons to number of unmethylated Cs values

destranded: logical value. If TRUE object is destranded, if FALSE it is not.

resolution: resolution of methylation information, allowed values: 'base' or 'region'

### **Details**

methylBase class extends data. frame class therefore providing novice and experienced R users with a data structure that is well known and ubiquitous in many R packages.

### **Subsetting**

In the following code snippets, x is a methylBase. Subsetting by x[i,] will produce a new object if subsetting is done on rows. Column subsetting is not directly allowed to prevent errors in the downstream analysis. see ?methylKit[.

#### Accessors

The following functions provides access to data slots of methylBase: getData,getAssembly, getContext

### Coercion

methylBase object can be coerced to GRanges object via as function.

38 methylBaseDB-class

#### **Examples**

```
data(methylKit)
library(GenomicRanges)
my.gr=as(methylBase.obj,"GRanges")
```

methylBase.obj

Example methylBase object.

### **Description**

methylBase, methylDiff and methylRawList. You can load the data using data(methylKit)

#### **Format**

methylBase.obj object stores the location and methylation information for bases that are covered in all samples. methylBase partially extends data. frame S3 class.

methylBaseDB-class

dbpath: path to flat file database(s)

dbtype: string for type of the flat file database, ex: tabix

An S4 class for storing methylation events sampled in multiple experiments as flat file database

# Description

This class is designed to contain methylation information such as coverage, number of methylated bases, etc... The class creates an object that holds methylation information and genomic location as flat file database. The object belonging to this class is produced by unite function.

### **Slots**

```
num.records: number of records (lines) in the object
sample.ids: character vector for ids of samples in the object
assembly: name of the genome assembly
context: context of methylation. Ex: CpG,CpH,CHH, etc
treatment: treatment vector denoting which samples are test and control
coverage.index: vector denoting which columns in the data correspond to coverage values
numCs.index: vector denoting which columns in the data correspond to number of methylatedCs
values
numTs.index: vector denoting which columns in the data correspond to number of unmethylated
Cs values
destranded: logical value. If TRUE object is destranded, if FALSE it is not.
resolution: resolution of methylation information, allowed values: 'base' or 'region'
```

methylDiff-class 39

#### **Details**

methylBaseDB class has the same functionality as methylBase class, but the data is saved in a flat database file and therefore allocates less space in memory.

### **Subsetting**

In the following code snippets, x is a methylBaseDB. Subsetting by x[i,] will produce a new methylBase object if subsetting is done on rows. Column subsetting is not directly allowed to prevent errors in the downstream analysis. see ?methylKit[.

#### Accessors

The following functions provides access to data slots of methylBaseDB: getData,getAssembly, getContext

# Coercion

methylBaseDB object can be coerced to: GRanges object via as function. methylBase object via as function.

### **Examples**

```
data(methylKit)
methylBaseDB.obj <- unite(methylRawList.obj,save.db=TRUE,dbdir="methylDB")
library(GenomicRanges)
my.gr=as(methylBaseDB.obj,"GRanges")

# remove Database again
rm(methylBaseDB.obj)
unlink("methylDB",recursive=TRUE)</pre>
```

methylDiff-class

An S4 class that holds differential methylation information

# Description

This class is designed to hold statistics and locations for differentially methylated regions/bases. It extends data.frame class. calculateDiffMeth function returns an object with methylDiff class.

### **Slots**

```
sample.ids ids/names of samples in a vector
assembly a name of genome assembly, such as :hg18,mm9, etc
context numeric vector identifying which samples are which group
treatment numeric vector identifying which samples are which group
destranded logical denoting if methylation inormation is destranded or not
resolution string either 'base' or 'region' defining the resolution of methylation information
.Data data.frame holding the locations and statistics
```

40 methylDiffDB-class

#### **Details**

methylDiff class extends data.frame class therefore providing novice and experienced R users with a data structure that is well known and ubiquitous in many R packages.

### **Subsetting**

In the following code snippets, x is a methylDiff. Subsetting by x[i,] will produce a new object if subsetting is done on rows. Column subsetting is not directly allowed to prevent errors in the downstream analysis. see ?methylKit[.

#### Coercion

methylDiff object can be coerced to GRanges object via as function.

#### Accessors

 $The following functions provides access to data slots of methyl Diff: {\tt getData,getAssembly,getContext} \\$ 

### **Examples**

```
data(methylKit)
library(GenomicRanges)
my.gr=as(methylDiff.obj,"GRanges")
```

methylDiff.obj

Example methylKit objects.

### **Description**

methylBase, methylDiff and methylRawList. You can load the data using data(methylKit)

#### **Format**

The Differential methylation information is stored in methylDiff.obj object. methylBase partially extends data.frame S3 class.

methylDiffDB-class

An S4 class that holds differential methylation information as flat file database

# Description

This class is designed to hold statistics and locations for differentially methylated regions/bases as flat file database. calculateDiffMeth function returns an object with methylDiffDB class.

methylDiffDB-class 41

#### **Slots**

```
dbpath: path to flat file database(s)
num.records: number of records (lines) in the object
sample.ids ids/names of samples in a vector
assembly a name of genome assembly, such as :hg18,mm9, etc
context numeric vector identifying which samples are which group
treatment numeric vector identifying which samples are which group
destranded logical denoting if methylation inormation is destranded or not
resolution string either 'base' or 'region' defining the resolution of methylation information
dbtype: string for type of the flat file database, ex: tabix
```

### **Details**

methylDiffDB class has the same functionality as methylDiff class, but the data is saved in a flat database file and therefore allocates less space in memory.

### **Subsetting**

In the following code snippets, x is a methylDiffDB. Subsetting by x[i, ] will produce a new object if subsetting is done on rows. Column subsetting is not directly allowed to prevent errors in the downstream analysis. see ?methylKit[.

#### Coercion

methylDiffDB object can be coerced to: GRanges object via as function. methylDiff object via as function.

### Accessors

The following functions provides access to data slots of methylDiffDB: getData,getAssembly, getContext

```
data(methylKit)
methylDiffDB.obj <- calculateDiffMeth(methylBase.obj,save.db=TRUE,dbdir="methylDB")
library(GenomicRanges)
my.gr=as(methylDiffDB.obj,"GRanges")
# remove Database again
rm(methylDiffDB.obj)
unlink("methylDB",recursive=TRUE)</pre>
```

42 methylRaw-class

methylRaw-class An S4 class for holding raw methylation data from an alignment pipeline.
--

### **Description**

This object stores the raw mehylation data that is read in through read function and extends data. frame. The raw methylation data is basically percent methylation values and read coverage values per genomic base/region.

### **Slots**

```
sample.id: string for an identifier of the sample
assembly: string for genome assembly, ex: hg18,hg19,mm9
context: methylation context string, ex: CpG,CpH,CHH, etc.
resolution: resolution of methylation information, 'base' or 'region'
```

#### **Details**

methylRaw class extends data.frame class therefore providing novice and experienced R users with a data structure that is well known and ubiquitous in many R packages.

### **Subsetting**

In the following code snippets, x is a methylRaw. Subsetting by x[i,] will produce a new object if subsetting is done on rows. Column subsetting is not directly allowed to prevent errors in the downstream analysis. see ?methylKit[.

#### Accessors

 $The following functions provides access to data slots of methylRaw: \verb|getData,getAssembly|, \verb|getContext| to the following functions provides access to data slots of methylRaw: \verb|getData,getAssembly|, \verb|getContext| to the following functions provides access to data slots of methylRaw: \verb|getData,getAssembly|, getContext| to the following functions provides access to data slots of methylRaw: \verb|getData,getAssembly|, getContext| to the following functions provides access to data slots of methylRaw: \verb|getData,getAssembly|, getContext| to the following functions provides access to data slots of methylRaw: \verb|getData,getAssembly|, getContext| to the following functions provides access to data slots of methylRaw: \verb|getData,getAssembly|, getContext| to the following functions access to data slots of methylRaw: \verb|getData,getAssembly|, getContext| to the following functions access to the$ 

### Coercion

methylRaw object can be coerced to GRanges object via as function.

methylRawDB-class 43

```
#coercing methylRaw object to GRanges object
my.gr=as(methylRawList.obj[[1]],"GRanges")
```

methylRawDB-class

An S4 class for storing raw methylation data as flat file database.

### **Description**

This object stores the raw mehylation data that is read in through read function as flat file database. The raw methylation data is basically percent methylation values and read coverage values per genomic base/region.

#### **Slots**

```
dbpath: path to flat file database
num.records: number of records (lines) in the object
sample.id: string for an identifier of the sample
assembly: string for genome assembly, ex: hg18,hg19,mm9
context: methylation context string, ex: CpG,CpH,CHH, etc.
resolution: resolution of methylation information, 'base' or 'region'
dbtype: string for type of the flat file database, ex: tabix
```

# Details

methylRawDB is created via read function and has the same functionality as methylRaw class, but the data is saved in a flat database file and therefore allocates less space in memory.

# **Subsetting**

In the following code snippets, x is a methylRawDB. Subsetting by x[i,] will produce a new methylRaw object if subsetting is done on rows. Column subsetting is not directly allowed to prevent errors in the downstream analysis. see ?methylKit[. x[] will return the methylRawDB object as new methylRaw object

#### Accessors

The following functions provides access to data slots of methylRawDB: getData,getAssembly, getContext

# Coercion

methylRawDB object can be coerced to: GRanges object via as function. methylRaw object via as function.

44 methylRawList-class

#### **Examples**

```
# example of a raw methylation data contained as a text file
read.table(system.file("extdata", "control1.myCpG.txt", package = "methylKit"),
header=TRUE,nrows=5)
methylRawDB.obj <- methRead(</pre>
system.file("extdata", "control1.myCpG.txt", package = "methylKit"),
                        sample.id = "ctrl1", assembly = "hg18",
                        dbtype = "tabix", dbdir = "methylDB")
# example of a methylRawDB object
methylRawDB.obj
str(methylRawDB.obj)
library(GenomicRanges)
#coercing methylRawDB object to GRanges object
my.gr=as(methylRawDB.obj,"GRanges")
#coercing methylRawDB object to methylRaw object
myRaw=as(methylRawDB.obj,"methylRaw")
# remove Database again
rm(methylRawDB.obj)
unlink("methylDB",recursive=TRUE)
```

methylRawList-class

An S4 class for holding a list of methylRaw objects.

### **Description**

This class stores the list of methylRaw objects. Functions such as lapply can be used on this list. It extends list class. This object is primarily produced by read function.

# Slots

```
treatment numeric vector denoting control and test samples .Data a list of methylRaw objects
```

```
data(methylKit)
#applying functions designed for methylRaw on methylRawList object
lapply(methylRawList.obj,"getAssembly")
```

methylRawList.obj 45

methylRawList.obj

Example methylRawList object.

### **Description**

methylBase, methylDiff and methylRawList. You can load the data using data(methylKit)

#### **Format**

Methylation data from multiple the samples regardless of common coverage are stored in methyl-RawList.obj object. methylRawList extends list S3 class

methylRawListDB-class An S4 class for holding a list of methylRawDB objects.

### **Description**

This class stores the list of methylRawDB objects. Functions such as lapply can be used on this list. It extends list class. This object is primarily produced by read function.

### **Slots**

treatment numeric vector denoting control and test samples
.Data a list of methylRawDB objects

46 normalizeCoverage

normalizeCoverage

normalize read coverage between samples

### **Description**

The function normalizes coverage values between samples using a scaling factor derived from differences between mean or median of coverage distributions

### Usage

```
normalizeCoverage(obj,method="median",chunk.size,save.db,...)
## S4 method for signature 'methylRawList'
normalizeCoverage(obj, method = "median",
    chunk.size = 1e+06, save.db = FALSE, ...)
## S4 method for signature 'methylRawListDB'
normalizeCoverage(obj, method = "median",
    chunk.size = 1e+06, save.db = TRUE, ...)
```

### **Arguments**

obj	methylRawList or methylRawListDB object
method	a string "mean" or "median" which denotes median or mean should be used to calculate scaling factor. (Default:median)
chunk.size	Number of rows to be taken as a chunk for processing the methylRawListDB objects. (Default: 1e6)
save.db	A Logical to decide whether the resulting object should be saved as flat file database or not, default: explained in Details sections
	optional Arguments used when save.db is TRUE suffix A character string to append to the name of the output flat file database, only used if save.db is true, default actions: append "_filtered" to current filename if database already exists or generate new file with filename "sampleID_filtered" dbdir The directory where flat file database(s) should be stored, defaults to getwd(), working directory for newly stored databases and to same directory for already existing database

### Value

a methylRawList or methylRawList object depending on class of input object

#### **Details**

The parameter chunk.size is only used when working with methylRawListDB objects, as they are read in chunk by chunk to enable processing large-sized objects which are stored as flat file database. Per default the chunk.size is set to 1M rows, which should work for most systems. If you encounter memory problems or have a high amount of memory available feel free to adjust the chunk.size.

PCASamples 47

The parameter save.db is per default TRUE for methylDB objects as methylRawListDB, while being per default FALSE for methylRawList. If you wish to save the result of an in-memory-calculation as flat file database or if the size of the database allows the calculation in-memory, then you might want to change the value of this parameter.

#### Author(s)

Altuna Akalin

### **Examples**

```
data(methylKit)
# normalize by the median coverage
newObj = normalizeCoverage(methylRawList.obj,method="median")
# normalize by mean coverage and save to database in folder methylDB
newDBObj = normalizeCoverage(methylRawList.obj,method="mean",
save.db=TRUE,dbdir="methylDB")
```

**PCASamples** 

Principal Components Analysis of Methylation data

### **Description**

The function does a PCA analysis using prcomp function using percent methylation matrix as an input.

# Usage

### **Arguments**

.Object a methylBase or methylBaseDB object

screeplot a logical value indicating whether to plot the variances against the number of

the principal component. (default: FALSE)

48 PCASamples

adj.lim	a vector indicating the propotional adjustment of xlim (adj.lim[1]) and ylim (adj.lim[2]). This is primarily used for adjusting the visibility of sample labels on the on the PCA plot. (default: $c(0.0004,0.1)$ )
scale	logical indicating if prcomp should scale the data to have unit variance or not (default: $TRUE$ )
center	logical indicating if prcomp should center the data or not (default: TRUE)
comp	vector of integers with 2 elements specifying which components to be plotted.
transpose	if TRUE (default) percent methylation matrix will be transposed, this is equivalent to doing PCA on variables that are regions/bases. The resulting plot will location of samples in the new coordinate system if FALSE the variables for the matrix will be samples and the resulting plot whill show how each sample (variable) contributes to the principle component.the samples that are highly correlated should have similar contributions to the principal components.
sd.filter	If TRUE, the bases/regions with low variation will be discarded prior to PCA (default:TRUE) $$
sd.threshold	A numeric value. If filterByQuantile is TRUE, the value should be between 0 and 1 and the features whose standard deviations is less than the quantile denoted by sd.threshold will be removed. If filterByQuantile is FALSE, then features whose standard deviations is less than the value of sd.threshold will be removed.(default:0.5)
filterByQuantil	e
	A logical determining if sd.threshold is to be interpreted as a quantile of all standard deviation values from bases/regions (the default), or as an absolute value

### Value

obj.return

chunk.size

The form of the value returned by PCASamples is the summary of principal component analysis by prcomp.

objects, default: 1e6

if the result of prcomp function should be returned or not. (Default:FALSE)

Number of rows to be taken as a chunk for processing the methylRawListDB

### **Details**

The parameter chunk.size is only used when working with methylBaseDB objects, as they are read in chunk by chunk to enable processing large-sized objects which are stored as flat file database. Per default the chunk.size is set to 1M rows, which should work for most systems. If you encounter memory problems or have a high amount of memory available feel free to adjust the chunk.size.

### Note

cor option is not in use anymore, since prcomp is used for PCA analysis instead of princomp

### **Examples**

```
data(methylKit)
```

# do PCA with filtering rows with low variation, filter rows with standard # deviation lower than the 50th percentile of Standard deviation distribution PCASamples(methylBase.obj,screeplot=FALSE, adj.lim=c(0.0004,0.1),

percMethylation 49

```
scale=TRUE,center=TRUE,comp=c(1,2),transpose=TRUE,sd.filter=TRUE,
sd.threshold=0.5,filterByQuantile=TRUE,obj.return=FALSE)
```

 $\begin{tabular}{ll} \it perc Methylation &\it get percent methylation scores from methylBase or methylBase DB object \\ \end{tabular}$ 

### **Description**

get percent methylation scores from methylBase or methylBaseDB object

#### Usage

```
percMethylation(methylBase.obj, rowids = FALSE, save.txt = FALSE,
    chunk.size = 1e+06)

## S4 method for signature 'methylBase'
percMethylation(methylBase.obj, rowids = FALSE)

## S4 method for signature 'methylBaseDB'
percMethylation(methylBase.obj, rowids = FALSE,
    save.txt = FALSE, chunk.size = 1e+06)
```

### Arguments

methylBase.obj a methylBase or methylBaseDB object

rowids if TRUE, matrix rownames have identifiers as base/region location (default:FALSE)

save.txt if TRUE, the matrix will be written to a text file, but only for methylBaseDB objects (default: FALSE)

chunk.size Number of rows to be taken as a chunk for processing the methylBaseDB objects (default: 1e6)

### Value

matrix with percent methylation values per base/region across all samples, row names would be base/region identifiers

#### **Details**

The parameter chunk.size is only used when working with methylBaseDB objects, as they are read in chunk by chunk to enable processing large-sized objects which are stored as flat file database. Per default the chunk.size is set to 1M rows, which should work for most systems. If you encounter memory problems or have a high amount of memory available feel free to adjust the chunk.size.

```
data(methylKit)
mat=percMethylation(methylBase.obj)
head(mat)
```

50 pool

pool

Pool replicates within groups to a single sample per group

# Description

The function sums up coverage, numCs and numTs values within each group so one representative sample for each group will be created in a new methylBase object

# Usage

```
pool(obj, sample.ids, chunk.size = 1e+06, save.db = FALSE, ...)
## S4 method for signature 'methylBase'
pool(obj, sample.ids, chunk.size = 1e+06,
    save.db = FALSE, ...)
## S4 method for signature 'methylBaseDB'
pool(obj, sample.ids, chunk.size = 1e+06,
    save.db = TRUE, ...)
```

# Arguments

obj	methylBase or methylBaseDB object with two groups or more and each group should have multiple samples
sample.ids	a character vector of new sample.ids ex:c("test","control"), should follow the same order as unique treatment vector, and should be equal to the length of the unique treatment vector
chunk.size	Number of rows to be taken as a chunk for processing the methylRawListDB objects, default: 1e6
save.db	A Logical to decide whether the resulting object should be saved as flat file database or not, default: explained in Details sections
	optional Arguments used when save.db is TRUE
	suffix A character string to append to the name of the output flat file database, only used if save.db is true, default actions: The default suffix is a 13-character random string appended to the fixed prefix "methylBase", e.g. "methylBase_16d3047c1a254.txt.bgz"
	dbdir The directory where flat file database(s) should be stored, defaults to getwd(), working directory for newly stored databases and to same directory for already existing database
	dbtype The type of the flat file database, currently only option is "tabix" (only

### Value

a methylBase or methylBaseDB object depending on class of input object

used for newly stored databases)

processBismarkAln 51

#### **Details**

The parameter chunk.size is only used when working with methylBaseDB objects, as they are read in chunk by chunk to enable processing large-sized objects which are stored as flat file database. Per default the chunk.size is set to 1M rows, which should work for most systems. If you encounter memory problems or have a high amount of memory available feel free to adjust the chunk.size.

The parameter save. db is per default TRUE for methylDB objects as methylBaseDB, while being per default FALSE for methylBase. If you wish to save the result of an in-memory-calculation as flat file database or if the size of the database allows the calculation in-memory, then you might want to change the value of this parameter.

### Author(s)

Altuna Akalin

### **Examples**

```
data(methylKit)

# methylBase.obj has two groups, each group has two samples,
# the following function will pool the samples in each group
# so that each group will be represented by one pooled sample
pooled.methylBase=pool(methylBase.obj,sample.ids=c("test","control"))
```

processBismarkAln

Get methylation percentage from sorted Bismark alignments

#### **Description**

The function calls methylation percentage per base from sorted Bismark SAM or BAM files and reads methylation information as methylKit objects. Bismark is a popular aligner for high-throughput bisulfite sequencing experiments and it outputs its results in SAM format by default, and can be converted to BAM. Bismark SAM/BAM format contains aligner specific tags which are absolutely necessary for methylation percentage calling using processBismarkAln. SAM/BAM files from other aligners will not work with this function.

#### Usage

```
processBismarkAln(location, sample.id, assembly, save.folder = NULL,
    save.context = c("CpG"), read.context = "CpG", nolap = FALSE,
    mincov = 10, minqual = 20, phred64 = FALSE, treatment,
    save.db = FALSE)

## S4 method for signature 'character, character, character'
processBismarkAln(location, sample.id,
    assembly, save.folder, save.context, read.context, nolap, mincov, minqual,
    phred64, save.db)

## S4 method for signature 'list,list,character'
processBismarkAln(location, sample.id, assembly,
```

52 processBismarkAln

```
save.folder = NULL, save.context = c("CpG"), read.context = "CpG",
nolap = FALSE, mincov = 10, minqual = 20, phred64 = FALSE, treatment,
save.db = FALSE)
```

# **Arguments**

location	location of sam or bam file(s). If multiple files are given this argument must be a list.
sample.id	the id(s) of samples in the same order as file. If multiple sam files are given this arugment must be a list.
assembly	string that determines the genome assembly. Ex: mm9,hg18 etc. This is just a string for book keeping. It can be any string. Although, when using multiple files from the same assembly, this string should be consistent in each object.
save.folder	The folder which will be used to save methylation call files, if set to NULL no methylation call file will be saved as a text file. The files saved can be read into R in less time using methRead function in methylKit
save.context	A character vector consisting following strings: "CpG", "CHG", "CHH". The methylation percentages for these methylation contexts will be saved to save.folder
read.context	One of the 'CpG','CHG','CHH' or 'none' strings. Determines what type of methylation context will be read-in to the memory which can be immediately used for analysis. If given as 'none', processBismarkAln will not return any object, but if a save.folder argument given it will save the methylation percentage call files.
nolap	if set to TRUE and the SAM/BAM file has paired-end reads, the one read of the overlapping paired-end read pair will be ignored for methylation calling.
mincov	minimum read coverage to call a methylation status for a base.
minqual	minimum phred quality score to call a methylation status for a base.
phred64	logical (default: FALSE) you will not need to set this TRUE, Currently bismark gives only phred33 scale
treatment	treatment vector only to be used when location and sample.id parameters are lists and you are trying to read-in multiple samples that are related to eachother in down-stream analysis.
save.db	A Logical to decide whether the resulting object should be saved as flat file database or not ( default: FALSE). With the default value, a text file containing methylation values will be saved. If TRUE, database will either be saved to location save.folder or if this is NULL, to a new folder in the current working directory named after this scheme: "methylDB <date> &lt;3randomlettersornumbers&gt;"</date>

### Value

methylRaw or methylRawList object

# Note

SAM files should be sorted with samtools sort or unix sort. Other sorting methods can alter the order of fields(columns) in the SAM file and that will result in an error when using processBismarkAln().

reconstruct 53

#### **Examples**

```
# reading one bismark file:
my.file=system.file("extdata", "test.fastq_bismark.sorted.min.sam",
                              package = "methylKit")
obj=processBismarkAln(my.file,"test",assembly="hg18",save.folder=NULL,
                 save.context="CpG",read.context="CpG")
# reading multiple files
file.list2=list(system.file("extdata", "test.fastq_bismark.sorted.min.sam",
package = "methylKit"),
               system.file("extdata", "test.fastq_bismark.sorted.min.sam",
               package = "methylKit"),
              system.file("extdata", "test.fastq_bismark.sorted.min.sam",
              package = "methylKit"),
               system.file("extdata", "test.fastq_bismark.sorted.min.sam",
                package = "methylKit"))
 objs=processBismarkAln(location=file.list2
             , sample.id=list("test1", "test2", "ctrl1", "ctrl1"), assembly="hg18",
             save.folder=NULL, save.context=NULL, read.context="CpG",
             nolap=FALSE,mincov=10,minqual=20,phred64=FALSE,
             treatment=c(1,1,0,0))
```

reconstruct

Reconstruct methylBase or methylBaseDB object based on a new methylation percentage matrix

### **Description**

The function reconstructs a new methylBase object from an input methylBase object and percent methylation matrix. Basically, it uses the read coverages in the input methylBase object and deduces new number of methylated Cs and unmethylated Cs based on the input percent methylation matrix. It is ideally to be used to reconstruct methylBase objects after batch correction on percent methylation values. The percent methylation matrix rows must match methylBase object rows in order ,and in addition column order (the order of samples) in input methylBase must match the order in percent methylation matrix.

# Usage

```
reconstruct(methMat, mBase, chunk.size = 1e+06, save.db = FALSE, ...)
## S4 method for signature 'ANY,methylBase'
reconstruct(methMat, mBase, chunk.size = 1e+06,
    save.db = FALSE, ...)
## S4 method for signature 'ANY,methylBaseDB'
reconstruct(methMat, mBase, chunk.size = 1e+06,
    save.db = TRUE, ...)
```

54 reconstruct

#### **Arguments**

methMat percent methylation matrix, row order and order of the samples same as the

methylBase object

mBase methylBase or methylBaseDB object to be reconstructed

chunk.size Number of rows to be taken as a chunk for processing the methylBaseDB objects

(default: 1e6)

save.db A Logical to decide whether the resulting object should be saved as flat file

database or not, default: explained in Details sections

.. optional Arguments used when save.db is TRUE

suffix A character string to append to the name of the output flat file database, only used if save.db is true, default actions: append "\_reconstructed" to current filename if database already exists or generate new file with filename "methyl-

Base\_reconstructed"

dbdir The directory where flat file database(s) should be stored, defaults to getwd(), working directory for newly stored databases and to same directory for

already existing database

dbtype The type of the flat file database, currently only option is "tabix" (only

used for newly stored databases)

### Value

new methylBase or methylBase object where methylation percentage matches input methMat and coverages matches input mBase

### **Details**

The parameter chunk.size is only used when working with methylBaseDB objects, as they are read in chunk by chunk to enable processing large-sized objects which are stored as flat file database. Per default the chunk.size is set to 1M rows, which should work for most systems. If you encounter memory problems or have a high amount of memory available feel free to adjust the chunk.size.

The parameter save.db is per default TRUE for methylDB objects as methylBaseDB, while being per default FALSE for methylBase. If you wish to save the result of an in-memory-calculation as flat file database or if the size of the database allows the calculation in-memory, then you might want to change the value of this parameter.

### Note

Batch effect correction (if any batch effect exists) is a tricky issue. We provide some simple ways to deal with it (see assocComp and removeComp ), But if you can find other ways to correct for batch effects and want to create a methylBase object with the corrected percent methylation values, you can use this function.

### Author(s)

Altuna Akalin

# Examples

```
data(methylKit)
```

# get percent methylation

regionCounts 55

```
mat=percMethylation(methylBase.obj)

# do some changes in the matrix
# this is just a toy example
# ideally you want to correct the matrix
# for batch effects
mat[mat==100]=80

# reconstruct the methylBase from the corrected matrix
newobj=reconstruct(mat,methylBase.obj)
```

regionCounts

Get regional counts for given GRanges or GRangesList object

### **Description**

Convert methylRaw, methylRawDB, methylRawList, methylRawListDB, methylBase or methylBaseDB object into regional counts for a given GRanges or GRangesList object.

# Usage

```
regionCounts(object,regions,cov.bases=0,strand.aware=FALSE,chunk.size,save.db,...)
## S4 method for signature 'methylRaw, GRanges'
regionCounts(object, regions, cov.bases = 0,
 strand.aware = FALSE, chunk.size = 1e+06, save.db = FALSE, ...)
## S4 method for signature 'methylBase, GRanges'
regionCounts(object, regions, cov.bases = 0,
 strand.aware = FALSE, chunk.size = 1e+06, save.db = FALSE, ...)
## S4 method for signature 'methylRaw, GRangesList'
regionCounts(object, regions, cov.bases = 0,
  strand.aware = FALSE, chunk.size = 1e+06, save.db = FALSE, ...)
## S4 method for signature 'methylBase,GRangesList'
regionCounts(object, regions,
 cov.bases = 0, strand.aware = FALSE, chunk.size = 1e+06,
  save.db = FALSE, ...)
## S4 method for signature 'methylRawList,GRanges'
regionCounts(object, regions, cov.bases = 0,
 strand.aware = FALSE, chunk.size = 1e+06, save.db = FALSE, ...)
## S4 method for signature 'methylRawList,GRangesList'
regionCounts(object, regions,
  cov.bases = 0, strand.aware = FALSE, chunk.size = 1e+06,
  save.db = FALSE, ...)
## S4 method for signature 'methylRawDB,GRanges'
regionCounts(object, regions, cov.bases = 0,
```

56 regionCounts

```
strand.aware = FALSE, chunk.size = 1e+06, save.db = TRUE, ...)
## S4 method for signature 'methylRawDB,GRangesList'
regionCounts(object, regions,
  cov.bases = 0, strand.aware = FALSE, chunk.size = 1e+06,
  save.db = TRUE, ...)
## S4 method for signature 'methylRawListDB,GRanges'
regionCounts(object, regions,
 cov.bases = 0, strand.aware = FALSE, chunk.size = 1e+06,
 save.db = TRUE, ...)
## S4 method for signature 'methylRawListDB,GRangesList'
regionCounts(object, regions,
 cov.bases = 0, strand.aware = FALSE, chunk.size = 1e+06,
 save.db = TRUE, ...)
## S4 method for signature 'methylBaseDB,GRanges'
regionCounts(object, regions, cov.bases = 0,
  strand.aware = FALSE, chunk.size = 1e+06, save.db = TRUE, ...)
## S4 method for signature 'methylBaseDB,GRangesList'
regionCounts(object, regions,
  cov.bases = 0, strand.aware = FALSE, chunk.size = 1e+06,
  save.db = TRUE, ...)
```

### Arguments

object a methylRaw, methylRawDB, methylRawList, methylRawListDB, methylBase

or methylBaseDB object

regions a GRanges or GRangesList object. Make sure that the GRanges objects are

unique in chr, start, end and strand columns. You can make them unique by using

unique() function.

cov.bases number minimum bases covered per region (Default:0). Only regions with base

coverage above this threshold are returned.

strand.aware if set to TRUE only CpGs that match the strand of the region will be summa-

rized. (default:FALSE)

chunk.size Number of rows to be taken as a chunk for processing the methylDB objects

(default: 1e6)

save.db A Logical to decide whether the resulting object should be saved as flat file

database or not, default: explained in Details sections

... optional Arguments used when save.db is TRUE

suffix A character string to append to the name of the output flat file database, only used if save.db is true, default actions: append "\_regions" to current filename if database already exists or generate new file with filename "sampleID\_regions"

or "methylBase\_filtered" dependent on input object

dbdir The directory where flat file database(s) should be stored, defaults to getwd(), working directory for newly stored databases and to same directory for

already existing database

dbtype The type of the flat file database, currently only option is "tabix" (only

used for newly stored databases)

removeComp 57

#### Value

a new methylRaw,methylBase or methylRawList object. If strand. aware is set to FALSE (default). Even though the resulting object will have the strand information of regions it will still contain methylation information from both strands.

#### **Details**

The parameter chunk.size is only used when working with methylRawDB, methylBaseDB or methylRawListDB objects, as they are read in chunk by chunk to enable processing large-sized objects which are stored as flat file database. Per default the chunk.size is set to 1M rows, which should work for most systems. If you encounter memory problems or have a high amount of memory available feel free to adjust the chunk.size.

The parameter save.db is per default TRUE for methylDB objects as methylRawDB, methylBaseDB or methylRawListDB, while being per default FALSE for methylRaw, methylBase or methylRawList. If you wish to save the result of an in-memory-calculation as flat file database or if the size of the database allows the calculation in-memory, then you might want to change the value of this parameter.

### **Examples**

```
data(methylKit)

# get the windows of interest as a GRanges object, this can be any set
# of genomic locations
library(GenomicRanges)
my.win=GRanges(seqnames="chr21",
ranges=IRanges(start=seq(from=9764513,by=10000,length.out=20),width=5000) )

# getting counts per region
regional.methylRaw=regionCounts(object=methylRawList.obj, regions=my.win,
cov.bases=0,strand.aware=FALSE)
```

removeComp

Remove principal components from a methylBase object

# **Description**

This function can remove a given principal componet from a given methylBase object. First, it calculates principal components from percent methylation matrix and removes the given component(s), reconstructs the methylation matrix then reconstructs number of methylated and unmethylated Cs per position based on the reconstructed percent methylation matrix, and finally returns a new methylBase object.

# Usage

```
removeComp(mBase, comp, chunk.size = 1e+06, save.db = FALSE, ...)
## S4 method for signature 'methylBase'
removeComp(mBase, comp, chunk.size = 1e+06,
    save.db = FALSE, ...)
```

58 removeComp

```
## S4 method for signature 'methylBaseDB'
removeComp(mBase, comp, chunk.size = 1e+06,
    save.db = TRUE, ...)
```

### **Arguments**

mBase methylBase or methylBaseDB object with no NA values, that means all bases

should be covered in all samples.

comp vector of component numbers to be removed

chunk.size Number of rows to be taken as a chunk for processing the methylBaseDB objects

(default: 1e6)

save.db A Logical to decide whether the resulting object should be saved as flat file

database or not, default: explained in Details sections

... optional Arguments used when save.db is TRUE

suffix A character string to append to the name of the output flat file database, only used if save.db is true, default actions: append "\_reconstructed" to current filename if database already exists or generate new file with filename "methyl-

Base\_reconstructed"

dbdir The directory where flat file database(s) should be stored, defaults to getwd(), working directory for newly stored databases and to same directory for

already existing database

dbtype The type of the flat file database, currently only option "tabix"

#### Value

new methylBase or methylBaseDB object

### **Details**

The parameter chunk.size is only used when working with methylBaseDB objects, as they are read in chunk by chunk to enable processing large-sized objects which are stored as flat file database. Per default the chunk.size is set to 1M rows, which should work for most systems. If you encounter memory problems or have a high amount of memory available feel free to adjust the chunk.size.

The parameter save.db is per default TRUE for methylDB objects as methylBaseDB, while being per default FALSE for methylBase. If you wish to save the result of an in-memory-calculation as flat file database or if the size of the database allows the calculation in-memory, then you might want to change the value of this parameter.

```
data(methylKit)
# remove 1st principal component
newObj=removeComp(methylBase.obj,comp=1)
# remove 3rd and 4th principal components
newObj=removeComp(methylBase.obj,comp=c(3,4))
```

reorganize 59

reorganize	Reorganize methylKit objects by creating new objects from subset of samples

### **Description**

The function creates a new methylRawList, methylRawListDB, methylBase or methylBaseDB object by selecting a subset of samples from the input object, which is a methylRawList or methylBase object. You can use the function to partition a large methylRawList or methylBase object to smaller object based on sample ids or when you want to reorder samples and/or give a new treatmet vector.

# Usage

```
reorganize(methylObj, sample.ids, treatment, chunk.size = 1e+06,
    save.db = FALSE, ...)

## S4 method for signature 'methylBase'
reorganize(methylObj, sample.ids, treatment,
    chunk.size = 1e+06, save.db = FALSE, ...)

## S4 method for signature 'methylRawList'
reorganize(methylObj, sample.ids, treatment,
    chunk.size = 1e+06, save.db = FALSE, ...)

## S4 method for signature 'methylRawListDB'
reorganize(methylObj, sample.ids, treatment,
    chunk.size = 1e+06, save.db = TRUE, ...)

## S4 method for signature 'methylBaseDB'
reorganize(methylObj, sample.ids, treatment,
    chunk.size = 1e+06, save.db = TRUE, ...)
```

# **Arguments**

methylObj	a methylRawList, methylRawListDB, methylBase or methylBaseDB object
sample.ids	a vector for sample.ids to be subset. Order is important and the order should be similar to treatment. sample.ids should be a subset or reordered version of sample ids in the input object.
treatment	treatment vector, should be same length as sample.ids vector
chunk.size	Number of rows to be taken as a chunk for processing the methylBaseDB or methylRawListDB objects, default: $1e6$
save.db	A Logical to decide whether the resulting object should be saved as flat file database or not, default: explained in Details sections
•••	optional Arguments used when save.db is TRUE suffix A character string to append to the name of the output flat file database, only used if save.db is true, default actions: For "methylBase": The default suffix is a 13-character random string appended to the fixed prefix "methylBase", e.g. "methylBase_16d3047c1a254.txt.bgz". For "methylRawList": ignored.

60 select

dbdir The directory where flat file database(s) should be stored, defaults to getwd(), working directory for newly stored databases and to same directory for already existing database

dbtype The type of the flat file database, currently only option "tabix"

#### Value

returns a methylRawList, methylRawListDB, methylBase or methylBaseDB object depending on the input object

#### **Details**

The parameter chunk.size is only used when working with methylBaseDB or methylRawListDB objects, as they are read in chunk by chunk to enable processing large-sized objects which are stored as flat file database. Per default the chunk.size is set to 1M rows, which should work for most systems. If you encounter memory problems or have a high amount of memory available feel free to adjust the chunk.size.

The parameter save. db is per default TRUE for methylDB objects as methylBaseDB and methylRawListDB, while being per default FALSE for methylBase and methylRawList. If you wish to save the result of an in-memory-calculation as flat file database or if the size of the database allows the calculation in-memory, then you might want to change the value of this parameter.

### **Examples**

select

selects rows from of methylKit objects

### Description

The function returns a subset of data contained in the methylKit objects.

select 61

#### Usage

```
select(x,i)
## S4 method for signature 'methylBase'
select(x, i)
## S4 method for signature 'methylRaw'
select(x, i)
## S4 method for signature 'methylDiff'
select(x, i)
## S4 method for signature 'methylRawDB'
select(x, i)
## S4 method for signature 'methylBaseDB'
select(x, i)
## S4 method for signature 'methylBaseDB'
select(x, i)
```

#### **Arguments**

x an methylBase,methylBaseDB, methylRaw,methylRawDB or methylDiff object

i a numeric or logical vector. This vector corresponds to bases or regions contained in methylKit objects. The vector is used to subset the data.

# Value

a methylBase,methylRaw or methylDiff object depending on the input object.

62 selectByOverlap

```
# remove Database again
rm(methylBaseDB.obj)
rm(methylRawDB.obj)
unlink("methylDB",recursive=TRUE)
```

selectByOverlap

selects records of methylDB objects lying inside a GRanges range

### **Description**

The function selects records from a methylBaseDB, methylRawDB or methylDiffDB object that lie inside the regions given by ranges of class GRanges and returns an object of class methylBase, methylRaw or methylDiff

### Usage

```
selectByOverlap(object,ranges)
## S4 method for signature 'methylRaw'
selectByOverlap(object, ranges)
## S4 method for signature 'methylRawList'
selectByOverlap(object, ranges)
## S4 method for signature 'methylBase'
selectByOverlap(object, ranges)
## S4 method for signature 'methylDiff'
selectByOverlap(object, ranges)
## S4 method for signature 'methylRawDB'
selectByOverlap(object, ranges)
## S4 method for signature 'methylRawListDB'
selectByOverlap(object, ranges)
## S4 method for signature 'methylBaseDB'
selectByOverlap(object, ranges)
## S4 method for signature 'methylDiffDB'
selectByOverlap(object, ranges)
```

### **Arguments**

object an methylBaseDB,methylRawDB or methylDiffDB object ranges a GRanges object specifying the regions of interest

#### Value

a methylBase,methylRaw or methylDiff object depending on the input object.

#### Author(s)

Alexander Gosdschan

### **Examples**

```
data(methylKit)
system.file("extdata", "control2.myCpG.txt", package = "methylKit") )
methylRawListDB.obj=methRead(file.list,
                        sample.id=list("test1","test2","ctrl1","ctrl2"),\\
                        assembly="hg18", treatment=c(1,1,0,0),
                        dbtype = "tabix",dbdir = "methylDB")
methylBaseDB.obj=unite(methylRawListDB.obj)
methylDiffDB.obj = calculateDiffMeth(methylBaseDB.obj)
# define the windows of interest as a GRanges object, this can be any set
# of genomic locations
library(GenomicRanges)
my.win=GRanges(seqnames="chr21",
ranges=IRanges(start=seq(from=9764513,by=10000,length.out=20),width=5000) )
# selects the records that lie inside the regions
myRaw <- selectByOverlap(methylRawListDB.obj[[1]],my.win)</pre>
# selects the records that lie inside the regions
myBase <- selectByOverlap(methylBaseDB.obj,my.win)</pre>
# selects the records that lie inside the regions
myDiff <- selectByOverlap(methylDiffDB.obj,my.win)</pre>
rm(methylRawListDB.obj)
rm(methylBaseDB.obj)
rm(methylDiffDB.obj)
unlink("methylDB",recursive=TRUE)
```

show, methylBase-method

show method for methylKit classes

### **Description**

 $The show method works for \verb|methylRaw|, methylRaw|DB|, \verb|methylRaw|List|DB|, methylBase|DB| and \verb|methylDiff| objects$ 

64 tileMethylCounts

### Usage

```
## S4 method for signature 'methylBase'
show(object)
## S4 method for signature 'methylRaw'
show(object)
## S4 method for signature 'methylRawList'
show(object)
## S4 method for signature 'methylDiff'
show(object)
## S4 method for signature 'methylRawDB'
show(object)
## S4 method for signature 'methylRawListDB'
show(object)
## S4 method for signature 'methylBaseDB'
show(object)
## S4 method for signature 'methylDiffDB'
show(object)
```

### **Arguments**

object any methylKit object

### **Examples**

```
data(methylKit)
methylDiff.obj
show(methylDiff.obj)
```

tileMethylCounts

Get methylated/unmethylated base counts for tilling windows

# Description

The function summarizes methylated/unmethylated base counts over tilling windows across genome. This function can be used when differential methylated analysis is preferable to tilling windows instead of base pairs.

tileMethylCounts 65

#### Usage

```
tileMethylCounts(object,win.size=1000,step.size=1000,cov.bases=0,mc.cores=1,save.db,...)
## S4 method for signature 'methylRaw'
tileMethylCounts(object, win.size = 1000,
  step.size = 1000, cov.bases = 0, mc.cores = 1, save.db = FALSE, ...)
## S4 method for signature 'methylRawList'
tileMethylCounts(object, win.size = 1000,
  step.size = 1000, cov.bases = 0, mc.cores = 1, save.db = FALSE, ...)
## S4 method for signature 'methylBase'
tileMethylCounts(object, win.size = 1000,
  step.size = 1000, cov.bases = 0, mc.cores = 1, save.db = FALSE, ...)
## S4 method for signature 'methylRawDB'
tileMethylCounts(object, win.size = 1000,
 step.size = 1000, cov.bases = 0, mc.cores = 1, save.db = TRUE, ...)
## S4 method for signature 'methylRawListDB'
tileMethylCounts(object, win.size = 1000,
  step.size = 1000, cov.bases = 0, mc.cores = 1, save.db = TRUE, ...)
## S4 method for signature 'methylBaseDB'
tileMethylCounts(object, win.size = 1000,
  step.size = 1000, cov.bases = 0, mc.cores = 1, save.db = TRUE, ...)
```

### **Arguments**

object	methylRaw, methylRawDB, methylRawList, methylRawListDB, methylBase or methylBaseDB object containing base pair resolution methylation information
win.size	an integer for the size of the tiling windows
step.size	an integer for the step size of tiling windows
cov.bases	minimum number of bases to be covered in a given window
mc.cores	number of cores to use when processing methylDB objects, default: 1, but always 1 for Windows)
save.db	A Logical to decide whether the resulting object should be saved as flat file database or not, default: explained in Details sections
	optional Arguments used when save.db is TRUE suffix A character string to append to the name of the output flat file database, only used if save.db is true, default actions: append "_tiled" to current filename if database already exists or generate new file with filename "sampleID_tiled" or "methylBase_tiled" dependent on input object dbdir The directory where flat file database(s) should be stored, defaults to getwd(), working directory for newly stored databases and to same directory for already existing database

### Value

methylRaw,methylBase or methylRawList object

66 unite

#### **Details**

The parameter chunk.size is only used when working with methylRawDB, methylBaseDB or methylRawListDB objects, as they are read in chunk by chunk to enable processing large-sized objects which are stored as flat file database. Per default the chunk.size is set to 1M rows, which should work for most systems. If you encounter memory problems or have a high amount of memory available feel free to adjust the chunk.size.

The parameter save.db is per default TRUE for methylDB objects as methylRawDB, methylBaseDB or methylRawListDB, while being per default FALSE for methylRaw, methylBase or methylRawList. If you wish to save the result of an in-memory-calculation as flat file database or if the size of the database allows the calculation in-memory, then you might want to change the value of this parameter

# **Examples**

unite

unite methylRawList to a single table

### **Description**

This functions unites methylRawList and methylRawListDB objects that only bases with coverage from all samples are retained. The resulting object is either of class methylBase or methylBaseDB depending on input.

#### Usage

```
unite(object, destrand = FALSE, min.per.group = NULL, chunk.size = 1e+06,
    mc.cores = 1, save.db = FALSE, ...)

## S4 method for signature 'methylRawList'
unite(object, destrand = FALSE,
    min.per.group = NULL, chunk.size = 1e+06, mc.cores = 1,
    save.db = FALSE, ...)

## S4 method for signature 'methylRawListDB'
unite(object, destrand = FALSE,
    min.per.group = NULL, chunk.size = 1e+06, mc.cores = 1,
    save.db = TRUE, ...)
```

### **Arguments**

object

a methylRawList or methylRawListDB object to be merged by common locations covered by reads

unite 67

destrand if TRUE, reads covering both strands of a CpG dinucleotide will be merged, do

not set to TRUE if not only interested in CpGs (default: FALSE). If the methyl-RawList object contains regions rather than bases setting destrand to TRUE will

have no effect.

min.per.group an integer denoting minimum number of samples per replicate needed to cover

a region/base. By default only regions/bases that are covered in all samples are united as methylBase object, however by supplying an integer for this argument users can control how many samples needed to cover region/base to be united as methylBase object. For example, if min.per.group set to 2 and there are 3 replicates per condition, the bases/regions that are covered in at least 2 replicates

will be united and missing data for uncovered bases/regions will appear as NAs.

chunk.size Number of rows to be taken as a chunk for processing the methylRawListDB

objects, default: 1e6

mc.cores number of cores to use when processing methylRawListDB objects, default: 1,

but always 1 for Windows)

A Logical to decide whether the resulting object should be saved as flat file save.db

database or not, default: explained in Details sections

optional Arguments used when save.db is TRUE

suffix A character string to append to the name of the output flat file database, only used if save.db is true, default actions: The default suffix is a 13-character

random string appended to the fixed prefix "methylBase", e.g. "methylBase\_16d3047c1a254.txt.bgz"

dbdir The directory where flat file database(s) should be stored, defaults to getwd(), working directory for newly stored databases and to same directory for

already existing database

dbtype The type of the flat file database, currently only option is "tabix" (only

used for newly stored databases)

#### Value

a methylBase or methylBaseDB object depending on input

### **Details**

The parameter chunk.size is only used when working with methylRawDB or methylRawListDB objects, as they are read in chunk by chunk to enable processing large-sized objects which are stored as flat file database. Per default the chunk.size is set to 1M rows, which should work for most systems. If you encounter memory problems or have a high amount of memory available feel free to adjust the chunk. size.

The parameter save.db is per default TRUE for methylDB objects as methylRawListDB, while being per default FALSE for methylRawList. If you wish to save the result of an in-memorycalculation as flat file database or if the size of the database allows the calculation in-memory, then you might change the value of this parameter.

```
data(methylKit)
## Following
my.methylBase=unite(methylRawList.obj)
my.methylBase=unite(methylRawList.obj,destrand=TRUE)
```

68 updateMethObject

updateMethObject $update\ methylKit\ objects\ The\ method\ updates\ object\ from\ ear\ sions\ (<\!v0.9.1)\ to\ latest\ object.$	lier ver-
---	-----------

# Description

update methylKit objects

The method updates object from earlier versions (<v0.9.1) to latest object.

# Usage

```
updateMethObject(object)
```

# Arguments

object a methylKit object: methylRaw, methylRawList, methylBase or methylDiff

### Value

```
methylRaw,methylDiff,methylBase or methylRawList object
@export @docType methods @rdname updateMethObject
```

# **Index**

```
[,methylBase,ANY,ANY,ANY-method
                                                  calculateDiffMeth,methylBase-method
        (extract), 15
                                                          (calculateDiffMeth), 7
[,methylBaseDB,ANY,ANY,ANY-method
                                                  calculateDiffMeth,methylBaseDB-method
                                                          (calculateDiffMeth), 7
        (extract), 15
[,methylDiff,ANY,ANY,ANY-method
                                                  calculateDiffMethDSS, 10
        (extract), 15
                                                  \verb|calculateDiffMethDSS|, methylBase-method|\\
                                                          (calculateDiffMethDSS), 10
[,methylDiffDB,ANY,ANY,ANY-method
        (extract), 15
                                                  clusterSamples, 11
                                                  clusterSamples, methylBase-method
[,methylRaw,ANY,ANY,ANY-method
                                                          (clusterSamples), 11
        (extract), 15
[,methylRawDB,ANY,ANY,ANY-method
                                                  clusterSamples, methylBaseDB-method
                                                          (clusterSamples), 11
        (extract), 15
adjust.methylC(adjustMethylC), 3
                                                  data.frame, 23, 37, 39, 40, 42
adjustMethylC, 3
                                                  dataSim, 9, 12
adjustMethylC, methylRaw, methylRaw-method
                                                  densityMclust, 35
        (adjustMethylC), 3
                                                  diffMethPerChr, 14
adjust \texttt{MethylC}, \texttt{methylRawDB}, \texttt{methylRawDB}-\texttt{method} \ diff \texttt{MethPerChr}, \texttt{methylDiff-method}
        (adjustMethylC), 3
                                                          (diffMethPerChr), 14
adjustMethylC,methylRawList,methylRawList-metdooffMethPerChr,methylDiffDB-method
        (adjustMethylC), 3
                                                          (diffMethPerChr), 14
\verb|adjustMethylC|, methylRawListDB|, methylRawListDB-method|
        (adjustMethylC), 3
                                                  extract, 15
annotate.WithFeature
                                                  fastseg, 35
        (genomation-deprecated), 18
                                                  filterByCoverage, 16
annotate.WithGenicParts
                                                  filterByCoverage, methylRaw-method
        (genomation-deprecated), 18
as, 37, 39–43
                                                          (filterByCoverage), 16
                                                  filterByCoverage,methylRawDB-method
assocComp, 5, 54
                                                          (filterByCoverage), 16
barplot, 14
                                                  filter By Coverage, methylRaw List-method
                                                          (filterByCoverage), 16
bedgraph, 5
bedgraph, methylDiff-method (bedgraph), 5
                                                  filter By Coverage, methyl Raw List DB-method\\
bedgraph, methylDiffDB-method
                                                          (filterByCoverage), 16
        (bedgraph), 5
                                                  genomation-deprecated, 18
bedgraph, methylRaw-method (bedgraph), 5
bedgraph, methylRawDB-method (bedgraph),
                                                  get.methylDiff (getMethylDiff), 26
                                                  getAssembly, 18, 37, 39-43
                                                  getAssembly,methylBase-method
bedgraph, methylRawList-method
                                                          (getAssembly), 18
        (bedgraph), 5
                                                  getAssembly,methylBaseDB-method
bedgraph, methylRawListDB-method
                                                          (getAssembly), 18
        (bedgraph), 5
                                                  getAssembly,methylDiff-method
calculateDiffMeth, 7, 39, 40
                                                          (getAssembly), 18
```

<pre>getAssembly,methylDiffDB-method</pre>	<pre>getMethylationStats,methylRaw-method</pre>
(getAssembly), 18	(getMethylationStats), 25
<pre>getAssembly,methylRaw-method</pre>	<pre>getMethylationStats,methylRawDB-method</pre>
(getAssembly), 18	(getMethylationStats), 25
<pre>getAssembly,methylRawDB-method</pre>	getMethylDiff, 26
(getAssembly), 18	<pre>getMethylDiff,methylDiff-method</pre>
getContext, 19, 37, 39–43	(getMethylDiff), 26
<pre>getContext,methylBase-method</pre>	<pre>getMethylDiff,methylDiffDB-method</pre>
(getContext), 19	(getMethylDiff), $26$
<pre>getContext,methylBaseDB-method</pre>	getSampleID, 27
(getContext), 19	getSampleID,methylBase-method
<pre>getContext,methylDiff-method</pre>	(getSampleID), 27
(getContext), 19	<pre>getSampleID,methylBaseDB-method</pre>
<pre>getContext,methylDiffDB-method</pre>	(getSampleID), 27
(getContext), 19	<pre>getSampleID, methylDiff-method</pre>
<pre>getContext,methylRaw-method</pre>	(getSampleID), 27
(getContext), 19	<pre>getSampleID,methylDiffDB-method</pre>
<pre>getContext,methylRawDB-method</pre>	(getSampleID), 27
(getContext), 19	getSampleID, methylRaw-method
getCorrelation, 20	(getSampleID), 27
<pre>getCorrelation,methylBase-method</pre>	getSampleID, methylRawDB-method
(getCorrelation), $20$	(getSampleID), 27
<pre>getCorrelation,methylBaseDB-method</pre>	<pre>getSampleID,methylRawList-method</pre>
(getCorrelation), $20$	(getSampleID), 27
getCoverageStats, 22	getSampleID,methylRawListDB-method
<pre>getCoverageStats,methylRaw-method</pre>	(getSampleID), 27
(getCoverageStats), 22	<pre>getSampleID&lt;- (getSampleID), 27</pre>
<pre>getCoverageStats,methylRawDB-method</pre>	<pre>getSampleID&lt;-,methylBase-method</pre>
(getCoverageStats), 22	(getSampleID), 27
getData, 15, 23, 37, 39–43	<pre>getSampleID&lt;-,methylBaseDB-method</pre>
<pre>getData, methylBase-method (getData), 23</pre>	(getSampleID), 27
<pre>getData,methylBaseDB-method(getData),</pre>	<pre>getSampleID&lt;-,methylDiff-method</pre>
23	(getSampleID), 27
<pre>getData, methylDiff-method (getData), 23</pre>	<pre>getSampleID&lt;-,methylDiffDB-method</pre>
<pre>getData,methylDiffDB-method(getData),</pre>	(getSampleID), 27
23	getSampleID<-,methylRaw-method
<pre>getData, methylRaw-method (getData), 23</pre>	(getSampleID), 27
<pre>getData,methylRawDB-method(getData), 23</pre>	getSampleID<-,methylRawDB-method
getDBPath, 24	(getSampleID), 27
getDBPath,methylBaseDB-method	<pre>getSampleID&lt;-,methylRawList-method</pre>
(getDBPath), 24	(getSampleID), 27
<pre>getDBPath,methylDiffDB-method</pre>	<pre>getSampleID&lt;-,methylRawListDB-method</pre>
(getDBPath), 24	(getSampleID), 27
<pre>getDBPath,methylRawDB-method</pre>	getTargetAnnotationStats
(getDBPath), 24	(genomation-deprecated), 18
<pre>getDBPath,methylRawListDB-method</pre>	getTreatment, 29
(getDBPath), 24	$\verb getTreatment,methylDiffDB-method \\$
getFeatsWithTargetsStats	(getTreatment), 29
(genomation-deprecated), 18	<pre>getTreatment,methylBase-method</pre>
<pre>getFlanks(genomation-deprecated), 18</pre>	(getTreatment), 29
${\tt getMembers} \; ({\tt genomation-deprecated}), \; 18$	<pre>getTreatment,methylBaseDB-method</pre>
getMethylationStats, 25	(getTreatment), 29

getTreatment,methylDiff-method	methSeg, 30, 35, 36, 37
(getTreatment), 29	methSeg2bed, <i>36</i> , <i>36</i>
getTreatment,methylDiffDB-method	methylBase, 5, 15, 18-20, 23, 27, 29, 31, 39,
(getTreatment), 29	<i>54–58, 61–63, 65, 68</i>
getTreatment,methylRawList-method	methylBase (methylBase-class), 37
(getTreatment), 29	methylBase-class, 37
getTreatment,methylRawListDB-method	methylBase.obj, 38
(getTreatment), 29	methylBaseDB, 5, 15, 18-20, 23, 24, 27-29,
<pre>getTreatment&lt;- (getTreatment), 29</pre>	31, 54–56, 58, 61, 62, 65
getTreatment<-,methylBase-method	methylBaseDB (methylBaseDB-class), 38
(getTreatment), 29	methylBaseDB-class,38
getTreatment<-,methylBaseDB-method	methylDiff, 5, 6, 14, 15, 18-20, 23, 26, 27,
(getTreatment), 29	29, 31, 35, 41, 61–63, 68
getTreatment<-,methylDiff-method	methylDiff(methylDiff-class), 39
(getTreatment), 29	methylDiff-class, 39
getTreatment<-,methylDiffDB-method	methylDiff.obj, 40
(getTreatment), 29	methylDiffDB, 5, 6, 24, 26–29, 31, 35, 62
getTreatment<-,methylRawList-method	methylDiffDB (methylDiffDB-class), 40
(getTreatment), 29	methylDiffDB-class, 40
<pre>getTreatment&lt;-,methylRawListDB-method</pre>	methylRaw, 5, 6, 15, 18–20, 23, 27, 31, 35, 43,
(getTreatment), 29	44, 55, 56, 61–63, 65, 68
GRanges, 30, 35–37, 39–43, 55	methylRaw (methylRaw-class), 42
GRangesList, 55	methylRaw-class, 42
	methylRawDB, 5, 6, 15, 18–20, 23, 24, 27, 31,
hclust, 11	35, 45, 55, 56, 61, 62, 65
hist, 22, 25	methylRawDB (methylRawDB-class), 43
	methylRawDB-class, 43
joinSegmentNeighbours, 30, 36	methylRawList, 5, 6, 27, 29, 31, 55, 56, 65, 68
	methylRawList (methylRawList-class), 44
list, <i>32</i> , <i>44</i> , <i>45</i>	methylRawList-class, 44
	methylRawList.obj, 45
makeMethylDB, 31	methylRawListDB, 5, 6, 24, 27–29, 31, 55, 56,
makeMethylDB,methylBase-method	65
(makeMethylDB), 31	methylRawListDB
makeMethylDB,methylBase-methods	(methylRawListDB-class), 45
(makeMethylDB), 31	
makeMethylDB,methylDiff-method	methylRawListDB-class, 45
(makeMethylDB), 31	1: 0 46
makeMethylDB,methylDiff-methods	normalizeCoverage, 46
(makeMethylDB), 31	normalizeCoverage, methylRawList-method
makeMethylDB,methylRaw-method	(normalizeCoverage), 46
(makeMethylDB), 31	normalizeCoverage, methylRawListDB-method
makeMethylDB,methylRaw-methods	(normalizeCoverage), 46
(makeMethylDB), 31	
makeMethylDB,methylRawList-method	p.adjust, <i>8</i> , <i>10</i>
(makeMethylDB), 31	PCASamples, 47
makeMethylDB,methylRawList-methods	PCASamples, methylBase-method
(makeMethylDB), 31	(PCASamples), 47
methRead, 32	PCASamples, methylBaseDB-method
methRead,character,character,character-me	
(methRead), 32	percMethylation, 49
methRead,list,list,character-method	percMethylation,methylBase-method
(methRead), 32	(percMethylation), 49

percMethylation, methylBaseDB-method	regionCounts, methylRawList, GRangesList-method
(percMethylation), 49	(regionCounts), 55
plotTargetAnnotation	regionCounts, methylRawListDB, GRanges-method
(genomation-deprecated), 18	(regionCounts), 55
pool, 9, 50	regionCounts, methylRawListDB, GRangesList-method
pool, methylBase-method (pool), 50	(regionCounts), 55
pool, methylBaseDB-method (pool), 50	removeComp, 54, 57
prcomp, 47	removeComp,methylBase-method
processBismarkAln, 51	(removeComp), 57
<pre>processBismarkAln,character,character,charac</pre>	
(processBismarkAln), 51	(removeComp), 57
$\verb processBismarkAln, list, list, character-method \\$	reorganize, $9,59$
(processBismarkAln), 51	reorganize,methylBase-method
	(reorganize), 59
qvalue, 8, 10	reorganize,methylBaseDB-method
1 42 45	(reorganize), 59
read, 43–45	reorganize, methylRawList-method
read (methRead), 32	(reorganize), 59
read.bed(genomation-deprecated), 18	reorganize, methylRawListDB-method
read.bismark(processBismarkAln), 51	(reorganize), 59
read.feature.flank	, , ,
(genomation-deprecated), 18	select, 60
read.table, 32	select, wo select, methylBase-method (select), 60
read.transcript.features	select, methylBaseDB-method (select), 60
(genomation-deprecated), 18	
reconstruct, 53	select, methylDiff-method (select), 60
reconstruct, ANY, methylBase-method	select, methylDiffDB-method (select), 60
(reconstruct), 53	select, methylRaw-method (select), 60
reconstruct, ANY, methylBaseDB-method	select, methylRawDB-method (select), 60
(reconstruct), 53	selectByOverlap, 62
reconstruct, methylBase-method	selectByOverlap,methylBase-method
(reconstruct), 53	(selectByOverlap), 62
reconstruct,methylBaseDB-method	selectByOverlap,methylBaseDB-method
(reconstruct), 53	(selectByOverlap), 62
regionCounts, 55	selectByOverlap,methylDiff-method
regionCounts, methylBase, GRanges-method	(selectByOverlap), 62
(regionCounts), 55	selectByOverlap,methylDiffDB-method
regionCounts, methylBase, GRangesList-method	(selectByOverlap), 62
(regionCounts), 55	selectByOverlap,methylRaw-method
regionCounts, methylBaseDB, GRanges-method	(selectByOverlap), 62
(regionCounts), 55	selectByOverlap,methylRawDB-method
regionCounts, methylBaseDB, GRangesList-method	(selectByOverlap), 62
(regionCounts), 55	selectByOverlap,methylRawList-method
regionCounts, methylRaw, GRanges-method	(selectByOverlap), 62
(regionCounts), 55	selectByOverlap,methylRawListDB-method
regionCounts, methylRaw, GRangesList-method	(selectByOverlap), 62
(regionCounts), 55	show, methylBase
regionCounts, methylRawDB, GRanges-method	(show, methylBase-method), 63
(regionCounts), 55	show, methylBase-method, 63
regionCounts, methylRawDB, GRangesList-method	show, methylBaseDB
(regionCounts), 55	(show, methylBase-method), 63
regionCounts, methylRawList, GRanges-method	show, methylBaseDB-method
(regionCounts), 55	(show, methylBase-method), 63
(1 05-0110041100), 55	(onon, me ony ibase me onou), oo

```
show, methylDiff
        (show, methylBase-method), 63
show, methylDiff-method
        (show, methylBase-method), 63
show, methylDiffDB
        (show, methylBase-method), 63
show, methylDiffDB-method
        (show, methylBase-method), 63
show, methylRaw
        (show, methylBase-method), 63
show, methylRaw-method
        (show, methylBase-method), 63
show, methylRawDB
        (show, methylBase-method), 63
show, methylRawDB-method
        (show, methylBase-method), 63
show, methylRawList
        (show, methylBase-method), 63
show, methylRawList-method
        (show, methylBase-method), 63
show, methylRawListDB
        ({\tt show,methylBase-method}),\,63
show, methylRawListDB-method
        (show, methylBase-method), 63
tileMethylCounts, 64
tileMethylCounts,methylBase-method
        (tileMethylCounts), 64
{\tt tileMethylCounts,methylBaseDB-method}
        (tileMethylCounts), 64
tileMethylCounts, methylRaw-method
        (tileMethylCounts), 64
{\tt tileMethylCounts,methylRawDB-method}
        (tileMethylCounts), 64
tileMethylCounts, methylRawList-method
        (tileMethylCounts), 64
tile {\tt MethylCounts,methylRawListDB-method}
        (tileMethylCounts), 64
unite, 37, 38, 66
unite, methylRawList-method (unite), 66
unite, methylRawListDB-method (unite), 66
updateMethObject, 68
```