ArrayExpress

October 25, 2011

ArrayExpress

R objects from ArrayExpress database

Description

ArrayExpress produces an AffyBatch, an ExpressionSet or a NChannelSet from a raw dataset from the ArrayExpress database. ArrayExpress needs an Internet connection.

Usage

```
ArrayExpress(input, path = getwd(), save = FALSE, rawcol = NULL)
```

Arguments

input	an ArrayExpress identifier. The function currently covers 96% of the raw data sets of the database.
path	the name of the directory in which the files downloaded on the ArrayExpress repository will be extracted. The default is the current directory.
save	if TRUE, the files downloaded from the database will not be deleted from path after executing the function.
rawcol	by default, for the raw data, the columns are automatically selected according to the scanner type. If the scanner is unknown or if the user wants to use different columns than the default, the argument 'rawcol' can be set. For two colour arrays it must be a list with the fields 'R', 'G', 'Rb' and 'Gb' giving the column names to be used for red and green foreground and background. For one colour arrays, it must be a character string with the column name to be used. These column names must correspond to existing column names of the expression files.

Value

The output is an object of class AffyBatch or ExpressionSet or NChannelSet with the raw expression values in the assayData of the object, the information contained in the .sdrf file in the phenoData, the adf file in the featureData and the idf file content in the experimentData.

If several array designs are used in the data set, the output is a list with an object for each array design.

2 extract.zip

Author(s)

Audrey Kauffmann Maintainer: <kauffmann@bergonie.org>

See Also

```
queryAE, getAE, magetab2bioc, getcolproc, procset
```

Examples

```
ETABM25.affybatch = ArrayExpress(input = "E-TABM-25")
print(ETABM25.affybatch)
sampleNames(ETABM25.affybatch)
colnames(pData(ETABM25.affybatch))
```

extract.zip

Unzip archives in a specified directory

Description

extract.zip extracts the files from a .zip archive in a specific directory.

Usage

```
extract.zip(file, extractpath = dirname(file)[1])
```

Arguments

file A file name.

extractpath A path to define where the files are to be extracted.

Value

```
Success is indicated by returning the directory in which the files have been extracted. If it fails, it returns an empty character string.
```

Author(s)

Audrey Kauffmann Maintainer: <kauffmann@bergonie.org>

getAE 3

getAE

Download MAGE-TAB files from ArrayExpress in a specified directory

Description

get AE downloads and extracts the MAGE-TAB files from an ArrayExpress dataset.

Usage

```
getAE(input, path = getwd(), type = "full", extract = TRUE)
```

Arguments

input is an ArrayExpress identifier. path is the name of the directory in which the files downloaded on the ArrayExpress repository will be extracted. can be 'raw' to download and extract only the raw data, 'processed' to download type and extract only the processed data or 'full' to have both raw and processed data. if FALSE, the files are not extracted from the zip archive.

Value

extract

A list with the names of the files that have been downloaded and extracted.

Author(s)

Audrey Kauffmann Maintainer: <kauffmann@bergonie.org>

See Also

ArrayExpress, magetab2bioc, getcolproc, procset

Examples

```
mexp1422 = getAE("E-MEXP-1422", type = "full")
## Build a an ExpressionSet from the raw data
MEXP1422raw = magetab2bioc(files = mexp1422)
## Build a an ExpressionSet from the processed data
cnames = getcolproc(mexp1422)
MEXP1422proc = procset(mexp1422, cnames[2])
```

4 getcolraw

getcolproc

Return the possible column names from processed MAGE-TAB files

Description

 ${\tt getcolproc\,extracts\,the\,column\,names\,from\,processed\,MAGE-TAB\,and\,return\,them.\,\,The\,output\,is\,needed\,to\,call\,the\,function\,\texttt{procset}.}$

Usage

```
getcolproc(files)
```

Arguments

files

A list as given from getAE function. Containing the following elements: **procfile** procfile is the name of the processed MAGE-TAB file to be read. **path** is the name of the directory where to find this file.

Author(s)

Audrey Kauffmann

Maintainer: <kauffmann@bergonie.org>

See Also

ArrayExpress, queryAE, getAE, procset

getcolraw

Return the possible column names from raw MAGE-TAB files

Description

getcolraw extracts the column names from raw MAGE-TAB and return them. The output can be use to set the argument 'rawcol' of the function magetab2bioc.

Usage

```
getcolraw(path, rawfiles)
```

Arguments

rawfiles rawfiles are the name of the raw MAGE-TAB files to be read.

path is the name of the directory where to find these files.

Author(s)

Audrey Kauffmann Maintainer: <kauffmann@bergonie.org>

See Also

```
ArrayExpress, queryAE, getAE
```

magetab2bioc 5

magetab2bioc

Convert MAGE-TAB files from raw data into a Bioconductor object

Description

magetab2bioc converts local MAGE-TAB files into a AffyBatch, an ExpressionSet or a NChannelSet.

Usage

```
magetab2bioc(files, rawcol = NULL, save = TRUE)
```

Arguments

files A list as given from getAE function. Containing the following elements:

rawfiles all the expression files to use to create the object. The content of the

raw.zip MAGE-TAB file.

sdrf the name of the sdrf file from MAGE-TAB.

idf the name of the idf file from MAGE-TAB.

adf the name of the adf file from MAGE-TAB.

path is the name of the directory containing these files.

rawcol by default, the columns are automatically selected according to the scanner type.

If the scanner is unknown or if the user wants to use different columns than the default, the argument 'rawcol' can be set. For two colour arrays it must be a list with the fields 'R', 'G', 'Rb' and 'Gb' giving the column names to be used for red and green foreground and background. For one colour arrays, it must be a character string with the column name to be used. These column names must

correspond to existing column names of the expression files.

save if TRUE, the files used to create the object will not be deleted from path after

executing the function.

Value

An object of class AffyBatch, ExpressionSet or NChannelSet with the raw expression values in the 'assayData' of the object, the information contained in the sdrf file in the 'phenoData', the adf file content in the 'featureData' and the idf file content in the 'experimentData'.

If several array designs are used in the dataset, the output is a list with an object for each array design.

Author(s)

Audrey Kauffmann Maintainer: <kauffmann@bergonie.org>

See Also

```
ArrayExpress, queryAE, getAE
```

Examples

An example can be found in the help of the getAE function.

6 queryAE

procset

Convert processed MAGE-TAB files into a Bioconductor object

Description

procset converts local MAGE-TAB files into an ExpressionSet.

Usage

```
procset(files, procol)
```

Arguments

files is the list with the names of the processed, the sdrf, the adf and the idf files and

the path of the data as given by getAE.

procol the name of the column to be extracted from the file. Obtained using getcolproc.

Author(s)

Audrey Kauffmann Maintainer: <kauffmann@bergonie.org>

See Also

```
queryAE, getAE, getcolproc
```

Examples

An example can be found in the help of the getAE function.

queryAE

XML query of the ArrayExpress repository

Description

queryAE queries the ArrayExpress database with keywords and give a dataframe with ArrayExpress identifiers and related information, as an output.

Usage

```
queryAE(keywords = NULL, species = NULL)
```

Arguments

keywords the keyword(s) of interest. To use several words, they must be separated by a

"+" as shown in the examples.

species the specie(s) of interest.

queryAE 7

Value

A dataframe with all the ArrayExpress dataset identifiers which correspond to the query in the first column. The following columns contain information about these datasets, such as the number of raw files, the number of data processed, the release date on the database, the pubmed ID, the species, the experiment design and the experimental factors.

Author(s)

Audrey Kauffmann

Maintainer: <kauffmann@bergonie.org>

See Also

ArrayExpress, getAE

Examples

```
## To retrieve all the identifiers of pneumonia data sets
pneumo = queryAE(keywords = "pneumonia")

## To retrieve all the identifiers of pneumonia data sets studied in human
pneumoHS = queryAE(keywords = "pneumonia", species = "homo+sapiens")
```

Index

```
*Topic datasets
    ArrayExpress, 1
    extract.zip, 2
    getAE, 3
    getcolproc, 4
    {\tt getcolraw}, {\tt 4}
    magetab2bioc, 5
    procset, 6
   queryAE, 6
AffyBatch, 1,5
ArrayExpress, 1, 3-5, 7
ExpressionSet, 1, 5, 6
extract.zip, 2
getAE, 2, 3, 4-7
getcolproc, 2, 3, 4, 6
getcolraw,4
magetab2bioc, 2, 3, 5
NChannelSet, 1, 5
procset, 2-4, 6
queryAE, 2, 4, 5, 6, 6
```